

Advasil Conform[®]

Silicone gel sheet for scar management

Product description


A sheet of Medical Grade self-adherent silicone gel for the management of scar tissue or recently healed wounds. Composed of a conformable layer of translucent Silfix[®] silicone gel and a smooth polyurethane outer layer to improve wear ability and patient compliance.

Silicone gel sheets are reported to help reduce the effects of scarring if worn for prolonged periods. Hypertrophic and keloid scarring may particularly benefit from Advasil Conform[®]. Advasil Conform[®] is made from medical grade silicone approved for skin contact. Advasil Conform[®] is supplied sterile for improved hygiene but once opened is non-sterile.

Use

Remove Advasil Conform[®] from the pack, remove the liner and place gel side to the scar. For smaller scars Advasil Conform[®] can be cut to size. For larger scars use additional sheets of Advasil Conform[®] and place side by side. Advasil Conform[®] is self adherent and reusable for up to 28 days.

When storing Advasil Conform[®] replace the liner and store in original packaging. It is important that the skin beneath the dressing is kept clean, and dry. Wash area daily. Wear time should be built up over a period of time to condition the skin after which Advasil Conform[®] can be comfortably worn for 24 hours.



Ordering information for Advasil Conform[®]

Size	Stock code	Pack size	NHS code	PIP code
10cm x 10cm	CR3845	5	EKH006	301-1640
10cm x 15cm	CR3832	5	ELY178	329-5979



Features

- ✓ Can be shaped to optimise coverage of scar area
- ✓ Conformable
- ✓ Washable
- ✓ Reusable for up to 28 days
- ✓ Self-securing with gentle release

Indications

- ✓ Hypertrophic & Keloid scars

Care and use symbols

	Caution - pay specific attention to any warnings on the IFU		CE marked medical device
	Latex free		Keep dry
	Single use only		Store out of direct sunlight
	Do not use if packaging is broken or damaged		Do not re-sterilize
	STERILE EO Sterilised using Ethylene Oxide		Read instructions for use

Contra-indications

Do not use if allergic to silicone.

Storage

Store at room temperature, out of direct sunlight and keep dry.

Advancis
Medical
www.advancis.co.uk

Another quality product selected and stocked by:

WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489
e enquiries@bamford.co.nz www.bamford.co.nz

When performance counts **BAMFORD** delivers

