

SURE-SCOPE™

Single-Use Laryngoscope with Built-In Light Source

*BRIGHTER
IS
BETTER*



OBP Medical's SURE-SCOPE is a single-use one piece laryngoscope with built-in, hands-free, LED light source. The ready-to-use device eliminates the time and expense of reprocessing and reduces the risk of cross-contamination in the medical setting by obsoleting reusable light sources, handles, and blades.

The SURE-SCOPE is designed to facilitate tracheal intubation during general anesthesia or cardiopulmonary resuscitation or for procedures on the larynx or other parts of the upper tracheobronchial tree.

OBP
MEDICAL

Simple. Bright. Ideas.

Another quality product selected and stocked by:
WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489
e enquiries@bamford.co.nz www.bamford.co.nz

when performance counts **BAMFORD** delivers



Stronger, brighter, safer than the rest.

Think your current laryngoscope is clean?

Think again. Recent studies demonstrate a clear potential for transmission of pathogens, and highlight a clear risk to both staff and patients.

- There are 4.5 hospital infections for every 100 patient admissions and nearly 100,000 deaths annually. (CDC)
- There are approximately 1.7 million hospital-acquired infections in US hospitals, costing over \$5 billion annually. (CDC)

How strong are SURE-SCOPE blades?

A study published in the *The Association of Anaesthetists of Gt Britain and Ireland* concluded that approximately 10lbs of force is needed to intubate the average patient. The SURE-SCOPE blade has been validated to support over 40lbs of blade force.

TO ORDER:

e orders@bamford.co.nz
p 0800-2263673

- Item # Mac 3: C090130
Mac 4: C090140
Miller 0: C090200
Miller 1: C090210
Miller 2: C090220
Miller 3: C090230
- Batteries are non-hazardous
- Batteries included in each unit
- Each unit comes packed sterile

Product Features:

- Single-use device is completely hygienic, eliminating the risk of cross-contamination from poorly cleaned reusable components
- Each SURE-SCOPE comes ready for use with a built-in single-use LED light source
- No assembly is required and no additional parts (replacement blades, light sources, chargers, batteries, sheaths, etc.) are needed
- High powered, battery-operated LED light source is ultra-bright, and can produce peak light for 30+ minutes
- LED light source generates absolutely no heat
- Comfortable ergonomic design
- One piece construction ensures blade integrity and eliminated risk of corrosion.
- Blunt tip helps reduce trauma.
- Ideal for OR, ICU, EMS, ER, Urgent Care, and Crash Carts.

About Us:

OBP Medical is the leading global developer of single-use, self-contained, illuminating medical devices. Our single-use lights are the brightest on the market, allowing for optimal visualization during procedures. Our innovative products are used in physician offices, surgery centers and more than 2,000 hospitals throughout the U.S., as well as healthcare facilities worldwide..

MedAssets Contract # MS02268
ISO13485 certified

OBP Medical Inc.
360 Merrimack St, Building 9
Lawrence, MA 01843
www.obpmedical.com

OBP
MEDICAL
Simple. Bright. Ideas.