


Actilite Protect[®]

The first time 3 dressings have been combined



Wound healing...
as nature intended

Another quality product selected and stocked by:
WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489
e enquiries@bamford.co.nz www.bamford.co.nz


when performance counts **BAMFORD** delivers

Actilite Protect®

The first time 3 dressings have been combined



At our state-of-the-art manufacturing facility in the UK, Advancis is constantly investing in the research and development of our range of advanced wound care products. This allows not only innovation, but enables us to deliver the unique combination of high quality products at competitive prices.



Actilite Protect® offers the benefits of 3 effective, well known, trusted products - combined into one new unique dressing:

Foam – a patient friendly absorbent dressing for light to medium exudate

Manuka Honey – nature's effective antimicrobial & wound healing creation

Silicone – an atraumatic wound contact layer for pain free removal

Benefits of using our 3-in-1 dressing

- ✓ Multifunctional – the dressing offers the key benefits of 3 dressings combined into 1
- ✓ The Activon Manuka Honey® is presented in a dry matrix, which gels when in contact with exudate. The resulting gel then interacts with the wound delivering the antimicrobial properties of the honey
- ✓ 10mm holes that allow the 2-way transfer of the benefits of Manuka Honey on to the wound and exudate back up into the foam layer
- ✓ 3.5mm foam layer providing increased absorbency and reducing frequency of dressing change
- ✓ Wide range of sizes – 10x10mm, 10x15mm, 10x25mm and 10x35mm
- ✓ Wider silicone adhesive borders providing secure dressing fixation, whilst being gentle on surrounding skin
- ✓ Clear backing allows easy assessment of wound exudate without requiring removal of the dressing
- ✓ Competitively priced in relation to its nearest competitors



The ideal choice to prevent infection

The dressing is an ideal choice to prevent infection but also to treat an existing wound infection. The dressing contains Activon Manuka Honey® and is an effective antimicrobial treatment. The versatility of this dressing means it is appropriate for every phase of wound healing and could be used for:

- ✓ Surgical wounds for preventing infection post operatively
- ✓ Caesarean section
- ✓ Leg ulcers in the final stages of healing
- ✓ Trauma wounds
- ✓ Superficial pressure ulcers
- ✓ Skin tears

Manuka honey contains the benefits of a broad spectrum antimicrobial, effective against *Staphylococcus aureus*, *Pseudomonas aeruginosa*, VRE and possibly the best known, MRSA - challenging the need for routine antibiotics to treat wound infection.

Manuka honey is a natural alternative to the more common wound therapies such as silver, with none of the perceived challenges that silver products face today, such as their cost, lack of efficacy and concerns around their safety and resistance¹.

¹Is Manuka Honey a Credible Alternative to Silver in Wound Care? by Cooper Rose; Gray David - Wounds UK 2012

Advancis Medical is a trading name of Brightwake Ltd