



Effective disinfection & lasting protection

- ✓ Advanced Active Technology™
- ✓ Kind to surfaces and your skin
- ✓ Broad spectrum antimicrobial efficacy



## Disinfectant wipes & sprays

### AZOACTIVE™

AZOACTIVE™ is based on a highly effective formulation of quaternary ammonium compounds (QACs) and polymeric biguanide (PHMB), working synergistically with other ingredients to deliver a potent and fast acting disinfectant which is kind and gentle to surfaces.

AZOACTIVE™ is designed for use in hospitals, nursing and care homes, laboratories, dental and GP surgeries. It has become increasingly popular in non medical areas like offices, beauty salons, nurseries and schools. It can be used on a wide range of surfaces from mattresses, bed frames, table tops, chairs, commodes, ambulance stretchers, cabinets, tile surfaces and many more. The range also includes wipes and sprays for disinfecting non invasive medical devices. AZOACTIVE™ offers additional benefits to alcohol based disinfectants like residual activity and broader antimicrobial efficacy including enveloped and non enveloped viruses e.g. Norovirus and Influenza H1N1.

### Instructions for use



Use AZOACTIVE wipe for disinfecting. Wipe over surface/equipment to be disinfected. Use more than one wipe if necessary so the surface stays visibly wet, then leave to dry.



Once used, discard the wipe into appropriate waste container.



Use AZOACTIVE spray to disinfect – spray onto the cleaned surface and wipe or leave to dry. To deodorise – spray onto contaminated surface and leave to dry.



Ensure the surface or a piece of equipment is clean prior to disinfection. Do not use one wipe on more than one piece of equipment. Do not mix with any other chemicals. Keep the canister/spray bottle/bucket closed. Store in a cool, dry place. Avoid direct contact with sunlight.

## Choose AZOACTIVE™

For effective disinfection and lasting protection

- **Long lasting, residual efficacy**  
keeps working for hours, even after it has dried out
- **Broad spectrum antimicrobial efficacy**  
effective against a wide range of micro-organisms
- **Chlorine, glutaraldehyde and alcohol free**  
non tainting & non corrosive
- **Kind to surfaces**  
reduces ongoing repair & replacement costs of equipment
- **Skin friendly**  
kind to your skin
- **Safe to environment**  
non hazardous & non toxic  
no specialist PPE or storage required
- **Pleasant fragrance**  
removes unpleasant odours
- **Multi-surface & medical devices wipes**  
wide range of applications
- **Available in a wide range of formats**  
pouches, canisters, buckets and sprays



Code	Product name	Product description	Size (mm)	Wipes qty	Case qty
AA050W	<b>AZOACTIVE™</b>	Multi - surface disinfectant wipes for universal use - pouch	330x220	50	20
AA100W	<b>AZOACTIVE™</b>	Multi - surface disinfectant wipes for universal use - canister	200x220	100	12
AA100NW	<b>AZOACTIVE™</b> <b>NEW</b>	Multi - surface disinfectant wipes for universal use - canister	130x180	100	12
AA200W	<b>AZOACTIVE™</b> <b>NEW</b>	Multi - surface disinfectant wipes for universal use - canister	180x200	200	12
AA300W	<b>AZOACTIVE™</b> <b>NEW</b>	Multi - surface disinfectant wipes for universal use - bucket	260x235	300	4
AA500S	<b>AZOACTIVE™</b> <b>SPRAY</b>	Multi - surface disinfectant spray for universal use - bottle	500ml	N/A	12
AA050WCE	<b>AZOACTIVE™</b> <b>NEW</b>	Disinfectant wipes for non-invasive medical devices, CE mark - pouch	300x220	50	30
AA100WCE	<b>AZOACTIVE™</b> <b>NEW</b>	Disinfectant wipes for non-invasive medical devices, CE mark - canister	130x180	100	12
AA200WCE	<b>AZOACTIVE™</b> <b>NEW</b>	Disinfectant wipes for non-invasive medical devices, CE mark - canister	180x200	200	12
AA500SCE	<b>AZOACTIVE™</b> <b>NEW</b> <b>SPRAY</b>	Disinfectant wipes for non-invasive medical devices, CE mark - bottle	500ml	N/A	12

Another quality product selected and stocked by:

WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489  
e enquiries@bamford.co.nz www.bamford.co.nz

When performance counts

