

# STERIL-PEEL SELF SEAL POUCHES

**Self-Seal pouches** provide a self contained method of sterile packaging for instruments without the use of heat sealing equipment. These steam sterilization pouches are ideal for table top sterilizers that are used in operating rooms, surgery centers, as well as medical and dental offices. Self seal pouches are also used as ethylene oxide pouches in hospitals.

**Pre-Folded** at the designated seal location, assures a quick, accurate seal without perforation holes.

**Triple seal** construction offers additional assurance of seal integrity.

## Features

- Tack points above chevron at corners reduce paper curl and protect top of pouch from airborne contaminants.
- Thumb notched for ease of loading and opening under sterile conditions.
- Unique blue film darkens to show seal strength, therefore unacceptable seals are clearly visible.
- 52 week shelf life study shows 100% sterile maintenance.
- All pouches have a class 1 process indicator.



<u>Code</u>	<u>Description</u>	<u>Packaging</u>
222100/705	Self Seal Pouch 89mm x 133mm	200/Box
222101/710	Self Seal Pouch 89mm x 228mm	200/Box
222102/715	Self Seal Pouch 133mm x 254mm	200/Box
222103/720	Self-Seal Pouch 191mm x 330mm	200/Box
222104/725	Self Seal Pouch 133mm x 381mm	200/Box
222218/730	Self Seal Pouch 304mm x 460mm	200/Box
222565/760	Self Seal Pouch 90mm x 520mm	200/Box



T 1800 074 085 Customer Support