

# SURGI STRIP



## SKIN CLOSURES

*Superior strength, exceptional porosity and the finest grade medical adhesive have been combined to produce a skin closure tape designed for maximum wound support and protection.*

High strength filament reinforcing provides greater wound support yet the extreme thinness of the backing ensures excellent conformability.

Rounded ends resist lifting and curling. SURGI STRIP remains clean and neat for longer.

SURGI STRIP skin closures allow earlier post-operative suture or staple removal to assist natural healing and minimize scarring.

A painless alternative for the closure of minor incisions, and for skin lacerations, enabling fast and easy application, with no need for anaesthetic, eliminating needle puncture marks.



*Colour coded packaging for quick size selection:*

C-01 Aqua 3mm x 75mm 50 pouches of 5  
C-02 Green 6mm x 75mm 50 pouches of 3  
C-03 Orange 6mm x 38mm 50 pouches of 6

C-04 Blue 6mm x 100mm 50 pouches of 10  
C-05 Red 13mm x 100mm 50 pouches of 6  
C-06 Purple 25mm x 125mm pouches of 4

**Also available in flesh colour non reinforced**

S-11 3mm x 75mm 50 pouches of 5  
S-13 6mm x 38mm 50 pouches of 6

S-12 6mm x 75mm 50 pouches of 3  
S-15 13mm x 100mm 50 pouches of 6

**Another quality product selected and stocked by:**

WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489  
e enquiries@bamford.co.nz [www.bamford.co.nz](http://www.bamford.co.nz)

**When performance counts BAMFORD delivers**

