

LIBERTY 7E-A PORTABLE SUCTION PUMP

Liberty[®]

The Liberty 7E-A Portable Suction Unit, is a new design, oil-free lubrication pump and based on extensive experience and development trends from domestic and overseas products.

The Liberty 7E-A Suction Pump is specially designed for absorbing thick liquid such as blood, phlegm, etc. Not only can this unit be used as an aspirator in the operating room or emergency room, but also as a health care product for the home.

Features

- Portable, compact small size
- Light-weight, easy to carry
- High negative pressure
- Large Flux
- Oil-free lubrication pump
- Square negative pressure dial
- Built-in bottle cavity
- Low noise



Technical Information

Voltage:	AC 220/240V 50HZ
Power:	≤ 90 VA
Max Negative Pressure Value:	≥ 0.075 MPa
Pumping Rate:	≥ 18 Litres/Minute
Reservoir Volume:	1000 mL, 1 Piece
Noise:	≤ 65 Db (A)

CODE	DESCRIPTION
SPB	LIBERTY 7E-A SUCTION PUMP 240V 18L

SPARES/ACCESSORIES

SPBF	FILTER
SPBC12	SUCTION CANNULA 12FG
SPBSB1L	SPARE BOTTLE 1L
SPBSB2L	SPARE BOTTLE 2L

E&OE 15OCT08



Square Negative Pressure Dial



Another quality product selected and stocked by:

WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489
e enquiries@bamford.co.nz www.bamford.co.nz

When performance counts

