

LIBERTY 7E-D MULTI-VOLTAGE PORTABLE SUCTION PUMP

Liberty[®]



Features

- Multi-Voltage
- Fully Portable
- Oil-free lubrication pump
- Square negative pressure dial
- Built-in bottle cavity
- Low noise
- Small size (H22CM x L29CM x D20CM)
- Easy to carry
- Fully portable unit incorporating three modes of electricity: 240V AC, 12V DC and Battery
- Suitable for vehicles, clinic, home and outdoor places

Technical Information

Voltage:	AC 220V/240V 50HZ 12VDC
Power:	90 VA
Max Negative Pressure Value:	≥ 0.075 MPa
Pumping Rate:	≥ 15 Litres/Minute
Reservoir Volume:	1,000 mL
Noise:	≤ 65 Db (A)

CODE	DESCRIPTION
SPBB	LIBERTY 7E-D SUCTION PUMP 240V/12V/BATTERY 15L

SPARES/ACCESSORIES

SPBF	FILTER
SPBC12	SUCTION CANNULA 12FG
SPBSB1L	SPARE BOTTLE 1L
SPBSB2L	SPARE BOTTLE 2L

E&OE 14OCT08



Another quality product selected and stocked by:

WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489
e enquiries@bamford.co.nz www.bamford.co.nz

When performance counts

