

PORTABLE SUCTION PUMP ZEINER RESCUE, MULTI-VOLTAGE 25 LITRE

Battery Operated Aspirator for First Aid Emergency



Features

- Ideal for ambulance and first-aid environments
- Internal Ni-MH battery pack, 45 minutes run time
- Car or ambulance cable and universal voltage (100 – 240VAC) adapter and internal rechargeable battery pack, 45 minutes run time
- Reliable, maintenance free and silent 12VDC double piston MEDEL pump
- Shock-resistant, easy-to-carry and easy-to-clean chassis
- Suction adjustment and vacuum gauge
- 800 ml canister with overflow valve
- Antibacterial filter

Technical Specifcation

Size:	240 x 125 x 216 mm (525 x 140 x 216 mm with double canisters brackets)
Weight:	2,8 Kg
Current:	3,5 A
Max vacuum:	-82 kPa; -610 mmHg
Max flow:	25 l/min
Battery	3300 mAh (Ni-MH battery cells)
Warranty:	2 years

CODE	DESCRIPTION
SPRESCUE	SUCTION PUMP ZEINER RESCUE 240V/12V/BATTERY 25LITRE

E&OE First Published: 22 March 2010
Last Modified: 16 April 2010
File Name: info_SPRESCUE_SUCTION-PUMP-RESCUE



Another quality product selected and stocked by:

WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489
e enquiries@bamford.co.nz www.bamford.co.nz

When performance counts

