



- Single use - reduces the risk of cross infection.
- High gloss, low-friction outer surface improves tube glide over bougie.
- Low point-of-pressure tip reduces risk of accidental trauma.
- Available in 600mm and 800mm length, each with or without angled (coudé) tip, allowing greater choice for clinicians.
- All available in 10ch and 15ch, suitable for use with all tubes from 4mm(ID) to 10mm(ID).
- Highly visible colour for improved positioning in poor light.
- Extra tip graduations allow for improved laryngeal inlet entry and exit guidance.
- Individually packed sterile, 10 units within highly visible outer tube for ease of product identification in emergency.
- Robust outer tube prevents damage in transit and storage.

- Barium tip for improved X-Ray guidance if required.



Code	Description	Size
S10600	Bougie, straight	10ch x 600mm
S15600	Bougie, straight	15ch x 600mm
S10800	Bougie, straight	10ch x 800mm
S15800	Bougie, straight	15ch x 800mm
A10600	Bougie, angled (coudé)	10ch x 600mm
A15600	Bougie, angled (coudé)	15ch x 600mm
A10800	Bougie, angled (coudé)	10ch x 800mm
A15800	Bougie, angled (coudé)	15ch x 800mm

P3 Medical Limited

Another quality product selected and stocked by:

WM Bamford & Co Ltd  
www.bamford.co.nz

p 04-576-2100  
f 04-569-6489  
e enquiries@bamford.co.nz

