

The Ultimate T-Bag™

Oxygen Enhancement Device



Recovery Room Management Made Easier

- Reliable FiO_2 concentrations of 70% at 6L/min and 50% at 3L/min.
- Holds approximately 300ml O_2 in its reservoir
- Lightweight - will not displace the airway device
- No significant rebreathing
- Confirmation of respiration through T-Bag movement
- Single use. Materials recyclable
- Available in boxes of 100 units



Single Use Only

Lightweight

**Confirms
Respiration**

Latex Free

Manufactured by

Ultimate Medical

Ultimate Medical Pty Ltd

484 Swan Street, Richmond, Australia 3121

Phone 61 3 9421 0500 Fax 61 3 9421 0555

**International Patent Pending.
Another innovation by Ultimate**