

# Mannual Resuscitator



**FORTUNE MEDICAL INSTRUMENT CORP.**

6Fl., No. 29, Sec 2, Jungjeng E. Rd. Danshuei Jen, Taipei, TAIWAN

Email: [info@fortunemed.com](mailto:info@fortunemed.com) [sales1@fortunemed.com](mailto:sales1@fortunemed.com)

TEL: (886)-2-2624 2233 FAX: (886)-2-2624-2266

[www.FortuneMed.com](http://www.FortuneMed.com)



ISO 9001 & CE Certified

## Fortune<sup>®</sup> Manual Resuscitator Features:

- Equipped with pressure relief valve at 40 cmH<sub>2</sub>O with overridden mechanism.
- 3 different materials of Silicone, Rubber and PVC for reuse or disposable selections.
- 3 sizes are available: infant, child and adult.
- ISO standard connectors.

Model	Infant	Child	Adult	Description
<b>Silicone Resuscitator</b>	SR-Infant	SR-Child	SR-Adult	Reusable
<b>Rubber Resuscitator</b>	RR-Infant	RR-Child	RR-Adult	Disposable
<b>PVC Resuscitator</b>	PR-Infant	PR-Child	PR-Adult	Disposable


### Silicone Resuscitator:

Product	REF. No.	Size	Description
<b>SR-Infant</b>	1610-0001	Infant	Silicone mask #0/ #1/ #2 Patient valve w/ pressure relief valve Silicone bag 240 ml Reservoir valve Oxygen reservoir 600 ml
<b>SR-Child</b>	1610-0002	Child	Silicone mask #3 Patient valve w/ pressure relief valve Silicone bag 500 ml Reservoir valve Oxygen reservoir 2600 ml
<b>SR-Adult</b>	1610-0003	Adult	Silicone mask #5 Patient valve w/ pressure relief valve Silicone bag 1600 ml Reservoir valve Oxygen reservoir 2600 ml



### Rubber Resuscitator:

Product	REF. No.	Size	Description
<b>RR-Infant</b>	1650-0001	Infant	Rubber mask #0/ #1/ #2 Patient valve w/ pressure relief valve Rubber bag 240 ml
<b>RR-Child</b>	1650-0002	Child	Rubber mask #3 Patient valve w/ pressure relief valve Rubber bag 500 ml
<b>RR-Adult</b>	1650-0003	Adult	Rubber mask #5 Patient valve w/ pressure relief valve Rubber bag 1600 ml



### PVC Resuscitator:

Product	REF. No.	Size	Description
<b>PR-Infant</b>	1690-0001	Infant	PVC Air-Cush mask #0/ #1/ #2 Patient valve w/ pressure relief valve PVC bag 240 ml Reservoir valve Oxygen reservoir 600 ml
<b>PR-Child</b>	1690-0002	Child	PVC Air-Cush mask #3 Patient valve w/ pressure relief valve PVC bag 500 ml Reservoir valve Oxygen reservoir 2600 ml
<b>PR-Adult</b>	1690-0003	Adult	PVC Air-Cush mask #5 Patient valve w/ pressure relief valve PVC bag 1600 ml Reservoir valve Oxygen reservoir 2600 ml



### Fortune<sup>®</sup> Air-Cush Mask Features:

- 100% medical grade PVC material.
- Clear shell provides inspection visibility
- Air cushioned dome provides a soft and friendly contact to patient's face.
- Equipped with inflatable check valve to adjust air amount for best cushion fitting outo patient's face.
- 6 sizes are available.
- Color coded.




### Fortune<sup>®</sup> Sil-Cush Mask Features:

- 100% medical grade silicone material for superior biocompatibility.
- Clear shell provides inspection visibility
- Air cushioned dome provides a soft and friendly contact to patient's face.
- Equipped with inflatable check valve to adjust air amount for best cushion hardness.
- Easy disassembling and assembling.
- Spare parts are available for replacement.
- 3 sizes are available: # 3, # 4 and # 5.
- Color coded.

### Mask:

Product	REF. No.	Size	Application	Description
Silicone Mask	1611-9000	# 0	Infant	Silicone Air cushioned Check valve
	1611-9001	# 1		
	1611-9002	# 2	Child	
	1611-9003	# 3		
	1611-9004	# 4	Adult	
1611-9005	# 5			
Sil-Cush Mask	1611-9093	# 3	Child	Silicone Air cushioned Check valve
	1611-9094	# 4	Adult	
	1611-9095	# 5		
Rubber Mask	1651-9000	# 0	Infant	Rubber
	1651-9001	# 1		
	1651-9002	# 2	Child	
	1651-9003	# 3		
	1651-9004	# 4	Adult	
1651-9005	# 5			
Air-Cush Mask	1691-9000	# 0	Infant	PVC Air cushioned Check valve
	1691-9001	# 1		
	1691-9002	# 2	Child	
	1691-9003	# 3		
	1691-9004	# 4	Adult	
1691-9005	# 5			



## Fortune<sup>®</sup> Airway Features:

- Proper rigidity to prevent deformation upon patient bite.
- Elastic plastic material to protect patient teeth.
- 2 models are available: Berman and Guedel.

### Airway:

Product	REF. No.	Size	Description
Berman Airway	1603-9040	40 mm	Infant
	1603-9050	50 mm	
	1603-9060	60 mm	
	1603-9070	70 mm	Child
	1603-9080	80 mm	
	1603-9090	90 mm	Adult
	1603-9100	100 mm	
Guedel Airway	1603-0040	40 mm	Infant
	1603-0050	50 mm	
	1603-0060	60 mm	
	1603-0070	70 mm	Child
	1603-0080	80 mm	
	1603-0090	90 mm	Adult
	1603-0100	100 mm	
	1603-0110	110 mm	



## Fortune<sup>®</sup> Carrying Bag Features:

- Bright orange Nylon bag with hanger for handy carrying purpose.
- Sufficient space to store pre-assembled resuscitator for emergency use.
- Clear accessories compartment provides storage and visibility for airways, teeth opener and necessities.
- Clear opening side with zipper provides content inspection visibility.
- Clear packet for placing manual or instruction.

### Accessories & Replacement Parts:

Product	REF. No.	Size	Description
Carrying Bag	1609-9001	Universal	Nylon Orange Clear side with zipper
Patient Valve	1601-9001	22 mm O.D.	Pressure relief valve @ 40 cmH <sub>2</sub> O
	1601-9002		Pressure relief valve @ 40 cmH <sub>2</sub> O w/ overridden lock
	1601-9003		w/o pressure relief valve
Reservoir Valve	1602-9001	22/25 mm	22/25 mm to oxygen reservoir 22 mm I.D. to intake valve
Oxygen Reservoir	1602-9060	600 ml	Infant
	1602-9260	2600 ml	Child & Adult
Oxygen Tubing	1604-9007	7 feet.	PVC
PEEP Valve Diverter	1605-8001	22/30 mm	Polysulfone
	1605-9001	22/30 mm	PC
PEEP Valve	1606-8010	10 cmH <sub>2</sub> O	Polysulfone
	1606-8020	20 cmH <sub>2</sub> O	
	1606-9010	10 cmH <sub>2</sub> O	PC
	1606-9020	20 cmH <sub>2</sub> O	
Teeth Opener	1607-9001	Universal	PC



Carrying Bag



PEEP Valve & Diverter



Teeth Opener