

## Writing Instruments - FIT for the PURPOSE

### WATERPROOF AUTOCLAVE PEN

**USE:** Fine writing on documentation which may be subjected to steam sterilization or for any documentation requiring exceptional long life legibility.

**POINT:** Strong, stainless steel roller ball, writing width 0.3.

**BODY:** Slim round white plastic incorporating stainless steel pocket clip and safety cap.

**INK:** Water-based, solvent free pigment ink, permanent.

**COLOURS AVAILABLE:**  
Black only.



Code: TTS209P

### AQUAMARC

**USE:** For permanent marking and writing on impervious surfaces which will be subjected to water, heat, chemicals and abrasion. eg. for recording items through a disinfectant/washer on waterproof paper.

**POINT:** Fine steel shielded microporous fibre, writing width 0.5.

**BODY:** Slim round black plastic pocket-sized incorporating safety cap with pocket clip.

**INK:** Permanent, longlife, water chemical and enzyme resistant - quick drying.

**COLOURS AVAILABLE:**  
Black only.



Code: WDM/003/00

### MEDIBAC ANTI-BACTERIAL PEN

**USE:** General purpose writing where cross contamination poses a high degree of risk.

**POINT:** Stainless steel ball point with reinforced ferrule.

**BODY:** Slimline clear plastic incorporating an anti-bacterial inorganic agent, White cap incorporating pocket clip and safety vent.

**INK:** Quick drying oil based suspension.

**COLOURS AVAILABLE:**  
Black only.



Code: AB/210

### MEDIMARC

**USE:** Permanent writing and marking on medical packaging which will not compromise the microbial integrity of a sterile pack. Bleed resistant.

**POINT:** Bullet tipped.

**BODY:** Tapered black plastic body incorporating colour coded safety cap with anti-roll feature.

**INK:** Permanent, non-solvent based ink, waterproof and fade resistant.

**COLOURS AVAILABLE:**  
Black, Blue, Green and Red.



Code: MM/02T

### STERITRAX MARKER

**USE:** Permanent marker that won't disappear under extreme conditions. Ideal for writing on surfaces that are exposed to water, heat, chemicals and abrasion. Suitable for most sterilization processes.

**POINT:** Ultra Fine Tip, Line Width: 0.4.

**BODY:** Slim round grey plastic incorporating pocket clip and safety cap.

**INK:** Permanent, longlife, water chemical and enzyme resistant - quick drying.

**COLOURS AVAILABLE:**  
Black only.



Code: SM04