



ENT  
CONSUMABLES


**Overview**

.....

CoreMed is an Australian based medical distribution company, specialising in the sales and marketing of ENT products.

We provide a professional, efficient and experienced approach towards the product and service requirements of hospital operating theatres, emergency departments and medical practices.

We are proud to be Spiggle & Theis and EndoClarix Exclusive Distributor in Australia and New Zealand.



**Contents**

.....



**Page**

1. Nasal Splints	1
2. Disposable Suction Tubes	2
3. Ventilation Tubes	3
4. Nasal Packings	4
5. Ear Wicks	4
6. Sterile Covers	5
7. Sheaths for Nasopharyngoscopes	5
8. Bandages	6–7

## Internal Nasal Splints

- Sterile, ready to use
- Punched suture holes
- Soft material for easy insertion and removal
- Surface with Lotus leaf effect (Silicone splints only)


Code	Feature	Packaging	Material
2030100	standard / with airway	5 pairs / box	silicone
2030110	standard / without airway	5 pairs / box	silicone
2042025	large / 0.25 mm thickness	5 pairs / box	fluoroplastic
2042050	large / 0.50 mm thickness	5 pairs / box	fluoroplastic



## Septal Button, Silicone

- Sterile, ready to use
- For non-surgical closure of septal perforation
- Soft material for easy insertion and removal
- Silicone prevents adhesion
- No scabbing, epistaxis or other septal perforation symptoms


Code	Size	Feature	Packaging
2031100	3.2 cm	one part	1 piece
2031110	3.2 cm	two parts	1 piece



## External Aluminium Nasal Splints

- Easily malleable
- Self-adhesive back, self-adhesive Dorsum foam pad included
- 10 pieces per box


Code	Dimensions	Size
2060200	56 x 38 x 30 mm	small
2060300	64 x 40 x 36 mm	medium
2060400	66 x 42 x 45 mm	large



## External Thermoplastic Nasal Splints

- Much less time and work required than with plaster of paris
- Includes three-part foam bridge
- Self-adhesive back
- 10 pieces per box

Code	A x B x C	Thickness	Size
2050100	56 x 38 x 30 mm	2.4 mm	small
2050200	64 x 40 x 36 mm	2.4 mm	medium
2050300	66 x 42 x 45 mm	2.4 mm	large
2051100	56 x 38 x 30 mm	1.6 mm	small
2051200	64 x 40 x 36 mm	1.6 mm	medium
2051300	66 x 42 x 45 mm	1.6 mm	large



## Disposable Suction Tubes - Micro

Code	Length	Diameter	Colour	Connection
301004	60 mm	0.4 mm	grey	Luer
301007	60 mm	0.7 mm	black	Luer
301010	60 mm	1.0 mm	white	Luer
301014	60 mm	1.4 mm	red	Luer
302004	70 mm	0.4 mm	grey	Luer
302007	70 mm	0.7 mm	black	Luer
302008	70 mm	0.8 mm	green	Luer
302010	70 mm	1.0 mm	white	Luer
302012	70 mm	1.2 mm	pink	Luer
302014	70 mm	1.4 mm	red	Luer
302016	70 mm	1.6 mm	white	Luer
301020	70 mm	2.0 mm	green	Luer
303004	80 mm	0.4 mm	grey	Luer
303007	80 mm	0.7 mm	black	Luer
303014	80 mm	1.4 mm	red	Luer
303020	80 mm	2.0 mm	green	Luer
303025	80 mm	2.5 mm	transparent	Luer-Lock
303030	80 mm	3.0 mm	transparent	Luer-Lock
303040	80 mm	4.0 mm	transparent	Luer-Lock
304010	100 mm	1.0 mm	white	Luer
304014	100 mm	1.4 mm	red	Luer
304020	100 mm	2.0 mm	green	Luer
304025	100 mm	2.5 mm	transparent	Luer-Lock
304030	100 mm	3.0 mm	transparent	Luer-Lock



- Angled cannula made of stainless steel
- Luer / Luer-Lock connection
- Atraumatic rounded cannula point
- Single sterile packaging
- Wide range of different sizes
- Packaging of 50 pieces per box

## Suction Handle

Code	Description	Package
301000	With tip control, Luer connection, sterile	50 pcs./box



## Disposable Suction Tubes

Code	Length	Diameter
3012020	120 mm	2.0 mm
3012025	120 mm	2.5 mm
3012030	120 mm	3.0 mm
3012035	120 mm	3.5 mm
3012040	120 mm	4.0 mm
3012050	120 mm	5.0 mm
3014020	140 mm	2.0 mm
3014025	140 mm	2.5 mm
3014030	140 mm	3.0 mm
3014035	140 mm	3.5 mm
3014040	140 mm	4.0 mm
3014050	140 mm	5.0 mm
3016020	160 mm	2.0 mm
3016025	160 mm	2.5 mm
3016030	160 mm	3.0 mm
3016035	160 mm	3.5 mm
3016040	160 mm	4.0 mm
3016050	160 mm	5.0 mm



- Tip control
- Blue tube connection
- Angled cannula made of stainless steel
- Atraumatic rounded cannula point
- Single sterile packaging
- Packaging of 50 pieces per box

### For microlaryngoscopy:



3025020	250 mm	2.0 mm	• Distal opening (adherence protection)	• Packaging of 10 pieces per box
3025030	250 mm	3.0 mm		

## Vent tubes


### Titanium

Code	Type	Inner Diameter	Feature	Packaging
10076	Collar Button	0.76 mm		10pcs / box
10100	Collar Button	1.00 mm		10pcs / box
10100F	Collar Button	1.00 mm	w/wire	10pcs / box
10125	Collar Button	1.25 mm		10pcs / box
10125F	Collar Button	1.25 mm	w/wire	10pcs / box
10150	Collar Button	1.50 mm		10pcs / box
10150F	Collar Button	1.50 mm	w/wire	10pcs / box
101100	Shepard	1.00 mm		10pcs / box
101100F	Shepard	1.00 mm	w/wire	10pcs / box
101125	Shepard	1.25 mm		10pcs / box
101125F	Shepard	1.25 mm	w/wire	10pcs / box



### Silicone

Code	Type	Inner Diameter	Feature	Packaging
103100	T-tube, 12 mm	1.10 mm		10pcs / box
103150	T-tube, 5 mm	1.10 mm		10pcs / box
103200	Grommet T-tube	1.27 mm	w/tab	10pcs / box
103300	Paparella I	1.10 mm	w/tab	10pcs / box
103310	Paparella II	1.50 mm	w/tab	10pcs / box



### Fluoroplastic

Code	Type	Inner Diameter	Feature	Packaging
104100	Collar Button	1.00 mm		10pcs / box
104100F	Collar Button	1.00 mm	w/wire	10pcs / box
104125	Collar Button	1.25 mm		10pcs / box
104125F	Collar Button	1.25 mm	w/wire	10pcs / box
104150	Collar Button	1.50 mm		10pcs / box
104150F	Collar Button	1.50 mm	w/wire	10pcs / box
104115	Shepard	1.15 mm		10pcs / box
104115F	Shepard	1.15 mm	w/wire	10pcs / box
104200	Shah	1.15 mm	w/tab	10pcs / box
104210	Shah	1.15 mm		10pcs / box



### Self-Cutting Titanium Ventilation Tube

- Time-saving, separate myringotomy blade is no longer needed
- Myringotomy blade is always sharp
- Easy positioning
- Excellent drainage and ventilation results
- Familiar collar button design
- Sits firmly
- To be used in combination with our 1.2 mm suction tube (Code 302012)



Code	Inner Diameter	Packaging
10125SC	1.25 mm	10 pieces / single, sterile



## Foam Nasal Packing with Latex cover

	Code	Size (mm)
• Application after rhino / septoplasty	35057	50 x 10 x 10
• X-ray contrast strip	35059	70 x 10 x 10
• Non-stick, with thread (30 cm)	35044	50 x 20 x 10
• Sterile packed in pairs of 25 per box	35043	70 x 20 x 10



## Foam Nasal Packing with Latex-free

	Code	Size (mm)
• Application after rhino / septoplasty		
• Water repellent		
• Non-adhesive with thread (25 cm)	35087	35 x 10 x 8
• X-ray contrast strip	35084	50 x 20 x 10
• Sterile packed in pairs of 25 per box	35112	70 x 15 x 10



## Expanding PVA Nasal Tampons

	Code	Size (mm)
• Application after rhinological procedures	66100	45
• Period of usage should not exceed 4 days	66120	45 (with ventilation)
• With 15 cm safety thread	66210	55
• 20 pieces per box	66160	75 (slender)
	66110	80
	66130	80 (with ventilation)
	66200	100



PVA - with airway



PVA - without airway

\* Other sizes of Nasal packings available on request

## PVA Ear Wicks

Post-operative tamponing of the external acoustic canal after tympanoplasty; treatment of inflammation of the external acoustic canal.

With use of a tweezer, place carefully in a dry and compressed state into the external acoustic canal. Afterwards apply thoroughly a water based solution to the tampon, so that it expands. Also available with a central lumen.

### IMPORTANT

In order to eliminate adhesion of the tampon to the tissue after tympanoplasty, it is absolutely necessary, to apply strips of silicone to the eardrum and the walls of the external acoustic canal.

Do not be used for implantation!

Code	Size (mm)	Material	Packaging
66400	9 x 15 (with lumen)	PVA	50 pieces
66410	12 x 24	PVA	20 pieces



## Sterile Covers


For Leica and Wild-Leitz microscopes:

Code	Dimensions	Shape	Lense	Feature	Packaging
2500110H	115 x 165 mm	D	70 mm	angled lense for Leica	10 pieces / box
2500120H	115 x 260 mm	D	70 mm	angled lense for Leica	10 pieces / box
2500130H	137 x 265 mm	D	70 mm	angled lense for Leica	5 pieces / box
2500140H	115 x 300 mm	D	70 mm	angled lense for Leica	10 pieces / box
2500140F	115 x 300 mm	D	70 mm	HQ lense for Leica	10 pieces / box



For Zeiss, Kaps or Moller-Wedel microscopes:

Code	Dimensions	Shape	Lense	Feature	Packaging
2548050F	115 x 180 cm	D	48 mm		10 pieces / box
2565070F	50 x 160 cm	A	65 mm		10 pieces / box
2565010F	125 x 67 x 155 cm	B	65 mm		10 pieces / box
2565080H	125 x 67 x 155 cm	B	65 mm		10 pieces / box
2565060F	135 x 80 x 155 cm	C	65 mm		10 pieces / box
2565020F	115 x 180 cm	D	65 mm		10 pieces / box
2565020H	115 x 180 cm	D	65 mm		10 pieces / box
2565030F	115 x 260 cm	D	65 mm		10 pieces / box
2565030H	115 x 260 cm	D	65 mm		10 pieces / box
2565035H	120 x 260 cm	D	65 mm	with handgrips	10 pieces / box
2565040F	115 x 300 cm	D	65 mm		10 pieces / box
2565040H	115 x 300 cm	D	65 mm		10 pieces / box
2565050F	115 x 210 cm	D	65 mm		10 pieces / box
2565050H	115 x 210 cm	D	65 mm		10 pieces / box
2565090H	137 x 260 cm	D	65 mm		5 pieces / box
2565120F	50 cm		65 mm	cap-shaped	10 pieces / box
2565120FXL	130 cm		65 mm	cap-shaped	10 pieces / box



## Sundry Covers

Code	Dimensions	Type	Packaging
2500010	universal	Microscope handgrips	20 pairs / box
2800070	8 x 150 cm	Drill cable cover	30 pieces / box

## Sheaths for flexible nasopharyngoscopes

- New lens profile for maximum optical integrity
- Innovative sheath material and protective tube facilitate smooth insertion and removal
- Modular connector adaptable to all flexible nasopharyngoscopes
- Tested according to FDA protocol and against BSE prion
- Ultimate patient comfort, atraumatic and easy to use
- Complete product range for both adult and paediatric nasopharyngoscopes
- Expertise and innovation providing custom made solutions
- Boxes of 50 units, sterile and single use, latex free

Code	Colour	Use with
15.1000	red	Pentax FNL 10S/10P2/10BS/10RP3/10RP2, Vision Sciences ENT 2000, 3000, 5000
15.1001	green	
15.1002	blue	Storz 11101 RP1 - Storz 11101 RP2
15.1003	black	Storz 11101 RP
15.1004	white	Olympus ENF P4/GP
15.1005	yellow	Pentax FNL 7RP3, Vision Sciences ENT, 1000 Olympus ENF XP
15.1006	light green	Storz 11101 SK1
15.1	clear	Olympus ENF P2/P3



### External nasal bandage

- For post-operative nasal surgery
- Absorbs nasal secretion
- Cotton compress contains a stretch net bandage
- Fits the nasal contour
- Non-sterile

Code	Size	Packaging
10904	45 mm	40 pieces / box



### Ear Bandages

- High quality material
- Ready for application
- Applied within seconds
- Optimal patient comfort
- Single, sterile packaging



### Ear bandage one-sided

- Easy and quick application
- Post-operative bandage
- One-sided
- Compress included
- Flexible and comfortable
- Single, sterile packaging

Code	Size	Packaging
35000	small	10 pieces / box
35001	medium	10 pieces / box
35002	large	10 pieces / box



### Ear bandage one-sided with Velcro

Code	Size	Packaging
35180	universal	10 pieces / box





Ear bandage two-sided

Code	Size	Packaging
35030	small	5 pieces / box
35029	medium	5 pieces / box
35028	large	5 pieces / box



Ear bandage for otoplasty

Code	Size	Packaging
35033	small	5 pieces / box
35032	medium	5 pieces / box
35031	large	5 pieces / box



Elastic head bandages

- For post-operative head surgery
- Elastic stretch bandage with face opening
- Simple application
- Non-sterile

Code	Size	Packaging
90006	universal	25 pieces / box



Elastic compression bandage

- Compression for incisions in the area below the ear
- With ear slit and fixation strap
- Velcro sealing strip
- Simple application
- Non-sterile

Code	Size	Packaging
90007	universal	5 pieces / box





Exclusive Distributor for Spiggle & Theis

PO Box 565, Balmain, NSW 2041, Australia P +61 02 9810 6947 F +61 02 8072 1556 E [admin@coremed.com.au](mailto:admin@coremed.com.au)  
W [www.coremed.com.au](http://www.coremed.com.au)