


Secu-Pack®

SENGEWALD 

STERIFLEX® Sterilisation Products

The Range of Accessories for Your Sterilisation Packaging.



Our Service

Full service from advice on the packaging and punctual delivery to comprehensive customer support on site by our sales representatives
Well-founded technical know-how

Customer-specific packaging solutions in individual sizes as required.

Please ask your sales representative!

STERIFLEX® products are the ideal supplement to Sengewald packaging and offer the highest level of sterilisation safety.

Indicators for steam formalene or EtO-gas are your visible control assurance.

SENGEWALD 
clinic innovations



Sengewald Klinikprodukte GmbH
Adlerstrasse 2
83101 Rohrdorf-Thansau
Germany
Telephone: +49-80 31-72 24-0
Telefax: +49-80 31-72 24-95
www.sengewald.de
infosengewald@pregis.com

more quality
more safety
more flexibility

STERIFLEX® Adhesive Tapes

- Adhesive tapes with and without printed indicators for steam and gas sterilisation.
- They assure safe fastening of the paper overwrap, i.e. they neither fall off nor tear, even at 134°C in the autoclave.
- After sterilisation it is easy to remove the tapes without residue to open the packets.

Adhesive Tapes

Code #	Size	Packaging Unit
92 19 01	19 mm x 55 m with steam indicator	48
92 19 02	19 mm x 55 m with gas indicator	48
92 19 03	19 mm x 50 m without indicator	48
92 19 04	19 mm x 50 m with steam indicator and stronger adhesive for nonwoven	48
92 19 05	25 mm x 50 m without indicator	36
92 19 06	25 mm x 50 m with steam indicator	36

On request tapes of 50 mm width are also available.

STERIFLEX® Test Cards

There are test cards in two different formats:

- as 3 x 20 cm test cards with triple indicators. They are suitable for sterilisation control and documentation with steam, gas and formalene sterilisation.
- as 4 x 6 cm container control cards with steam indicator for the documentation and control of steam sterilisation.

Test cards with Triple Indicators for Steam, Gas and Formalene Sterilisation

Code #	Size	Packaging Unit
96 03 20	3 x 20 cm	500

Container Control Cards with Steam Indicator

96 04 06	4 x 6 cm	1250
----------	----------	------

STERIFLEX® Semi Crepe Wrapping Papers

- Semi crepe wrapping papers in sheets and on rolls, in white and green, 60 g/m².
- Reliable bacterial impermeability
- Effective steam and gas permeability
- Good foldability and drapeability for creped paper
- High tear resistance
- Good water resistance

Semi creped wrapping papers are suitable as interior and exterior wraps for OR laundry, instrument sieves, OR packs etc. The colours can be used to recognise the interior or exterior overwraps or to recognise use-by dates.

Semi Crepe Wrapping Papers in Sheets and on Rolls

Code #	Size	Packaging Unit
90 40 40	40 x 40 cm	500
90 50 50	50 x 50 cm	500
90 60 60	60 x 60 cm	500
90 75 75	75 x 75 cm	250
90 90 90	90 x 90 cm	250
90 10 10	100 x 100 cm	250
90 12 12	120 x 120 cm	100
91 75 10	roll 75 cm x 100 m	3

STERIFLEX®

High Quality Wrapping Materials

Super 66 Soft Reinforced Paper

Has the following advantages over sterilisation paper:

- It is softer and therefore has far greater drapeability. This makes it easier to work with and allows cleaner and denser packets.
- Memory free, i.e. provides a functional sterile field for sterile work after unpacking.
- Higher stability, i.e. safe protection of the sterile contents and less re-packing.

Super 66 Soft Reinforced Paper

Code #	Size	Packaging Unit
98 40 40	40 x 40 cm	500
98 50 50	50 x 50 cm	500
98 60 60	60 x 60 cm	500
98 75 75	75 x 75 cm	250
98 90 90	90 x 90 cm	250
98 10 10	100 x 100 cm	250
98 12 12	120 x 120 cm	100

Super 77 Nonwoven

The advantages over textiles as a wrapping material are:

- Bacterial impermeability in wet and dry environments
- Water repellent
- Lintfree
- Alcohol repellent
- No reprocessing such as washing and repairs, therefore saving staff
- Permanently available
- Low space requirement

Super 77 Nonwoven

Code #	Size	Packaging Unit
87 50 50	50 x 50 cm	500
87 60 60	60 x 60 cm	500
87 75 75	75 x 75 cm	250
87 90 90	90 x 90 cm	250
87 10 10	100 x 100 cm	250
87 12 12	120 x 120 cm	100
87 13 15	130 x 150 cm	100

WM Bamford & Co Ltd

12 Victoria Street
Private Bag 31346
Lower Hutt 5040, New Zealand

p 04-576-2100
f 04-569-6489
e enquiries@bamford.co.nz
www.bamford.co.nz

