



**Contisure**<sup>®</sup>

# DISPOSABLES

CONTINENCE CARE from Synergy Healthcare



**Contisure**

new and improved range!

for Stress  to Severe  incontinence

# CONTISURE PREMIUM ELASTICATED SHAPED PADS

Synergy Healthcare's comprehensive range of Contisure disposable continence care products are designed to meet the needs of the NHS and the private healthcare sector.

Not only have we revised and improved our range, but have also introduced new products and developed old favourites.

Synergy customers will continue to enjoy first class support from our team of nurse advisors.

As always, a range of clear and concise educational literature is available for your patients.

Our range of Sahara washable continence care products, Conti patient cleansing wipes and Conti urology products mean that Synergy Healthcare really can offer you complete continence care.

**Contisure Premium is a range of elasticated shaped pads that are soft and fluffy for excellent comfort, fit and absorbency.**

**They are suitable for men and women with stress to severe incontinence.**









### The benefits

- 2 strands of Lycra® elastics in the crotch area for an extremely comfortable fit
- Soft, water-repellent standing cuffs which act as a leakage barrier (on Normal, Extra, Super and Super Plus only)
- Integrated wetness indicator
- Anatomically-shaped crotch area to fit the contours of the body
- Rustle-free plastic backing sheet
- Profiled core with Super Absorbent Polymers (SAP) for absorbency where it's needed most
- Stay-dry topsheet to reduce wetback and promote dry skin
- Packed in 28's



PADS MOST EFFECTIVE with Stretch Pants See page 6 for details

### CONTISURE PREMIUM ELASTICATED SHAPED PADS

Product	Extra Light	Light	Normal	Extra	Super	Super Plus
Symbol						
Level of incontinence	Stress	Light	Moderate	Heavy	Heavy	Severe
Colour	White	White	Blue	Yellow	Green	Lilac
Absorbency	500ml	800ml	1300ml	1750ml	2100ml	2350ml
Standing Cuffs	-	-	✓	✓	✓	✓
Case Quantity	8 x 28	6 x 28	4 x 28	4 x 28	3 x 28	3 x 28
Synergy Code	CSP160	CSP170	CSP400	CSP410	CSP420	CSP425
NPC Code	CFP412	CFP348	CFP349	CFP376	CFP377	CFP414

# CONTISURE PREMIUM MICRO SHAPED PADS

Contisure Premium Micro Shaped Pads are small, slim and discreet, offering leakage protection for stress to light incontinence.

With an adhesive backing strip, they offer the ideal solution for women and men with active lifestyles.








## The benefits

- Elasticated for an improved fit
- Anatomically-shaped for improved comfort
- Adhesive backing strip to secure pad in place
- White rustle-free plastic backing sheet for discreetness
- Compact packaging to reduce space and for easier handling
- Mini, Micro, Micro Plus and Micro Super now available in packs of 28
- Mini, Micro, Micro Plus individually wrapped



## CONTISURE PREMIUM MICRO SHAPED PADS

Product	Mini	Micro	Micro Plus	Micro Super	Male Pouch
Symbol					
Level of incontinence	Stress	Stress	Light	Light	Stress
Colour	White	White	White	White	White
Absorbency	200ml	340ml	600ml	930ml	265ml
Case Quantity	12 x 28	12 x 28	12 x 28	12 x 28	12 x 14
Synergy Code	CSP110N1	CSP120N1	CSP130N1	CSP140N1	168600
NPC Code	CFP410	CFP411	CFP208	CFP413	-



for Stress  to Heavy  incontinence

and **Containment** for faecal incontinence

# SELECT

## CONTISURE SELECT NON-ELASTICATED SHAPED PADS

Suitable for men and women with stress to heavy urinary incontinence, the Contisure Select range of non-elasticated shaped pads are high on performance yet low on cost.

The range now includes pads for faecal incontinence, available in two sizes.



### The benefits






- Improved pulp mat with Super Absorbent Polymers (SAP) is now even stronger and more absorbent, and gives better all-round wetness protection (not faecal pad)
- 2 new pads to tackle faecal incontinence
- Rustle-free plastic backing sheet
- Integrated wetness indicator
- Contoured crotch area for improved fit
- Compact packaging
- Packed in 28's

PADS MOST EFFECTIVE  
with  
**Stretch  
Pants**  
See page 6 for details





**New  
faecal  
pads!**

### CONTISURE SELECT NON-ELASTICATED SHAPED PADS

Product	Extra Light	Light	Normal	Extra	Super
<b>Symbol</b>					
<b>Level of incontinence</b>	Stress	Light	Moderate	Heavy	Heavy
<b>Colour</b>	White	White	Blue	Yellow	Green
<b>Absorbency</b>	500ml	800ml	1300ml	1750ml	2100ml
<b>Case Quantity</b>	6 x 28	4 x 28	4 x 28	4 x 28	3 x 28
<b>Synergy Code</b>	755060	755070	CSS300	CSS310	CSS320

### FAECAL PADS

Small	Large
	
Faecal Containment	Faecal Containment
White	White
870ml	1000ml
4 x 28	4 x 28
<b>755070F</b>	<b>CSS330F</b>

All absorbencies quoted are Total Absorbency as determined by a Total Saturation Test over 30 minutes. International Standard Method ISO 11948-1

# CONTISURE RECTANGULAR PADS

With 5 products in the range, Contisure Rectangular Pads are designed for men and women with stress to moderate incontinence.

Slim and discreet, they are comfortable to wear and contain Super Absorbent Polymers (SAP) for improved absorbency and drier skin.



## The benefits






- Rustle-free waterproof backing for security
- All rectangular pads have an adhesive backing strip to secure pad in place (excluding insert)
- Maxi pad can be used for faecal incontinence
- Soft surface ensures patient comfort
- Now available in packs of 28



PADS MOST EFFECTIVE  
with  
**Stretch  
Pants**  
See page 6 for details



## CONTISURE RECTANGULAR PADS

Product	Mini (with adhesive)	Midi (with adhesive)	Midi Plus (with adhesive)	Maxi (with adhesive)	Insert
Symbol					
Level of incontinence	Stress	Light	Light	Moderate	Light
Size	11 x 34cm	11 x 43cm	15 x 46cm	15 x 60cm	12 x 35cm
Colour	Blue	Blue	Blue	Blue	White
Absorbency	360ml	550ml	800ml	1100ml	620ml
Case Quantity	12 x 28	8 x 28	8 x 28	4 x 28	8 x 28
Synergy Code	BN021	BN121	BN131	BN221	BN010



# CONTISURE STRETCH PANTS

Contisure disposable pads are most effective when worn as a 2-piece system with Contisure Stretch Pants.

Manufactured in breathable fabric, they are designed to hold the pads firmly in place to reduce the risk of leakage and to ensure the patient feels confident and comfortable.








## The benefits

- Comfortable for men, women and children
- Longer leg increases patient comfort and security and helps secure the pad in place
- Colour-coded waist and leg band for easy identification and selection
  - Air-permeable, breathable fabric
  - No side seams provide improved comfort for sensitive skin
  - Completely washable
  - XXL now has even more stretch
- Cost-effective packs of 10



## CONTISURE STRETCH PANTS

Product	Small	Medium	Large	Extra Large	XX Large
Symbol					
Colour	Red	Blue	Brown	Green	Blue/Green
Waist (inches) (centimetres)	Up to 31" 80cm	23 - 35" 60 - 90cm	31 - 43" 80 - 110cm	39 - 55" 100 - 140cm	55"+ 140cm+
Case Quantity	5 x 10	5 x 10	5 x 10	5 x 10	5 x 10
Synergy Code	SPS10	SPM10	SPL10	SPX10	SPXX10N
NPC Code	CFP417	CFP418	CFP419	CFP415	CFP416

# PULL-UPS

Pull-Ups are designed to be pulled up and down like normal underwear and are ideal for men and women with incontinence who are both mobile and independent.

Available in small, medium and large they are the ultimate disposable pant range.






## The benefits

- Can be pulled up and down like normal underwear without any adhesive
- Tear-away side seams allow quick and easy removal
- Designed to fit most body shapes
- Soft and rustle-free with textile-feel fabric
- Standing cuffs and elastic at crotch area ensure leakage protection



## PULL-UPS

Product	Small	Medium	Large
Symbol			
Level of incontinence	Moderate	Moderate	Moderate
Waist (inches) (centimetres)	23 - 35" 60 - 90cm	31 - 47" 80 - 120cm	39 - 57" 100 - 145cm
Colour	White	White	White
Absorbency	1320ml	1460ml	1590ml
Case Quantity	8 x 14	8 x 14	8 x 14
Synergy Code	1301	1302	1303

All absorbencies quoted are Total Absorbency as determined by a Total Saturation Test over 30 minutes. International Standard Method ISO 11948-1

# CONTIFIT ALL-IN-ONES

This improved range of all-in-ones features an Advanced Odour Reduction System and is designed for men and women with moderate to very severe incontinence.

The 3-layer retention core is now slimmer, with improved absorption speed and better liquid distribution.





A wide range of sizes are available, from extra small to extra large and can easily be identified using our new colour coding.







## The benefits

- 3-layer core now slimmer and more discreet
- Advanced odour reduction
- 4 re-sealable tapes provide security and easy fastening
- Elasticated waist offers a comfortable fit
- Skin-friendly top sheet
- Standing cuffs provide extra security against leaks
- Integrated wetness indicator
- Now available in packs of 28 and 14



CONTIFIT ALL-IN-ONES				
Product	Extra Small	Small	Medium	Large
<b>Symbol</b>				
<b>Level of incontinence</b>	Moderate	Heavy	Heavy	Severe
<b>Colour</b>	White	Yellow	Blue	Green
<b>Absorbency</b>	1500ml	2000ml	2150ml	2400ml
<b>Waist (inches)</b>	20 - 26"	24 - 35"	28 - 51"	39 - 59"
<b>(centimetres)</b>	50 - 65cm	60 - 90cm	70 - 130cm	100 - 150cm
<b>Case Quantity</b>	4 x 28	8 x 14	4 x 14	4 x 14
<b>Synergy Code</b>	BN055	EBN0601	EBN07014	EBN17014

CONTIFIT PLUS ALL-IN-ONES				
	Small	Medium	Large	Extra Large
<b>Symbol</b>				
<b>Level of incontinence</b>	Severe	Very Severe	Very Severe	Very Severe
<b>Colour</b>	Yellow	Blue	Green	White
<b>Absorbency</b>	2350ml	2900ml	3200ml	3700ml
<b>Waist (inches)</b>	24 - 35"	28 - 51"	39 - 59"	51 - 67"
<b>(centimetres)</b>	60 - 90cm	70 - 130cm	100 - 150cm	130 - 170cm
<b>Case Quantity</b>	8 x 14	4 x 14	4 x 14	4 x 14
<b>Synergy Code</b>	EBN060P1	EBN070P14	EBN170P14	14348140



# CONTISURE UNDERPADS & SUPERDRY PROTECTION SHEET

## CONTISURE UNDERPADS





Available in 4 sizes to cover a variety of uses, the Contisure Underpads range is made with virgin fluff pulp. They are fully sealed with no plastic edges exposed to the skin to ensure optimum patient comfort and maximum security.

### The benefits

- Diamond embossing for improved wicking of fluid
- Compact packaging for easier storage
- 4 sizes



### CONTISURE UNDERPADS

Size	60 x 40cm	60 x 60cm	60 x 75cm	60 x 90cm
Symbol				
Level of incontinence	Light	Moderate	Moderate	Heavy
Colour	Blue	Blue	Blue	Blue
Absorbency	950ml	1400ml	1600ml	1930ml
Case Quantity	4 x 60	4 x 25	4 x 30	4 x 25
Synergy Code	BN111D	BN211D	BN311D	BN411D

## CONTISURE SUPERDRY PROTECTION SHEET





In a choice of 4 sizes, our unique Contisure SuperDry Protection Sheet combines 100% virgin fluff pulp and Airlaid SuperCore® technology with Super Absorbent Polymers (SAP) to offer outstanding benefits for a wide range of applications.

### The benefits

- High volumes can be absorbed and wicked away quickly
- Skin stays dry and near to normal pH maintaining skin integrity
- Helps reduce risk of bacteria colonisation
- No puddling reducing change frequency
- Packed in 28's



### CONTISURE SUPERDRY PROTECTION SHEET

Size	60 x 40cm	60 x 60cm	60 x 75cm	60 x 90cm
Symbol				
Level of incontinence	Moderate	Moderate	Heavy	Heavy
Colour	Blue	Blue	Blue	Blue
Absorbency	1045ml	1540ml	1760ml	2123ml
Case Quantity	8 x 28	4 x 28	4 x 28	4 x 28
Synergy Code	CUP604	CUP606	CUP607	CUP609

All absorbencies quoted are Total Absorbency as determined by a Total Saturation Test over 30 minutes. International Standard Method ISO 11948-1



For more information simply contact Synergy Healthcare on **0161 624 5641** or visit our website at **[www.synergyhealthcare.plc.uk](http://www.synergyhealthcare.plc.uk)**

Your Synergy Approved Supplier

Synergy™  Healthcare

Synergy Healthcare (UK) Ltd, Patient Care Division, Lion Mill, Fitton Street, Royton, Oldham OL2 5JX. England.  
Tel: 0161 624 5641 Fax: 0161 627 0902 Email: [patientcare@synergyhealthcare.plc.uk](mailto:patientcare@synergyhealthcare.plc.uk)

**[www.synergyhealthcare.plc.uk](http://www.synergyhealthcare.plc.uk)**

Sahara® is a Registered Trade Mark of Synergy Healthcare plc. Contisure, Clinisan, Conticath and Contibag are Trade Marks of Synergy Healthcare plc.  
Lycra® is a Registered Trade Mark of Dupont. SuperCore® is a Registered Trade Mark of McAirmaid's.



Certificate No: 096/93