



REFERENCE	DESCRIPTION	PACKAGING
HPF 13040	1 Breast shield with FTC, 1 feeding bottle 120 ml with cap, 100 cm soft tubing	60 pcs/box
HPF 13086	1 Breast shield 22 mm, with valve, 1 feeding bottle 50 ml with cap, 100 cm soft tubing with connector	60 pcs/box
HPF 13087	1 Breast shield 22 mm, with valve, 1 feeding bottle 120 ml with cap, 100 cm soft tubing with connector	60 pcs/box
HPF 13070	1 Breast shield 25 mm, with valve, 1 feeding bottle 50 ml with cap, 100 cm soft tubing with connector	60 pcs/box
HPF 13071	1 Breast shield 25 mm, with valve, 1 feeding bottle 120 ml with cap, 100 cm soft tubing with connector	60 pcs/box
HPF 13089	1 Breast shield 30 mm, with valve, 1 feeding bottle 50 ml with cap, 100 cm soft tubing with connector	60 pcs/box
HPF 13088	1 Breast shield 30 mm, with valve, 1 feeding bottle 120 ml with cap, 100 cm soft tubing with connector	60 pcs/box
HPF 13051	2 Breast shields with FTC, 2 feeding bottles 120 ml with cap, 60 cm with Y-connector and 2 x 40 cm soft tubing	30 pcs/box
HPF 7042	Y-Connector (for double pumping) for the connection of 2 breast milk collector sets to one milk pump	50 pcs/box



BELDICO

... We simply love babies



The disposable breast milk collector sets

'A breast shield with a proven efficiency'



BELDICO S.A./N.V.

rue Andre Feher 1 • Zoning de Aye
6900 Marche-en-Famenne • Belgium
Tel.: +32 (0)84 31 37 80 • Fax: +32 (0)84 31 37 87
E-mail: info@beldico.be • Website: www.beldico.be

Member of the:
Beldico
group



INTENTIONS OF USE

The Beldico disposable breast shield allows the easy collection of the mother's milk.

The mother will benefit from the advantages of a disposable hygienic breast shield. The baby will fully profit from the benefits of the mother's milk.

Single breast shield 30 mm with valve



Beldico breast shields



Single breast shields with FTC or valve and feeding bottle



USER FRIENDLY

The milk is collected directly into the Beldico disposable feeding bottles of 50, 120 or 240 ml.

These collection bottles can be stored in the refrigerator or the freezer.

Following the mother's needs, the disposable breast milk collector sets are available in different sizes: 22 mm, 25 mm or 30 mm. Our sets are also suitable for double pumping.



HYGIENIC

For hygiene reasons, the disposable breast shields are presented as single, ready to use products.

The Beldico breast shields are supplied individually wrapped and ETO sterilised.



COMFORTABLE

The form and sizes suit all types of breasts, even in case of inverted nipples.

This disposable breast shields will help the mother to pump smoothly.

It allows an efficient expression of the mother's milk.

Thanks to the ergonomic design and the choice between different sizes, the breast shields assure a perfect and comfortable fit to the mother.



PRACTICAL USE

The breast milk collector sets are compatible with most of the milk pumps. One set is available for use with vacuum pump systems in the hospitals.

The disposable breast milk collector sets with valve are available in different sizes. This range of sets guarantees for each mother a smooth contact with the breast shield and a good aspiration with all the milk pumps available in the hospital.

One of our disposable breast shields has finger tip control, which makes the suction duration very flexible with a vacuum pump system.

The disposable breast milk collector sets are delivered with soft pre-connected tubing and connectors which fit on the different milk pumps.

