


# REFERENCE

# DESCRIPTION


# PACKAGING



37041

Pacifiers for newborns

8 x 50 pcs/box



37040

Pacifiers for premature babies

8 x 50 pcs/box

*Created by*



**BELDICO S.A.**

Rue Andre Feher 1 - Zoning de Aye

6900 Marche-en-Famenne - BELGIUM

Tel. +32 (0) 84 31 37 80

Fax +32 (0) 84 31 37 87

[www.beldico.be](http://www.beldico.be)

*Distributed by*



**BELDICO** ... *We simply love babies*



**The disposable pacifiers for  
neonates and paediatrics**

*'A new generation of disposable pacifiers'*





## INTENTIONS OF USE

- **To satisfy the “non-nutritive” suction needs** of newborn, pre born infants (neonate) and children with health problems and / or colic disturbances.

- **To stimulate the natural suck-reflex** (decrease the use of fingers and thumbs).

- **To facilitate the weaning process** (process by which one changes from one method of feeding to another).

- Beldico pacifiers can be used **for neonates and paediatrics who receive CPAP ventilation**.

- Beldico pacifiers can be used **for analgesic treatment** (30 % glucose or sucrose solutions in order to decrease pain, crying and stress).

- Beldico pacifiers can be used **with oral and nasal tubes**.



*Pacifier for premature babies (from 1 to 2 kg)*

### ODOUR-FREE AND FLAVOURLESS

The Beldico pacifiers are made of synthetic thermoplastic rubber and are therefore completely odour-free and flavourless.

### COMFORTABLE

This new generation of disposable pacifiers, created by Beldico, are made in one piece and thanks to the anatomical shape, they assure a perfect and comfortable position in the baby's mouth.

### HYGIENIC

In order to reduce the risk of contamination and infection in the hospital environment, the disposable pacifiers are presented as single use products. The Beldico pacifiers are delivered individually wrapped and gamma treated.



*Pacifier for newborn (> 2 kg and ≤ 6 kg)*

### SAFE

For the production of the pacifiers, Beldico uses a raw material that does not contain toxic derivatives. The Beldico pacifiers are conform EN1400-part3 of the European Regulation for soothers.