

Cost effective tubular bandaging



comfi



Comfiast®



A range of colour coded tubular bandages for use in dressing retention.

 **comfi**
elasticated viscose stockinette **fast**®

Available
on Drug Tariff
and NHS
Contract!



Also for wet & dry wrapping in the treatment of atopic eczema

Comfinette®, Comfigauz® and Comfitube®



**comfi
nette**[®]
plain viscose
stockinette

**comfi
gauz**[®]
plain cotton
stockinette

**comfi
tube**[®]
plain cellulose
stockinette

A cost effective range
of tubular bandages
used for dressing
retention



Available
on Drug Tariff
and NHS
Contract!

Comfigrip®



**comfi
grip**[®]
elasticated tubular
bandage B.P.

A range of tubular
bandages providing
light support for muscle
sprains, strains and
joints

Available
on NHS
Contract!

Ideal for dressing retention and for wet & dry wrapping in the treatment of atopic eczema.

Comfifast® is a colour coded range of tubular bandages for dressing retention and for wet wrapping in the treatment of atopic eczema.

Comfifast® can also be used with topical treatments to prevent staining of clothes and retaining stump dressings.

Made from a light breathable fabric, Comfifast® allows freedom of movement and helps reduce skin irritation.

- **Simple and easy to apply**
- **5 colour-coded sizes**
- **Double stripe for easy identification**
- **Cost-effective - washable and re-usable**
- **Comfortable, light breathable fabric allows complete freedom of movement**
- **A range of support materials available.**

Comfifast®: Elasticated Viscose Stockinette
Material: 92% viscose, 5% elastane, 3% polyamide



Ordering Details

DT - on Drug Tariff (Drug Tariff Specification 46)
NHS - on NHS Contract

Width (Unstretched)	Application	Limb Circumf.	1m REF (PIP)	3m REF (PIP)	5m REF (PIP)	10m REF (NPC)
Red Line 3.5cm	Small limbs	8 - 15cm	F12 ^{DT} 289-1745			F15 ^{NHS} EGP005
Green Line 5.0cm	Small/medium limbs	10 - 25cm	F22 ^{DT} 289-1729	F23 ^{DT} 289-1794	F24 ^{DT} 289-1828	F25 ^{NHS} EGP006
Blue Line 7.5cm	Large limbs	20 - 45cm	F32 ^{DT} 289-1760	F33 ^{DT} 289-1802	F34 ^{DT} 289-1836	F35 ^{NHS} EGP007
Yellow Line 10.75cm	Extra large limbs, heads, children's trunks	35 - 65cm	F42 ^{DT} 289-1778	F43 ^{DT} 289-1810	F44 ^{DT} 289-1844	F45 ^{NHS} EGP008
Beige Line 17.5cm	Adult trunks	50 - 120cm	F52 ^{DT} 289-1786			F55 ^{NHS} EGP009

+ Healthcare
comfi fast
Elasticated Viscose
Stockinette

Tubular Bandage



Synergy + Healthcare
comfi fast
Elasticated Viscose Stockinette

Tubular Bandage

for Dressing Retention and Wet Wrapping

Beige Lines - 17.5cm x 10 metre
for adult trunks

Synergy + Healthcare
comfi fast
Elasticated Viscose Stockinette

Tubular Bandage

for Dressing Retention

Yellow Lines - 10.25cm
for extra large limbs and arms

Synergy + Healthcare
comfi fast
Elasticated Viscose Stockinette

Tubular Bandage

for Dressing Retention and Wet Wrapping

Red Lines - 3.5cm x 1 metre
for small limbs

Light support for muscle strains, sprains and joints.

A simple to use elasticated tubular bandage providing light support for muscle sprains, strains and also joints.



- Wide range of sizes in natural shade designed to fit all areas of the body
- Remains secure without using pins or tape
- Cost-effective - washable and re-usable
- Comfortable for patients - allows complete freedom of movement

Comfigrip®: Elasticated Tubular Bandage B.P.
Material: 83% cotton, 9% elastodiene, 8% polyamide

Ordering Details

DT - on Drug Tariff (Drug Tariff Specification 46)
NHS - on NHS Contract

Size	Width (Unstretched)	Application	Limb Circumf.	0.5m Natural REF (PIP)	1m Natural REF (PIP)	10m Natural REF (NPC)
A	4.50cm	Infant feet or arms	10 - 12.4cm			GAN5 ^{NHS} EGA020
B	6.25cm	Small hands and arms	12.5 - 14.5cm	GBN1 ^{DT} 289-1851	GBN2 ^{DT} 289-1927	GBN5 ^{NHS} EGA022
C	6.75cm	Medium arms, small ankles	15 - 24cm	GCN1 ^{DT} 289-1968	GCN2 ^{DT} 289-1919	GCN5 ^{NHS} EGA024
D	7.50cm	Large arms, medium ankles, small knees	25 - 35cm	GDN1 ^{DT} 289-1950	GDN2 ^{DT} 289-1901	GDN5 ^{NHS} EGA005
E	8.75cm	Large ankles, medium knees, small thighs	36 - 44.5cm	GEN1 ^{DT} 289-1943	GEN2 ^{DT} 289-1893	GEN5 ^{NHS} EGA006
F	10.00cm	Large knees, medium thighs	45 - 50cm	GFN1 ^{DT} 289-1935	GFN2 ^{DT} 289-1877	GFN5 ^{NHS} EGA007
G	12.00cm	Large thighs	51 - 60cm	GGN1 ^{DT} 289-1869	GGN2 ^{DT} 289-1885	GGN5 ^{NHS} EGA008
J	17.50cm	Small trunks	61 - 70cm			GJN5 ^{NHS} EGA026
K	21.50cm	Medium trunks	71 - 80cm			GKN5 ^{NHS} EGA028
L	32.50cm	Large trunks	81 - 90cm			GLN5 ^{NHS} EGA039
M	36.50cm	Ex large trunks	91 - 110cm			GMN5

Comfinette®, Comfigauz® and Comfitube®

A cost effective range of tubular bandages for dressing retention

Comfinette®, Comfigauz® and Comfitube® are a lower cost range of tubular bandages used for dressing retention.

Ideal as a skin protector prior to wool band padding on leg ulcers etc.

Also for prevention of soiling to clothes when using topical creams and for total body coverage for psoriasis treatment and eczema.

- **Comfinette®** - cost-effective for frequent dressing changes
- **Comfigauz®** - ideal for use where natural fibres are preferred
- **Comfitube®** - water-repellent containing natural oils
- Particularly ideal for fingers and toes
- Each available in various widths to fit all areas of the body



Comfinette®: Plain Viscose Stockinette
Material: 100% viscose

Comfigauz®: Plain Cotton Stockinette
Material: 100% bleached cotton

Comfitube®: Plain Cellulose Stockinette
Material: 100% cotton

Ordering Details

NHS - on NHS Contract

Size	Application	Comfinette® 20m REF (NPC)	Comfigauz® 20m REF (NPC)	Comfitube® 20m REF (NPC)
00	Toes	N006 NHS EGJ045	Z006 NHS EGJ051	T006
01	Fingers and toes	N016 NHS EGJ041	Z016 NHS EGJ035	T016 NHS EGJ055
12	Bulky dressings for fingers and toes	N126 NHS EGJ042	Z126 NHS EGJ036	T126 NHS EGJ056
34	Children's limbs. Adult hands	N346 NHS EGJ046	Z346 NHS EGJ052	T346 NHS EGJ057
56	Adult limbs	N566 NHS EGJ043	Z566 NHS EGJ037	T566 NHS EGJ058
78	Large adult limbs	N786 NHS EGJ044	Z786 NHS EGJ038	T786 NHS EGJ039
T1	Small body bandage. Adult head bandage	NT16 NHS EGJ047	ZT16 NHS EGJ053	TT16 NHS EGJ059
T2	Large body bandage	NT26 NHS EGJ048	ZT26 NHS EGJ054	TT26 NHS EGJ060

Also from the Comfi range...

Comfi[®] Easywrap Suits

Top to toe convenience for wrapping



 **comfi**
easywrap™ suits **fast**[®]



Available
on Drug Tariff
and NHS
Contract!

A range of
convenient ready to
wear clothing for
wet & dry wrapping

For further information contact Synergy Healthcare on **0161 624 5641**

Synergy™  Healthcare

Synergy Healthcare (UK) Limited, Patient Care Division
Lion Mill, Fitton Street, Royton, Oldham OL2 5JX, England
Tel: 0161 624 5641 Fax: 0161 627 0902
Email: patientcare@synergyhealthcare.plc.uk

www.synergyhealthcare.plc.uk

Comfi[®], Comfigrip, Comfinette, Comfigauz and Comfitube are registered Trade Marks of Synergy Healthcare plc.
The Comfi Range is manufactured in England by Synergy Healthcare (UK) Limited.



Certificate No: 096/93