



Product Catalogue

Effective March 2014

www.defries.com.au



Custom Procedure Packs
Surgical Drapes & Gowns
Single Use Medical Consumables

www.defries.com.au

Contents

| | | |
|--|----------------------------------|-----|
| | Introduction | 2 |
| | Protective Apparel | 6 |
| | Medical Towels | 10 |
| | Surgical and Equipment Drapes | 13 |
| | Custom Procedure Packs | 17 |
| | Disposable Bowls and Plasticware | 21 |
| | Anaesthetics | 27 |
| | Neuro | 32 |
| | Catheter Packs | 35 |
| | Theatre Consumables | 40 |
| | Bandages | 48 |
| | Protective Covers and Bags | 51 |
| | Surgical Tape | 54 |
| | Silicone Drains | 57 |
| | Silicone Loops and Boots | 60 |
| | Silicone Tubing Range | 62 |
| | Ward Consumables | 65 |
| | Custom Printed Labels | 71 |
| | Intensive Care | 75 |
| | Instruments | 79 |
| | Instrument Packs | 87 |
| | Medical Imaging | 91 |
| | Maternity Products | 95 |
| | ENT | 101 |
| | Index | 104 |





Defries Industries is a 100% Australian family owned medical business. As a world class manufacturer, our vision has remained clear over the years: We are committed to improving patient care by offering safe and effective single-use, sterile medical products and customised packs to healthcare professionals.

Quality Products

Defries Industries offer innovative products, designed in conjunction with healthcare professionals. We have steadily expanded our product range, in response to healthcare demand of hospitals.

Defries Industries supplies around 3000 products, including a broad range of single use products, from theatre and ward sterile consumables to single-use instruments, and custom-designed theatre procedure packs.





Procedure Packs can incorporate small packs including Catheter and Anaesthetic Packs, to Ward Packs such as Lumbar Puncture Packs and Thyroid Packs, and Theatre Customised Packs including Orthopaedic, Cardiac, Caesarean, Minor and Major Bundles, Neuro, Transplant Packs and in addition Medical Imaging Custom Procedure Packs including Cath Lab Packs and Angiography Packs.

Our worldwide suppliers are scrutinised carefully for the quality of their products, their communication skills and their delivery management. The TGA plays a large role in assuring suppliers chosen are also certified approved manufacturers.

Specialist Manufacturers

Defries Industries manufacture medical devices in an ISO 8 clean room facility. Our clean rooms are equipped with the latest technology, including electronically interlocked doors, clean room tables, ergonomically designed chairs and the latest extrusion systems. The cleanroom airlocks are fitted with the latest sensitised taps and automatic hand dryers.

Our automatic thermoforming packaging machines use top quality medical grade paper and tyvek with a 5-layer co-extruded film to protect the integrity of our products.

As a therapeutic manufacturer in Australia, we are obliged to maintain a TGA licence: please visit our website to view our current certification (www.defries.com.au).





Customer Service

The Customer Service Team at Defries Industries have excellent communication skills and strive to keep you informed about the delivery turnaround time and product availability.

Our Customer Service Manager is a trained nurse and can assist with product queries and medical questions.



Research & Development Program

By recognising international trends in the healthcare sector, Defries Industries designs new product.

We constantly revise and update our practices to be sensitive to environmental impacts. Please refer to our Environmental Impact Policy on www.defries.com.au.





Distribution Centre

Our Warehouse Management Team is committed to storing and transporting your products under ideal conditions. Our tracking system reliably traces your order from dispatch to delivery. This efficient turnaround ensures quick and reliable delivery.

All our products are double boxed or wrapped in plastic sheaths for added protection.

Our warehouse complies with all codes for the storage of sterile devices. Our medical goods are delivered by trained staff.

We also offer a stock impress service for Custom Packs. Please contact us for further information.

Online Service

Our website presents comprehensive product information. You are welcome to email questions and requests from our website. Please visit www.defries.com.au

Service Commitment

Defries Industries has a reputation for delivering to hospitals efficiently and on time. We understand how important turnaround delivery time is.

For Custom Packs we hold a three month supply of stock at Defries Warehouses. For other consumables we aim to have stock available for immediate delivery. For certain specialty products there is a manufacturing turnaround of up to 3 weeks.




> PROTECTIVE APPAREL



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

**PROTECTIVE APPAREL - Surgical Gowns**

- All gowns are designed from a unique material which offers different levels of protection. All fabric's are soft, light weight, comfortable and low linting.
- Adjustable neckline with Velcro Strips. Soft edge Velcro sewn on outer neck to reduce irritation.
- Gown seams are ultrasonically bonded, providing additional strength.
- Gown ties are attached securely.
- Cuffs are manufactured with a soft high level quality cotton.
- All gowns tested successfully for fire retardant CPSC 16 CFR Part 1610.
- All surgical gowns meet ASTM F 1670 and AAMI PB 70.
- All Defries surgical gowns meet the test requirements set by the Advancement of Medical Instrumentation - AAMI. Please contact your Defries sales representative for test results.

CLASSIC STANDARD SURGICAL GOWN AAMI LEVEL 3


- Soft 50gsm Gown wrapped with 1 x Huck Towel 40cm x 68cm and in a SMS Field 91cm x 91cm.
- Suitable for majority of Surgical Procedures.

| | | |
|---------|--------------|----|
| DEF5000 | Medium | 25 |
| DEF5001 | Large | 25 |
| DEF5002 | XLarge | 25 |
| DEF5003 | XXLarge | 25 |
| DEF5004 | XLarge XLong | 25 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

CLASSIC STANDARD SURGICAL GOWN WITH 2 X WHITE TOWELS 32cm x 32cm

| | | |
|---------|--------------|----|
| DEF6000 | Medium | 25 |
| DEF6001 | Large | 25 |
| DEF6002 | XLarge | 25 |
| DEF6003 | XXLarge | 25 |
| DEF6004 | XLarge XLong | 25 |

CLASSIC PLUS REINFORCED SURGICAL GOWN AAMI LEVEL 3

- 63gsm Gown wrapped with 1 x Huck Towel 40cm x 68cm and in SMS Field 91cm x 91cm.
- These gowns are suitable for surgical cases when an increased level of protection is required.
- Reinforced panel on front of gown and sleeves.

| | | |
|---------|--------------|----|
| DEF5005 | Medium | 25 |
| DEF5006 | Large | 25 |
| DEF5007 | XLarge | 25 |
| DEF5008 | XXLarge | 25 |
| DEF5009 | XLarge XLong | 25 |

CLASSIC PLUS REINFORCED SURGICAL GOWN WITH 2 X WHITE TOWELS 32cm x 32cm


| | | |
|---------|--------------|----|
| DEF6005 | Medium | 25 |
| DEF6005 | Large | 25 |
| DEF6005 | XLarge | 25 |
| DEF6005 | XXLarge | 25 |
| DEF6005 | XLarge XLong | 25 |

Code

Description

Box Qty.

COOLMAX SURGICAL GOWN AAMI LEVEL 4



- This gown is wrapped with 1 x HUCK Towel 40cm x 68cm and SMS Field 91cm x 91cm.
- Gowns are a 64 gsm fabric and include an advanced trilaminate fabric which is breathable and cool.

| | | |
|---------|--------------|----|
| DEF5011 | Large | 25 |
| DEF5012 | XLarge | 25 |
| DEF5013 | XLarge XLong | 25 |
| DEF5014 | XXLarge | 25 |

COOLMAX SURGICAL GOWN WITH 2 X WHITE TOWELS 32cm x 32cm

| | | |
|---------|--------------|----|
| DEF6010 | Medium | 25 |
| DEF6011 | Large | 25 |
| DEF6012 | XLarge | 25 |
| DEF6014 | XXLarge | 25 |
| DEF6013 | XLarge XLong | 25 |

MEDICAL SINGLE USE PATIENT GOWN



- Comfortable and breathable SMS three layer fabric which is lightweight and generous to fit.
- Includes secure tie.
- Ideal for patients in Hospital, Day Procedure, GP clinics, Radiology or any setting where a Patient is required to wear a disposable gown.
- Versatile sizing to accomodate all body shapes.

| | | |
|---------------------|--|----|
| PATIENT GOWN SS MED | Patient Gown Short Sleeve Medium Size | 50 |
| PATIENT GOWN SS LRG | Patient Gown Short Sleeve Large Size | 50 |
| PATIENT GOWN SS XL | Patient Gown Short Sleeve X Large Size | 50 |

Code

Description

Box Qty.

MEDICAL SINGLE USE SCRUB SUITS



- Comfortable and breathable fabric which is lightweight and generous to fit.
- Guaranteed performance compared to reusable scrub suits.
- Improvement in infection control as only worn once and low linting fabric
- Ideal for Nurses in the Operating room, ER departments, Doctors and GP's, Vets, Dentists and Students in the Medical setting.
- Versatile sizing range to accomodate all body shapes, total 3 deep pockets on top.

| | | |
|--------------------|---|-----|
| SCRUB SUIT SML TOP | Single Use Disposable Scrub Suit Small Top | 100 |
| SCRUB SUIT SML TOP | Single Use Disposable Scrub Suit Small Pants | 100 |
| SCRUB SUIT MED TOP | Single Use Disposable Scrub Suit Medium Top | 100 |
| SCRUB SUIT MED TOP | Single Use Disposable Scrub Suit Medium Pants | 100 |
| SCRUB SUIT LRG TOP | Single Use Disposable Scrub Suit Large Top | 100 |
| SCRUB SUIT LRG TOP | Single Use Disposable Scrub Suit Large Pants | 100 |
| SCRUB SUIT XL TOP | Single Use Disposable Scrub Suit X Large Top | 100 |
| SCRUB SUIT XL TOP | Single Use Disposable Scrub Suit X Large Pants | 100 |
| SCRUB SUIT XXL TOP | Single Use Disposable Scrub Suit XX Large Top | 100 |
| SCRUB SUIT XXL TOP | Single Use Disposable Scrub Suit XX Large Pants | 100 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

DISPOSABLE GLASSES

| | | |
|---------|-----------------------------------|---|
| DEF1755 | CLEAR (POLYCARBONATE) NON STERILE | 1 |
|---------|-----------------------------------|---|

**FACIAL PROTECTION - NON STERILE**

- Anti-fog both sides with wrap around frame.
- Protects from splashes and sprays from the top down.
- Fits over prescription glasses.

| | | |
|--------|-----------------------------------|----|
| 23-002 | FACE MASK FULL WITH ELASTIC STRAP | 50 |
|--------|-----------------------------------|----|



| | | |
|--------|----------------|----|
| 23-006 | HALF FACE MASK | 50 |
|--------|----------------|----|

COSMIC NEE-NOGGS

| | | |
|---------|---|----|
| 23-001C | FULL FACE SHIELD WITH VENTED FOAM & VELCRO STRAP 20CM | 50 |
|---------|---|----|



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|--------|--------------------------------------|-----|
| 23-501 | EYE SHIELD REPLACEMENT LENS (SPeyes) | 100 |
|--------|--------------------------------------|-----|

| | | |
|--------|--------------------------------|----|
| 23-502 | EYE SHIELD FRAME ONLY (SPeyes) | 25 |
|--------|--------------------------------|----|



| | | |
|--------|------------------------------------|-----|
| 23-505 | EYE SHIELDS PRE-ASSEMBLED (Speyes) | 100 |
|--------|------------------------------------|-----|

**DISPOSABLE BERET CAPS**

| | | |
|---------|-----------------------------|-----|
| DEF1142 | 53cm (21") BLUE NON STERILE | 250 |
|---------|-----------------------------|-----|

**DISPOSABLE OVERSHOES**

| | | |
|---------|---|-----|
| DEF1145 | NON SKID UNIVERSAL OVERSHOES POLYPROPOLYENE | 100 |
|---------|---|-----|



> MEDICAL TOWELS



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

HAND DRYING & STERILE FIELDS**HUCK TOWEL BLUE 40CM X 68CM**
High absorbancy cotton towels.

| | | |
|---------|-----------------|-----|
| DEF1300 | BLUE HUCK TOWEL | 100 |
|---------|-----------------|-----|



| | | |
|---------|-----------------------------------|----|
| DEF2018 | BLUE HUCK TOWEL 5's | 40 |
| DEF1976 | BLUE HUCK TOWEL 1's | 30 |
| DEF2219 | BLUE HUCK TOWEL WITH XRAY ELEMENT | 30 |

HUCK TOWEL WHITE 40CM X 68CM

| | | |
|---------|------------------|----|
| DEF2753 | WHITE HUCK TOWEL | 30 |
|---------|------------------|----|

HUCK TOWEL GREEN 40CM X 68CM

| | | |
|---------|----------------------|----|
| DEF2140 | GREEN HUCK TOWEL 5's | 40 |
|---------|----------------------|----|



| | | |
|---------|----------------------|-----|
| DEF2619 | GREEN HUCK TOWEL 2's | 100 |
| DEF2501 | GREEN HUCK TOWEL 4's | 50 |

HUCK TOWEL PURE COTTON

| | | |
|---------|--|-----|
| DEF2745 | PURE COTTON HUCK TOWEL 40cm x 68cm 1's | 30 |
| DEF2905 | BLUE PURE COTTON HUCK TOWEL 40cm x 68cm NON STERILE | 100 |

DRESSING TOWEL

| | | |
|--------|----------------------------|----|
| DEF470 | DRESSING TOWEL 40cm x 42cm | 50 |
|--------|----------------------------|----|



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

DRESSING TOWEL - T SMALL 32CM X 32CM

| | | |
|--------|--------------------|------|
| DEF865 | DRESSING TOWEL - T | 100 |
| DEF243 | DRESSING TOWEL - T | 1000 |



| | | |
|--------|------------------------|----|
| DEF862 | DRESSING TOWEL - T 2's | 50 |
|--------|------------------------|----|

DRESSING TOWEL - T LARGE 32CM X 60CM

| | | |
|--------|------------------------|-----|
| DEF879 | DRESSING TOWEL - T | 100 |
| DEF857 | DRESSING TOWEL - T 2's | 50 |

DRESSING TOWEL - D SMALL 32CM X 34CM

| | | |
|--------|--------------------|-----|
| DEF471 | DRESSING TOWEL - D | 500 |
|--------|--------------------|-----|



| | | |
|--------|-------------------------|-----|
| DEF656 | DRESSING TOWEL - D 2's | 50 |
| DEF831 | DRESSING TOWEL - D 2's | 300 |
| DEF832 | DRESSING TOWEL - D 10's | 100 |

DRESSING TOWEL - D LARGE 34CM X 60CM

| | | |
|--------|------------------------|----|
| DEF894 | DRESSING TOWEL - D | 50 |
| DEF912 | DRESSING TOWEL - D 2's | 50 |
| DEF910 | DRESSING TOWEL - D 5's | 25 |

PUR COTTON TOWEL - 30CM X 32CM

| | | |
|---------|-------------------------|-----|
| DEF2927 | PURE COTTON TOWEL | 50 |
| DEF2928 | PURE COTTON TOWEL - 2's | 100 |

PUR COTTON TOWEL - 30CM X 60CM

| | | |
|---------|-------------------------|-----|
| DEF2929 | PURE COTTON TOWEL | 100 |
| DEF2930 | PURE COTTON TOWEL - 2's | 100 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

HEALTHCARE WIPES

- Soft and absorbent wipe, which are individually folded in a dispenser box.
- Strong wet strength of wipe.
- Easy to store and use as dispenser boxes.

| | | |
|---------|---|-----|
| DEF3927 | SMALL HEALTHCARE WIPE 32.5cm x 30cm WHITE IN DISPENSER BOX | 100 |
|---------|---|-----|



| | | |
|---------|---|-----|
| DEF3928 | LARGE HEALTHCARE WIPE 32.5cm x 60cm WHITE IN DISPENSER BOX | 100 |
|---------|---|-----|



To order, please contact our
Customer Service Team on

1300 550 278

or see our website for details
www.defries.com.au



> SURGICAL AND EQUIPMENT DRAPES



Code Description Box Qty.

EQUIPMENT DRAPES

SUPREME BACK TABLE COVER

| | | |
|---------|---------------------------------------|----|
| DEF5015 | 112cm x 228cm (44" x 90") | 40 |
| DEF5016 | 152cm x 228cm (60" x 90") | 40 |
| DEF5017 | HEAVY DUTY 152cm x 228cm (60" x 90") | 40 |
| DEF5018 | HEAVY DUTY 178cm x 280cm (70" x 110") | 40 |

MAYO STAND COVERS

| | | |
|---------|--|-----|
| DEF5019 | SUPREME REINFORCED 127cm x 58.5cm, with Adsorbant pad 50cm x 65cm with Cuff and sealed end | 50 |
| DEF761 | SEALED ENDS 62cm x 91cm Plastic | 100 |
| DEF775 | BLUE PLASTIC WITH CUFF 70cm x 137cm | 50 |



| | | |
|---------|---------------------------------|----|
| DEF5022 | HEAVY DUTY X LARGE 76cm x 145cm | 40 |
|---------|---------------------------------|----|

CLEAR PLASTIC DRAPE



| | | |
|--------|-----------------------------|-----|
| DEF742 | 45cm X 61cm | 200 |
| DEF741 | 61cm X 91cm | 25 |
| DEF079 | 61cm X 91cm | 200 |
| DEF749 | 70cm X 120cm | 100 |
| DEF745 | 76cm X 106cm | 100 |
| DEF746 | 91cm X 150cm | 100 |
| DEF752 | 120cm X 140cm DOUBLE PACKED | 100 |
| DEF748 | 120cm X 140cm | 100 |

Code Description Box Qty.

POLYBACKED DRAPES

An absorbant cotton top layer, with a waterproof plastic backing. Soft and low linting
Ideal for catheter & minor procedures


| | | |
|---------|--|-----|
| DEF833 | 40cm x 42cm no Fenestration | 750 |
| DEF1287 | 40cm x 60cm no Fenestration | 100 |
| DEF1288 | 40cm x 60cm with 10cm Circular Fenestration | 100 |
| DEF1289 | 40cm x 60cm with 10cm x 15cm Fenestration and Adhesive on each corner of drape 2.5cm x 2.5cm | 100 |



| | | |
|---------|---|-----|
| DEF1706 | 50cm x 60cm with 10cm Circular Fenestration | 100 |
| DEF2072 | 75cm x 105cm with 10cm x 15cm Fenestration and 5cm Adhesive length along 10cm sides of Fenestration | 25 |

AZURE MINOR PROCEDURE DRAPES

- New viscose fabric combined with a PE Laminate to create a drape capable of delivering a high level of performance to most types of surgery.



| | | |
|---------|--|----|
| DEF5090 | 50cm x 60cm with 6.5 Circular Fenestration with Full Adhesive | 25 |
| DEF5094 | CIRCUMCISION DRAPE 38cm x 45cm with 3cm Circular Fenestration | 50 |
| DEF5091 | 90cm x 90cm with 10cm x 15cm Fenestration with Full Adhesive | 40 |
| DEF5092 | 200cm x 300cm with 15cm x 25cm Fenestration with Full Adhesive | 12 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

CLEAR PLASTIC EPIDURAL DRAPES WITH 10CM X 15CM FENESTRATION AND 2.5CM ADHESIVE AROUND FENESTRATION

| | | |
|---------|--------------|-----|
| DEF1481 | 76cm x 106cm | 100 |
|---------|--------------|-----|



| | | |
|------------|---------------|----|
| DEF748FA-S | 120cm x 140cm | 10 |
|------------|---------------|----|

| | | |
|---------|---------------|-----|
| DEF3394 | 120cm x 140cm | 100 |
|---------|---------------|-----|

CLEAR PLASTIC PATIENT DRAPE

| | | |
|---------|---|----|
| DEF2430 | 91cm x 115cm with 2.5cm x 115cm Adhesive down length of drape | 50 |
|---------|---|----|

PATIENT SURGICAL DRAPES CLASSIC STOCKINETTE

| | | |
|---------|---------------------|----|
| DEF5060 | MEDIUM 20cm x 122cm | 25 |
|---------|---------------------|----|



| | | |
|---------|--------------------|----|
| DEF5061 | LARGE 23cm x 122cm | 25 |
|---------|--------------------|----|

| | | |
|---------|---------------------|----|
| DEF5062 | XLARGE 30cm x 122cm | 25 |
|---------|---------------------|----|

SUPREME DRAPE SHEETS

- A Hydrophilic Trilaminar Fabric specifically designed for top level performance coupled with great absorbency, low lint and patient comfort.
- This fabric is suitable for all levels of surgery.

| | | |
|---------|--------------------------|----|
| DEF5024 | HALF DRAPE 112cm x 145cm | 20 |
|---------|--------------------------|----|

| | | |
|---------|----------------------------|----|
| DEF5025 | MEDIUM DRAPE 112cm x 180cm | 20 |
|---------|----------------------------|----|

| | | |
|---------|-------------------------|----|
| DEF5026 | 3/4 DRAPE 142cm x 193cm | 20 |
|---------|-------------------------|----|

| | | |
|---------|---------------------------|----|
| DEF5028 | LARGE DRAPE 183cm x 215cm | 20 |
|---------|---------------------------|----|

| | | |
|---------|-----------------------------|----|
| DEF5029 | X LARGE DRAPE 178cm x 254cm | 12 |
|---------|-----------------------------|----|

| | | |
|---------|--|----|
| DEF5030 | TOP/HEAD DRAPE 250cm x 135cm with Reinforcement and Adhesive | 15 |
|---------|--|----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|---------|--|----|
| DEF5031 | SIDE DRAPE 110cm x 100cm with Reinforcement and Adhesive | 20 |
|---------|--|----|

| | | |
|---------|---|----|
| DEF5032 | BOTTOM/FOOT DRAPE 190cm x 185cm with Reinforcement and Adhesive | 20 |
|---------|---|----|

SUPREME IMPERVIOUS SPLITSHEET DRAPE

| | | |
|---------|-----------------------------------|----|
| DEF5049 | 152cm x 213cm, Split 76cm x 20cm. | 25 |
|---------|-----------------------------------|----|

| | | |
|---------|----------------------------------|----|
| DEF5050 | 152cm x 183cm, Split 52cm x 15cm | 25 |
|---------|----------------------------------|----|

| | | |
|---------|-----------------------|---|
| DEF5034 | UNIVERSAL PACK | 6 |
|---------|-----------------------|---|

1 x Top Drape 250cm x 135cm.
Reinforcement and Adhesive.

2 x Side Drapes 110cm x 100cm.
Reinforcement and Adhesive.

1 x Bottom Drape 190cm x 185cm.
Reinforcement and Adhesive

1 x Mayo Stand Cover 58cm x 141cm

1 x Backtable Cover 152cm x 225cm

4 x Utility Drapes

| | | |
|---------|-------------------------------------|---|
| DEF5111 | UNIVERSAL PACK - SHORT SIDES | 6 |
|---------|-------------------------------------|---|

1 x Top Drape 250cm x 135cm.
Reinforcement and Adhesive.

2 x Side Drape 110cm x 100cm
Reinforcement and Adhesive.

1 x Bottom Drape 190cm x 185cm
Reinforcement and Adhesive.

1 x Back Table Cover 152cm x 228cm

| | | |
|----------|--|---|
| DEF5110A | VASCULAR HD UNIVERSAL PACK - LONG SIDES | 6 |
|----------|--|---|

1 x Top Drape 250cm x 135cm.
Reinforcement and Adhesive.

2 x Side Drapes 110cm x 183cm
Reinforcement and Adhesive.

1 x Bottom Drape 190cm x 185cm
Reinforcement and Adhesive.

1 x Back Table Cover 152cm x 228cm

Code

Description

Box Qty.

Defries Industries will design and customize your procedural drape. Below are some sample drapes from our range. Please contact your Sales Representative to discuss your drape requirements.



SPECIALTY PROCEDURE DRAPES

ORTHOPAEDIC DRAPES

| | | |
|----------|--|----|
| DEF5035A | ARTHROSCOPY DRAPE 305cm x 450cm with Armboard Covers & Pouch with Port | 6 |
| DEF5036B | HIP DRAPE 295cm x 368cm with Tube Holders and Bilateral Pouches | 6 |
| DEF5037 | LARGE SPLIT U DRAPE 200cm x 300cm with Reinforcement & Adhesive, Split 16cm x 118cm | 12 |
| DEF5038A | EXTREMITY DRAPE WITH ARMBOARDS, 315cm x 365cm, Reinforcement & Adhesive with 6cm Conformable Fenestration | 6 |
| DEF5039A | HAND DRAPE 290cm x 370cm with Velcro Tab | 6 |
| DEF5040 | SHOULDER ARTHROSCOPY DRAPE 254cm x 160cm Teardrop Fenestration 17.5cm x 28cm, Fluid Pouch 115mm x 55cm, Tube holders | 8 |
| DEF5041 | BODY SPLIT DRAPE 196cm x 295cm SPLIT 25cm x 65cm | 12 |

GYNAECOLOGY AND OBSTETRIC DRAPES

| | | |
|---------|--|----|
| DEF5042 | CAESAREAN DRAPE (LSCS) 249cm x 312cm with Pouch Cut Out Incise | 8 |
| DEF5043 | LITHOTOMY T DRAPE - 140cm x 195cm x 74cm, 7.5cm x 15cm Perineal Fenestration | 16 |
| DEF5044 | LAPAROSCOPY DRAPE - 235cm x 250cm. 18cm x 18cm x 23cm Trapezoidal Abdominal Fenestration. 18cm x 7.5cm Perineal Fenestration | 12 |
| DEF5048 | UNDER BUTTOCK DRAPE 90cm x 124cm with Pouch | 20 |
| DEF5047 | UNDER BUTTOCKS DRAPE 90cm x 124cm | 20 |

GENERAL SURGICAL DRAPES

| | | |
|---------|--|----|
| DEF5051 | LAPAROTOMY DRAPE 260cm x 312cm with Fenestration 10cm x 30.5cm and Built-In Incise | 12 |
| DEF5052 | ENDOSCOPY DRAPE 260cm x 312cm x 200cm with Fenestration 28cm x 25.5cm, Pockets & Armboard Covers | 12 |
| DEF5188 | BARIATRIC DRAPE 280cm x 280cm SPLIT 135cm with Reinforcement and Adhesive | 8 |

Code

Description

Box Qty.

CARDIAC AND ANGIOGRAPHY DRAPES

| | | |
|---------|--|----|
| DEF2847 | ANGIOGRAPHY PROCEDURE DRAPE 205cm x 400cm with Formers | 10 |
| DEF5107 | CARDIAC DRAPE 275cm x 360cm Extra Long without Pouches Split 18cm x 250cm | 6 |
| DEF5053 | CV SPLIT DRAPE 254cm x 213cm, SPLIT 18cm x 250cm | 6 |
| DEF5085 | CARDIAC SADDLE 24cm x 110cm with Pocket | 40 |
| DEF5163 | CARDIAC SPLIT DRAPE 275cm x 360cm, SPLIT 18cm x 250cm, Extra Large Reinforced Area | 6 |
| DEF5118 | CV GROIN DRAPE 208cm x 191cm with 2 x SPLITS 50cm x 16cm with Adhesive | 8 |

EENT DRAPES

| | | |
|---------|---|----|
| DEF5057 | EENT SPLIT 196cm x 280cm, Split 7cm x 65cm with Reinforcement & Adhesive Full Length of Split | 12 |
| DEF5058 | HEAD DRAPE TURBAN STYLE 109cm x 85cm with Adhesive & Turban | 12 |
| DEF5120 | TRANSPHENOIDAL DRAPE 198cm x 367cm Diamond Shaped Fenestration 13cm x 13cm | 12 |
| DEF5120 | TRANSPHENOIDAL DRAPE 198cm x 367cm Diamond Shaped Fenestration 13cm x 13cm | 12 |

OPHTHALMIC DRAPES

| | | |
|---------|---|---|
| DEF5123 | OPHTHALMIC DRAPE 198cm x 281cm Clear Incise. Bilateral Pouches. Former and Spare Incise | 6 |
| DEF5065 | OPHTHALMIC DRAPE 166cm x 254cm. 10cm x 15.5cm Fenestration | 6 |

UROLOGY DRAPES

| | | |
|---------|---|---|
| DEF5067 | TUR DRAPE 154cm x 165cm. 6cm Circular Abdominal Built-In Leggings. Pouch with Sieve. Latex Free Finger Cot. | 6 |
|---------|---|---|

NEURO DRAPES

| | | |
|---------|--|----|
| DEF5081 | NEURO DRAPE 256cm x 345cm with Pouch, Reinforcement and Oval Fenestration with Built-In Incise | 8 |
| DEF5120 | TRANSPHENOIDAL DRAPE 198cm x 367cm Diamond Shaped Fenestration 13cm x 13cm | 12 |
| DEF5105 | UNIVERSAL SPINE DRAPE 290cm x 355cm 10cm x 30cm Fenestration. Reinforcement and Adhesive. 2 x Instrument Pouches | 8 |

PAEDIATRIC DRAPES

| | | |
|---------|--|----|
| DEF5122 | NEO NATAL DRAPE 244cm x 170cm Clear Anaesthetic Screen. Circular Fenestration 7cm | 12 |
| DEF5137 | PAEDIATRIC LAPAROTOMY DRAPE (WITH ARMBOARDS) 260cm x 280cm x 200cm Fenestration 15cm x 15cm Reinforcement and Adhesive | 12 |
| DEF5138 | PAEDIATRIC DRAPE (NO ARMBOARDS) 200cm x 280cm. Fenestration 15cm x 15cm Reinforcement and Adhesive | 12 |

> CUSTOM PROCEDURE PACKS



Code Description Box Qty.

CUSTOM PROCEDURE PACKS

- A Custom Procedure pack is designed to contain all the consumables used during a procedure.
- The Individual consumables can range from Single Use Drapes and Gowns, to Diathermy, Suction Tubing, Packs and Gauze.
- Defries assembles these custom packs to your hospitals' specifications.



Clinical Cost Savings


- Custom Packs save setup time for all procedures.
- They ensure all components are available for every procedure.
- A Customised Pack will reduce packaging waste.



- Infection control is maximised by using single use Drapes and Gowns. A reduction in linen ensures minimal lint.
- Custom Packs reduce both inventory and time required to order individual components.
- Emergency cases can be attended to more quickly with a Custom Procedure Pack.
- Our warehouse will hold 3 months stock of Custom Packs, therefore ensuring quick delivery when required.
- Delivery can be to stores or direct to relevant Departments.
- Flexible delivery options can be daily, biweekly or weekly to reduce storage requirements.

Code Description Box Qty.

- Flexibility to change components of Custom Packs.



- Easy traceability of all components in pack.



- Easy accountability of costs per procedure.

Implementation

- Contact your Defries Representative or Customer Service 1300 550 278. We will work with you to discuss the Custom Procedure Pack and components you require.



- Our relationship with our suppliers enables us to competitively quote.
- We now have our own range of Single Use Gowns, Drapes and Plasticware ensuring a complete Defries Custom Pack.
- We offer non sterile or sterile samples, so early changes can be identified.
- Ongoing education and followup is vital with Custom Packs, our Representatives are always available to discuss your specific needs.

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

CUSTOM THEATRE PACKS ORTHOPAEDIC PACKS

| |
|------------------------|
| |
| Ankle Arthroscopy |
| Knee Arthroscopy |
| Shoulder Arthroscopy |
| Arthroscopy |
| Knee Pack |
| Hip Pack |
| Shoulder Pack |
| Bilateral Knee |
| Major Orthopaedic Pack |
| Minor Orthopaedic Pack |
| Wrist/Hand Pack |
| Carpal Tunnel Pack |
| Hand Pack |
| ACL Pack |
| Foot/Ankle Pack |
| Extremity Pack |

LAMINECTOMY PACK

SPINAL CUSTOM PACK

| |
|----------------------|
| Spinal Pack |
| Cervical Fusion Pack |

GENERAL SURGERY

| |
|----------------------------|
| Colposcopy Pack |
| Major Pack |
| Minor Pack |
| Thyroidectomy Pack |
| Laparoscopy Pack |
| Transverse Laparotomy Pack |
| D&C/Hemorrhoidectomy Pack |

CAESAREAN PACKS

| |
|-----------------|
| Caesarean Pack |
| Delivery Bundle |

CARDIAC PACKS

| |
|----------------------|
| Radial Feeder Pack |
| CABG Pack |
| Valve Pack |
| Donor Retrieval Pack |
| Thoracic Pack |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

NEURO PACKS

| |
|-----------------|
| Craniotomy Pack |
| VP Shunt Pack |
| Spinal Pack |



GYNAEOLOGICAL PACK

| |
|----------------------|
| Gynaecological Pack |
| PV Set |
| LAVH Pack |
| Laparoscopy Pack |
| D & C Procedure Pack |

VASCULAR PACKS

| |
|------------------------|
| Femoral/Popliteal Pack |
| Phlebectomy Pack |
| Vascular Pack |
| Varicose Veins Pack |

TRANSPLANT CUSTOM PACK

| |
|------------------------|
| Kidney Transplant Pack |
|------------------------|

LUMBAR PUNCTURE PACK

PLASTICS CUSTOM PACK

| |
|-----------------------|
| Face Lift Pack |
| Breast Reduction Pack |
| Lipectomy Pack |

UNIVERSAL CUSTOM PACK

UROLOGY CUSTOM PACK

| |
|------------------------------|
| Transurethral Resection Pack |
| Urology Prep Set |

Code

Description

Box Qty.

SPECIALIZED CUSTOM PACK

- Tunnel Venus Cath
- Venous Access Pack
- Cervical Pack Spinal and Neuro
- CJD Minor Pack
- TVOR Pack
- Vertebroplasty pack

EENT PROCEDURE PACKS

- Oral Nasal Pack
- Oral Pack
- Excision Set
- Transhenoidal Pack
- Myringotomy Procedure Pack
- Dental Tray
- Dental Bowl Set
- Dental Pack
- Nasal Custom Pack
- Major Jaw Pack

OPHTHALMIC PACKS

- Eye Pack
- Eye Procedure Pack
- Intra-Vitreal Pack
- Extra Ocular pack
- Vitrectomy Pack

MEDICAL IMAGING/CATH LAB CUSTOM PACKS

- Joint Injection Pack
- PTCA Procedure Pack
- Medium Dressing Tray
- Pacemaker Procedure Pack
- Pacemaker Tray



Code

Description

Box Qty.

ANGIOGRAPHY PACKS

- Angio Procedure Pack
- Peripheral Angiography Pack
- CCL Angio Pack
- Coronary Angio Pack



RADIOLOGY PACKS

- Radiology Pack
- Radial Feeder Pack

PICC PACKS

- PICC Line Procedure Pack

PACEMAKER PACKS

- Pacemaker Pack



To order, please contact our Customer Service Team on

1300 550 278

or see our website for details

www.defries.com.au

> DISPOSABLE BOWLS & PLASTICWARE



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

MEASURING CUP

| | | |
|--------|------|-----|
| DEF459 | 60ml | 50 |
| DEF460 | 60ml | 250 |

**non sterile**

| | | |
|---------|------|-----|
| DEF1526 | 60ml | 100 |
|---------|------|-----|

GALLIPOTS

| | | |
|--------|-----------------------------------|----|
| DEF867 | 60ml - 2cm Height, 7.5cm Diameter | 50 |
|--------|-----------------------------------|----|



| | | |
|--------|------|-----|
| DEF860 | 60ml | 250 |
|--------|------|-----|

| | | |
|---------|------------------------------------|----|
| DEF1147 | 150ml - 5cm Height, 7.5cm Diameter | 25 |
|---------|------------------------------------|----|



| | | |
|--------|-------|-----|
| DEF888 | 150ml | 150 |
|--------|-------|-----|

| | | |
|---------|---------------------|----|
| DEF2792 | 150ml PACKED IN 3'S | 25 |
|---------|---------------------|----|

non sterile


| | | |
|---------|-------|-----|
| DEF1114 | 150ml | 500 |
|---------|-------|-----|

| | | |
|---------|-------|-----|
| DEF2362 | 150ml | 200 |
|---------|-------|-----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SPECIMEN JARS

| | | |
|--------|-----------------|----|
| DEF758 | 70ml YELLOW LID | 50 |
|--------|-----------------|----|




| | | |
|--------|-----------------|-----|
| DEF843 | 70ml YELLOW LID | 150 |
|--------|-----------------|-----|

| | | |
|--------|-----------------------------------|----|
| DEF278 | 70ml DOUBLE PACKED 4'S YELLOW LID | 40 |
|--------|-----------------------------------|----|

| | | |
|--------|------------------|-----|
| DEF842 | 120ml YELLOW LID | 150 |
|--------|------------------|-----|



| | | |
|--------|------------------|----|
| DEF840 | 250ml YELLOW LID | 80 |
|--------|------------------|----|



| | | |
|---------|------------------|----|
| DEF1109 | 250ml YELLOW LID | 25 |
|---------|------------------|----|

| | | |
|---------|-----------------|----|
| DEF1139 | 500ml WHITE LID | 36 |
|---------|-----------------|----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SCALPEL HOLDER

| | | |
|---------|--|----|
| DEF2820 | SCALPEL HOLDER (SCALPEL BLADE AND HANDLE NOT INCLUDED) | 50 |
|---------|--|----|

**BLADESAFE SCALPEL TRANSFER DEVICE**

| | | |
|---------|---|----|
| DEF2631 | BLADESAFE SCALPEL TRANSFER DEVICE (SCALPEL BLADE AND HANDLE NOT INCLUDED) | 24 |
|---------|---|----|

**DISPOSABLE JUGS**

| | | |
|---------|---|----|
| DEF2092 | 1.2 LITRE DOUBLE PACKED - LEFT Increments in 100ml to 1000ml | 12 |
| DEF2773 | 1 LITRE TALL - RIGHT Increments in 100ml to 1000ml | 12 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

DISPOSABLE SQUARE TRAYS

| | | |
|---------|--|----|
| DEF2546 | SQUARE FULLY DEEP TRAY 27cm x 25cm x 4cm Deep | 20 |
| DEF321 | DRUG TRAY - DRESSING TRAY 3 COMPARTMENT 23cm x 9cm x 3cm Deep | 80 |



| | | |
|---------|--|----|
| DEF2690 | DRESSING TRAY - TRAY 3 COMPARTMENT 25.5cm x 20.5cm x 2.5cm Deep | 36 |
|---------|--|----|

**MIXING BOWL AND SPATULA**

| | | |
|---------|--|----|
| DEF1530 | MIXING BOWL 600ml & SPATULA DOUBLE PACKED | 10 |
|---------|--|----|



| | | |
|---------|------------------------------------|----|
| DEF1532 | PLASTIC SPATULA 3cm x 17cm x 2.5cm | 10 |
|---------|------------------------------------|----|



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

DISPOSABLE BOWLS AND RECEIVER SETS

| | | |
|---------|-------------------------|----|
| DEF2901 | 500ml BLUE | 20 |
| DEF2902 | 1 LITRE BLUE | 20 |
| DEF2489 | 1.5 LITRES BLUE D/P | 20 |
| DEF2550 | 1.5 LITRE BLUE | 20 |
| DEF2553 | 2.5 LITRE DOUBLE PACKED | 10 |

BOWL DRAPE

- Flexible Bowl Drape fits over re-usable Splash Bowl.
- Tough plastic material that will not tear or break.
- Size 138cm x 138cm, bowl size reinforcement 36cm x 36cm, wrapped in laminated field 50 x 50cm
- Sterile and Latex free.

| | | |
|-------|---|----|
| 28868 | BOWL DRAPE STERILE 138 x 138cm wrapped in a laminated field 50 x 50cm | 24 |
|-------|---|----|

**SPLASH BOWLS**

| | | |
|----------|---|---|
| DEF2650A | DISPOSABLE SPLASH BOWL SET 6 LITRE WRAPPED IN 1 x SMS FIELD 130cm x 150cm | 6 |
| DEF2434A | DISPOSABLE SPLASH BOWL SET 6 LITRE WRAPPED IN 2 x SMS FIELD 130cm x 150cm | 6 |

GUIDEWIRE BOWL

| | | |
|---------|----------------|----|
| DEF2579 | GUIDEWIRE BOWL | 10 |
|---------|----------------|----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

KIDNEY DISHES

| | | |
|---------|------------|----|
| DEF2791 | 700ml BLUE | 35 |
|---------|------------|----|



| | | |
|---------|-------------|----|
| DEF1174 | 700ml CLEAR | 35 |
|---------|-------------|----|



| | | |
|---------|--------------|----|
| DEF2109 | 700ml YELLOW | 35 |
|---------|--------------|----|



| | | |
|---------|--------------|----|
| DEF2551 | 1 LITRE BLUE | 20 |
|---------|--------------|----|

non sterile

| | | |
|---------|-------------|-----|
| DEF1661 | 400ml GREEN | 640 |
|---------|-------------|-----|

| | | |
|---------|--------------------------------------|----|
| DEF409A | KIDNEY DISH AND GALLIPOT PACK | 35 |
|---------|--------------------------------------|----|

1 x Kidney Dish 700ml Clear
1 x Gallipot 150ml
Field Blue Plastic 75cm x 75cm

| | | |
|---------|------------------------|----|
| DEF2837 | KIDNEY DISH SET | 20 |
|---------|------------------------|----|

2 x Gallipot 150ml
2 x Kidney Dish 700ml Clear
Field Blue Plastic 50cm x 50cm

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

BOWL SETS**DEF2600D MINOR BOWL SET 15**

- 2 x Kidney Dish Clear 700ml
- 1 x Gallipot 150ml
- 1 x 500ml Blue Bowl
- 4 x Sharps Tray Yellow Rectangular 9.5cm x 26cm x 4cm deep
- SMS Field 75cm x 75cm

**DEF2601E MAJOR BOWL SET 6**

- 2 x Kidney Dish Clear 700ml
- 1 x Gallipot 150ml
- 2 x 500ml Blue Bowl
- 1 x 6L Blue Bowl
- 4 x Sharps Tray Yellow Rectangular
- SMS Field 130cm x 150cm

**DEF2787 GENERAL BOWL SET 4**

- 1 x 6L Blue Bowl
- 1 x 2.5L Blue Bowl
- 2 x 500ml Blue Bowl
- 5 x Gallipot 150ml
- 1 x Kidney Dish Yellow 700ml
- 2 x Kidney Dish Blue 700ml
- 2 x Kidney Dish Green 700ml
- Double Wrapped Field SMS 130cm x 150cm

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

DEF2541A PERIOPERATIVE KIDNEY DISH SET 20

- 1 x Kidney Dish Clear 700ml
- 2 x Kidney Dish Green 700ml
- 3 x Gallipot 150ml
- Field SMS 75cm x 75cm

DEF2607 BOWL SET 6

- 2 x Kidney Dish Blue 700ml
- 2 x 600ml Bowl Green
- 1 x 2.5L Blue Bowl
- 1 x Bladesafe Scalpel Transfer Device
- 1 x Kidney Dish Yellow 700ml
- 3 x Gallipot 150ml
- Field SMS 91cm x 91cm

DEF2652 SPLASH BOWL SET 6

- 1 x 2.5L Bowl
- 1 x 6L Bowl
- SMS Field 130cm x 150cm

DEF2093 SUTURE BOWL SET 36

- 1 x Dressing Tray 3 Cavity- 25.5cm x 20.5cm x 2.5cm
- 1 x Towel Blue Cotton 40cm x 68cm
- 10 x Gauze Swabs 10cm x 10cm Plain
- Field Plastic 91cm x 115cm

DEF2014A RECEIVER SET 6

- 1 x 2.5 Litre Bowl
- 1 x 1.0 Litre Bowl
- 2 x Kidney Dishes 700ml Clear
- 4 x Gallipots 150ml
- Field SMS 91.5cm x 91.5cm

DEF2433C DISPOSABLE BOWL SET 5

- 1 x Bowl 600ml
- 1 x Bowl 1.5L
- 1 x Bowl 6L
- 4 x Gallipot 150ml
- 1 x Kidney Dish Yellow 700ml
- 3 x Kidney Dish Green 700ml
- 1 x Label "Has time out been done"
- 2 x Field SMS 130cm x 150cm

DEF2549 MAJOR BOWL SET 5

- 1 x Bowl 6L
- 1 x Bowl 2.5L
- 2 x Kidney Dishes 700ml Clear
- 2 x Gallipots 150ml
- 1 x Tray 2 Compartments
- Field SMS 130cm x 150cm

| Code | Description | Box Qty. |
|----------|--|----------|
| DEF2215B | RECEIVER SET | 6 |
| | 1 x 2.5 Litre Bowl | |
| | 5 x Kidney Dishes 700ml Clear | |
| | 2 x Gallipots 150ml | |
| | 1 x Kidney Dish 700ml Yellow | |
| | 1 x Light Handle Cover | |
| | Field Backtable Cover 60" x 90" (152cm x 228cm) | |
| DEF2949 | MAJOR BOWL SET | 6 |
| | 2 x Blue Bowl 500ml | |
| | 1 x Blue Bowl 6 Litre | |
| | 1 x Kidney Dish 700ml Blue | |
| | 1 x Gallipot 150ml | |
| | 2 x Kidney Dish 700ml Clear | |
| | 1 x Kidney Dish 700ml Yellow | |
| | Field SMS (150cm x 200cm) | |
| DEF3180A | STANDARD BOWL SET | 5 |
| | 2 x Blue Bowl 250ml | |
| | 1 x Blue Bowl 6 Litre | |
| | 1 x Kidney Dish 700ml Clear | |
| | 1 x Kidney Dish 1000ml Blue | |
| | 1 x Kidney Dish 700ml Yellow | |
| | Field SMS 150cm x 150cm | |
| DEF3132A | BOWL SET | 10 |
| | 2 x Blue Bowl 500ml | |
| | 1 x Blue Bowl 1.5 Litre | |
| | 2 x Gallipot 150ml | |
| | 2 x Kidney Dish 700ml Clear | |
| | 1 x Kidney Dish 700ml Yellow | |
| | Field SMS 130cm x 150cm | |

| Code | Description | Box Qty. |
|----------|--|----------|
| DEF2066D | MAJOR BOWL SET | 6 |
| | 1 x Suction Tubing 5.5MM X 8.4MM X 3M | |
| | 5 x Abdominal sponges X-RAY 30CM X 48CM | |
| | 1 x Blue Bowl 1 Litre | |
| | 2 x Blue Bowl 2.5 Litre | |
| | 2 x Kidney Dish 700ml Blue | |
| | Foam pad 50mm X 60mm X 25mm | |
| | 1 x Gallipot 150ml | |
| | 1 x Yankauer Sucker without Vent | |
| | 5 x Gauze Swabs 10cm X 10cm X/RAY | |
| | 1 x Needle Card Double magnetic | |
| | 1 x Bladesafe Scalpel Transfer Device | |
| | 1 x Back Table Cover 60" X 90" | |
| DEF3179 | BOWL SET | 16 |
| | 2 x Blue Bowl 250ml | |
| | 1 x Kidney Dish 700ml Clear | |
| | 1 x Blue Bowl 1 Litre | |
| | 1 x Kidney Dish 700ml Yellow | |
| | Field SMS 91cm x 91cm | |



To order, please contact our
Customer Service Team on

1300 550 278

or see our website for details
www.defries.com.au


> ANAESTHETICS



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

GREEN GAUZE SWABS/SPONGES

| | | |
|--------|-----------------|-----|
| DEF190 | 10cm x 10cm 2's | 100 |
|--------|-----------------|-----|



| | | |
|---------|------------------|-----|
| DEF192 | 10cm x 10cm 3's | 100 |
| DEF207 | 10cm x 10cm 5's | 100 |
| DEF225 | 10cm x 10cm 20's | 50 |
| DEF1178 | 10cm x 20cm 10's | 50 |

BITE BLOCK

| | | |
|---------|----------------------------------|-----|
| DEF1828 | BITE BLOCK 7.5cm X-RAY WITH TAIL | 100 |
|---------|----------------------------------|-----|



| | | |
|---------|---|-----|
| DEF1246 | BITE BLOCK 7.5cm X-RAY WITH TAIL NON STERILE | 100 |
|---------|---|-----|

SILTAPPE

Soft silicone perforated tape roll used for taping eyelids in theatre



| | | |
|--------|------------|---|
| CR3938 | 2cm x 3m | 1 |
| CR3939 | 4cm x 1.5m | 1 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

DRUG TRAY

| | | |
|--------|---|----|
| DEF321 | DRUG TRAY - DRESSING TRAY WITH 3 COMPARTMENTS 23cm x 9cm x 3cm Deep | 80 |
|--------|---|----|

**ANAESTHETIC/REGIONAL/EPIDURAL PACKS**

| | | |
|---------|-------------------------|----|
| DEF1524 | ANAESTHETIC TRAY | 80 |
|---------|-------------------------|----|

2 x Green Swabs 10cm x 10cm
1 x Dressing Towel 40cm x 42cm
1 x Long Dressing Tray 23cm x 9cm x 3cm

| | | |
|--------|-------------------------|----|
| DEF406 | ANAESTHETIC PACK | 60 |
|--------|-------------------------|----|

1 x Long Dressing Tray
23cm x 9cm x 3cm Deep
5 x Green Swabs 10cm x 10cm
1 x Polybacked Towel 40cm x 42cm

| | | |
|--------|----------------------|----|
| DEF811 | REGIONAL TRAY | 60 |
|--------|----------------------|----|

1 x Epidural Clear Plastic Drape
76cm x 106cm with Full 2.5cm Adhesive
around 10cm x 15cm Fenestration and
Adhesive on Corner tabs
1 x Long Dressing Tray
23cm x 9cm x 3cm Deep
6 x Green Gauze Swabs 10cm x 10cm
1 x Green Disposable Forcep
Field Blue Plastic 75cm x 75cm



| Code | Description | Box Qty. |
|---------|---|----------|
| DEF1796 | EPIDURAL PACK | 50 |
| | 1 x Epidural Clear Plastic Drape 76cm x 106cm with Full 2.5cm Adhesive around 10cm x 15cm Fenestration and Adhesive on Corner tabs | |
| | 2 x Long Dressing Tray 23cm x 9cm x 3cm Deep | |
| | 5 x Green Swabs 10cm x 10cm | |
| | 1 x Dressing Towel - T Large 32cm x 60cm | |
| | 2 x Green Plastic Disposable Forceps | |
| | Field Blue Plastic 91cm x 115cm | |
| DEF1479 | REGIONAL TRAY WITH LARGE DRAPE | 60 |
| | 1 x Epidural Clear Drape 120cm x 140cm with 2.5cm Full Adhesive around 10cm x 15cm Fenestration and Adhesive on Corner Tabs | |
| | 1 x Long Dressing Tray 23cm x 9cm x 3cm Deep | |
| | 4 x Green Gauze Swabs 10cm x 10cm | |
| | 2 x Forceps Plastic Green | |
| | 1 x Towel White 2 Ply | |
| | Field Blue Plastic 91cm x 115cm | |
| DEF1618 | ANAESTHETIC PACK | 60 |
| | 1 x Blue Plastic Drape 75cm x 75cm with 10cm x 15cm Fenestration with Adhesive tabs on Corners of Fenestration and Corners of Drape | |
| | 1 x Long Dressing Tray 23cm x 9cm x 3cm Deep | |
| | 5 x Green Swabs | |
| | 1 x Green Disposable Forcep | |
| | 1 x 5ml L/S Syringe | |
| | 1 x 3ml L/S Syringe | |
| | 1 x Dressing Towel - T Small 30cm x 32cm | |
| | 1 x 19g x 11/2" Needle | |
| | 1 x 18g x 11/2" Blunt Drawing up Needle | |
| | 1 x 25g x 1" Needle | |
| | Field Blue Plastic 75cm x 75cm | |



| Code | Description | Box Qty. |
|----------|---|----------|
| DEF2255A | ANAESTHETIC PROCEDURE PACK | 60 |
| | 1 x Epidural Clear Plastic Drape 76cm x 106cm with Full 2.5cm Adhesive around 10cm x 15cm Fenestration and Adhesive on Corner tabs | |
| | 1 x Long Dressing Tray 23cm x 9cm x 3cm Deep | |
| | 1 x Syringe 3ml L/Lock | |
| | 1 x Syringe 5ml L/Lock | |
| | 2 x Plastic Pearman Forceps | |
| | 1 x Hypodermic Needle 18g x 1.5" Blunt | |
| | 1 x Hypodermic Needle 25g x 1.5" | |
| | 2 x Non Woven Balls | |
| | 5 x Gauze Sponges Plain 10cm x 10cm | |
| | 1 x Gallipot 60ml | |
| | 1 x Particle Filter 5 Micron | |
| | 1 x Plastic Drape 91cm x 115cm | |
| | Field Blue Plastic 91cm x 115cm | |
| DEF1480 | COMPLETE EPIDURAL PACK | 60 |
| | 1 x Epidural Clear Plastic Drape 76cm x 106cm with Full 2.5cm Adhesive around 10cm x 15cm Fenestration and Adhesive on Corner tabs | |
| | 1 x Long Dressing Tray 23cm x 9cm x 3cm Deep | |
| | 5 x Green Sponges | |
| | 1 x Green Disposable Forcep | |
| | 1 x 5ml L/Slip Syringe | |
| | 1 x 3ml L/Slip Syringe | |
| | 1 x 19g x 1.5" Hypodermic Needle | |
| | 1 x 25g x 1" Hypodermic Needle | |
| | 1 x 18g x 1.5" Hypodermic Needle Blunt | |
| | 1 x Dressing Towel - T Small 32cm x 32.5cm | |
| | Field Blue Plastic 75cm x 75cm | |



| Code | Description | Box Qty. |
|----------|---|----------|
| DEF1794A | EPIDURAL PACK | 15 |
| | 1 x Epidural Clear Plastic Drape 76cm x 106cm with Full 2.5cm Adhesive around 10cm x 15cm Fenestration and Adhesive on Corner tabs | |
| | 1 x Long Dressing Tray 23cm x 9cm x 3cm Deep | |
| | 3 x Foam Prep Stick Blue 6" (15cm) | |
| | 5 x Green Swabs 10cm x 10cm | |
| | 1 x Syringe 10ml L/Slip | |
| | 1 x Syringe 5ml L/Slip | |
| | 1 x Syringe 3ml L/Slip | |
| | 1 x Hypodermic Needle 18g x 1.5" Blunt | |
| | 1 x Hypodermic Needle 23g x 1.25" | |
| | 1 x Gallipot 60ml | |
| | 1 x Huck Towel 40cm x 68cm | |
| | 1 x Surgical Gown X/Large Unreinforced | |
| | Field SMS 91cm x 91cm | |

| | | |
|----------|--|----|
| DEF1988B | EPIDURAL PACK | 50 |
| | 1 x Epidural Clear Drape 120cm x 140cm with 2.5cm Full Adhesive around 10cm x 15cm Fenestration and Adhesive on Corner Tabs | |
| | 2 x Long Dressing Tray 23cm x 9cm x 3cm Deep | |
| | 5 x Green Swabs 10cm x 10cm | |
| | 1 x Syringe 5ml L/Slip | |
| | 1 x Syringe 3ml L/Slip | |
| | 2 x Hypodermic Needle 18g x 1.5" Blunt | |
| | 1 x Hypodermic Needle 25g x 1.5" | |
| | 1 x Hypodermic Needle 25g x 1" | |
| | 1 x Dressing Towel - T Small 32 x 32.5cm | |
| | 3 x Absorbent Foam Applicator 20cm | |
| | Field Blue Plastic 75cm x 75cm | |



| Code | Description | Box Qty. |
|---------|---|----------|
| DEF1823 | REGIONAL PACK - DOUBLE PACKED | 60 |
| | 1 x Epidural Clear Plastic Drape 76cm x 106cm with Full 2.5cm Adhesive around 10cm x 15cm Fenestration and Adhesive on Corner tabs | |
| | 2 x Long Dressing Tray 23cm x 9cm x 3cm Deep | |
| | 2 x Green Disposable Plastic Forceps | |
| | 3 x Plain Green Sponges | |
| | 1 x White Towel 2 Ply | |
| | 1 x Hypodermic Needle 19g x 1.5" | |
| | 1 x Hypodermic Needle 25g x 1" | |
| | 1 x Hypodermic Needle 18g x 1.5" Blunt | |
| | 1 x Syringe 5ml L/Slip | |
| | 1 x Syringe 3ml L/Slip | |
| | Field Blue Plastic 75cm x 75cm | |

| | | |
|----------|--|----|
| DEF2345A | EPIDURAL PACK/SPINAL PACK | 24 |
| | 1 x Epidural Clear Plastic Drape 120cm x 140cm with Full 2.5cm Adhesive around 10cm x 15cm Fenestration and Adhesive on Corner tabs | |
| | 1 x Hypodermic Needle 18g x 1.5" Blunt | |
| | 1 x Hypodermic Needle 25g x 1.5" | |
| | 1 x Particle Filter 5 Micron | |
| | 5 x Green Swabs 10cm x 10cm | |
| | 1 x Gallipot 60ml | |
| | 1 x Foam Prep Stick Blue 6" (15cm) | |
| | 2 x Dressing Towel - T Small 32 x 32.5cm | |
| | 1 x Dressing Tray - 3 Cavity 23cm x 9cm x 3cm Deep | |
| | 1 x Syringe 3ml L/Slip | |
| | 1 x Syringe 5ml L/Lock | |
| | 1 x Syringe 10ml L/Lock | |
| | Field Blue Plastic 91cm x 115cm | |



| Code | Description | Box Qty. |
|---------|---|----------|
| DEF3085 | EPIDURAL PACK/SPINAL PACK | 24 |
| | 1 x Dressing Tray with 3 compartments 23cm x 9cm x 3cm | |
| | 5 x Green Gauze Swabs 10cm x 10cm | |
| | 1 x Syringe 5ml Luer Slip | |
| | 1 x Syringe 3ml Luer Slip | |
| | 1 x Hypodermic Needle 25g x 1.5" | |
| | 1 x Hypodermic Needle 18g x 1.5" | |
| | 1 x Hypodermic Needle 25g x 1" | |
| | 1 x Dressing Towel Small T 32cm x 32cm | |
| | 1 x Epidural Drape 120cm x 140cm with 2.5cm full adhesive around 10cm x 15cm Fenestration and adhesive on corner tabs | |



To order, please contact our
Customer Service Team on


1300 550 278

or see our website for details

www.defries.com.au



> NEURO



| Code | Description | Box Qty. |
|--------------------|-----------------|----------|
| COTTON WOOL | | |
| DEF859 | 10cm x 10cm | 25 |
| DEF615 | 10cm x 15cm D/P | 100 |

FOLDED GAUZE

| | | |
|--------|----------------------|----|
| DEF347 | 7cm x 21cm 5's X-RAY | 50 |
|--------|----------------------|----|



| | | |
|---------|-----------------------|----|
| DEF1074 | 10cm x 20cm X-RAY 5's | 25 |
|---------|-----------------------|----|


| | | |
|--------|------------------------------|----|
| DEF611 | HEAD DRESSINGS D/P | 35 |
| | 8 x Gauze Swabs 10cm x 20cm | |
| | 20 x Gauze Swabs 10cm x 10cm | |

| | | |
|---------|-----------------------------|---|
| DEF2789 | NEURO BOWL SET | 8 |
| | 1 x Jug Tall 1 Litre | |
| | 2 x Bowl 500ml | |
| | 1 x Kidney Dish Clear 700ml | |
| | 1 x Kidney Dish Green 700ml | |

| | | |
|----------|------------------------------------|---|
| DEF2490A | NEURO RECEIVER SET | 8 |
| | 8 x Gallipots 150ml | |
| | 1 x Jug 1.2 Litre | |
| | 1 x Label "Has Time out been Done" | |
| | 1 x Plastic Bag | |

**NEURO SKIN HOOKS**

| | | |
|--------|----------------------------|----|
| DEF405 | 5 ON A TRAY Simmons Needle | 20 |
| DEF336 | 5 ON A TRAY D/P Small Hook | 20 |



| | | |
|---------|----------------------------------|----|
| DEF1032 | SIZE: 3/0 5 ON A TRAY Large Hook | 20 |
|---------|----------------------------------|----|

NEURO 1/2 (53CM) YARD STRIPS

| | | |
|--------|-----------|----|
| DEF858 | X-RAY 5'S | 50 |
|--------|-----------|----|

Code

Description

Box Qty.

NEURO SUCTION TUBING SILICONE
WITH MEDIUM CONNECTOR

| | | |
|---------|------------------|----|
| DEF1186 | 5mm x 9mm 30cm | 50 |
| DEF306 | 5mm x 9mm x 50cm | 50 |



| | | |
|---------|------------------|----|
| DEF1337 | 5mm x 9mm x 1m | 50 |
| DEF1712 | 5mm x 9mm x 60cm | 40 |

Code

Description

Box Qty.



To order, please contact our
Customer Service Team on
1300 550 278
or see our website for details
www.defries.com.au

> CATHETER PACKS



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

CATHETER PACKS XRAY DETECTABLE

| Code | Description | Box Qty. |
|---------|--|----------|
| DEF1601 | CATHETER PACK | 60 |
| | 1 x Drape Polybacked 50cm x 60cm with 10cm Circular Fenestration | |
| | 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm Deep | |
| | 5 x Gauze Swabs X-ray 10cm x 10cm | |
| | 1 x Syringe 10ml Luer Slip | |
| | 2 x Green Disposable Forcep | |
| | 1 x Blue Huck towel 40cm x 68cm | |
| | Field Blue Plastic 75cm x 75cm | |



| Code | Description | Box Qty. |
|----------|---|----------|
| DEF1999A | CATHETER PACK - THEATRE | 24 |
| | 1 x Drape SMS 91cm x 91cm with 10cm x 15cm Fenestration | |
| | 5 x Gauze Swabs 10cm x 10cm X-Ray | |
| | 1 x Clear Kidney Dish 700ml | |
| | 1 x Blue Huck towel 40cm x 68cm | |
| | 1 x 24cm Plastic Sponge Holder | |
| | 2 x Gallipot 150ml | |
| | 1 x 10ml Luer Slip Syringe | |
| | Field Blue Plastic 75cm x 75cm | |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| Code | Description | Box Qty. |
|---------|---|----------|
| DEF2074 | CATHETER PACK | 24 |
| | 1 x Polybacked Towel 75cm x 105cm with 15cm x 15cm Fenestration | |
| | 1 x Mediwipe Towel 35cm x 60cm | |
| | 2 x Plastic Pearman Forceps | |
| | 1 x 10ml Luer Slip Syringe | |
| | 5 x Gauze Swabs 10cm x 10cm X-Ray | |
| | 2 x Clear Kidney Dishes 700ml | |
| | 1 x Gallipot 150ml | |
| | Field Plastic Blue 91cm x 115cm | |



| Code | Description | Box Qty. |
|----------|--|----------|
| DEF2061A | CATHETER PACK | 24 |
| | 1 x Polybacked Towel 75 x 105cm with 5cm x 15cm Fenestration | |
| | 1 x 10ml Syringe L/Slip | |
| | 2 x Gallipot 150ml | |
| | 1 x Plastic Pearman Forcep | |
| | 1 x Clear Kidney Dish 700ml | |
| | 1 x Green Kidney Dish 700ml | |
| | 5 x Gauze Swabs 10 x 10cm X-Ray | |
| | 1 x Dressing Towel - T Small 32cm x 32.5cm | |
| | 1 x 19cm Plastic Sponge Holder | |
| | Field Blue Plastic 91cm x 115cm | |



| Code | Description | Box Qty. |
|---------|--|----------|
| DEF1262 | CATHETER PACK | 30 |
| | 1 x Drape 60cm x 60cm with 10cm x 10cm Fenestration | |
| | 1 x Dressing Tray - with 3 Compartments 23cm x 9cm x 3cm | |
| | 1 x Gallipot 60ml | |
| | 1 x Kidney Dish Clear 700ml | |
| | 1 x Plastic Forcep Pearman | |
| | 1 x Syringe with Leur Slip (10ml) | |
| | 1 x Dressing Towel - T Small 32cm x 32.5cm | |
| | 1 x Tube Clamp with Lock | |
| | 5 x Gauze Sponges 10cm x 10cm X-ray Detectable | |
| | Field Blue Plastic 75cm x 75cm | |

| | | |
|----------|--|----|
| DEF2649A | CATHETER PACK | 18 |
| | 1 x Polybacked Towel 75cm x 105cm with 5cm x 15cm Fenestration | |
| | 2 x Clear Kidney Dishes 700ml | |
| | 5 x Gauze Swabs 10cm x 10cm X-Ray | |
| | 1 x 10ml Luer Slip Syringe | |
| | 1 x 24cm Plastic Sponge Holder | |
| | 2 x Gallipots 150ml | |
| | Field SMS 130cm x 150cm | |

CATHETER PACKS WITH PLAIN GAUZE

| | | |
|---------|--|----|
| DEF1522 | CATHETER PACK | 60 |
| | 3 x Poly Backed Towel 40cm x 42cm | |
| | 2 x Dressing Tray - with 3 Compartments 23cm x 9cm x 3cm | |
| | 5 x Gauze Sponges Plain 10cm x 10cm | |
| | 2 x Green Disposable Forceps | |
| | 1 x 10ml Leur Slip Syringe | |
| | Field Blue Plastic 75cm x 75cm | |



| | | |
|----------|---|----|
| DEF1650A | CATHETER PACK | 60 |
| | 1 x Minor Procedure Drape 50cm x 60cm with 6cm Circular Fenestration with Full Adhesive | |
| | 2 x Dressing Tray - with 3 Compartments 23cm x 9cm x 3cm | |
| | 2 x Green Disposable Plastic Forceps | |
| | 1 x 2 Ply White Towel 40cm x 42cm | |
| | 1 x 10ml Luer Slip Syringe | |
| | Field Blue Plastic 75cm x 75cm | |

| | | |
|----------|--|----|
| DEF2253B | CATHETER PACK | 60 |
| | 1 x Drape Polybacked 50cm x 60cm with 10cm Circular Fenestration | |
| | 1 x Dressing Tray - with 3 Compartments 23cm x 9cm x 3cm | |
| | 1 x Dressing Towel - T Small 32cm x 32cm | |
| | 2 x Green Disposable Forceps | |
| | 2 x Nitrile Gloves - Large | |
| | 2 x 10ml Luer Slip Syringe | |
| | 1 x Lube Jelly sachet 2.7g | |
| | 5 x Cotton Balls - Large | |
| | Field Blue Plastic 75cm x 75cm | |

| | | |
|---------|---|----|
| DEF1922 | CATHETER PACK | 50 |
| | 1 x SMS Drape 91cm x 91cm with 10cm Circular Fenestration | |
| | 2 x Dressing Tray - with 3 Compartments 23cm x 9cm x 3cm | |
| | 1 x 19cm Plastic Sponge Holder | |
| | 1 x Blue Cotton Huck Towel 40cm x 68cm | |
| | Field Blue Plastic 91cm x 115cm | |

| | | |
|---------|-----------------------------------|----|
| DEF1881 | CATHETER PACK | 30 |
| | 2 x Kidney Dish 700ml Clear | |
| | 1 x Gallipot 150ml | |
| | 7 x Non Woven Balls | |
| | 2 x Non Woven Swabs 7.5cm x 7.5cm | |
| | 2 x Plastic Pearman Forceps | |
| | 3 x Blue Huck Towels 40cm x 68cm | |
| | Field Blue Plastic 75cm x 75cm | |

| | | |
|---------|---|----|
| DEF1036 | CATHETER PACK | 30 |
| | 1 x SMS Drape 91cm x 91cm with 10cm Circular Fenestration | |
| | 2 x Clear Kidney Dish 700ml | |
| | 2 x Gallipot 150ml | |
| | 2 x Green Disposable Forceps | |
| | 1 x Dressing Towel - T Large 32cm x 60cm | |
| | 6 x Cotton Balls - Large | |
| | Field Blue Plastic 75cm x 75cm | |

| Code | Description | Box Qty. |
|----------|--|----------|
| DEF1038 | CATHETER PACK | 30 |
| | 1 x SMS Drape 91cm x 91cm with 10cm Circular Fenestration | |
| | 1 x Clear Kidney Dish 700ml | |
| | 1 x Gallipot 150ml | |
| | 1 x Tube Clamp with Lock | |
| | Field Blue Plastic 75cm x 75cm | |
| DEF2465 | CATHETER PACK | 24 |
| | 1 x SMS Drape 91cm x 91cm with 10cm Circular Fenestration | |
| | 1 x Clear Kidney Dish 700ml | |
| | 1 x Dressing Tray - with 3 Compartments 23cm x 9cm x 3cm | |
| | 5 x Gauze Swabs Plain 10cm x 10cm | |
| | 2 x Green Disposable Forceps | |
| | 1 x Gallipot 60ml | |
| | 1 x 10ml Luer Slip Syringe | |
| | 2 x Dressing Towel - T Small 32cm x 32cm | |
| | Field Blue Plastic 91cm x 115cm | |
| DEF1335 | CATHETER PACK | 50 |
| | 1 x SMS Drape 91cm x 91cm with 10cm Circular Fenestration | |
| | 2 x Dressing Tray - with 3 Compartments 23cm x 9cm x 3cm | |
| | 5 x Cotton Balls - Large | |
| | 2 x Green Disposable Forceps | |
| | 1 x Dressing Towel - T Small 32cm x 32cm | |
| | Field Plastic Drape 91cm x 115cm | |
| DEF1333T | CATHETER PACK | 60 |
| | 1 x SMS Drape 91cm x 91cm with 10cm Circular Fenestration | |
| | 1 x Dressing Tray - with 3 Compartments 23cm x 9cm x 3cm | |
| | 1 x 10ml Luer Slip Syringe | |
| | 2 x Green Disposable Forceps | |
| | 1 x Dressing Towel - T Small 32cm x 32cm | |
| | Field Blue Plastic 75cm x 75cm | |

| Code | Description | Box Qty. |
|----------|---|----------|
| DEF3009C | CATHETER PACK | 10 |
| | 7 x Non Woven Swabs 5cm X 5cm | |
| | 1 x Forcep Green Plastic Disposable | |
| | 1 x Forcep Crile Straight 14cm with Green Handle | |
| | 1 x Forcep 13cm Blunt with Green Handle | |
| | 1 x Gallipot 150ml | |
| | 2 x Gallipot 60ml | |
| | 2 x Kidney Dish 1000ml Blue | |
| | 1 x SMS Drape 130CM x 150CM with 14cm slit fenestration | |
| | 1 x SMS 91CM x 91CM | |
| | 2 x Towel Clamp U Shaped White | |
| | 2 x Blue Huck Towel 40cm x 68cm | |
| | 1 x Field SMS 130cm x 150cm | |
| DEF1478 | PROCEDURE PACK | 50 |
| | 2 x Green Gauze Swabs 10cm x 10cm | |
| | 1 x Dressing Tray with 3 compartments 23cm x 9cm x 3cm | |
| | 2 x Green Disposable Plastic Forceps | |
| | 4 x Cotton Balls Large | |
| | 1 x Gallipot 150ml | |
| | 1 x Dressing Towel Small T 32cm x 32cm | |
| | 1 x Field Blue Plastic 91cm x 115cm | |
| DEF2202 | DELIVERY SUITE CATHETER PACK | 50 |
| | 10 x Cotton Balls Large | |
| | 1 x Dressing Tray with 3 compartments 23cm x 9cm x 3cm | |
| | 1 x Green Disposable Plastic Forceps | |
| | 1 x Gallipot 150ml | |
| | 1 x Field Blue Plastic 50cm x 50cm | |
| DEF3920 | CATHETER PACK | 30 |
| | 5 x Dextex Gauze Swab 10cm x 10cm X-ray | |
| | 1 x Dressing Tray with 3 compartments 23cm x 9cm x 3cm | |
| | 1 x Green Disposable Plastic Forceps | |
| | 1 x Kidney Dish 700ml Clear | |
| | 1 x Sponge Holder 24cm Plastic | |
| | 1 x Syringe 10ml Luer Slip Concentric Tip | |
| | 1 x Blue Huck Towel 40cm x 68cm | |
| | 1 x PolyBacked Drape 50cm x 60cm with 10cm circular Fenestration | |
| | 1 x Field Blue Plastic 75cm x 75cm | |

| Code | Description | Box Qty. |
|----------|---|----------|
| DEF1650A | CATHETER PACK | 60 |
| | 2 x Dressing Tray with 3 compartments 23cm x 9cm x 3cm | |
| | 2 x Green Disposable Plastic Forceps | |
| | 1 x Azure Minor Procedure Drape 50cm x 60cm with 6.5cm circular Fenestration and Adhesive around fenestration | |
| | 1 x Syringe 10ml Luer Slip Concentric Tip | |
| | 1 x Dressing Towel white 40cm x 42cm | |
| | 1 x Field Blue Plastic 75cm x 75cm | |
| DEF1932 | THEATRE CATHETER PACK | 30 |
| | 2 x Gallipot 150ml | |
| | 1 x Kidney Dish 700ml Clear | |
| | 5 x Gauze Swab 10cm x 10cm X-ray | |
| | 1 x Sponge Holder 24cm Plastic | |
| | 1 x Syringe 30ml Eccentric Tip | |
| | 1 x Towel D 32cm x 34cm | |
| | 1 x PolyBacked Drape 50cm x 60cm with 10cm circular Fenestration | |
| | 1 x Field Blue Plastic 91cm x 115cm | |
| DEF2649A | CATHETER PACK | 30 |
| | 1 x Polybacked Towel 75cm x 105cm with 5cm x 15cm Fenestration | |
| | 1 x Clear Kidney Dishes 700ml | |
| | 5 x Gauze Swabs 10cm x 10cm X-Ray | |



To order, please contact our
Customer Service Team on

1300 550 278

or see our website for details
www.defries.com.au

> THEATRE CONSUMABLES



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

BLADESAFE SCALPEL TRANSFER DEVICE

- Ideal for the safe transfer of scalpel blade & handle between scrub nurse and surgeon.

| | | |
|---------|--|----|
| DEF2631 | BLADESAFE SCALPEL TRANSFER DEVICE | 24 |
| | Scalpel Blade and Handle not included. | |

**BURNS PRODUCTS**

- Non adherant, highly absorbant wound dressing.



| | | |
|--------|---------------------|----|
| DEF987 | 50cm x 40cm MELOLIN | 50 |
| DEF679 | 50cm x 70cm D/P | 35 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

BURNS DRESSINGS GAUZE

- Soft roll suitable for burns and wounds.

| | | |
|---------|---------------------------------|-----|
| DEF342 | SMALL 90cm x 2m 1's | 50 |
| DEF344 | LARGE 90cm x 2m 4's | 25 |
| DEF1662 | 90cm x 5m | 10 |
| DEF650 | BURNS PLASTIC BAG 390mm x 305mm | 100 |

DIATHERMY PENCIL

- Comfortable, ergonomic design,
- Full range of accessories, including Teflon coated flat blade.

| | | |
|---------|------------------------------|----|
| DEF2863 | ROCKER SWITCH AND TEFLON TIP | 40 |
|---------|------------------------------|----|



| | | |
|---------|---------------------------------------|----|
| DEF2864 | ROCKER SWITCH AND STAINLESS STEEL TIP | 40 |
| DEF2865 | FINE TOUCH AND TEFLON TIP | 40 |
| DEF2866 | FINE TOUCH AND STAINLESS STEEL TIP | 40 |

| | | |
|---------|-----------------------|----|
| DEF2863 | DIATHERMY PACK | 10 |
|---------|-----------------------|----|

1 x Diathermy Pencil with Rocker Switch & Teflon Blade

1 x Disposable Quiver 5cm x 25cm with metal ring

| | | |
|---------|-----------------------|----|
| DEF2864 | DIATHERMY PACK | 10 |
|---------|-----------------------|----|

1 x Diathermy Pencil with Rocker Switch & Stainless Steel Blade

1 x Diathermy Scratch Pad Tip Cleaner

1 x Quiver 5cm x 25cm with metal ring

| | | |
|---------|-----------------------|----|
| DEF2865 | DIATHERMY PACK | 10 |
|---------|-----------------------|----|

1 x Diathermy Pencil with Fine Touch & Teflon Blade

1 x Quiver 5cm x 25cm with metal ring

Code

Description

Box Qty.

DIATHERMY QUIVER

- Single use Polyethylene Quiver.

DEF2772

DISPOSABLE QUIVER 5cm x 25cm LENGTH WITH METAL RING

90



DEF3068

DISPOSABLE QUIVER 5cm x 25cm WITH CABLE TIE

90

DEF3291

LAP INSTRUMENT QUIVER 9cm x 35cm WITH METAL RING

40

DEF3296

LAP INSTRUMENT QUIVER 9cm x 35cm WITH CABLE TIE

40

DRESSING NON WOVEN COMBINE

- Combine dressings can be used either in direct contact for medium to heavy exuding wounds or as a secondary dressing.



DEF711

9cm x 20cm

25

DEF719

9cm x 20cm 5's

100

DEF720

9cm x 20cm 10's

50

DEF715

20cm x 20cm

25

DEF80

20cm x 60cm

100

DEF1160

20cm x 80cm

50

DRESSING NON WOVEN COMBINE ROLL

DEF714

9cm x 10m D/P

16

DEF1004

20cm x 2.5m

25

Code

Description

Box Qty.

DRESSING NON WOVEN AMPUTATION

DEF804

NON WOVEN 70cm x 40cm

50



DEF2289

MELOLIN T-SHAPED 40cm x 20cm

35

FOAM PADS - POLYURETHANE LATEX FREE

DEF792

100mm x 100mm x 13mm THICK

100

DEF769

100mm x 100mm x 25mm THICK

50



DEF783

100mm x 180mm x 13mm THICK

50

DEF784

100mm x 180mm x 25mm THICK

40

DEF738

150mm x 150mm x 13mm THICK

50

DEF914

PLASTICS FOAM 150mm x 200mm x 25mm

25

I.V. PRESSURE PADS - LATEX FREE

- Unique application - Allows you to apply the pressure pad without contaminating the pad or adhesive. Hypoallergenic adhesive. Adequate compression by 3mm thick pad.

DEF929

36mm WHITE

100



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

INSTRUMENT POUCH

- Plastic pouch used to organise instruments in the sterile field.
- Adhesive strips on side of pouch to hold into place on drape.

| | | |
|---------|----------|----|
| DEF2406 | 2 POCKET | 25 |
|---------|----------|----|

**LIGHT HANDLE COVER**

- Gives control of light direction to surgeon and / or scrub personnel.
- Ability to re-cover if sterility is compromised.
- Compatible Handle Adaptors.

| | | |
|---------|----------------------|----|
| DEF2082 | GREEN - FRICTION FIT | 50 |
|---------|----------------------|----|



| | | |
|---------|---------------------|----|
| DEF2919 | GREEN 2's | 40 |
| DEF2557 | FOAM BLUE | 50 |
| DEF2675 | BLUE - CLICK ON | 50 |
| DEF3170 | BLUE - CLICK ON 2'S | 40 |



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

FORCEP PLASTIC

| | | |
|--------|-------|-----|
| DEF686 | GREEN | 100 |
|--------|-------|-----|

**MARKING PENS**

- For smooth clean line marking on patient's skin prior to surgical procedure. Non Toxic, Non tattooing, quick drying.

| | | |
|--------|-------------------------------|----|
| 26-002 | FINE POINT WITH 15cm/5" RULER | 50 |
|--------|-------------------------------|----|



| | | |
|--------|-----------------------------|----|
| 26-001 | STANDARD WITH 15cm/5" RULER | 50 |
|--------|-----------------------------|----|



| | | |
|--------|--------------------------|----|
| DEF464 | PREPACKED MARKING PENCIL | 24 |
|--------|--------------------------|----|

non sterile

| | | |
|---------|-------------------------------|-----|
| DEF2183 | SKIN MARKERS STANDARD TIP | 50 |
| DEF2669 | MINI PRE-OP SKIN MARKERS BULK | 100 |

Code Description Box Qty.

NEEDLE CASE

- Easy to read numbers & description facilitate correct count of blades & needles. Red colour signifies Biohazard waste.

25-1701 NEEDLE COUNTER RED, 15 COUNT WITH 50
BLADE REMOVER FOAM STRIPS/MAG



25-1702 NEEDLE COUNTER RED, 30 COUNT WITH 27
BLADE REMOVER FOAM STRIPS/MAG



25-1401 NEEDLE COUNTER RED, 20 COUNT WITH 20
BLADE REMOVER FOAM BLOCK/MAG



Code Description Box Qty.

25-1401Y NEEDLE COUNTER YELLOW, 20
20 COUNT WITH BLADE REMOVER
FOAM BLOCK/MAG



25-0702 NEEDLE COUNTER RED, 30 COUNT 27
FOAM STRIPS/MAG



25-1108 NEEDLE COUNTER RED, 60 COUNT WITH 27
BLADE REMOVER DOUBLE MAGNETIC



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|---------|---------------------------|-----|
| 25-0900 | SIMPLICITY NEEDLE COUNTER | 108 |
|---------|---------------------------|-----|




| | | |
|---------|-------------------------------------|-----|
| 25-0901 | SIMPLICITY NEEDLE COUNTER WITH FOAM | 108 |
|---------|-------------------------------------|-----|

**PREP APPLICATORS**

- All in one applicator for preping.

| | | |
|---------|----------------------------------|-----|
| DEF921 | 50mm x 60mm x 25mm | 100 |
| DEF1847 | WITH HANDLE 150mm BLUE PACK OF 2 | 100 |



| | | |
|---------|-------------------------|-----|
| DEF2035 | WITH HANDLE 20cm ORANGE | 100 |
|---------|-------------------------|-----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

PREP SWAB BALLS X/Ray

| | | |
|--------|-----------|-----|
| DEF255 | LARGE 1's | 100 |
| DEF251 | LARGE 2's | 100 |



| | | |
|--------|---------------|-----|
| DEF254 | LARGE 3's | 100 |
| DEF264 | LARGE D/P 5's | 50 |

RUBBER BANDS

- Latex bands.

| | | |
|---------|----------|-----|
| DEF1590 | #16 2's | 100 |
| DEF340 | #28 5's | 100 |
| DEF1405 | #31 2's | 100 |
| DEF625 | #31 10's | 100 |



| | | |
|--------|-------------|-----|
| DEF669 | #33 5'S D/P | 25 |
| DEF619 | #34 5's | 25 |
| DEF600 | #34 10's | 25 |
| DEF624 | #64 | 100 |
| DEF626 | #64 5's D/P | 100 |

SAFETY PINS

| | | |
|----------|-----------|----|
| DEF071SP | #2 (38mm) | 50 |
| DEF072SP | #3 (50mm) | 50 |
| DEF073SP | #4 (57mm) | 50 |

**non sterile**

| | | |
|--------|-----------|-----|
| DEF635 | #3 (50mm) | 144 |
| DEF636 | #4 (57mm) | 144 |

Code

Description

Box Qty.

SPONGES ABDOMINAL X-RAY

- High quality 6 Ply Packs.

| | | |
|---------|-----------------|-----|
| DEF1350 | 23cm x 23cm 5's | 100 |
| DEF1351 | 30cm x 30cm 5's | 75 |



| | | |
|---------|-----------------|----|
| DEF1352 | 30cm x 48cm 5's | 60 |
|---------|-----------------|----|

ABDOMINAL SPONGES WITH TAPE

- High quality 6 Ply Packs.

| | | |
|--------|---------------------------|-----|
| DEF162 | 23cm x 23cm 1's WITH TAPE | 100 |
| DEF160 | 23cm x 23cm 3's WITH TAPE | 100 |
| DEF168 | 23cm x 23cm 5's WITH TAPE | 100 |



SPONGES PEANUT

- Small, firm, quality gauze sponges used in delicate surgical procedures.

| | | |
|--------|--|----|
| DEF316 | INSTAPACK LARGE 5's IN AN ADHESIVE FOAM CASE FOR EASE OF COUNT | 50 |
|--------|--|----|



Code

Description

Box Qty.

| | | |
|--------|--|----|
| DEF315 | INSTAPACK (SMALL) 5's IN AN ADHESIVE FOAM CASE FOR EASE OF COUNT | 50 |
| DEF046 | 5 ON A SAFETY PIN X-RAY | 50 |



SPONGES GAUZE X-RAY

- 8 Ply xray swabs with 24 x 20 thread count.


| | | |
|--------|-----------------------|-----|
| DEF202 | 7.5cm x 7.5cm 5's D/P | 100 |
| DEF044 | 10cm x 10cm 5's D/P | 100 |



| | | |
|---------|---------------------|-----|
| DEF165 | 10cm x 10cm 5's D/P | 500 |
| DEF1689 | 10cm x 10cm 10's | 250 |

SPONGES GAUZE - SPECIALIZED

| | | |
|--------|--|-----|
| DEF986 | DONOR SITE GAUZE DRESSING 5cm x 10cm D/P | 100 |
| DEF169 | RECOVERY GAUZE 5cm x 5cm PLAIN 32Ply | 100 |



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SILICONE TOURNIQUET - STERILE

- Soft stretchy neoprene material. Latex free.

| | | |
|---------|----------------------|----|
| DEF1006 | 6.5cm x 2m BLUE D/P | 10 |
| DEF2764 | 10cm x 2.7m BLUE D/P | 20 |




| | | |
|---------|----------------------|----|
| DEF1081 | 10cm x 3m BLUE | 5 |
| DEF2668 | 15cm x 2.7m BLUE D/P | 20 |

TIP CLEANER X-RAY

- Effective removal of debris from non teflon coated electrodes
- Radio-opaque strip ensures xray detection.

| | | |
|---------|--|----|
| DEF1886 | 5cm x 5cm WITH ADHESIVE AND RADIO-OPAQUE STRIP | 50 |
|---------|--|----|

**PVC TUBING SNUGGERS**

| | | |
|---------|-----------------------|----|
| DEF3051 | 16Fr x 17cm X-RAY 6's | 25 |
| DEF3052 | 16Fr x 22cm X-RAY 6's | 25 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

TUBE CLAMPS

- Ideal for drape clamping management.

| | | |
|---------|-----------------------------|----|
| DEF1575 | WITH LOCK BLUE PLASTIC 12cm | 50 |
| DEF2075 | TOWEL CLAMP, U-SHAPED WHITE | 50 |



| | | |
|---------|---------------------------------|----|
| DEF2190 | TOWEL CLIP GREEN PLASTIC | 25 |
| DEF2488 | TOWEL CLIP STAINLESS STEEL 11CM | 10 |
| DEF2298 | TUBE CLAMPS 13CM | 10 |

TUBE HOLDER

- Appropriate for managing tubing with easy to place and remove velcro.

| | | |
|---------|-----------------|----|
| DEF1889 | VELCRO ADHESIVE | 25 |
|---------|-----------------|----|

**DRAPE CLIP - NON STERILE - MIXED COLOURS**

1



> BANDAGES



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

CREPE BANDAGE**medium quality**

| | | |
|--------|------------------|----|
| DEF265 | 2.5cm x 1.6m S/P | 48 |
| DEF266 | 5cm x 1.6m S/P | 48 |
| DEF267 | 7.5cm x 1.6m S/P | 48 |
| DEF268 | 10cm x 1.6m S/P | 48 |



| | | |
|-------------------------|------------------|----|
| DEF269 | 15cm x 1.6m S/P | 48 |
| hospital quality | | |
| DEF819 | 2.5cm x 1.5m D/P | 48 |
| DEF821 | 7.5cm x 1.5m D/P | 48 |
| DEF820 | 5cm x 1.5m D/P | 48 |
| DEF822 | 10cm x 1.5m D/P | 48 |



| | | |
|--------|-----------------|----|
| DEF823 | 15cm x 1.5m D/P | 36 |
|--------|-----------------|----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

heavy duty quality

| | | |
|--------|------------------|----|
| DEF270 | 5cm x 2.3m D/P | 48 |
| DEF271 | 7.5cm x 2.3m D/P | 48 |
| DEF272 | 10cm x 2.3m D/P | 48 |



| | | |
|--------|-------------|----|
| DEF273 | 15cm x 2.3m | 24 |
|--------|-------------|----|

COHESIVE BANDAGE

| | | |
|---------|-----------------|----|
| DEF1056 | 6cm x 1.5m D/P | 25 |
| DEF257 | 10cm x 1.5m D/P | 24 |

CONFORMING BANDAGE

| | | |
|---------|------------------|----|
| DEF1140 | 2.5cm x 1.5m D/P | 48 |
| DEF1100 | 5cm x 3.5m D/P | 48 |



| | | |
|---------|---------------|----|
| DEF1222 | 10cm x 3m D/P | 18 |
|---------|---------------|----|

H CONFORMING

| | | |
|--------|------------------|----|
| DEF258 | 5cm x 1.5m D/P | 48 |
| DEF259 | 7.5cm x 1.5m D/P | 48 |
| DEF260 | 10cm x 1.5m D/P | 48 |
| DEF261 | 15cm x 1.5m D/P | 48 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

UNDERCAST BANDAGES**natural**

| | | |
|---------|------------------|----|
| DEF830N | 7.5cm x 2.7m D/P | 48 |
| DEF824N | 10cm x 2.7m D/P | 48 |



| | | |
|---------|-----------------|----|
| DEF825N | 15cm x 2.7m D/P | 24 |
|---------|-----------------|----|

S Undercast

| | | |
|--------|------------------|----|
| DEF829 | 5cm x 2.7m D/P | 48 |
| DEF830 | 7.5cm x 2.7m D/P | 48 |
| DEF824 | 10cm x 2.7m D/P | 48 |



| | | |
|--------|-----------------|----|
| DEF825 | 15cm x 2.7m D/P | 24 |
|--------|-----------------|----|

V Undercast

| | | |
|--------|-------------------|----|
| DEF826 | 10cm x 2.75m D/P | 48 |
| DEF827 | 15cm x 2.75m D/P | 24 |
| DEF828 | 7.5cm x 2.75m D/P | 48 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

W Undercast

| | | |
|--------|----------------|----|
| DEF673 | 5cm x 3.6m D/P | 48 |
| DEF674 | 7.5cm x 3.6m | 48 |
| DEF675 | 10cm x 3.6m | 48 |



| | | |
|--------|-------------|----|
| DEF676 | 15cm x 3.6m | 24 |
|--------|-------------|----|

VIC BANDAGES

| | | |
|--------|-----------------|----|
| DEF517 | 10cm x 2.2m D/P | 48 |
|--------|-----------------|----|



| | | |
|--------|--------------------------|----|
| DEF518 | 15cm x 2.2m D/P WITH PIN | 24 |
| DEF516 | 7.5cm x 2.2m D/P | 48 |

> PROTECTIVE COVERS AND BAGS



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

CAMERA AND CORD COVERS

| | | |
|---------|--------------|----|
| DEF1164 | 15cm x 230cm | 50 |
|---------|--------------|----|



| | | |
|---------|----------------------------------|----|
| DEF1201 | D/P (CUFF) 15cm X 2.3m S/E & T/P | 50 |
|---------|----------------------------------|----|

CHAIR COVER

| | | |
|--------|--------------------------------------|-----|
| DEF772 | 575mm x 620mm & #33 RUBBER BANDS 2's | 100 |
|--------|--------------------------------------|-----|

CLIP SEAL BAG


| | | |
|---------|---------------|-----|
| DEF1703 | 280mm x 380mm | 25 |
| DEF640 | 200mm x 250mm | 100 |

DRAWSTRING BAG

| | | |
|--------|---------------|-----|
| DEF861 | 445mm x 465mm | 100 |
| DEF864 | 445mm x 465mm | 25 |

DRILL SLEEVE COVER

| | | |
|---------|------------------------------------|-----|
| DEF1224 | 5cm x 80cm WITH 3 LABELS | 100 |
| DEF1368 | 5cm x 110cm WITH TAPE | 100 |
| DEF618 | 10cm x 152cm WITH TAPE D/P | 25 |
| DEF1516 | 7.5cm x 1m WITH ADHESIVE BOTH ENDS | 100 |
| DEF1154 | 7.5cm x 1.5m WITH R/BAND | 100 |



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

DONOR ORGAN BAG

| | | |
|--------|--------------------------|----|
| DEF847 | 850mm x 500mm DRAWSTRING | 25 |
|--------|--------------------------|----|

**FREEZER BAG**

| | | |
|-------|-----------------------------------|-----|
| DEF66 | 300mm x 450mm & 2-#34 R.BANDS D/P | 100 |
|-------|-----------------------------------|-----|

**LEAD SCREEN COVERS**

| | | |
|---------|--------------------------|----|
| DEF766T | 575mm x 600mm WITH TAPES | 50 |
|---------|--------------------------|----|

PLASTIC SLEEVE COVER

| | | |
|--------|---------------------------------|-----|
| DEF454 | 15cm x 1m & 2 - #31 R/BANDS D/P | 100 |
|--------|---------------------------------|-----|



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

PROBE COVERS

| | | |
|---------|----------------------|-----|
| DEF1482 | 23cm x 1m SEALED END | 100 |
|---------|----------------------|-----|



| | | |
|---------|--------------------------|----|
| DEF648 | DRAWSTRING 17.5cm x 1m | 25 |
| DEF2630 | DRAWSTRING 17.5cm x 1.5m | 25 |

PLASTIC BAGS

| | | |
|---------|--------------------------------|-----|
| DEF642 | WITH COTTON TAPE 455mm x 510mm | 100 |
| DEF861 | 445mm x 465mm DRAW STRING | 100 |
| DEF1878 | 140mm x 300mm | 100 |



| | | |
|---------|---|-----|
| DEF1276 | 140mm x 300mm WITH BLANK ADHESIVE LABEL | 100 |
| DEF650 | BURNS PLASTIC BAG 305mm x 390mm | 100 |
| DEF647 | 200mm x 200mm | 100 |
| DEF612 | 455mm x 510mm 5's | 100 |
| DEF638 | 280mm x 325mm | 100 |

PROTECTIVE COVER

| | | |
|---------|--------------|----|
| DEF613 | 15cm x 1m | 50 |
| DEF763 | 104cm x 1.5m | 40 |
| DEF1573 | 208cm x 1.5m | 50 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SPECIMEN PLASTIC BAG

| | | |
|--------|--------------------------|-----|
| DEF808 | 575mm x 620mm HEAVY DUTY | 100 |
|--------|--------------------------|-----|

TROLLEY BAG COVER

| | | |
|--------|----------------|----|
| DEF641 | 610cm x 1020cm | 50 |
|--------|----------------|----|



| | | |
|---------|---------------|----|
| DEF1239 | 200cm x 300cm | 25 |
|---------|---------------|----|

ULTRASOUND PROBE COVER

| | | |
|---------|------------------------|-----|
| DEF1572 | SEALED END 17cm x 50cm | 100 |
| DEF1088 | SEALED END 23cm x 80cm | 100 |



| | | |
|---------|---|-----|
| DEF643 | SEALED POINTED END 23cm x 1.5m | 25 |
| DEF2409 | SEALED END 25cm x 1.5m WITH #64 RUBBER BAND | 100 |

> SURGICAL TAPE



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

COTTON TAPE

| | | |
|--------|------------------|-----|
| DEF173 | 1.2cm x 60cm D/P | 100 |
|--------|------------------|-----|



| | | |
|---------|---------------|-----|
| DEF2106 | 1.25cm x 1.5m | 100 |
|---------|---------------|-----|

non sterile roll

| | | |
|--------|--------------|---|
| DEF918 | 1.2cm x 500m | 1 |
|--------|--------------|---|

| | | |
|--------|--------------|---|
| DEF919 | 2.5cm x 250m | 1 |
|--------|--------------|---|

SURGICAL TAPE - F

| | | |
|---------|-------------|-----|
| DEF1504 | 10cm x 12cm | 100 |
|---------|-------------|-----|



| | | |
|---------|-------------|-----|
| DEF1505 | 10cm x 20cm | 100 |
|---------|-------------|-----|

| | | |
|---------|-------------|-----|
| DEF1506 | 10cm x 45cm | 100 |
|---------|-------------|-----|

| | | |
|---------|-------------|----|
| DEF1614 | 15cm x 70cm | 50 |
|---------|-------------|----|

| | | |
|---------|-------------|-----|
| DEF1507 | 20cm x 20cm | 100 |
|---------|-------------|-----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SURGICAL TAPE - H

| | | |
|---------|--------------|----|
| DEF1623 | 2.5cm x 30cm | 50 |
|---------|--------------|----|

| | | |
|---------|-----------|----|
| DEF1647 | 5cm x 8cm | 50 |
|---------|-----------|----|

| | | |
|---------|------------|----|
| DEF1995 | 5cm x 15cm | 66 |
|---------|------------|----|

| | | |
|---------|------------|----|
| DEF1599 | 5cm x 25cm | 50 |
|---------|------------|----|

| | | |
|---------|------------|----|
| DEF1624 | 5cm x 30cm | 50 |
|---------|------------|----|

| | | |
|---------|------------|----|
| DEF1797 | 5cm x 50cm | 60 |
|---------|------------|----|

| | | |
|---------|-------------|-----|
| DEF1044 | 10cm x 10cm | 100 |
|---------|-------------|-----|



| | | |
|--------|-------------|----|
| DEF944 | 10cm x 25cm | 50 |
|--------|-------------|----|

| | | |
|---------|---------------|-----|
| DEF1057 | 10cm x 36.5cm | 100 |
|---------|---------------|-----|

| | | |
|---------|-------------|----|
| DEF1798 | 10cm x 50cm | 60 |
|---------|-------------|----|

| | | |
|---------|-------------|----|
| DEF1600 | 15cm x 30cm | 50 |
|---------|-------------|----|

| | | |
|---------|-------------|-----|
| DEF1026 | 15cm x 60cm | 100 |
|---------|-------------|-----|

| | | |
|---------|--------------|-----|
| DEF1025 | 15cm x 100cm | 100 |
|---------|--------------|-----|

| | | |
|--------|-------------|-----|
| DEF678 | 20cm x 30cm | 100 |
|--------|-------------|-----|

| | | |
|--------|-----------------|----|
| DEF475 | 20cm x 75cm D/P | 50 |
|--------|-----------------|----|


| | | |
|---------|------------|-----|
| DEF2160 | 8cm x 20cm | 125 |
|---------|------------|-----|

| | | |
|--------|-------------|----|
| DEF943 | 10cm x 20cm | 50 |
|--------|-------------|----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SURGICAL TAPE NICHIMESH


| | | |
|--------|---------------------|----|
| DEF947 | 5cm x 30cm 10's D/P | 50 |
|--------|---------------------|----|

**non sterile roll**

| | | |
|--------|-------------|---|
| DEF925 | 3.8cm x 10m | 1 |
| DEF927 | 10cm x 10m | 1 |

SURGICAL RETRACTION TAPE R

| | | |
|--------|------------|-----|
| DEF167 | 2.5cm x 2m | 100 |
|--------|------------|-----|



To order, please contact our
Customer Service Team on

1300 550 278

or see our website for details
www.defries.com.au

> SILICONE DRAINS



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SILICONE DRAINS DOUBLE PACKED

- Soft, flexible and inert for patient comfort and faster healing.
- Closed Wound Drain Systems are used to remove blood and excess fluid from wounds, so as to prevent or minimise haematomas and eliminate dead space.
- Promotes the approximation of skin flaps and reduces the potential of post-op infections.
- Closed Wound Drain Systems are composed of 2 components - Surgical Drain and Evacuators/Reservoirs.
- Trocars can be included or can be attached during the surgical procedure.

SILICONE BULB

- Silicone drainage Bulb is compatible with the following drains below: Flat Drains, Fluted Drains, Round Drains.
- Squeeze the Bulb and then attach to the drain for suction.
- As the Silicone Bulb expands, it creates a vacuum to a maximum of 125mmHg pressure, and slowly sucks fluid and blood from the wound site, preventing Haematomas forming.
- The Bulb is usually removed 1 - 2 days after the operation.

| Code | Description | Box Qty. |
|-------------|-------------------------|----------|
| 313300380DS | EVACUATOR BULB 100ml/cc | 10 |

**FLAT SILICONE DRAINS**

- Hubless drains
- Soft and comfortable silicone flat drains with multiple drainage holes 36 x 2 Holes, to enhance the fluid removal in a wound in the patient.
- Ergonomic Flat shape enabling easier drain placement.
- Indicated for Abdominal Surgery, and suitable for use as a Pleural Space Drain.

| Code | Description | Box Qty. |
|---------------|-----------------------|----------|
| 11100200780DS | WHITE X/R 7mm x 190mm | 10 |



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|---------------|------------------------|----|
| 11100201080DS | WHITE X/R 10mm x 190mm | 10 |
|---------------|------------------------|----|

**FLUTED SILICONE DRAINS**

- Flexible silicone drains with four channel unique H shape.
- This drain shape allows for greater multi-channel drainage, and also ensures ease of removal due to no solid silicone outer drain perimeter.
- Inner core of H design offers drain strength and flexibility.
- Suitable for Plastic Surgery particularly Breast, Abdominoplasty and Rhytidectomies. Also used extensively in Cardiac Surgery. Useful as an adjunct to the following surgical procedures: Orthopedic, Neurological, Gynecological, Head and Neck.

**round**

| | | |
|---------------|------------------------------------|----|
| 11003131080DS | CLEAR FULL FLUTED 10FR with Trocar | 10 |
|---------------|------------------------------------|----|

| | | |
|---------------|------------------------|----|
| 11003021080DS | CLEAR FULL FLUTED 10FR | 10 |
|---------------|------------------------|----|



| Code | Description | Box Qty. |
|---------------|------------------------------------|----------|
| 11003831080DS | WHITE WITH HUB 10FR | 10 |
| 11003131580DS | CLEAR FULL FLUTED 15FR with Trocar | 10 |
| 11003131980DS | CLEAR FULL FLUTED 19FR with Trocar | 10 |



| | | |
|---------------|---------------------|----|
| 11003811980DS | WHITE 19FR with Hub | 10 |
|---------------|---------------------|----|



| | | |
|---------------|--|----|
| 11003122480DS | CLEAR FULL FLUTED 24FR Hubless flat white | 10 |
| 11101310780DS | 3/4 FLUTED 7mm with Trocar, with Hub | 10 |
| 11101311080DS | 3/4 FLUTED 10mm with Trocar, with Hub | 10 |
| 11101410780DS | FULL FLUTED 7mm with Trocar, with Hub | 10 |
| 11101411080DS | FULL FLUTED 10mm with Trocar, with Hub | 10 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SILICONE ROUND CLEAR DRAINS X-RAY

- Drain consists of 5 side holes with a drain length of 50cm.
- Used for General Surgery, Chest and Neck Surgery.



| | | |
|---------------|---------------------|----|
| 11000301080DS | 10 FR 1.6mm x 3.2mm | 10 |
| 11000301280DS | 12 FR 2mm x 4mm | 10 |
| 11000301580DS | 15 FR 3mm x 5mm | 10 |
| 11000301880DS | 18FR 4mm x 6mm | 10 |
| 11000302180DS | 21FR 4.5mm x 7mm | 10 |
| 11000302480DS | 24 FR 5mm x 8mm | 10 |
| 11000302780DS | 27FR 6mm x 9mm | 10 |
| 11000303080DS | 30 FR 7mm x 10mm | 10 |
| 11000303680DS | 36 FR 8mm x 12mm | 10 |

CAPILLARY ACTION DRAINS - PASSIVE WOUND DRAINAGE MULTILUMEN DRAINS X-RAY

- Clear Drains with varying lumens according to drainage required.
- Lumens can be separated out and positioned as necessary. The drainage takes place both through, under and over the lumens.

| | | |
|---------|-------------------|----|
| DEF1110 | 2 LUMENS x 12.5cm | 10 |
| DEF1076 | 5 LUMENS x 25cm | 10 |




| | | |
|---------|-----------------------|----|
| DEF1788 | 9 LUMENS 2.5cm x 20cm | 10 |
| DEF1024 | 9 LUMENS 2.5cm x 25cm | 10 |
| DEF1750 | 16 LUMEN x 20cm | 10 |
| DEF1751 | 16 LUMEN x 40cm | 10 |

EASY FLOW DRAINS -

- 300mm Drain Length Radio-opaque Striped Clear drain.
- Drain incorporates inner "teeth" to ensure open drainage.

| | | |
|----------|-----------|----|
| 21740632 | FLAT 40mm | 10 |
| 21706632 | ROUND 6mm | 10 |

> SILICONE LOOPS & BOOTS




| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

IDENTILOOPS

- Silicone radio-opaque Identi Loops are used in surgical situations requiring retraction, identification or occlusion of arteries, veins, nerves, ureters or tendons.
- They are oval in cross-section, allowing for even pressure distribution, minimizing damage to tissues.

| | | |
|---------|--|----|
| DEF1318 | 2 BLUE, 1 RED, 1 YELLOW, 1 WMAXI-WHITE 40cm 5's | 10 |
|---------|--|----|



| | | |
|-----------|--------------------------|----|
| 2200620-2 | MINI-BLUE 40cm 2's | 10 |
| 2202620-2 | MAXI-BLUE 40cm 2's | 10 |
| DEF1282 | MAXI-BLUE 40cm 2's D/P | 50 |
| DEF1193 | SUPER MAXI BLUE 40cm 2's | 10 |



| | | |
|-----------|--------------------------|----|
| DEF2182 | SUPER MAXI BLUE 60cm 2's | 10 |
| 2200220-2 | MINI-RED 40cm 2's | 10 |
| 2200220-4 | MINI-RED 40cm 4's | 10 |
| 2202220-2 | MAXI-RED 40cm 2's | 10 |
| 2202220-4 | MAXI-RED 40cm 4's | 10 |
| 2200420-2 | MINI-YELLOW 40cm 2's | 10 |
| 2202420-2 | MAXI-YELLOW 40cm 2's | 10 |
| 2200920-2 | MINI-WHITE 40cm 2's | 10 |
| 2202920-2 | MAXI-WHITE 40cm 2's | 10 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SILCLAMPS

- The Silclamps are covers for surgical clamps, hemostats and forceps.
- They are soft and pliable, providing even pressure distribution, minimizing tissue trauma.

| | | |
|---------|------------------------|-----|
| DEF1628 | BLUE 3/10" x 10cm 2's | 50 |
| DEF1189 | RED 1/10" x 10cm SMALL | 100 |

**SUTURE BOOTS SILICONE WITH ADHESIVE TRAY DP**

- Suture Boot Covers offer a non-damaging way to firmly hold delicate sutures without slippage.
- Radio-opaque.

| | | |
|---------|-----|----|
| DEF1093 | 2's | 25 |
|---------|-----|----|



| | | |
|---------|------|----|
| DEF1094 | 4's | 25 |
| 30-702 | 6's | 25 |
| DEF1096 | 8's | 25 |
| 30-701 | 10's | 25 |
| DEF2284 | 12's | 25 |

> SILICONE TUBING RANGE



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SILICONE TUBING

- 100% Medical Grade, clear Silicone Tubing available in both Metric and Inch Range.
- Please contact Customer Service if you prefer a different size.

SILICONE TUBING METRIC

| | | |
|---------|-----------------------------------|-----|
| DEF1564 | 1.57mm x 3.18mm x 20cm | 25 |
| DEF2617 | 1mm x 3mm x 50cm | 50 |
| DEF473 | 1mm x 3mm x 50cm 4's | 20 |
| DEF2522 | 1mm x 2mm x 50cm | 50 |
| DEF1428 | 1.57mm x 2.41mm x 15cm D/P | 100 |
| DEF1944 | 3mm x 5mm x 165cm | 15 |
| DEF2314 | 3.18mm x 6.35mm x 25cm | 60 |
| DEF2038 | 4mm x 6mm x 45cm | 55 |
| DEF2786 | 6.4mm x 9.5mm x 2cm PACKED IN 2's | 375 |
| DEF2060 | 5mm x 8mm x 3.1m | 8 |
| DEF1268 | 5mm x 7mm x 42cm | 25 |
| DEF481 | 5mm x 8mm x 45cm | 50 |
| DEF1261 | 5mm x 9mm x 1.5m | 50 |
| DEF2032 | 7mm x 10mm x 3m | 8 |
| DEF482 | 7mm x 12mm x 85mm 2's D/P | 100 |
| DEF1204 | 8mm x 14mm x 3m | 5 |
| DEF2468 | 10mm x 16mm x 150cm | 16 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SILICONE TUBING INCH RANGE

| | | |
|---------|------------------------------|----|
| DEF477 | 0.62" x 0.95" x 30cm 2'S D/P | 25 |
| DEF2307 | 0.012" x 0.025" x 20cm | 75 |
| DEF1545 | 0.062" x 0.125" x 1m | 15 |
| DEF1540 | 0.078" x 0.125" x 21cm | 70 |

TUBING CONNECTOR

| | | |
|--------|-------|-----|
| DEF465 | SMALL | 100 |
|--------|-------|-----|

**SILICONE SNUGGERS CARDIAC**

- Used during Cardiac Bypass Surgery

| | | |
|---------|-------------------------------|----|
| DEF1226 | 4.77mm x 9.5mm x 15cm 2's D/P | 50 |
| DEF1665 | 5mm x 7mm x 20cm 2's D/P | 25 |
| DEF2271 | 18FR x 15cm 6's D/P | 30 |

SILICONE TUBING WITH MEDIUM CONNECTOR ATTACHED

- Flexible Tubing for suctioning

| | | |
|---------|----------------------|----|
| DEF735 | 5mm x 8mm x 1m | 50 |
| DEF1186 | 5mm x 9mm x 30cm D/P | 50 |
| DEF1712 | 5mm x 9mm x 60cm | 40 |



| | | |
|----------|---------------------|----|
| DEF1337 | 5mm x 9mm x 1m | 50 |
| DEF1231A | 5mm x 9mm x 1.5m | 30 |
| DEF305 | 7mm x 11mm x 1m D/P | 50 |
| DEF1253 | 5mm x 8mm x 1.6m | 50 |
| DEF306 | 5mm x 9mm x 50cm | 50 |

| Code | Description | Box Qty. |
|---|---|----------|
| SILICONE TUBING NON STERILE | | |
| • Please contact Customer Service if you prefer a different size. | | |
| 203005009 | 0.020" x .037" / 0.51mm x 0.94mm 15m | 1 |
| 203095159 | 0.375" x .625" / 9.50mm x 15.90mm 15m | 1 |
| 203158222 | 0.625" x .875" / 15.87mm x 22.23mm 15m | 1 |
| 201005010 | 0.50mm ID x 1mm OD 25m | 1 |
| 201020030 | 2mm (ID) x 3mm (OD) 25m | 1 |
| 201020040 | 2mm (ID) x 4mm (OD) 25m | 1 |
| 201020060 | 2mm (ID) x 6mm (OD) 25m | 1 |
| 201030070 | 3mm (ID) x 7mm (OD) 25m | 1 |
| 201040080 | 4mm (ID) x 8mm (OD) 25m | 1 |
| 201050100 | 5mm (ID) x 10mm (OD) 25m | 1 |
| 201060080 | 6mm (ID) x 8mm (OD) 25m | 1 |
| 201060100 | 6mm (ID) x 10mm (OD) 25m | 1 |
| 201060110 | 6mm (ID) x 11mm (OD) 25m | 1 |
| 201070120 | 7mm (ID) x 12mm (OD) 25m | 1 |
| 201080140 | 8mm (ID) x 14mm (OD) 25m | 1 |
| 201080160 | 8mm (ID) x 16mm (OD) 25m | 1 |
| 201090130 | 9mm (ID) x 13mm (OD) 25m | 1 |
| 201100120 | 10mm (ID) x 12mm (OD) 25m | 1 |
| 201100160 | 10mm (ID) x 16mm (OD) 25m | 1 |
| 201100180 | 10mm (ID) x 18mm (OD) 25m | 1 |
| 201120170 | 12mm (ID) x 17mm (OD) 25m | 1 |
| 201150210 | 15mm (ID) x 21mm (OD) 25m | 1 |
| 201160220 | 16mm (ID) x 22mm (OD) 25m | 1 |
| 201160240 | 16mm (ID) x 24mm (OD) 25m | 1 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|



To order, please contact our
Customer Service Team on
1300 550 278
or see our website for details
www.defries.com.au

> WARD CONSUMABLES



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

COTTON BALLS

| | | |
|--------|-----------|-----|
| DEF725 | LARGE 4'S | 100 |
| DEF074 | LARGE 5'S | 100 |

**COTTON WOOL ROLL**

| | | |
|--------|------------------|----|
| DEF874 | 15cm x 55cm D/P | 50 |
| DEF875 | 30cm x 95cm D/P | 25 |
| DEF880 | 38cm x 165cm D/P | 15 |

I.V. PRODUCTS**I.V. PRESSURE PADS**

- Unique application - Allows you to apply the pressure pad without contaminating the pad or adhesive.
- Hypoallergenic adhesive. Latex free.

| | | |
|--------|---|-----|
| DEF929 | 36mm Total size, Pad size is 16mm x 3mm white | 100 |
|--------|---|-----|



| | | |
|---------|-----------------------------|----|
| DEF1251 | 15mm x 27mm x 6mm PAD WHITE | 50 |
| DEF1252 | 15mm x 27mm x 9mm PAD WHITE | 30 |



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

IV INSERTION PACKS

- For the insertion of IV therapy.

| | | |
|--------|-----------------------|-----|
| DEF920 | BRONNY IV PACK | 200 |
|--------|-----------------------|-----|

2 x Gauze Plain Swab 5cm x 5cm

1 x Plastic Drape 60cm x 45cm



| | | |
|---------|---------------------------------------|-----|
| DEF1676 | BRONNY IV PACK WITH 4 SWABS DP | 200 |
|---------|---------------------------------------|-----|

4 x Gauze Plain Swab 5cm x 5cm

1 x Plastic Drape 60cm x 45cm

**FORCEP PLASTIC**

| | | |
|--------|-------|-----|
| DEF686 | GREEN | 100 |
|--------|-------|-----|



| Code | Description | Box Qty. |
|----------|---|----------|
| DEF219 | PLASTICS PACK | 60 |
| | 1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm Deep | |
| | 1 x Gallipot 60ml | |
| | 12 x Gauze Swabs 10cm x 10cm | |
| | 2 x Towel White 2 Ply | |
| | 2 x Swab Sticks 15cm with Wooden Shaft | |
| | Field Blue Plastic 50cm x 50cm | |
| DEF2234A | LARGE DRESSING PACK | 10 |
| | 1 x Dressing Tray 23.5cm x 29.5cm x 5cm | |
| | 1 x Kidney Dish Clear 700ml | |
| | 2 x Gallipot 150ml | |
| | 2 x Dressing Towel - T Small 32cm x 32cm | |
| | 10 x Non Woven Balls | |
| | 3 x Green Disposable Forceps | |
| | 10 x Non Woven Swab 7.5cm x 7.5cm | |
| | Field Blue Plastic 91cm x 115cm | |
| DEF2746B | LUMBAR PUNCTURE PACK | 24 |
| | 1 x Dressing Tray - 3 Cavity 25.5cm x 20.5cm x 2.5cm | |
| | 5 x Gauze Sponges 5cm x 5cm | |
| | 2 x Green Disposable Forceps | |
| | 1 x Spinal Manometer 0-34 | |
| | 1 x Hypodermic Needle 25g x 5/8" | |
| | 1 x Syringe 5ml Luer Lock | |
| | 1 x IV 3000 Dressing 6cm x 8.5cm | |
| | 1 x Specimen Jar 5ml | |
| | 3 x Specimen tube 10ml with "V" Bottom | |
| | 1 x SMS Drape 61cm x 61cm with 14cm Circular Fenestration and with Corner Tabs | |
| | 1 x Dressing Towel - T Small 32cm x 32cm | |
| | 1 x Label DEF4618, DEF4619, DEF4620, DEF4621 | |
| | Field Blue Plastic 75cm x 75cm | |

| Code | Description | Box Qty. |
|---------|---|----------|
| DEF2765 | LUMBAR PUNCTURE PACK | 15 |
| | 1 x Set of 3 Labelled Specimen Tubes 10ml "V" Bottom | |
| | 1 x Plastic Bag 126mm x 152mm (Specimen Tubes Enclosed) | |
| | 2 x Gallipot 60ml | |
| | 2 x Green Disposable Forceps | |
| | 1 x Spinal Manometer 0-34 | |
| | 5 x Green Sponges 10cm x 10cm | |
| | 1 x Tegaderm Dressing 4.4cm x 4.4cm | |
| | 1 x Syringe 5ml Luer Slip | |
| | 1 x Hypodermic Needle 23g x 1.25" | |
| | 1 x Stopcock 3 Way low Pressure-Off | |
| | 1 of 4 Lumbar Puncture Information Sheets | |
| | 1 x Dressing Tray - 3 Cavity - 25.5cm x 20.5cm x 2.5cm | |
| | 1 x Dressing Towel - T Small 32cm x 32cm | |
| | 1 x Clear Plastic Drape 76cm x 106cm with Full 2.5cm Adhesive around 10cm x 15cm Fenestration and Adhesive on Corner tabs | |
| | Field Blue Plastic 75cm x 75cm | |

JUMBO MOUTH & THROAT SWAB STICKS PAPERSHAFT

- 100% Paper shaft, 100% Cotton head. Head securely attached.

| | | |
|--------|----------|-----|
| DEF400 | 15cm 2's | 150 |
| DEF401 | 15cm 3's | 150 |
| DEF402 | 15cm 5's | 150 |
| DEF455 | 20cm 2's | 150 |
| DEF407 | 20cm 3's | 150 |



| | | |
|---------|--------------------|-----|
| DEF1402 | 20cm 4's | 100 |
| | non sterile | |
| DEF055 | 15cm | 200 |
| DEF054 | 20cm | 200 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

MOUTH TOILET PACK

| | | |
|--------|---|----|
| DEF399 | 10 x 15cm M&T Swab Stick & 1 x 150ml Gallipot | 50 |
|--------|---|----|

| | | |
|---------|-----------------|----|
| DEF1866 | MSU PACK | 50 |
|---------|-----------------|----|

1 x Kidney Dish 400ml
 1 x Gallipot 60ml
 3 x Non-Woven Swabs 7.5cm x 7.5cm
 7 x Non Woven Balls
 2 x Vinyl Gloves Large
 Field Blue Plastic 50cm x 50cm

**NASAL BOLSTER**

- Non woven cover, Extremely absorbant and Comfortable

| | | |
|----------|---------|-----|
| DEF1381 | BOLSTER | 20 |
| DEF 1382 | BOLSTER | 100 |

**non sterile**

| | | |
|---------|---------|----|
| DEF1380 | BOLSTER | 20 |
|---------|---------|----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|---------|----------------|----|
| DEF1407 | PV TRAY | 60 |
|---------|----------------|----|

1 x Dressing Tray with 3 compartments
 23cm x 9cm x 3cm
 1 x Prep Swabs Small
 1 x Prep Swabs Large
 1 x Dressing Towel - T Small 32cm x 32cm
 Field Blue Plastic 75cm x 75cm

**RENAL PRODUCTS**

| | | |
|---------|----------------------------|----|
| DEF1525 | RENAL DIALYSIS PACK | 80 |
|---------|----------------------------|----|

1 x Dressing Tray 12.5cm x 10.8cm x 2.5cm
 2 x Green Disposable Forceps
 12 x Cotton Balls Large
 2 x Towels White 2ply 40 x 42cm
 8 x Gauzes Sponges Plain 10cm x 10cm
 1 x Towel Polybacked 40cm x 60cm
 Field Blue Plastic 75cm x 75cm

| | | |
|---------|--|----|
| DEF1250 | CHUSHAVEN PRESSURE PAD 27mm x 27mm x 8 mm THICK | 60 |
|---------|--|----|

| | | |
|---------|--|----|
| DEF1251 | STEPTY PRESSURE PAD 15mm x 27mm x 6mm THICK | 50 |
|---------|--|----|

| | | |
|---------|--|----|
| DEF1252 | STEPTY "P" PRESSURE PAD 15mm x 27mm x 9mm THICK | 30 |
|---------|--|----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

GAUZE SPONGES PLAIN 8 PLY

| | | |
|--------|--------------------|-----|
| DEF183 | 7.5cm x 7.5cm 3's | 100 |
| DEF181 | 7.5cm x 7.5cm 5's | 100 |
| DEF371 | 7.5cm x 7.5cm 20's | 50 |
| DEF189 | 10cm x 10cm 2's | 100 |
| DEF186 | 10cm x 10cm 3's | 100 |
| DEF045 | 10cm x 10cm 5's | 100 |



| | | |
|--------------------|-------------------------------|-----|
| DEF178 | 10cm x 10cm 20's | 25 |
| DEF1033 | 10cm x 20cm 4's | 100 |
| DEF201 | 10cm x 20cm 10's | 250 |
| DEF195 | 10cm x 20cm 10's | 50 |
| DEF853 | 10cm x 80cm 4's 16 PLY FOLDED | 50 |
| DEF213 | 10cm x 10cm 10's | 300 |
| DEF208 | 10cm x 10cm 10's | 50 |
| DEF177 | 10cm x 10cm 25's | 25 |
| DEF214 | 10cm x 10cm 25's | 250 |
| DEF1696 | 10cm x 10cm 50's | 50 |
| DEF198 | 10cm x 20cm 5's | 250 |
| non sterile | | |
| DEF136 | 10cm x 10cm PLAIN BULK | 100 |

SPONGES NON WOVEN SWABS

| | | |
|---------|---------------|-----|
| DEF1131 | 5cm x 5cm 5's | 100 |
|---------|---------------|-----|

SAFETY SCALPEL AND HANDLE

| | | |
|---------|-----|-----|
| DEF2700 | #10 | 100 |
| DEF2701 | #11 | 100 |



| | | |
|---------|-----|-----|
| DEF2702 | #15 | 100 |
|---------|-----|-----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SANITARY BAG & DISPENSER**non sterile**

| | | |
|--------|--|----|
| DEF692 | BAG DISPENSER | 1 |
| DEF709 | SANITARY DISPOSAL BAG 100's FLORAL DESIGN | 20 |

STAPLE REMOVER

| | | |
|---------|----------------|----|
| DEF1826 | PLASTIC HANDLE | 50 |
|---------|----------------|----|

**non sterile**

| | | |
|---------|----------------|----|
| DEF1825 | PLASTIC HANDLE | 20 |
|---------|----------------|----|

7.5cm SWAB STICKS WITH WOODEN SHAFT

| | | |
|--------|---------|-----|
| DEF057 | S/E 3's | 100 |
|--------|---------|-----|



| | | |
|--------|---------|-----|
| DEF442 | S/E 5's | 100 |
|--------|---------|-----|

15cm SWAB STICKS WITH WOODEN SHAFT

| | | |
|--------|---------|-----|
| DEF439 | S/E 2's | 100 |
|--------|---------|-----|



| | | |
|--------|---------|-----|
| DEF059 | S/E 3's | 100 |
| DEF443 | S/E 5's | 100 |

non sterile

| | | |
|--------|----------|------|
| DEF062 | S/E BULK | 1000 |
|--------|----------|------|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

7.5cm SWAB STICKS WITH PAPER SHAFT

| | | |
|--------|----------|-----|
| DEF424 | D/E 5's | 100 |
| DEF421 | D/E 10's | 100 |



| | | |
|--------|------------------------------------|------|
| DEF104 | D/E 50's D/P non sterile | 80 |
| DEF061 | D/E BULK | 1250 |

7.5cm SWAB STICK WITH PLASTIC SHAFT

| | | |
|--------|---------|-----|
| DEF058 | D/E 3's | 100 |
| DEF447 | D/E 3's | 200 |

**SWAB STICK PROCTOSCOPE**

| | | |
|--------|----------------|----|
| DEF954 | 40cm PAPER S/E | 50 |
|--------|----------------|----|



To order, please contact our
Customer Service Team on

1300 550 278

or see our website for details
www.defries.com.au



> CUSTOM PRINTED LABELS



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

CUSTOMISED LABELS

- Please contact Customer Service to discuss your label requirements.

**NON STERILE LABELS**

White Background / Blue, Black, Red, Magenta, Yellow & Green Writing only in Rolls of 100 or 1000

- Size: 15mm x 30mm (Anaesthetic Size)
- Size: 21mm x 61mm
- Size: 32mm x 50mm
- Size: 32mm x 70mm (Standard Size)
- Size: 50mm x 70mm
- Size: 40mm x 102mm
- Size: 40mm x 170mm
- Size: 80mm x 102mm (2 Rolls of 500)

Full Colour Background (Any colour combination can be made either for the background or writing) in Rolls of 100 or 1000

- Size: 15mm x 30mm (Anaesthetic Size)
- Size: 21mm x 61mm
- Size: 32mm x 50mm
- Size: 32mm x 70mm (Standard Size)
- Size: 50mm x 70mm
- Size: 40mm x 102mm
- Size: 40mm x 170mm
- Size: 80mm x 102mm (2 Rolls of 500)

Fluorescent Labels in Orange, Yellow, Red, Green, Purple, Blue or Pink with Black writing in Rolls of 100 or 1000

- Size: 32mm x 50mm - ONLY Fluorescent Red
- Size: 30mm x 70mm
- Size: 50mm x 70mm - ONLY Fluorescent Green
- Size: 80mm x 102mm



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

STERILE LABELS

White Background / Blue, Black, Red, Magenta, Yellow or Green Writing only in Boxes 100 (15mm x 30mm)

- Packed in 1's
- Packed in 2's
- Packed in 3's
- Packed in 6's
- Packed in 12's
- Packed in 18's
- Packed in 24's



Full Colour Background (Any colour combination can be made either for the background or writing) in Boxes 100 (15mm x 30mm)

- Packed in 1's
- Packed in 2's
- Packed in 3's
- Packed in 6's
- Packed in 12's
- Packed in 18's
- Packed in 24's

PIGGY BACK LABELS

PBLABELS - B (Blank) in Rolls of 100 or 1000

- 25mm x 50mm
- 32mm x 50mm
- 36mm x 105mm

PBLABELS - P (Printed) in Rolls of 100 or 1000

- 25mm x 50mm
- 32mm x 50mm
- 36mm x 105mm

PBLABELS - S (Sterrad) in Rolls of 100 or 1000

- 25mm x 50mm
- 32mm x 50mm
- 50mm x 70mm

Code Description Box Qty.


STERILE PERI-OP LABELS

- Labels printed, folded and packaged in a peel down pouch, then sterilised.
- Label sheet is folded into a 3 fold to minimise packaging size.
- Labels can be customised according to the National Recommendations for user-applied Labelling.
- Labels are waterproof and can be written on by surgical markers: DEF3219 Surgical Marker Permanent Fine Tip.
- LAB9009 and LAB9010 are currently used in hospitals.
- Cost effective price and ready to use for every procedure.
- Sterile and Latex free.

LAB9009 Sterile Medication Label Sheet – specialty General, Ortho, Neuro, Cardiac, Vascular, Urology, Obstetrics & Gynaecology, Breast & Endocrine specific. 100



LAB9010 Sterile Medication Label Sheet – specialty ENT, Max/Facial, Plastic, Ophthalmology and Neuro Surgery. 100



Code Description Box Qty.

DRUG LABELS - AUSTRALIAN STANDARDS non sterile

| | | |
|------------|--|-----|
| AL/AE | ANTI-EMETIC - SALMON | 500 |
| AL/AE-5 | ONDANSETRON HCl(ZOFRAN) - SALMON | 500 |
| AL/BW | BLANK NO TEXT - WHITE | 500 |
| AL/HA-6 | CLONIDINE - VIOLET STRIPE | 500 |
| AL/HA-7 | GTN - VIOLET/ WHITE | 500 |
| AL/HA-8 | VERAPAMIL - VIOLET | 500 |
| AL/HA-9 | PHENOXYBENZAMINE HCl - VIOLET/WHITE DIAG | 500 |
| AL/IA-2 | PROPOFOL mg/ml -YELLOW | 500 |
| AL/IA-6 | TRAMADOL - BLACK ON YELLOW | 500 |
| AL/LA-2 | PROCAINE - GREY | 500 |
| AL/LAB | BLANK - GREY | 500 |
| AL/MES-18 | CHANGE ON - WHITE | 500 |
| AL/MES-21 | ALLERGY - RED | 500 |
| AL/MES-34 | CVP BLACK ON BLUE | 500 |
| AL/MES-35 | REGIONAL INF N/S | 500 |
| AL/MES-9 | PULMONARY ARTERY - BLACK ON YEL | 500 |
| AL/MIS-27 | GENTAMYCIN - BLACK ON WHITE | 500 |
| AL/MIS-30 | THROMBIN - BLACK ON WHITE | 500 |
| AL/MIS-37 | SODIUM BIOC. NON STERILE - WHITE | 500 |
| AL/MIS-39 | METHADONE - BLACK ON WHITE | 500 |
| AL/MIS-47 | ABCIXIMAB - WHITE | 500 |
| AL/MIS-50 | KETOROLAC - BLACK ON WHITE B/GROUND | 500 |
| AL/MIS-62 | TRAMADOL - BLACK ON WHITE | 500 |
| AL/MIS-70 | POTASSIUM | 500 |
| AL/MR-1 | SUXAMETHONIUM mg/ml - RED | 500 |
| AL/N | NARCOTIC - BLUE | 500 |
| AL/N-2 | FENTANYL mg/ml - BLUE | 500 |
| AL/NB | BLANK - BLUE | 500 |
| AL/PATHC04 | PATHOLOGY LABELS - LOWER LOBE - CARINE BK/ORANGE | 500 |
| AL/VB | BLANK NO TEXT - VIOLET | 500 |

| Code | Description | Box Qty. |
|---------------|---|----------|
| LABELS | | |
| ME1 | DUST COVER LABEL WITH USE BY | 400 |
| ME15 | LOAN SET | 500 |
| ME16 | ETO GAS STERILIZE 2cm x 6cm - FLURO ORANGE | 1 ROLL |
| ME2 | DECONTAMINATED | 400 |
| ME21 | LOAN ITEM - BLACK ON PINK | 800 |
| ME23 | CAUTION - PATIENT WITH SAME OR SIMILAR NAME | 500 |
| ME25 | ADDITIVE DRUG LABEL ADDITIVE - RED/WHT BORDER | 500 |
| ME6G | STERILE IF PACK IS UNOPENED - GREEN | 400 |
| ME6P | STERILE IF PACK IS UNOPENED - PINK | 400 |
| ME6Y | STERILE IF PACK IS UNOPENED - YELLOW | 400 |
| ME7 | URGENT E STEAM STERILSABLE - RED ON WHITE | 800 |
| ME9 | TRAY INCOMPLETE LABEL - ORANGE | 400 |
| MEP11 | PLAIN WHITE POLYPROPYLENE | 500 |
| CP1/W | STERILIZED BATCH LABEL - WHITE PIGGY BACK | 500 |
| PBLA-BELS-P | PIGGY BACK LABELS (PRINTED) | 1000 |
| PBLA-BELS-S | PIGGY BACK LABELS "STERRAD" | 1000 |

POLY INSTRUMENT REPAIR TAG non sterile

- Please contact Customer Service to have Customized Poly Instrument Repair Tags.

| | | |
|---------|--|-----|
| DEF1293 | PAPER INSTRUMENT REPAIR TAG WITH NYLON FASTENERS | 500 |
| DEF1468 | POLY INSTRUMENT REPAIR TAG WITH NYLON FASTENER | 100 |



| | | |
|-------------|--|------|
| DEF1469 | SECURI-TIE NON STERILE | 5000 |
| DEF1469-100 | SECURI-TIE NON STERILE | 100 |
| DEF1725 | POLY LOAN EQUIP PROCESSING TAG WITH NYLON FASTENER | 100 |
| DEF2010 | POLY TAG-INSTRUMENT TRAY 150 x 120 BIG TAG W/TIES | 100 |

| Code | Description | Box Qty. |
|--------------|---|----------|
| DEF2011 | POLY TAG-LOAN SET 50 x 90mm W/TIES | 100 |
| DEF2021 | BME REPAIR TAG WITH TIE - CABRINI | 100 |
| DEF2031 | POLY TAG - 32 x 60mm INSTRUMENT NAME W/TIES | 100 |
| DEF2065 | INSTRUMENT REPAIR TAG - BLANK WITH NYLON TIE | 100 |
| PT-109G | POLY TAG 50mm x 90mm - LOAN SET - GREEN | 100 |
| PT001 | FAST TRACK POLY TAG 5cm x 9cm - RED ON WHITE | 100 |
| PT013/BLUE | POLY TAGS - NO TEXT W/NYLON TIE - BLUE | 100 |
| PT013/GREEN | POLY TAGS - NO TEXT W/NYLON TIE - GREEN | 100 |
| PT013/ORANGE | POLY TAGS - NO TEXT W/NYLON TIE - ORANGE | 100 |
| PT013/PINK | POLY TAGS - NO TEXT W/NYLON TIE - PINK | 100 |
| PT013/RED | POLY TAGS - NO TEXT W/NYLON TIE - RED | 100 |
| PT013/VIOLET | POLY TAGS - NO TEXT W/NYLON TIE - VIOLET | 100 |
| PT013/YELLOW | POLY TAGS - NO TEXT W/NYLON TIE - YELLOW | 100 |
| PT015 | INSTRUMENT REPAIR TAG - POLY - TRAY OR SPECIALTY | 100 |
| PT076-D | PRIORITY 7cm x 9cm - RED ON WHITE | 100 |
| PT107 | FAST TRACK POLY TAG 5cm x 7cm | 100 |
| PT108 | INSTRUMENT REPAIR TAG 5cm x 7cm | 100 |
| PT153/01G | THEATRE NUMBER 1 - GREEN | 100 |
| PT153/02O | THEATRE NUMBER 2 - ORANGE | 100 |
| PT153/03B | POLY TAG W/NYLON TIE THEATRE NUMBER 3 - BLUE | 100 |
| PT153/04Y | POLY TAG W/NYLON TIE THEATRE NUMBER 4 - YELLOW | 100 |
| PT153/05P | POLY TAG W/NYLON TIE THEATRE NUMBER 5 - PINK | 100 |
| PT153/06R | POLY TAG W/NYLON TIE THEATRE NUMBER 6 - RED | 100 |
| PT153/07S | POLY TAG W/NYLON TIE THEATRE NUMBER 7 - GOLD | 100 |
| PT152/08 PUR | POLY TAG W/NYLON TIE THEATRE NUMBER 8 - PURPLE | 100 |
| PT166 5045MM | INSTRUMENT LABELS - NO TIES 50mm x 45mm WITH A HOLE | |

PERSPEX DISPENSER

| | | |
|---------|--------------------------|---|
| PD10 | HOLDS 10 BOXES OF LABELS | 1 |
| PD15 | HOLDS 15 BOXES OF LABELS | 1 |
| PD15-2T | HOLDS 30 BOXES OF LABELS | 1 |
| PD20 | HOLDS 20 BOXES OF LABELS | 1 |

> INTENSIVE CARE



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|----------|----------------------|---|
| DEF2161A | ARTERIAL TRAY | 5 |
|----------|----------------------|---|

1 x Polybacked Towel 40cm x 42cm with Fenestration 8cm x 8cm square
 1 x Dressing Tray 3 Cavity 25.5cm x 20.5cm x 2.5cm
 6 x Gauze Sponges 7.5cm x 7.5cm Plain
 2 x Huck Towels 48cm x 60cm
 1 x Scissors 12.5cm Sharp-Blunt Straight
 1 x Scissors Iris 11.5cm
 1 x Mayo Hagar Needle Holder 13cm
 2 x Forcep Artery Mosquito 12.5cm Straight
 2 x Forcep Artery Mosquito 12.5cm Curved
 1 x Iris Forcep Plain 10cm
 1 x Iris Forcep Toothed 10cm
 1 x Lacrimal Probe Size 1-2
 1 x Lacrimal Probe Size 3-4
 Field Clear Plastic 75cm x 75cm



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|----------|---------------------------|----|
| DEF2393A | CVC INSERTION PACK | 25 |
|----------|---------------------------|----|

1 x Scissors 12.5cm Sharp-Blunt Straight GH
 5 x Gauze Sponges 10cm x 10cm Plain
 2 x Huck Towel 48cm x 60cm
 1 x Kidney Dish Clear 700ml
 2 x Gallipot 150ml
 5 x Non Woven balls
 1 x Needle Holder 13cm GH
 1 x Forcep Adson Micro Toothed GH
 2 x Forceps Spencer Wells Straight 13cm GH
 Field Blue Plastic 91cm x 115cm



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| Code | Description | Box Qty. |
|----------|--------------------------|----------|
| DEF2766B | CENTRAL LINE TRAY | 10 |

- 1 x Dressing Tray 24cm x 15cm x 4cm
- 1 x Gown Unreinforced XL
- 1 x Central Line Drape 160cm x 245cm
- 2 x Dressing Towel - T Small 32cm x 32cm
- 1 x Iris Scissor Sharp Sharp with Green Handle
- 1 x Needle Holder 16cm Green Handle
- 1 x Forcep Adson Plain Green Handle
- 2 x Gallipot 150ml
- 10 x Green Gauze Sponges 10cm x 10cm
- 1 x Absorbant Foam Prep Applicator 20cm
- 1 x SMS Field 130cm x 150cm



| Code | Description | Box Qty. |
|---------|--|----------|
| DEF1609 | INTERCOSTAL CATHETER INSERTION TRAY | 24 |

- 2 x Kidney Dish Clear 700ml
- 2 x Gallipot 150ml
- 6 x Cotton balls Large
- 1 x Needle Holder 16cm
- 2 x Forceps Spencer Wells Straight 13cm
- 2 x Forceps Spencer Wells Curved 13cm
- 4 x Gauze Sponges 10cm x 10cm Plain
- 4 x Dressing Towel - T Small 32cm x 32cm
- 1 x Scissors 12.5cm Sharp-Sharp Straight
- 1 x Scalpel Handle #4
- 1 x Forcep Gillies Toothed
- Field SMS 91cm x 91cm

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| Code | Description | Box Qty. |
|---------|-------------------|----------|
| DEF2275 | TUBE CLAMP | 5 |

- | | | |
|---------|-----------------------|---|
| DEF2275 | 18cm GREEN HANDLE 2's | 5 |
|---------|-----------------------|---|



- | | | |
|---------|------------------------------|----|
| DEF1781 | 18cm GREEN HANDLE | 1 |
| DEF2775 | TUBE CLAMP 18cm GREEN HANDLE | 10 |

| Code | Description | Box Qty. |
|---------|--------------------------|----------|
| DEF2373 | TRACHEAL DILATORS | 1 |

- | | | |
|---------|---------|---|
| DEF2373 | 3 PRONG | 1 |
|---------|---------|---|



- | | | |
|---------|-------------------|---|
| DEF2585 | TROUSSEAU 2 PRONG | 5 |
| DEF1467 | TROUSSEAU 2 PRONG | 1 |



- | | | |
|---------|--|---|
| DEF2445 | TRACHY PROCEDURE PACK | 1 |
| | 1 x Tracheal Dilator Trousseau 2 prong | |
| | 1 x Retractor Berman 16.5cm | |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|----------|-------------------------------------|---|
| DEF2396A | POST OP TRACHY 3 PRONG PACK | 5 |
| | 1 x Kidney Dish Clear 700ml | |
| | 2 x Gallipot 150ml | |
| | 3 x Blue Huck Towel | |
| | 2 x Retractor Berman 16.5cm | |
| | 2 x Forcep Spencer Wells | |
| | 13cm Curved Green Handle | |
| | 1 x Tracheal Dilator 3 prong | |
| | 1 x Field 91cm x 115cm Plastic Blue | |

TRACHEOSTOMY SPONGE

| | | |
|---------|----------------------------------|-----|
| DEF2163 | 5cm x 5cm Y CUT NON WOVEN 2's | 100 |
| DEF2164 | 7.5 x 7.5cm Y CUT NON WOVEN 2's | 100 |
| DEF2165 | 10 x 10cm Y CUT NON WOVEN 2's | 100 |
| DEF803 | 10cm x 10cm 2 TAPES 12mm x 600mm | 100 |



| | | |
|--------|---|----|
| CR3900 | ADVAZORB T LITE ULTRA-THIN HYDROPHILLIC FOAM DRESSING 8cm x 8cm | 20 |
| CR4212 | ADVAZORB PLUS T STANDARD HYDROPHILLIC FOAM DRESSING 8cm x 8cm | 20 |



| | | |
|--------|-------------------|-----|
| DEF053 | PLAIN 10cm x 10cm | 100 |
|--------|-------------------|-----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

TRACHEOSTOMY TAPE

| | | |
|--------|-----------------|-----|
| DEF541 | 12mm x 90cm 2's | 100 |
| DEF066 | 12mm x 1m D/P | 25 |

> INSTRUMENTS



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

DENTAL INSTRUMENTS

| | | |
|---------|---------------------------------|-----|
| DEF1450 | 1 x MIRROR | 100 |
| DEF1457 | 1 x MIRROR 1 x PROBE 1 x FORCEP | 100 |



- Each instrument is engraved with Defries Industries and Single-Use, green dipped handles visually identifies single use
- Instruments are ultrasonically cleaned and each instrument is then checked at Defries by the Quality Assurance Team, packed and sterilized locally
- Single-use Instruments are a medical grade stainless steel, with an exacting mechanical finish
- Technical design of each instrument modelled on reusable instruments
- Benefit of no repeat sterilisation costs for hospitals
- By using single use, this removes the risk of cross infection
- Can be recycled into meltdown metal, or incinerated to become non-toxic file ash
- Instrument Repair Tags.

FORCEP ADSON

| | | |
|---------|---------------------------|-----|
| DEF2147 | PLAIN 12.5cm GREEN HANDLE | 10 |
| DEF2148 | PLAIN 12.5cm GREEN HANDLE | 100 |



| | | |
|---------|-----------------------------|----|
| DEF2142 | TOOTHED 12.5cm GREEN HANDLE | 10 |
|---------|-----------------------------|----|



| | | |
|---------|--|-----|
| DEF2143 | TOOTHED 12.5cm GREEN HANDLE | 100 |
| DEF2346 | TOOTHED SUPER FINE 12.5cm GREEN HANDLE | 10 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

FORCEP MICRO ADSON GREEN HANDLE

| | | |
|------------|----------------|-----|
| DEF2135-10 | PLAIN 12.5cm | 10 |
| DEF2135 | PLAIN 12.5cm | 100 |
| DEF2134-10 | TOOTHED 12.5cm | 10 |
| DEF2134 | TOOTHED 12.5cm | 100 |

**FORCEP ROCHESTER PEAN GREEN HANDLE**

| | | |
|---------|-------------|----|
| DEF2920 | 20cm CURVED | 10 |
|---------|-------------|----|



| | | |
|----------|---------------|---|
| DEF1782A | 22CM STRAIGHT | 1 |
|----------|---------------|---|

**FORCEP MOSQUITO ARTERY GREEN HANDLE**

| | | |
|------------|-----------------|-----|
| DEF2413-10 | STRAIGHT 12.5cm | 10 |
| DEF2413 | STRAIGHT 12.5cm | 100 |
| DEF2414-10 | CURVED 12.5cm | 10 |
| DEF2414 | CURVED 12.5cm | 100 |



| | | |
|---------|--------------|----|
| DEF2108 | MICRO CURVED | 10 |
|---------|--------------|----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

FORCEP SPENCER WELLS ARTERY GREEN HANDLE

| | | |
|------------|---------------|-----|
| DEF2277-10 | STRAIGHT 13cm | 10 |
| DEF2277 | STRAIGHT 13cm | 100 |



| | | |
|------------|-------------|----|
| DEF2415-10 | CURVED 13cm | 10 |
|------------|-------------|----|



| | | |
|------------|---------------|-----|
| DEF2415 | CURVED 13cm | 100 |
| DEF2417-10 | STRAIGHT 16cm | 10 |
| DEF2417 | STRAIGHT 16cm | 100 |



| | | |
|------------|-------------|-----|
| DEF2416-10 | CURVED 16cm | 10 |
| DEF2416 | CURVED 16cm | 100 |



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

FORCEP CRILE ARTERY GREEN HANDLE

| | | |
|---------|---------------|-----|
| DEF2441 | STRAIGHT 14cm | 10 |
| DEF1073 | STRAIGHT 14cm | 100 |



| | | |
|---------|-------------|----|
| DEF2442 | CURVED 14cm | 10 |
|---------|-------------|----|

**FORCEP DRESSING**

| | | |
|---------|-------------------------|-----|
| DEF2153 | 13cm BLUNT GREEN HANDLE | 100 |
|---------|-------------------------|-----|



| | | |
|---------|-------------------------|-----|
| DEF2152 | 13cm BLUNT GREEN HANDLE | 10 |
| DEF1072 | 16cm BLUNT GREEN HANDLE | 100 |



| | | |
|---------|---------------------------|-----|
| DEF1264 | 13cm POINTED GREEN HANDLE | 100 |
|---------|---------------------------|-----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

FORCEP GILLIES 15.5CM GREEN HANDLE

| | | |
|------------|-----------|-----|
| DEF2418-10 | 1/2 TOOTH | 10 |
| DEF2418 | 1/2 TOOTH | 100 |

**FORCEP MCINDOE 15.5CM**

| | | |
|---------|--------------------------|----|
| DEF2677 | NON TOOTHED GREEN HANDLE | 10 |
|---------|--------------------------|----|



| | | |
|---------|--------------------------|-----|
| DEF2948 | NON TOOTHED GREEN HANDLE | 100 |
|---------|--------------------------|-----|

FORCEP MCGILL

| | | |
|---------|-------------------------|----|
| DEF1807 | 20cm CHILD PLAIN HANDLE | 10 |
| DEF2083 | 25cm ADULT PLAIN HANDLE | 10 |

FORCEP ADDITIONAL

| | | |
|--------|----------------------|-----|
| DEF686 | FORCEP PLASTIC GREEN | 100 |
|--------|----------------------|-----|



| | | |
|---------|---------------------------------------|----|
| DEF1578 | WILDE FORCEP 15cm ANGLED PLAIN HANDLE | 10 |
|---------|---------------------------------------|----|



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|---------|--------------------------------------|----|
| DEF2578 | KELMAN-McPHERSON FORCEP GREEN HANDLE | 10 |
|---------|--------------------------------------|----|



| | | |
|---------|----------------------------------|----|
| DEF2443 | LISTER SINUS FORCEP GREEN HANDLE | 10 |
|---------|----------------------------------|----|



| | | |
|---------|--|---|
| DEF1579 | BOZMAN FORCEP STRAIGHT 25cm PLAIN HANDLE | 5 |
|---------|--|---|



| | | |
|---------|--|----|
| DEF2174 | SPLINTER FORCEP FEICHENFELDT 11cm PLAIN HANDLE | 10 |
|---------|--|----|

| | | |
|---------|---|----|
| DEF2750 | HARTMANN CROC-ACTION FORCEP 14cm GREEN HANDLE | 10 |
|---------|---|----|



| | | |
|---------|--------------------------------------|----|
| DEF3141 | FINE CROC-ACTION FORCEP GREEN HANDLE | 10 |
|---------|--------------------------------------|----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|---------|-------------------------------|----|
| DEF2695 | IRIS FORCEP 10CM PLAIN HANDLE | 10 |
|---------|-------------------------------|----|




| | | |
|---------|--------------------------------|----|
| DEF2683 | TILLEY EAR FORCEP PLAIN HANDLE | 10 |
|---------|--------------------------------|----|

**HOOKS**

| | | |
|---------|---------------------------------------|----|
| DEF1577 | GILLIES SKIN HOOK 16.5cm PLAIN HANDLE | 10 |
|---------|---------------------------------------|----|

| | | |
|--------|----------------------------------|----|
| DEF336 | NEURO FISH HOOKS 5 ON A TRAY D/P | 20 |
|--------|----------------------------------|----|



| | | |
|---------|--------------------------------|----|
| DEF1032 | NEURO FISH HOOKS SIZE: 3/0 5's | 20 |
|---------|--------------------------------|----|

| | | |
|--------|------------------------------|----|
| DEF405 | NEURO SKIN HOOKS 5 ON A TRAY | 20 |
|--------|------------------------------|----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

NASAL / EAR INSTRUMENTS

| | | |
|---------|---------------------------|---|
| DEF2240 | EAR SPECULUM GREEN HANDLE | 5 |
|---------|---------------------------|---|

NASAL SPECULUM THUDICHUM PLAIN HANDLE

| | | |
|---------|----------------------------|----|
| DEF1356 | #1 - THUDICHUM 20mm x 70mm | 25 |
|---------|----------------------------|----|

| | | |
|---------|----------------------------|----|
| DEF1357 | #2 - THUDICHUM 22mm x 70mm | 25 |
|---------|----------------------------|----|

| | | |
|---------|----------------------------|----|
| DEF1358 | #3 - THUDICHUM 25mm x 70mm | 25 |
|---------|----------------------------|----|

**NASAL SPECULUM MOD WEIN PLAIN HANDLE**

| | | |
|---------|-------|----|
| DEF2292 | SMALL | 10 |
|---------|-------|----|



| | | |
|---------|------|----|
| DEF2293 | MED. | 10 |
|---------|------|----|

| | | |
|---------|-------|----|
| DEF2294 | LARGE | 10 |
|---------|-------|----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

NEEDLE HOLDERS

| | | |
|---------|------|----|
| DEF2149 | 13cm | 10 |
|---------|------|----|



| | | |
|---------|------|-----|
| DEF2136 | 13cm | 100 |
|---------|------|-----|

| | | |
|---------|------|----|
| DEF2150 | 16cm | 10 |
|---------|------|----|



| | | |
|---------|------|-----|
| DEF2151 | 16cm | 100 |
|---------|------|-----|

| | | |
|---------|---------|----|
| DEF2707 | NIEVERT | 10 |
|---------|---------|----|

PROBE

| | | |
|---------|-------------------------------|----|
| DEF1580 | 14cm STAINLESS STEEL WITH EYE | 20 |
|---------|-------------------------------|----|

RETRACTORS

| | | |
|---------|------------------------|----|
| DEF2402 | CATSPAW RETRACTORS 2's | 10 |
|---------|------------------------|----|



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

STAPLE REMOVERS

| | | |
|---------|---------------------------|----|
| DEF1826 | PLASTIC HANDLE DISPOSABLE | 50 |
|---------|---------------------------|----|

**SCISSORS DRESSING SHARP/BLUNT 12.5CM**

| | | |
|---------|-------------|----|
| DEF1871 | SHARP/BLUNT | 10 |
|---------|-------------|----|



| | | |
|---------|-------------|-----|
| DEF1870 | SHARP/BLUNT | 100 |
|---------|-------------|-----|

non sterile

| | | |
|---------|-------------|----|
| DEF1280 | SHARP/BLUNT | 12 |
|---------|-------------|----|

SCISSORS SHARP/SHARP 12.5CM

| | | |
|---------|-------------|----|
| DEF1311 | SHARP/SHARP | 10 |
|---------|-------------|----|

| | | |
|---------|-------------|-----|
| DEF2276 | SHARP/SHARP | 100 |
|---------|-------------|-----|

**SCISSORS BLUNT/BLUNT 12.5cm**

| | | |
|---------|--------------------|-----|
| DEF1637 | BLUNT/BLUNT 12.5cm | 100 |
|---------|--------------------|-----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SCISSORS IRIS

| | | |
|---------|--------------------|----|
| DEF2144 | 7.5cm SHARP/SHARP | 10 |
| DEF2141 | 11.5cm SHARP/SHARP | 10 |



| | | |
|---------|---------------------------|-----|
| DEF2133 | 11.5cm SHARP/SHARP | 100 |
| DEF2460 | 11.5cm SHARP/SHARP CURVED | 10 |

SCISSORS MAYO 17CM

| | | |
|---------|----------|----|
| DEF2117 | CURVED | 10 |
| DEF2408 | STRAIGHT | 10 |

**SCISSORS METZENBAUM**

| | | |
|---------|-------------|----|
| DEF2511 | 16cm CURVED | 10 |
|---------|-------------|----|



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SCISSORS UNIVERSAL

| | | |
|---------|-------------------|----|
| DEF2851 | BLACK HANDLE 18CM | 25 |
|---------|-------------------|----|

**SCALPEL HANDLE**

| | | |
|---------|-----------------|----|
| DEF2145 | #3 GREEN HANDLE | 10 |
|---------|-----------------|----|



| | | |
|---------|-----------------|-----|
| DEF2137 | #3 GREEN HANDLE | 50 |
| DEF2146 | #3-GREEN HANDLE | 100 |

SAFETY SCALPEL AND HANDLE

| | | |
|---------|-----|-----|
| DEF2700 | #10 | 100 |
| DEF2701 | #11 | 100 |



| | | |
|---------|-----|-----|
| DEF2702 | #15 | 100 |
|---------|-----|-----|

SPONGE HOLDERS

| | | |
|---------|--|----|
| DEF2609 | RAMPLEY SPONGE HOLDER 24cm PLAIN HANDLE STAINLESS STEEL | 10 |
| DEF1923 | PLASTIC 19cm BLUE | 10 |
| DEF2071 | PLASTIC 24cm BLUE | 10 |



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

SUCKERS GREEN HANDLE

| | | |
|---------|---------------------|----|
| DEF2329 | YANKUER SUCKER 30CM | 10 |
|---------|---------------------|----|



| | | |
|---------|------------------|----|
| DEF2554 | POOL SUCKER 25CM | 10 |
|---------|------------------|----|

| | | |
|---------|------------------------------|----|
| DEF2498 | PEDIATRIC SUCKER MEDIUM 29CM | 10 |
|---------|------------------------------|----|



| | | |
|---------|-----------------------------|----|
| DEF2497 | PEDIATRIC SUCKER SMALL 27CM | 10 |
|---------|-----------------------------|----|

**FRAZIER SUCKER 17CM**

| | | |
|---------|-----|----|
| DEF2325 | 6FR | 10 |
|---------|-----|----|

| | | |
|---------|-----|----|
| DEF2326 | 8FR | 10 |
|---------|-----|----|

| | | |
|---------|------|----|
| DEF2327 | 10FR | 10 |
|---------|------|----|

| | | |
|---------|------|----|
| DEF2328 | 12FR | 10 |
|---------|------|----|



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

TUBE CLAMPS

| | | |
|---------|------|----|
| DEF2298 | 13CM | 10 |
|---------|------|----|

| | | |
|---------|-----------------|-----|
| DEF2304 | 13CM SMOOTH JAW | 100 |
|---------|-----------------|-----|

| | | |
|---------|------|---|
| DEF1781 | 18CM | 1 |
|---------|------|---|


| | | |
|---------|-----------------|----|
| DEF2775 | TUBE CLAMP 18cm | 10 |
|---------|-----------------|----|



| | | |
|---------|-----------------|---|
| DEF2275 | TUBE CLAMP 18cm | 5 |
|---------|-----------------|---|



| | | |
|---------|-----------------------------|----|
| DEF1575 | WITH LOCK BLUE PLASTIC 12cm | 50 |
|---------|-----------------------------|----|

**WIRE CUTTER**

| | | |
|---------|----------------|----|
| DEF1800 | UNIVERSAL 12cm | 50 |
|---------|----------------|----|

> INSTRUMENT PACKS



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

DISPOSABLE SUTURE SET

| | | |
|----------|---------------------------------------|----|
| DEF1328A | PLASTIC SET WITH GREEN HANDLES | 10 |
|----------|---------------------------------------|----|

1 x Iris Scissor 11.5cm
1 x McIndoe Forcep 15.5cm Plain

| | | |
|----------|--|----|
| DEF 1486 | SUTURE REMOVAL SET WITH GREEN HANDLES | 10 |
|----------|--|----|

1 x Forcep Adson Toothed 13cm
1 x Scissor Iris 11.5cm

| | | |
|----------|---|----|
| DEF 1656 | FINE SUTURE SET WITH GREEN HANDLES | 25 |
|----------|---|----|

1 x Scissor Iris 11.5cm
1 x Forcep Adson Micro Plain

| | | |
|---------|--|----|
| DEF1560 | SUTURE SET SMALL WITH GREEN HANDLES | 10 |
|---------|--|----|

1 x Adson Toothed Forcep
1 x Iris Scissor 11.5cm Sharp/Sharp
1 x Needle Holder 13cm



| | | |
|---------|--|----|
| DEF1561 | SUTURE SET LARGE WITH GREEN HANDLES | 10 |
|---------|--|----|

1 x Adson Toothed Forcep
1 x Needle Holder 16cm
1 x Iris Scissor 11.5cm Sharp/Sharp



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|---------|---|----|
| DEF1494 | DISPOSABLE SUTURE SET WITH GREEN HANDLES | 25 |
|---------|---|----|

1 x Iris Scissor 11.5cm Sharp/Sharp
1 x Forcep Blunt 13cm
1 x Needle Holder 16cm

| | | |
|-------------|--------------------------------------|----|
| DEF 1339-10 | SUTURE SET WITH GREEN HANDLES | 10 |
| DEF 1339-50 | | 50 |

1 x Scissor Sharp/Blunt Straight
1 x Forcep Adson Toothed 13cm
1 x Needle Holder 16cm

| | | |
|-------------|--|-----|
| DEF1636-10 | FINE FACIAL MICRO SUTURE SET WITH GREEN HANDLES | 10 |
| DEF1636-25 | | 25 |
| DEF1636-100 | | 100 |

1 x Forcep Adson Micro Toothed
1 x Iris Scissor 11.5cm Sharp/Sharp
1 x Needle Holder 13cm



| | | |
|----------|--|----|
| DEF1671A | SUTURE KIT WITH COTTON BALLS WITH GREEN HANDLES | 10 |
|----------|--|----|

1 x Dressing Tray - Long Catheter
23cm x 9cm x 3cm
3 x Gauze Sponges Plain 10cm x 10cm
1 x Green Disposable Forcep
1 x Forcep 13cm Blunt
1 x Forcep Adson Toothed 13cm
1 x Iris Scissor 11.5cm Sharp/Sharp
1 x Scalpel Handles #3
1 x Towels White 2 Ply 400mm x 420mm
1 x Dressing Tray Field 50cm x 50cm
5 x Cotton Balls Large
1 x Needle Holder 16cm

| | | |
|---------|--------------------------------------|----|
| DEF1811 | SUTURE KIT WITH GREEN HANDLES | 25 |
|---------|--------------------------------------|----|

1 x Scalpel Handle and no. 11 Blade
1 x Needle Holder 16cm

| Code | Description | Box Qty. |
|----------|--|----------|
| DEF2214 | SUTURE REMOVAL SET WITH GREEN HANDLES | 10 |
| | 1 x Forcep Adson Plain Micro | |
| | 1 x Iris Scissor 11.5cm Sharp/Sharp | |
| DEF2226 | FINE SUTURE SET WITH GREEN HANDLES | 10 |
| | 1 x Dressing Tray - Long Catheter 23cm x 9cm x 3cm | |
| | 10 x Gauze Sponges Plain 10cm x 10cm | |
| | 1 x Dressing Towel - T Small 32cm x 32cm | |
| | 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration & Adhesive Tabs in Drape Corners | |
| | 1 x Green Disposable Forcep | |
| | 1 x Needle Holder 13cm | |
| | 1 x Iris Scissor 11.5cm Sharp/Sharp | |
| | 1 x Forcep Adson Micro Toothed | |
| | 1 x Forcep Artery - Mosquito 12.5cm Curved | |
| | 1 x Syringe 5mL L/Slip | |
| | 1 x Hypodermic Needle 25g x 1" | |
| | Field Blue Plastic 75cm x 75cm | |
| DEF2209 | SUTURE INCISION SET WITH GREEN HANDLES | 10 |
| | 1 x Dressing Tray - Long Catheter 23cm x 9cm x 3cm | |
| | 2 x Galli Pot 60mL | |
| | 1 x Iris Scissor 11.5cm Sharp/Sharp | |
| | 1 x Scalpel Handles no. 3 | |
| | 1 x Forcep Adson Micro Toothed | |
| | 5 x Cotton Balls Large | |
| | 4 x Gauze Sponges Plain 10cm x 10cm | |
| | 1 x Dressing Towel - T Small 32cm x 32cm | |
| | 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration | |
| | 1 x Needle Holder Collier 12.5cm | |
| | 1 x Forcep 13cm Blunt | |
| | Field Blue Plastic 50cm x 50cm | |
| DEF2050A | SUTURE TRAY WITH GREEN HANDLES | 10 |
| | 1 x Dressing Tray - Long Catheter 23cm x 9cm x 3cm | |
| | 1 x Scissor 12.5cm SH/BL Straight | |
| | 5 x Gauze Sponges Plain 10cm x 10cm | |
| | 2 x Non Woven Swabs 7.5cm x 7.5cm | |
| | 1 x Towels Blue Cotton 40cm x 68cm | |
| | 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration | |
| | 1 x Forcep Adson Plain | |
| | 1 x Iris Scissor 11.5cm Sharp/Sharp | |
| | 1 x Needle Holder 13cm | |
| | Field Blue Plastic 50cm x 50cm | |

| Code | Description | Box Qty. |
|----------|--|----------|
| DEF2051 | FINE PLASTIC SUTURE TRAY WITH GREEN HANDLES | 10 |
| | 1 x Dressing Tray - Long Catheter 23cm x 9cm x 3cm | |
| | 1 x Needle Holder 13cm | |
| | 5 x Gauze Sponges Plain 10cm x 10cm | |
| | 2 x Non Woven Swabs 7.5cm x 7.5cm | |
| | 1 x Towels Blue Cotton 40cm x 68cm | |
| | 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration | |
| | 1 x Scalpel Handle #3 | |
| | 1 x Iris Scissor 11.5cm Sharp/Sharp | |
| | 1 x Forcep Adson Micro Plain | |
| | 1 x Forcep Adson Micro Toothed | |
| | Field Blue Plastic 50cm x 50cm | |
| DEF2179 | SUTURE KIT WITH GREEN HANDLES | 50 |
| | 1 x Dressing Tray - Long Catheter 23cm x 9cm x 3cm | |
| | 5 x Gauze Sponges Plain 10cm x 10cm | |
| | 1 x Disposable Green Forcep | |
| | 1 x Gillies Forcep Toothed | |
| | 1 x Iris Scissor 11.5cm Sharp/Sharp | |
| | 1 x Scissor 12.5cm SH/SH Straight | |
| | 1 x Needle Holder 16cm | |
| | 2 x Towels White 2 Ply 400mm x 420mm | |
| | 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration | |
| | Field Blue Plastic 75cm x 75cm | |
| DEF1913B | INSTRUMENT PACK WITH GREEN HANDLES | 10 |
| | 1 x Needle Holder 13cm | |
| | 1 x Forcep Adson Toothed | |
| | 1 x Forcep Adson Plain | |
| | 1 x Iris Scissor 11.5cm Sharp/Sharp | |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|---------|--|----|
| DEF1677 | INSTRUMENT DRESSING PACK WITH GREEN HANDLES | 60 |
|---------|--|----|

1 x Long Dressing Tray 23cm x 9cm x 3cm
 3 x Gauze Sponges Plain
 1 x Disposable Green Forcep
 1 x Forcep 13cm Blunt
 1 x Forcep 13cm Pointed
 1 x Scissor Sharp/Blunt
 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration
 Field Blue Plastic 75cm x 75cm



| | | |
|----------|---|----|
| DEF2453A | FINE PROCEDURE PACK WITH GREEN HANDLES | 30 |
|----------|---|----|

1 x Iris Scissor 11.5cm Sharp/Sharp
 1 x Needle Holder 13cm
 1 x Scalpel Handle # 3
 1 x Forcep Mosquito Artery Curved 12.5cm
 1 x Gallipot 60ml
 1 x Kidney Dish Clear 700ml
 1 x Forcep Adson Micro Toothed
 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration
 4 x Gauze Sponges Plain 7.5cm x 7.5cm
 Field Blue Plastic 75cm x 75cm

| | | |
|---------|--|----|
| DEF2486 | PHARMACY INSTRUMENT PACK WITH GREEN HANDLES | 10 |
|---------|--|----|

1 x Kidney Dish Clear 700ml
 1 x Forcep Crile 14cm
 1 x Forcep Blunt 13cm
 1 x Scissor Sharp/Blunt
 Field Blue Plastic 75cm x 75cm

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|---------|--|----|
| DEF2502 | MEDIUM DRESSING TRAY WITH GREEN HANDLES | 30 |
|---------|--|----|

1 x Towel Huck Blue
 1 x Kidney Dish Clear 700ML
 1 x Forcep Adson Plain
 1 x Gallipot 150ml
 1 x Forcep Spencer Wells 13cm Curved
 1 x Forcep Spencer Wells 13cm Straight
 14 x Gauze Sponges Plain 7.5cm x 7.5cm
 1 x Dressing Tray 23cm x 9cm x 3cm
 Field Blue Plastic 91cm x 115cm

| | | |
|----------|--|----|
| DEF1834A | THYROID PACK WITH PLAIN HANDLES | 20 |
|----------|--|----|

2 x Forcep Spencer Wells Straight 13cm
 10 x Gauze Sponges Plain 10cm x 10cm
 1 x Towel Blue Huck 40 cm x 68cm
 1 x Scissor Sharp/Blunt 12.5cm
 1 x Staple Remover Scissor Action with Plastic Handle
 Field SMS 61cm x 61cm



> MEDICAL IMAGING



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

PROTECTIVE GLASSES

| | | |
|---------|---------------------|---|
| DEF1755 | CLEAR POLYCARBONATE | 1 |
|---------|---------------------|---|

**SPLASH PROTECTION FOR EYES - SPEYES NON STERILE**

| | | |
|--------|-----------------------------------|----|
| 23-002 | FACE MASK FULL WITH ELASTIC STRAP | 50 |
|--------|-----------------------------------|----|



| | | |
|--------|--------------------------------|-----|
| 23-006 | HALF FACE MASK | 50 |
| 23-501 | EYE SHIELD REPLACEMENT LENS | 100 |
| 23-502 | EYE SHIELD FRAME ONLY (Speyes) | 25 |



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

COSMIC NEE NOGGS

| | | |
|---------|---|----|
| 23-0016 | FULL FACE SHIELD WITH VENTED FOAM & VELCRO STRAP 20CM | 50 |
|---------|---|----|



| | | |
|--------|---------------------------|-----|
| 23-505 | EYE SHIELDS PRE-ASSEMBLED | 100 |
|--------|---------------------------|-----|

**SURGICAL MARKING PEN**

| | | |
|---------|-------------------------------|----|
| DEF1888 | FINE POINT WITH 15cm/5" RULER | 25 |
| 26-001 | STANDARD WITH 15cm/5" RULER | 25 |




| | | |
|---------|---------------------|-----|
| DEF2669 | MINI PRE-OP MARKERS | 100 |
|---------|---------------------|-----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

DEFOGGER™

- Specially formulated solution that safely prevents lens fogging
- Will not damage scope
- Each bottle comes with non-abrasive, x-ray detectable foam pad
- Place a few drops on pad for fog-free surgery
- Adhesive-backed pad allows universal placement
- Sterile
- Latex Free

| | | |
|---------|-----------------------|----|
| 28-0101 | DEFOGGER™ OVAL BOTTLE | 20 |
|---------|-----------------------|----|

**GUIDEWIRE BOWL**

| | | |
|---------|----------------|----|
| DEF2579 | GUIDEWIRE BOWL | 10 |
|---------|----------------|----|

**HUCKL TOWEL**

| | | |
|---------|-----------------|----|
| DEF2018 | BLUE COTTON 5's | 40 |
|---------|-----------------|----|



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

HANDLE SLEEVE COVERS

| | | |
|---------|--|-----|
| DEF1102 | 14cm x 30cm WITH ELASTIC OPENING | 50 |
| DEF1053 | 14cm x 30cm WITH ADHESIVE ALONG TOP OPEN END OF BAG 2.5cm x 10cm | 100 |

**I.I. COVER**

| | | |
|---------|--|----|
| DEF1045 | 1 SECURING TAPE, CUFF & LABELS D/P | 15 |
| DEF1348 | 2 SECURING TAPES CUFF & LABELS STERILE D/P | 15 |
| DEF764 | CAP-SMALL 61cm x 91cm D/P | 40 |
| DEF767 | MEDIUM 91cm x 115cm | 50 |
| DEF774 | LARGE 120cm x 140cm | 40 |
| DEF770 | SMALL 76cm x 106cm D/P | 50 |



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

LEAD SCREEN COVER

| | | |
|---------|----------------------------------|----|
| DEF766T | WITH TAPES 57.5cm x 60cm | 50 |
| DEF1245 | 80cm x 80cm WITH ELASTIC OPENING | 50 |
| DEF1762 | 104cm x 105cm | 50 |

**ULTRASOUND PROBE COVER**

| | | |
|---------|--------------------------|-----|
| DEF643 | 23cm x 1.5m | 25 |
| DEF2409 | S/END 23cm x 1.5m | 100 |
| DEF1572 | SEALED END 17.5cm x 50cm | 100 |
| DEF1088 | SEALED END 23cm x 80cm | 100 |



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

DEF2485 X-RAY INSTRUMENT PACK WITH GREEN HANDLES 10

- 1 x Forcep Gillies Toothed
- 1 x Forcep Crile 14cm Straight
- 1 x Forcep Artery Mosquito 12.5cm Curved
- 1 x Forcep McIndoe 15.5cm Plain
- 1 x Scissor Sharp/Blunt
- 1 x Needle Holder 16cm

DEF2768 ANGIO INSTRUMENT SET WITH GREEN HANDLES 10

- 1 x Sponge Holder Plastic 24cm
- 1 x Scissor Sharp/Blunt
- 1 x Forcep Mosquito Curved 12.5cm
- 1 x Forcep Mosquito Straight 12.5cm



To order, please contact our
Customer Service Team on

1300 550 278

or see our website for details
www.defries.com.au

> MATERNITY PRODUCTS



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

BABY CUP WITH LID

| | | |
|---------|-------------------|-----|
| DEF2129 | BABY CUP WITH LID | 100 |
|---------|-------------------|-----|



BABY TEATS

| | | |
|---------|-----------------------|-----|
| DEF1059 | TEAT CHERRY SLOW FLOW | 100 |
|---------|-----------------------|-----|




| | | |
|---------|--------------------------------|-----|
| DEF1349 | TEAT CHERRY SLOW FLOW AND RING | 100 |
|---------|--------------------------------|-----|



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

TEAT BAGS

| | | |
|--------|----------------------|-----|
| DEF299 | SIZE: #21 PACK OF 10 | 100 |
|--------|----------------------|-----|



non sterile

| | | |
|--------|--------------------------------|---|
| DEF410 | NIPPLE SHIELDS GENTLE FEED 2'S | 2 |
|--------|--------------------------------|---|



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

CORD CLAMP CUTTING DEVICES

| | | |
|---------|---|----|
| DEF2535 | UMBILICAL CORD CUTTING SCISSORS non sterile | 10 |
| DEF2168 | UMBILICAL CORD CLAMP CUTTER PLASTIC HANDLE | 10 |

**MATERNITY PADS**

| | | |
|--------|------------|----|
| DEF721 | 9cm x 30cm | 50 |
|--------|------------|----|



| | | |
|---------|---|-----|
| DEF1163 | 9cm x 30cm 5's non sterile | 50 |
| DEF700 | BULK 10cm x 30cm | 300 |
| DEF238 | 9cm x 30cm 10's PACKED IN CLIPSEAL BAG | 25 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

PLACENTA BASIN

| | | |
|---------|--|----|
| DEF2822 | PLACENTA BASIN 26cm x 18.5cm x 7cm DEEP | 12 |
|---------|--|----|



| | | |
|----------|--|----|
| DEF2131A | MATERNITY SUTURE SET WITH GREEN HANDLES | 10 |
|----------|--|----|

Contents:

- 1 x Long Dressing Tray 23cm x 9cm x 3cm
- 1 x Gallipot 60ml
- 1 x Needle Holder 16cm
- 1 x Forcep Gillies Toothed 15.5cm
- 1 x Forcep Mosquito Curved 12.5cm
- 1 x Forcep Mosquito Straight 12.5cm
- 1 x Ins Scissor 11.5cm
- 2 x Packing Gauze 10cm x 2metre Xray
- 1 x Field 75cm x 75cm

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|---------|---|----|
| DEF2118 | DELIVERY SUTURE SET WITH GREEN HANDLES | 25 |
|---------|---|----|

Contents

- 1 x Scissor Mayo 17cm Straight
- 1 x Needle Holder 16cm
- 2 x Crile Forcep curved 14cm
- 1 x Forcep Bonney 18cm 1/2 Tooth



| | | |
|---------|---------------------|----|
| DEF2752 | DELIVERY SET | 15 |
|---------|---------------------|----|

Contents

- 2 x Forcep Rochester Pean 22cm Straight Plain Handle
- 1 x Mayo Scissor 17cm Curved Green Handle
- 2 x Cord Clamps
- 1 x Baby Receiver
- 1 x Maternity Pad
- 20 x Non Woven Swabs 10cm x 10cm
- 3 x Towels White 2 Ply
- 1 x Kidney Dish Clear 700ml
- 1 x Bowl 500ml
- 1 x Drape 76cm x 106cm SMS
- Field SMS 91cm x 91cm

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|----------|---|----|
| DEF2671A | EPISIOTOMY SUTURE SET D/P WITH GREEN HANDLES | 10 |
|----------|---|----|

- 1 x Needle Holder 16cm
- 1 x Forcep Gillies Toothed 15.5cm
- 1 x Iris Scissor 11.5cm

| | | |
|----------|------------------------|----|
| DEF2130A | DELIVERY BUNDLE | 10 |
|----------|------------------------|----|

- 2 x Rochester Pean Forcep 22cm Straight
- 1 x Mayo Scissor Straight 17cm Green Handle
- 1 x Mayo Scissor Curved 17cm Green Handle
- 1 x Dressing Towel - T Small 32cm x 32cm
- Field SMS 91cm x 91cm

| | | |
|---------|---------------------------------------|----|
| DEF2526 | LABOUR PACK WITH PLAIN HANDLES | 10 |
|---------|---------------------------------------|----|

- 2 x Rochester Pean Forcep 22cm Straight
- 1 x Mayo Scissor Straight 23cm Green Handle
- 1 x Bowl 2.5 Litre
- 1 x Bowl 600ml
- 1 x Gallipot 150ml
- Field SMS 91cm x 91cm

| | | |
|----------|--|----|
| DEF2670A | BIRTHROOM BUNDLE WITH PLAIN HANDLES | 10 |
|----------|--|----|

- 1 x Placenta Basin 26cm x 18.5cm x 7cm Deep
- 1 x Bowl 500ml
- 1 x Scissors Sharp/Blunt 13cm
- 1 x Rochester Pean Forcep 22cm Straight
- Field SMS 91cm x 91cm

CLEAR PLASTIC EPIDURAL DRAPES WITH 10CM X 15CM FENESTRATION AND 2.5CM ADHESIVE AROUND FENESTRATION

| | | |
|---------|--------------|-----|
| DEF1481 | 76cm x 106cm | 100 |
|---------|--------------|-----|



| | | |
|------------|---------------|----|
| DEF748FA-S | 120cm x 140cm | 10 |
|------------|---------------|----|

| | | |
|---------|---------------|-----|
| DEF3394 | 120cm x 140cm | 100 |
|---------|---------------|-----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

I.V. PRODUCTS

- Unique application - Allows you to apply the pressure pad without contaminating the pad or adhesive. Hypoallergenic acrylic adhesive. Latex Free.
- Adequate compression by 3mm thick pad.

| | | |
|--------|------------|-----|
| DEF929 | 36mm WHITE | 100 |
|--------|------------|-----|

**GYNAEOLOGICAL SINGLE USE INSTRUMENTS**

| | | |
|---------|--|----|
| DEF2248 | ROCHESTER PEAN FORCEP CURVED SIZE: 20cm | 10 |
|---------|--|----|



| | | |
|----------|---|----|
| DEF1782A | ROCHESTER PEAN FORCEPS 22CM STRAIGHT | 10 |
|----------|---|----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|---------|------------------------------------|----|
| DEF2400 | SIMMS VAGINAL SPECULUM SIZE: LARGE | 20 |
|---------|------------------------------------|----|



| | | |
|---------|-----------------------------|----|
| DEF2399 | SIMMS SPECULUM SIZE: MEDIUM | 20 |
|---------|-----------------------------|----|

| | | |
|---------|----------------------------|----|
| DEF2398 | SIMMS SPECULUM SIZE: SMALL | 20 |
|---------|----------------------------|----|

VAGINAL SPECULUMS

- Built-in battery operated LED light source capable of producing peak light for 20 minutes.
- LED light source generates NO heat.
- Single-handled ratcheted adjusters for both vertical and blade opening.
- Comfortable ergonomic design for improved balance and easy maneuverability during exams.
- Batteries used in ER-SPEC light source are landfill friendly and require no special disposal instructions.
- For use in Hospital Emergency, IVF centre's, Women's Clinic.

Sterile Speculums with Light Source

| | | |
|---------|--------|----|
| DEF3140 | SMALL | 18 |
| DEF3139 | MEDIUM | 18 |
| DEF3138 | LARGE | 16 |

Non Sterile Speculums with Light Source

| | | |
|----------------|--------|----|
| ER-SPEC-SMALL | SMALL | 18 |
| ER-SPEC-MEDIUM | MEDIUM | 18 |
| ER-SPEC-LARGE | LARGE | 16 |

Sterile Speculums without Light Source

| | | |
|---------|--------|----|
| DEF3157 | SMALL | 20 |
| DEF3158 | MEDIUM | 20 |
| DEF3159 | LARGE | 20 |

Non Sterile Speculums without Light Source

| | | |
|-----------|--------|----|
| DEF3157NS | SMALL | 20 |
| DEF3158NS | MEDIUM | 20 |
| DEF3159NS | LARGE | 20 |

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

IUD INSERTION KIT SINGLE-USE & STERILE

- Kit provided sterile, easy to use with quality consumables and medical grade stainless steel instruments. Saves setup time and sterilisation time for all IUD Insertions.
- Time saving on ordering and ensures all components are available for every insertion. Reduction on packaging waste.
- Benefit of no repeat sterilisation costs and by using single use, this removes the risk of cross infection.

| | | |
|---------|--------------------------|----|
| DEF2752 | IUD INSERTION KIT | 10 |
|---------|--------------------------|----|

Contents

2 x Swab Sticks 20 cm Paper Shaft

2 x Gallipots 150ml

1 x Kidney Dish 700ml clear

1 x Vaginal Speculum Large

1 x Mayo Scissor 23cm Straight with Green Handles

1 x Tenaculum Forcep 25cm with Green Handles

1 x Sponge Holder 19cm Plastic

1 x Uterine Sound (Hysterometer CH14) 4.67mm with Nozzle

Wrapped in a blue field 75cm x 75cm



To order, please contact our
Customer Service Team on

1300 550 278

or see our website for details

www.defries.com.au



> ENT PRODUCTS



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

DENTAL PELLETS

| | | |
|---------|------------------|-----|
| DEF1180 | SIZE: #0 5's D/P | 100 |
|---------|------------------|-----|

**DENTAL PRESSURE PACK**

| | | |
|--------|------------------|-----|
| DEF799 | WITH COTTON TAIL | 100 |
|--------|------------------|-----|



| | | |
|---------|---------------|-----|
| DEF1063 | WITH TAPE 2's | 100 |
|---------|---------------|-----|



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

| | | |
|--------|-----------------|----|
| DEF144 | X-RAY SMALL 2's | 50 |
|--------|-----------------|----|



| | | |
|--------|-----------------|----|
| DEF145 | X-RAY LARGE 2's | 50 |
|--------|-----------------|----|

DENTAL ROLLS

| | | |
|---------|-----|-----|
| DEF1209 | 5's | 100 |
|---------|-----|-----|

EYE PADS

| | | |
|--------|----------|-----|
| DEF069 | EYE PADS | 100 |
|--------|----------|-----|

**EYE SPEARS****Highly Absorbant**

| | | |
|--------|-----|-----|
| DEF484 | 2's | 100 |
|--------|-----|-----|



| | | |
|--------|-----|----|
| DEF485 | 5's | 20 |
|--------|-----|----|


| | | |
|--------|-----|----|
| DEF486 | 6's | 20 |
|--------|-----|----|

| | | |
|--------|------|----|
| DEF487 | 10's | 20 |
|--------|------|----|

| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

EYE STROLLS**Highly Absorbant**

| | | |
|--------|-----|----|
| DEF112 | 5's | 50 |
|--------|-----|----|



| | | |
|--------|------|----|
| DEF081 | 10's | 50 |
|--------|------|----|

| | | |
|--------|------|-----|
| DEF704 | BULK | 500 |
|--------|------|-----|

TOOTHBRUSH

| | | |
|--------|------------|----|
| DEF690 | TOOTHBRUSH | 25 |
|--------|------------|----|



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|



> INDEX



A

| | |
|--|-------|
| ABDOMINAL SPONGES | 45-46 |
| ANAESTHETIC/REGIONAL/EPIDURAL PACKS | 28-31 |
| ANGIO INSTRUMENT SET - DISPOSABLE STAINLESS STEEL | 94 |
| ARTERIAL LINE TRAY | 76 |

B

| | |
|--|--------|
| BABY CUP WITH LID | 96 |
| BAGS: | |
| CLIP SEAL | 52 |
| DONOR ORGAN | 52 |
| DRAWSTRING | 52 |
| FREEZER | 52 |
| PLASTIC | 53 |
| SPECIMEN | 53 |
| BANDAGES: | |
| CREPE | 49 |
| COHESIVE | 49 |
| CONFORMING | 49 |
| H CONFORMING | 49 |
| UNDERCAST | 50 |
| VIC | 50 |
| BIRTHROOM BUNDLE | 98 |
| BITE BLOCK | 28 |
| BLADESAFE SCALPEL TRANSFER DEVICE | 23, 41 |
| BOWLS PLASTIC: | 24 |
| BONE | 24 |
| GUIDEWIRE | 24-93 |
| SPLASH | 24 |
| BOWL SETS PLASTIC | 25 |
| RECEIVER SET | 25 |
| BRONNY IV PACKS | 66 |
| BURNS PRODUCTS | 41 |
| BURNS PAD | 41 |
| BURNS PLASTIC BAG | 41 |
| GAUZE | 41 |

C

| | |
|----------------------------------|-------|
| CAMERA AND CORD COVERS | 52 |
| CAPS - BERET | 9 |
| CATHETER PACKS: | 36-39 |
| PLAIN | 37-39 |
| XRAY | 36-37 |
| CENTRAL LINE TRAY | 77 |
| CHAIR COVER | 52 |
| CLIP SEAL BAG | 52 |
| COMBINE DRESSING | 42 |
| CORD CLAMP CUTTING DEVICE | 97 |
| COTTON BALLS | 66 |
| COTTON TAPE | 55 |
| COTTON WOOL ROLL | 66 |

COVERS:

| | |
|-------------------------------|--------|
| CAMERA | 52 |
| CHAIR | 52 |
| CORD | 52 |
| DRILL SLEEVE | 52 |
| HANDLE SLEEVE | 53, 93 |
| IMAGE INTENSIFIER | 93 |
| LEAD SCREEN | 52, 94 |
| PROBE | 53 |
| PROTECTIVE | 53 |
| TROLLEY BAG | 53 |
| ULTRASOUND | 53 |
| CUPS MEASURING | 22 |
| CUSTOM PROCEDURE PACKS | 17 |
| CVC INSERTION PACK | 76 |

D

| | |
|--|-------|
| DELIVERY SUTURE SET | 98 |
| DELIVERY BUNDLE | 98 |
| DELIVERY SET | 98 |
| DEFOGGER | 93 |
| DENTAL INSTRUMENTS - PLASTIC | 80 |
| DENTAL PELLETS | 102 |
| DENTAL PRESSURE PACK | 102 |
| DENTAL ROLLS | 102 |
| DIATHERMY: | 41 |
| PENCIL | 41 |
| QUIVER | 42 |
| DONOR ORGAN BAG | 52 |
| DRAINS: | 58-59 |
| CAPILLARY ACTION DRAINS - PASSIVE WOUND DRAINAGE | 59 |
| EASY FLOW DRAINS | 59 |
| FLAT SILICONE DRAINS | 58 |
| FLAT WHITE | 59 |
| FLUTED SILICONE DRAINS HUBLESS | 58 |
| MULTILUMEN DRAINS X-RAY | 59 |
| ROUND | 58-59 |
| SILICONE ROUND DRAINS X-RAY | 59 |
| DRAPE CLIP | 47 |
| DRAPES EQUIPMENT: | 14 |
| BACK TABLE COVER | 14 |
| CLEAR PLASTIC | 14 |
| MAYO STAND COVERS | 14 |
| DRAPES PATIENT: | |
| CLEAR PLASTIC EPIDURAL | 15 |
| CLEAR PLASTIC PATIENT | 15 |
| DRAPE SHEETS | 15 |
| MINOR PROCEDURE | 14 |
| POLYBACKED | 14 |
| SPLITSHEET DRAPE | 15 |
| UNIVERSAL PACK | 15 |

DRAPES SPECIALITY:

| | |
|---|----|
| CARDIAC AND ANGIOGRAPHY | 16 |
| EENT | 16 |
| GENERAL SURGICAL | 16 |
| GYNAECOLOGY AND OBSTETRIC | 16 |
| OPHTHALMIC | 16 |
| ORTHOPAEDIC | 16 |
| UROLOGY | 16 |
| DRAWSTRING BAG | 52 |
| DRESSING PACK - LARGE | 67 |
| DRESSINGS: | 42 |
| AMPUTATION | 42 |
| COMBINE NON WOVEN | 42 |
| COMBINE ROLL NON WOVEN | 42 |
| DRILL SLEEVE COVER | 52 |
| DRUG LABELS - AUSTRALIAN STANDARDS | 73 |
| DRUG TRAY PLASTIC | 28 |

E

| | |
|--|-------|
| EAR SPECULUM | 83 |
| EPIDURAL/REGIONAL/ANAESTHETIC PACKS | 28-31 |
| EPIDURAL DRAPES | 98 |
| EYE PROTECTION: | 9, 92 |
| EYE SHIELD FRAME ONLY (SPeyes) | 9, 92 |
| EYE SHIELD REPLACEMENT LENS (SPeyes) | 9, 92 |
| EYE SHIELDS PRE-ASSEMBLED (Speyes) | 9, 92 |
| EYE PADS | 102 |
| EYE SPEARS | 102 |
| EYE STROLLS | 103 |

F

| | |
|---|-------|
| FACE SHIELD: | 9, 92 |
| FACE MASK FULL WITH ELASTIC STRAP | 9, 92 |
| FACE MASK HALF | 9, 92 |
| FINE PROCEDURE PACK - DISPOSABLE STAINLESS STEEL | 90 |
| FOAM PADS | 42 |
| FORCEP: | 43 |
| ADSON | 80 |
| BOZMAN | 82 |
| CRILE ARTERY GREEN HANDLE | 81 |
| DRESSING | 81 |
| GILLIES 15.5CM | 82 |
| HARTMANN CROC-ACTION | 82 |
| IRIS FORCEP PLAIN HANDLE | 83 |
| KELMAN-McPHERSON | 82 |
| LISTER SINUS | 82 |
| MCGILL | 82 |
| MCINDOE 15.5CM | 82 |
| MICRO ADSON GREEN HANDLE | 80 |
| MOSQUITO ARTERY GREEN HANDLE | 80 |
| PLASTIC GREEN | 66 |
| ROCHESTER PEAN GREEN HANDLE | 80 |

| | |
|-----------------------------------|----|
| SPENCER WELLS ARTERY GREEN HANDLE | 81 |
| SPLINTER FEICHENFELDT | 82 |
| TILLEY EAR | 83 |
| WILDE | 82 |
| FRAZIER SUCKER | 86 |
| FREEZER BAG | 52 |

G

| | |
|--|-------|
| GALLIPOTS | 22 |
| GAUZE SWABS/SPONGES: | |
| BURNS | 41 |
| GREEN | 28 |
| FOLDED | 33 |
| PEANUT | 46 |
| PLAIN | 69 |
| RECOVERY | 46 |
| SPECIALIZED | 46 |
| X-RAY | 46 |
| GAUZE SPONGES PLAIN | 69 |
| GLASSES SAFETY | 9, 92 |
| GUIDEWIRE BOWL | 93 |
| GYNAEOLOGICAL SINGLE USE INSTRUMENTS: | 100 |
| ROCHESTER PEAN FORCEP | 99 |
| SIMMS SPECULUM SIZE: LARGE | 99 |

H

| | |
|----------------------------|--------|
| HANDLE SLEEVE COVER | 93 |
| HOOKS: | 83 |
| GILLIES | 83 |
| NEURO FISH | 31, 83 |

I

| | |
|---|------------|
| I.I. COVER | 93 |
| IDENTILOOPS | 61 |
| IMAGE INTENSIFIER COVER | 93 |
| IUD INSERTION KIT | 100 |
| INSTRUMENT PACK - DISPOSABLE STAINLESS STEEL | 90 |
| INTERCOSTAL CATHETER INSERTION TRAY | 77 |
| IV INSERTION PACKS | 66 |
| I.V. PRESSURE PADS | 42, 66, 99 |
| JUGS PLASTIC | 23 |

K

| | |
|------------------------------|----|
| KIDNEY DISHES PLASTIC | 24 |
|------------------------------|----|

L

| | |
|-----------------------------|--------|
| LABELS - CUSTOMISED | 72 |
| LABELS - STERILE | 72 |
| LABELS - PIGGYBACKED | 72 |
| LABOUR PACK | 98 |
| LEAD SCREEN COVERS | 52, 94 |
| LIGHT HANDLE COVER | 43 |
| LUMBAR PUNCTURE PACK | 67 |

M

| | |
|--|-------|
| MARKING PENS | 43 |
| MATERNITY PADS | 97 |
| MATERNITY SUTURE SET | 97-99 |
| MEDIUM DRESSING TRAY - DISPOSABLE STAINLESS STEEL | 90 |
| MOUTH & THROAT SWAB STICKS PAPERSHAFT JUMBO | 67 |
| MOUTH TOILET PACK | 68 |
| MSU PACK | 68 |

N

| | |
|--------------------------|----|
| NASAL BOLSTER | 68 |
| NASAL SPECULUM THUDICHUM | 83 |
| NASAL SPECULUM MOD WEIN | 83 |
| NEEDLE CASE/COUNTERS | 44 |
| NEEDLE HOLDERS | 84 |
| NEEDLE HOLDER NIEVERT | 84 |
| NEURO: | |
| 1/2 YARD STRIPS | 33 |
| BOWL SET | 33 |
| COTTON WOOL | 33 |
| GAUZE | 33 |
| HEAD DRESSINGS | 33 |
| RECEIVER SET | 33 |
| SKIN HOOKS | 33 |
| SUCTION TUBING | 34 |
| NIPPLE SHIELDS | 96 |

O

| | |
|-----------|---|
| OVERSHOES | 9 |
|-----------|---|

P

| | |
|--|-------|
| PACKS - ABDOMINAL | 45-46 |
| PERSPEX DISPENSER | 74 |
| PHARMACY INSTRUMENT PACK - DISPOSABLE STAINLESS STEEL | 90 |
| PLACENTA BASIN | 97 |
| PLASTIC: | |
| BONE BOWL | 24 |
| BOWLS | 24 |
| JUGS | 23 |
| KIDNEY | 24 |
| SPATULA | 23 |
| TRAYS | 23 |
| PLASTIC BAGS | 53 |
| PLASTICS PACK | 67 |
| POLY INSTRUMENT REPAIR TAG | 74 |
| POST OP TRACHY 3 PRONG PACK | 78 |
| PREP APPLICATORS | 45 |
| PREP SWAB BALLS X/Ray | 45 |
| PROBE COVERS | 53 |
| PROBE | 84 |

| | |
|------------------|----|
| PROTECTIVE COVER | 53 |
| PV TRAY | 68 |

R

| | |
|-------------------------------------|-------|
| REGIONAL/EPIDURAL/ANAESTHETIC PACKS | 28-31 |
| RENAL PRODUCTS | 68 |
| RENAL DIALYSIS PACK | 68 |
| PRESSURE PAD - CHUSHAVEN & STEPTY | 68 |
| RETRACTORS - CATSPA W | 84 |
| RUBBER BANDS | 45 |

S

| | |
|--|-------|
| SAFETY PINS | 45 |
| SAFETY SCALPEL AND HANDLE | 69 |
| SANITARY BAG & DISPENSER | 69 |
| SCALPEL HOLDER | 23 |
| SCALPEL HANDLE | 85 |
| SCISSORS: | 84-85 |
| BLUNT/BLUNT | 84 |
| DRESSING SHARP/BLUNT | 84 |
| IRIS | 85 |
| MAYO 17CM | 85 |
| METZENBAUM | 85 |
| SHARP/SHARP | 84 |
| UNIVERSAL | 85 |
| SCRATCHIES TIP CLEANER | 47 |
| SILCLAMPS | 61 |
| SILICONE BULB 100ml/cc | 58 |
| SILICONE SNUGGERS CARDIAC | 63 |
| SILICONE TUBING METRIC | 63 |
| SILICONE TUBING INCH RANGE | 63 |
| SILICONE TUBING WITH CONNECTOR ATTACHED | 63 |
| SILICONE TUBING NON STERILE | 64 |
| SIL TAPE | 28 |
| SKIN MARKERS | 43 |
| SLEEVE COVERS PLASTIC | 52 |
| SPATULA PLASTIC | 23 |
| SPECIMEN JARS | 22 |
| SPECIMEN PLASTIC BAG | 53 |
| SPONGE HOLDERS | 85 |
| SPONGES: | 45-46 |
| ABDOMINAL | 45-46 |
| BURNS | 41 |
| GREEN | 28 |
| NEURO | 33 |
| PEANUT | 46 |
| PLAIN | 69 |
| RECOVERY | 46 |
| SPECIALIZED | 46 |
| X-RAY | 46 |
| SNUGGERS | 47 |

| | |
|--|-------|
| SPONGES NON WOVEN SWABS | 69 |
| STAPLE REMOVER PLASTIC HANDLE | 69 |
| STAPLE REMOVERS | 84 |
| STOCKINETTE | 15 |
| SUCKERS: | 86 |
| FRAZIER SUCKER | 86 |
| PEDIATRIC SUCKER | 86 |
| POOL SUCKER | 86 |
| YANKUER SUCKER | 86 |
| SURGICAL GOWNS: | 7 |
| CLASSIC PLUS COOLMAX SURGICAL GOWN | 8 |
| CLASSIC PLUS REINFORCED SURGICAL GOWN | 7 |
| CLASSIC STANDARD SURGICAL GOWN | 7 |
| SURGICAL MARKING PEN | 92 |
| SURGICAL TAPE - F | 55 |
| SURGICAL TAPE - H | 55 |
| SURGICAL TAPE NICHIMESH | 56 |
| SURGICAL RETRACTION TAPE R | 56 |
| SUTURE BOOTS SILICONE | 61 |
| SUTURE SET - DISPOSABLE STAINLESS STEEL | 88-90 |
| SUTURE SET - EPISIOTOMY | 98 |
| SUTURE SET MATERNITY | 98-99 |
| SWAB STICKS WITH WOODEN SHAFT | 69 |
| SWAB STICKS WITH PAPER SHAFT | 70 |
| SWAB STICK WITH PLASTIC SHAFT | 70 |
| SWAB STICK PROCTOSCOPE | 70 |

T

| | |
|--|-------|
| TAPE: | 55-56 |
| COTTON | 55 |
| NICHIMESH | 56 |
| RETRACTION R | 56 |
| SURGICAL TAPE - F | 55 |
| SURGICAL TAPE - H | 55 |
| TEATS | 96 |
| TEAT BAGS | 96 |
| THYROID PACK - DISPOSABLE STAINLESS STEEL | 90 |
| TIP CLEANER X-RAY | 47 |
| TOOTHBRUSH | 103 |
| TOURNIQUET SILICONE | 47 |
| TOWEL - HAND: | |
| TOWEL BLUE HUCK | 11 |
| TOWEL DRESSING | 11 |
| TOWEL GREEN HUCK | 11 |
| TOWEL PURE COTTON | 11 |
| TOWEL WHITE HUCK | 11 |
| TRAYS PLASTIC | 23 |
| TRACHEAL DILATORS | 77 |
| TRACHEOSTOMY PRODUCTS | 78 |
| TRACHEOSTOMY SPONGE | |

| | |
|------------------------------|--------|
| TRACHEOSTOMY TAPE | 78 |
| TRACHY PROCEDURE PACK | 77 |
| TROLLEY BAG COVER | 53 |
| TUBE CLAMP | 77 |
| TUBE CLAMPS | 47, 86 |
| TUBE HOLDER | 47 |
| TUBING CONNECTOR | 63 |

U

| | |
|--|--------|
| ULTRASOUND PROBE COVER | 53, 94 |
| UMBILICAL CORD CUTTING SCISSORS | 97 |

V

| | |
|-------------------------|----|
| VAGINAL SPECULUM | 99 |
|-------------------------|----|

W

| | |
|--------------------|----|
| WIRE CUTTER | 86 |
|--------------------|----|

X

| | |
|------------------------------|----|
| X-RAY INSTRUMENT PACK | 94 |
|------------------------------|----|



Hospitals Trust Defries

Custom Procedure Packs
Surgical Drapes & Gowns
Single Use Medical Consumables



Head Office

18-24 Villas Road
Dandenong South VIC 3175

t: 61 3 9706 3600

f: 61 3 9706 3619

freecall: 1300 550 278

e: sales@defries.com.au

www.defries.com.au

Manufacturing Site

17-20 & 22 Kempson Court
Keysborough VIC 3173