


> PROTECTIVE APPAREL



Code Description Box Qty.

**PROTECTIVE APPAREL - Surgical Gowns**


- All gowns are designed from a unique material which offers different levels of protection. All fabric's are soft, light weight, comfortable and low linting.
- Adjustable neckline with Velcro Strips. Soft edge Velcro sewn on outer neck to reduce irritation.
- Gown seams are ultrasonically bonded, providing additional strength.
- Gown ties are attached securely.
- Cuffs are manufactured with a soft high level quality cotton.
- All gowns tested successfully for fire retardant CPSC 16 CFR Part 1610.
- All surgical gowns pass ASTM F 1670.

CLASSIC STANDARD SURGICAL GOWN

- Soft 50gsm Gown wrapped with 2 x Hand Towels 32cm x 60cm Pur Cotton and in a SMS Field 91cm x 91cm.
- Suitable for majority of Surgical Procedures.


DEF5000	Medium	25
DEF5001	Large	25
DEF5002	XLarge	25
DEF5003	XXLarge	25
DEF5004	XLarge XLong	25

Code Description Box Qty.

CLASSIC PLUS REINFORCED SURGICAL GOWN

- 63gsm Gown wrapped with 2 x Hand Towels 32cm x 60cm Pur Cotton and in SMS Field 91cm x 91cm.
- These gowns are suitable for surgical cases when an increased level of protection is required.
- Reinforced panel on front of gown and sleeves.

DEF5005	Medium	25
DEF5006	Large	25
DEF5007	XLarge	25
DEF5008	XXLarge	25
DEF5009	XLarge XLong	25

COOLMAX SURGICAL GOWN

- This gown is wrapped with 2 x Hand Towels 32cm x 60cm Pur Cotton and SMS Field 91cm x 91cm.
- The gowns are a 64 gsm fabric and include an advanced trilaminate fabric meeting AAMI level 4 protection on front panel and sleeves whilst still offering a breathable and cool level of comfort.
- Coolmax gowns pass ASTM F 1671 on front panel and sleeve reinforcement.

DEF5011	Large	25
DEF5012	XLarge	25
DEF5013	XLarge XLong	25
DEF5014	XXLarge	25

Code Description Box Qty.

DISPOSABLE GLASSES

DEF1755 CLEAR (POLYCARBONATE) 1



FACIAL PROTECTION - NON STERILE

- Anti-fog both sides with wrap around frame.
- Protects from splashes and sprays from the top down.
- Fits over prescription glasses.

23-002 FACE MASK FULL WITH ELASTIC STRAP 50



23-006 HALF FACE MASK 50

23-501 EYE SHIELD REPLACEMENT LENS (SPeyes) 100

23-502 EYE SHIELD FRAME ONLY (SPeyes) 25



Code Description Box Qty.

23-503 NECK CORDS FOR EYE SHIELDS 25



23-505 EYE SHIELDS PRE-ASSEMBLED (SPeyes) 100



DISPOSABLE BERET CAPS

DEF1142 53cm (21") BLUE NON STERILE 250



DISPOSABLE OVERSHOES

DEF1145 NON SKID UNIVERSAL OVERSHOES POLYPROPOLYENE 100

