


Fine Surgical  
Accessories from  
PMT.

## SURGICAL SUCTION AND IRRIGATION

# MICROMAT™

The MicroMat™ is designed with special features that provide the surgeon with a clear field of view at the surgical site.



### *Benefits...*

#### **MicroMat™**

- The silicone based background material features a unique MicroRib surface.
- Allows tissue elevation for suturing purposes. The MicroRib's 1mm channel spacing aides in accurate measurement capabilities.
- Disposable, prepackaged sterile and can be trimmed to a desired size or shape.
- Available in blue or yellow backgrounds for improved contrast.
- Available with or without suction.

**PMT®**

Quality, Innovation, Customer Driven

# Product Information

# MICROMAT™

## *Description:*

**PMT MicroMat™ Microsurgical Background Material.**

**Model 3103-3G with suction, 3FR**


**Model 3103-5G with suction, 5FR**

**Model 3103-7G with suction, 7FR**

**Model 3105-G without suction, MicroRib design**

**Model 3105-S without suction, smooth surface**

Opaque blue color provides superior contrast  
(also available in yellow).



## *Our Commitment...*

### *Customer Satisfaction*

PMT is committed to providing you the highest quality products and the finest customer service. We are available to assist you Toll-Free 24 hours a day, 7 days a week, 365 days a year. If it's more convenient, you can even fax us your requirements.

All of our products are thoroughly tested and only the highest quality materials are used in their construction. Innovative thought and unique products have proven us to be a leader in the medical manufacturing industry. Continued commitment to customer support allows us to retain this position both now and in the future.

▶ **Local sales representatives are available to serve you.  
Contact PMT Corporation for more information.**

These products are protected by one or more of the following U.S. patents:  
4,533,352

PMT and MANKIND are registered trademarks and MICROMAT is a trademark of PMT Corporation  
© Copyright 2004 PMT Corporation. All Rights Reserved

**PMT®**  
*Quality, Innovation, Customer Driven*

PMT CORPORATION  
1500 Park Road  
Chanhassen, MN USA 55317

Call Toll Free... 24 Hours a Day/365 Days a Year  
**1-800-MANKIND®** (626-5463)  
Phone 952-470-0866 Fax 952-470-0865