



FLASH STERILISER module




features : process to patient tracking

components

The MEDITRAX® Flash Steriliser Module includes :

Qty	Description	Code
1	MEDITRAX® flash label applicator	MHG1400
1	flash steriliser module log book	FAM003
1	ink roller	IRT1400
5 rolls	flash cycle labels	FAM402
1000	flash instrument labels	FAM01
1	patient record of flash sterilised instrument	FAM004
1	steriliser/applicator ID labels	FAM009
400	flash integrator labels	FAM002

manual recording



accountability & traceability

The MEDITRAX® Flash Steriliser Module is a two step system, where the first step covers the traceability and logging system for items processed by Flash Sterilisers. The system tracks the staff , items, patient, cycle and steriliser information, with notes on quality control checks and whether cycle passed or failed.

The second step covers traceability from process-to-patient and patient-to-process, with accurate records cross referencing patient and steriliser details.

implementation

MEDITRAX® Flash Steriliser Module includes at no extra cost ongoing :-

- staff training
- auditing of system
- quality assurance consultation
- maintenance program

summary

The MEDITRAX® Flash Steriliser Module benefits

- Customised for Flash Steriliser
- Steriliser records reference patient
- Patient records reference steriliser
- Staff data linked to the above
- All materials supplied
- Easy to learn

The MEDITRAX® Flash Steriliser Log Book is manufactured from Sterisyn waterproof synthetic paper for essential record protection in harsh wet environments.

quick reference

MEDITRAX® Flash Steriliser Module includes ;



BATCH RECORD



ITEM RECORD



STERILISER
DATA



STAFF DATA
RECORDED



PATIENT RECORD

related MEDITRAX® products :










- Total Traceability System
- Batch Recall Plus II module
- Cold Sterilisation module
- Steris module





FLASH STERILISER module

consumables for Meditrax Flash Steriliser Module :

Code		Description	Packaging	Order Qty
MHG1400		MEDITRAX® Flash label applicator	1 each	
IRT1400		ink roller	1 each	
FAM002		flash integrator labels	400/box	
FAM01		flash steriliser instrument labels	1000/packet	
FAM402		flash cycle labels	1000 /roll 5 rolls/box	
FAM003		flash steriliser log book	1 each	
FAM004		patient record of flash sterilised instruments	50 sheets/pad 5 pads /packet	

T 1800 074 085 customer support

accountability & traceability

training


Austmel provides MEDITRAX® System Inservice at the point of implementation. Courses are run on a person-to-person basis ensuring rapid absorption. Advanced courses are designed to train the trainer where required. MEDITRAX® training systems comply with the AS/NZS 4187:2003 standards and documentation.

one-on-one MEDITRAX® customer support

Austmel Pty Ltd is able to offer one-on-one support during office hours on our 1800 074 085 service. This is a valuable service to healthcare personnel who require immediate technical support to use the system. This service is a complete one-on-one person to person level. Additionally IT Remote Access Tech-Support provides quick on-line diagnosis and repair of systems to clients requiring rapid response under maintenance contracts.

documentation

MEDITRAX® documentation is built on the philosophy of paper waste minimisation, taking into account policies and procedures relating to the processing of items being sterilised in healthcare facilities. All aspects of the complex requirements are clearly and easily followed.



summary

The MEDITRAX® Flash Steriliser Module benefits

- Customised for Flash Steriliser
- Steriliser records reference patient
- Patient records reference steriliser
- Staff data linked to the above
- All materials supplied
- Easy to learn

The MEDITRAX® Flash Steriliser Log Book is manufactured from Sterisyn waterproof synthetic paper for essential record protection in harsh wet environments.

quick reference

MEDITRAX® Flash Steriliser Module includes ;



background

Austmel Pty. Ltd is an Australian owned company established in 1974. Austmel have been directly involved in the supply and development of Sterility Assurance products for over 30 years. Our unique products and systems extend into both the Health Care and Industrials markets, including :

- MEDITRAX® Accountability & Traceability Systems
- Biological and Chemical Sterilization Indicators
- Bowie-Dick Tests
- Sterilization Packaging
- Quality Assurance Consultation
- TMI Electronic Data Loggers
- Surgical Disposables

related MEDITRAX® products :

- Total Traceability System
- Batch Recall Plus II module
- Cold Sterilisation module
- Steris module

