

Foetal Monitoring Belt



Bamford Foetal Monitor Belts are manufactured from superior quality lightweight elasticised fabric. Button holes evenly spaced every 1.5cm allow easy securing of the Foetal Monitor to the patient for contraction and heartbeat monitoring.

Length 1.5m
Width 5cm
Code: BHE50

Another quality product selected and stocked by:

WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489
e enquiries@bamford.co.nz www.bamford.co.nz

When performance counts

