





total traceability system
GAS PLASMA starter pack



total traceability system

The MEDITRAX® TTS Gas Plasma module offers all the same benefits of the Total Traceability System but with the use of synthetic non cellulose based labels and documentation with specific process indicators suitable for the Gas Plasma process.

Unique identification of all items processed with strong process control provides a high degree of accountability ensuring all items are processed according to policy and procedure.



features : plasma sterilization tracking

standard recording

components

The MEDITRAX® TTS Gas Plasma starter pack includes :

Qty	Description	Code
1	MEDITRAX® label applicator	MHG1800
1	non gassing ink roller	IRT/NG05
5 rolls	plasma batch labels	MHI104P
1000	MEDITRAX® plasma labels	TTS4P
250	MEDITRAX® plasma batch cycle cards	TTS202P
1	sterilizer ID labels	MHI113/L
1	applicator ID labels	MHI111/S
50	MEDITRAX® patient record sheets	TTS201
6	MEDITRAX® binders	TTS/S01
50	daily audit sheets	TTS/AM/11
50	weekly audit sheets	TTS/AM/05
50	batch clips	TTS/S20
1	clipboard for batch clips	TTS/AM/007
1	SteriTrax™ marker	SM/04

implementation

MEDITRAX® TTS Gas Plasma starter pack includes ongoing :-

- staff training
- auditing of system
- quality assurance consultation
- maintenance program at no extra cost.

summary

The MEDITRAX® TTS Gas Plasma starter pack benefits

- Full audit trail traceability
- Batch record data captured
- Patient record data captured
- Instrument data captured
- Sterilizer data captured
- Nonconformance monitored
- Staff data linked to the above
- Closed loop accountability
- All materials supplied
- Easy to learn

quick reference

MEDITRAX® TTS Gas Plasma starter pack module includes :



BATCH RECORD



ITEM RECORD



PATIENT RECORD



STAFF DATA RECORDED

related MEDITRAX® products :

- Archive module
- Total Traceability System
- TTS Non Conformance module
- TTS Multi Batching module



Austmel Pty Ltd, P PO Box 2345 Head Office Claude House T 1800 074 085 F 07 5597 5244 A 16 Jade Drive, Nerang QLD 4211, AUSTRALIA W www.austmel.com E info@austmel.com

Austmel Pty Ltd, P PO Box 2345 Head Office Claude House T 1800 074 085 F 07 5597 5244 A 16 Jade Drive, Nerang QLD 4211, AUSTRALIA W www.austmel.com E info@austmel.com



- Archive module
- Total Traceability System
- TTS Non Conformance module
- TTS Multi Batching module

related MEDITRAX® products :



total traceability system
GAS PLASMA starter pack

total traceability system

consumables for MEDITRAX® TTS Gas Plasma starter pack :

Code		Description	Packaging	Order Qty
MHG1800		MEDITRAX® label applicator	1 each	
IRT/NG05		non gassing ink roller	1 each	
873170		standard ink roller	1 each	
MHI04P		plasma batch labels	1000/roll	
TTS4P		MEDITRAX® plasma labels	1000/roll	
TTS202P		MEDITRAX® plasma batch cycle cards	250/packet	
TTS201		MEDITRAX® patient record sheets	5 pads/packet 50 sheets/pad	
TTS/S01		MEDITRAX® binders (specify Sterilizer numbers)	6/box	
TTS/AM/11		daily audit sheets	50/pad	
TTS/S20		batch clips	50/packet	
TTS/AM/007		clipboard for batch cards	1 each	
SM04		SteriTrax™ marker	25/packet	
TTS/AM/05		weekly audit sheets	50/packet	

training

Austmel provides MEDITRAX® System Inservice at the point of implementation. Courses are run on a person-to-person basis ensuring rapid absorption. Advanced courses are designed to train the trainer where required. MEDITRAX® training systems comply with the AS/NZS 4187 standards and documentation.

one-on-one MEDITRAX® customer support

Austmel Pty Ltd is able to offer one-on-one support during office hours on our 1800 074 085 service. This is a valuable service to healthcare personnel who require immediate technical support to use the system. This service is a complete one-on-one person to person level. Additionally IT Remote Access Tech-Support provides quick on-line diagnosis and repair of systems to clients requiring rapid response under maintenance contracts.

documentation

MEDITRAX® documentation is built on the philosophy of paper waste minimisation, taking into account policies and procedures relating to the processing of items being sterilized in healthcare facilities. All aspects of the complex requirements are clearly and easily followed.

T 1800 074 085 **customer support**



summary

The MEDITRAX® TTS Gas Plasma starter pack benefits

- Full audit trail traceability
- Batch record data captured
- Patient record data captured
- Instrument data captured
- Sterilizer data captured
- Nonconformance monitored
- Staff data linked to the above
- Closed loop accountability
- All materials supplied
- Easy to learn

quick reference

MEDITRAX® TTS Gas Plasma starter pack includes :



background

Austmel Pty. Ltd is an Australian owned company established in 1974. Austmel have been directly involved in the supply and development of Sterility Assurance products for over 30 years. Our unique products and systems extend into both the Health Care and Industrials markets, including:

- Meditrax® Accountability & Traceability Systems
- Biological and Chemical Sterilization Indicators
- Bowie-Dick Tests
- Sterilization Packaging
- Quality Assurance Consultation
- TMI Electronic Data Loggers
- Surgical Disposables