



UltraEnhancer®

for permanent makeup and medical tattooing

The first permanent makeup device to obtain both the CE and UL marks for medical equipment safety and quality.

 Listed to UL medical device standards for USA and Canada and CE medical device standards for Europe. 

1-year manufacturer's warranty.

Powerful enough for all micropigmentation procedures including nipple-areola tattoo, scar and hair camouflage and permanent makeup.

The equipment comes complete with:

Power Control Unit

- Hospital Grade medical casing
- International voltage switch
- Speed selection knob (0-8000 rpm)
- Manual or Foot Pedal Operation
- Hospital grade power cord

Precision Handpiece

- Accommodates all Permark® needles (12 sizes) and pigments (53 shades)
- Straight design for maximum visibility and comfort
- Reciprocating needle mechanism
- Nose cones for accurate depth control
- Can be sterilized by gas, autoclave or Sterrad®

Foot Pedal (Variable or fixed speed)



*Made with unsurpassed
craftsmanship in the USA*

PMT® / PERMARK®

1500 Park Road
Chanhassen, MN 55317
sales@pmtcorp.com

PMT® / PERMARK®
1-800-PERMARK www.permark.com

©2005 PMT/Permark. All Rights Reserved. PMT Corporation.