

**MEDI PLAST<sup>®</sup>**

# Products



Mediplast AB  
P.O. Box 9504, SE-200 39 Malmö, Sweden  
Visiting address: Stenyxegatan 25 C, 213 76 Malmö  
Phone: +46 40 671 23 00 Fax: +46 40 671 23 10



## METAL SUCTION TIPS

- Single use, for surgery and examination of the ear
- Sterile, packed in peel pack
- ETO sterilised
- Special angle on the tip for ease of sight during examination and operations with microscope
- Quick connector system – it is easy to change tip during ongoing examination or surgery

### New Metal Suction Tip Stapes

- Tube in tube for rigidity and fine tip end lumen
- For stapedotomy and stapedectomy procedures

Material: PP/ Medical grade stainless steel  
(virgin steel)  
MDD Rule: Rule 6  
MDD Class: IIa  
CE: Acc. existing directive



<b>Art No</b>	<b>Size mm</b>		<b>Ward box / Trsp box</b>	<b>Peel pack</b>
6066500040	0,4 x 70	lilac	60/900	1 pc
6066507040	0,4/0,7 x 77	lilac	10/140	1 pc
6066500070	0,7 x 70	black	60/900	1 pc
6066500100	1,0 x 70	brown	60/900	1 pc
6066500140	1,4 x 70	white	60/900	1 pc
6066500140BP	1,4 x 70	white	60/720	2 x 30 pcs
6066500208	2,0 x 70	green	60/900	1 pc
6066500208BP	2,0 x 70	green	60/720	2 x 30 pcs
6066500045	0,4 x 80	lilac	60/900	1 pc
6066500072	0,7 x 80	black	60/900	1 pc
6066500102	1,0 x 80	brown	60/900	1 pc
6066500142	1,4 x 80	white	60/900	1 pc
6066500142BP	1,4 x 80	white	60/720	2 x 30 pcs
6066500200	2,0 x 80	green	60/900	1 pc
6066500200BP	2,0 x 80	green	60/720	2 x 30 pcs
6066500250	2,5 x 80	orange	30/900	1 pc
6066500250BP	2,5 x 80	orange	30/720	2 x 30 pcs
6066500300	3,0 x 80	blue	30/900	1 pc
6066500300BP	3,0 x 80	blue	30/720	2 x 30 pcs
6066500400	4,0 x 80	yellow	30/900	1 pc
6066500400BP	4,0 x 80	yellow	30/720	2 x 30 pcs
6066509100	1,0 x 100	brown	30/900	1 pc
6066509140	1,4 x 100	white	30/900	1 pc
6066509140BP	1,4 x 100	white	30/720	2 x 30 pcs
6066509200	2,0 x 100	green	30/900	1 pc
6066509200BP	2,0 x 100	green	30/720	2 x 30 pcs
6066509250	2,5 x 100	orange	30/900	1 pc
6066509250BP	2,5 x 100	orange	30/720	2 x 30 pcs
6066509300	3,0 x 100	blue	30/900	1 pc
6066509300BP	3,0 x 100	blue	30/720	2 x 30 pcs

ward box = min qty / order

### Metal Suction Tips - BigPack

- Sterile BigPack for out patient clinical use. 2 x 30 pcs/pack



6066500208BP



## REGULATOR - SUCTION TIP HANDLE

- The Suction Tip Handle, with or without suction control, should be used as a connector between the Suction Tip and the Mediplast Suction Tube
- Single use, sterile, single packed in peel-pack
- White ABS-plastic

<u>Art.no.</u>		<u>ward.box/trsp.box</u>
6066500992	with control	75/1125
6066505009	without control	75/1125
6066500990	with control	25/25

Ward.box = min. qty/order

Material: ABS  
MDD Rule: Rule 1  
MDD Class: I sterile  
CE: Acc to regulation



## MEDIBENDER

A bending instrument for Mediplast Single Use  
Metal Suction Tips:  
2,0 mm - 1,4 mm - 1,0 mm diameters

- Stainless Steel
- Non sterile
- Can be autoclaved

How to use it:

1. Grip the instrument with the thumb on the Mediplast engraving and with the holes visible.
2. Take the wanted size of metal suction tip.
3. Put the metal suction tip through the appropriate sized hole:  
left = 2,0mm , middle = 1,4mm , right = 1,0mm diameters
4. Bend the metal suction tip gently in the length direction of the instrument to the angle wanted.

**Art No**

6066500

**ward box/trsp. box**

1/1

Material: Stainless Steel  
MDD Rule: Rule 1  
MDD Class: I non-sterile  
CE: Acc to regulation



## AURAL SPECULUMS

- New oval design of the cone and side cut end for improved line of sight
- New, wider, rimb on the grip side for better stability when hand held
- For examination procedures
- Fits Peter's or Siegler's pneumatic glass
- Hygienically packed
- Latex and PVC-free

<b>Art No</b>	<b>mm</b>	<b>Ward. Box/Trsp. Box</b>
606625	2,5	100/2800
606640	4	100/2800
606650	5	100/2800
606660	6	100/2800
606670	7	100/2800
606680	8	100/2800

Material: PP  
MDD Rule: Rule 5 - 1  
MDD Class: I  
CE: According to existing directives

ward. box = min qty / order



## MYRINGOTOMY SCALPEL

- Single use Myringotomy Scalpel, triple cut edges
- Used for Paracenticis
- Sterile, single packed in peel-pack

<b>Art.no.</b>	<b>mm</b>	<b>ward.box/trsp.box</b>
6066121370	1,3 x 70	30/450

ward.box = min.qty/order

Material: PP/ Medical Grade Stainless Steel  
(virgin steel)  
MDD Rule: Rule 6  
MDD Class: IIa  
CE: Acc to regulation



## MYRINGOTOMY KIT

- A ready made Kit for Paracentesis Surgery
- Single use
- Sterile, single packed in peel-pack



<b><u>Art.no.</u></b>	<b><u>ward.box/trsp.box</u></b>
<b>6066143300</b>	<b>30/180</b>

Contents:

- Metal Suction Tip
- 1,4 x 70 mm, 2,0 x 70 mm, 2,5 x 80 mm
- Myringotomy Scalpel
- Suction Tip Handle
- Aural Speculas (2,5 mm, 4 mm, 5 mm)
- 2 Non woven gauzes
- Sterile field

<b>6066143400</b>	<b>30/180</b>
-------------------	---------------

Contents:

- Metal Suction Tip
- 1,4 x 70 mm, 2,0 x 70 mm, 4,0 x 80 mm
- Myringotomy Scalpel
- Suction Tip Handle
- Aural Speculas (5 mm, 6 mm, 7 mm)
- 2 Non woven gauzes
- Sterile field

ward.box = min.qty/order

Material: PP/ ABS/ Medical Grade Stainless Steel  
(virgin steel)  
MDD Rule: Rule 6  
MDD Class: IIa  
CE: Acc to regulation



## MIDDLE EAR SURGERY KIT

- A ready made Kit for Middle Ear Surgery
- Single use, sterile, single packed in peel-pack

<b>Art.no.</b>	<b>ward.box/trsp.box</b>
----------------	--------------------------

6066502070	30/30
------------	-------

Contents:

Metal Suction Tip

1,0 x 70 mm (2 pcs)

1,4 x 70 mm (2 pcs)

2,0 x 80 mm (2 pcs)

Suction Tip Handle

2 m Medioplast Suction Tube

6066532070	30/30
------------	-------

Contents:

Metal Suction Tip

1,0 x 70 mm (2 pcs)

1,4 x 70 mm (2 pcs)

2,0 x 80 mm (2 pcs)

Suction Tip Handle

3,5 m Medioplast Suction Tube

ward.box = min.qty/order

Material: PP/ PVC/ Medical Grade Stainless Steel  
(virgin steel)  
MDD Rule: Rule 6  
MDD Class: IIa  
CE: Acc to regulation



## NASAL SUCTION TIP

- For surgical operations and rhinological examinations
- Single use, sterile, single packed in peel-pack
- Medical grade stainless steel, luer slip handle (PP)

<b>Art.no.</b>	<b>mm</b>		<b>ward.box/trsp.box</b>
6066516020	2,0 x 150	green	40/600
6066516030	3,0 x 150	blue	40/600
6066516040	4,0 x 150	yellow	40/600

ward.box = min.qty/order

Material: PP/ Medical Grade Stainless Steel (virgin steel)  
MDD Rule: Rule 6  
MDD Class: IIa  
CE: Acc to regulation



## NASAL SUCTION TIP WITH TRANSPARENT HANDLE

- For surgical operations and rhinological examinations
- Single use, sterile, single packed in peel-pack
- Medical grade stainless steel, transparent handle (PMMA)

<b>Art.no.</b>	<b>mm</b>	<b>ward.box/trsp.box</b>
6066516120	2,0 x 150	25/375
6066516130	3,0 x 150	25/375
6066516140	4,0 x 150	25/375

ward.box = min.qty/order

Material: PMMA/ Medical Grade Stainless Steel  
(virgin steel)  
MDD Rule: Rule 6  
MDD Class: IIa  
CE: Acc to regulation



## NASAL SUCTION TIP TYPE FERG

- For surgical operations and rhinological examinations
- Single use, sterile, single packed in peel-pack
- Medical grade stainless steel, transparent handle (PMMA)

<b>Art.no.</b>	<b>mm</b>	<b>ward.box/trsp.box</b>
6066501020	2,0 x 150	25/375
6066501030	3,0 x 150	25/375
6066501040	4,0 x 150	25/375

ward.box = min.qty/order

Material: PMMA/ Medical Grade Stainless Steel  
(virgin steel)  
MDD Rule: Rule 6  
MDD Class: IIa  
CE: Acc to regulation



## NASAL SUCTION TIP TYPE FERG WITH SIDEHOLE

- For delicate surgery procedures in Rhinology
- 1 side hole at tip end, to avoid biopsy effect (tissue damage)
- Single use, sterile, packed in peel-pack
- Medical grade stainless steel, transparent handle of PMMA with built in vacuum control

<b>Art.No.</b>	<b>mm</b>	<b>Ward box/trsp.box</b>
6066511020	2,0 x 150	25/375
6066511030	3,0 x 150	25/375
6066511040	4,0 x 150	25/375

ward box. = minimum order qty

Material: PMMA/ Medical grade stainless steel ("virgin steel")  
MDD Rule: Rule 6  
MDD Class: IIa  
CE: Acc. to existing directive



## ENDONASAL SUCTION TIP WITH TRANSPARENT HANDLE

- For rhinological surgery, FESS operations
- Olive shaped end for non-traumatic entrance via ostium
- Single use, sterile, single packed in peel-pack
- Medical grade stainless steel, transparent handle (PMMA)

<b>Art.no.</b>	<b>mm</b>	<b>ward.box/trsp.box</b>
6066503020	2,0 x 150	25/375
6066503025	2,5 x 150	25/375
6066503030	3,0 x 150	25/375
6066503035	3,5 x 150	25/375
6066503040	4,0 x 150	25/375

ward.box = min qty/ order

Material: PMMA/ Medical Grade Stainless Steel  
(virgin steel)  
MDD Rule: Rule 6  
MDD Class: IIa  
CE: Acc to regulation






## NEURO/OTO/RHINO SUCTION TIPS

- For delicate surgery procedures in Neuro/Oto/Rhinology
- 2 side holes at tip end, to avoid biopsy effect (tissue damage)
- Single use, sterile, packed in peel-pack
- Medical grade stainless steel, transparent handle of PMMA with built in vacuum control

<b>Art. No.</b>	<b>mm</b>	<b>Ward box/trsp.box</b>
6066509020	2,0 x 110	25/375
6066509025	2,5 x 110	25/375
6066509030	3,0 x 110	25/375
6066509035	3,5 x 110	25/375

ward box = minimum order qty

Material: PMMA/ Medical grade stainless steel  
("virgin steel")  
MDD Rule: Rule 6  
MDD Class: IIa  
CE: Acc. to existing directive



## NASAL SPECULA

- Single use, hygienically and single packed in peel pouch
- White ABS plastic, spring made of stainless steel
- Reversed spring effect allows Physician to use both hands

<b>Art.no.</b>	<b>ward.box/trsp.box</b>
6066516100	20/200

ward.box = min.qty/order

Material: ABS / Stainless Steel  
MDD Rule: Rule 5, subrule 1  
MDD Class: I non sterile  
CE: Acc to regulation



## PLASTIC SUCTION TIPS

- Plastic Suction Tips with smoothly rounded tips
- Funnel end for direct connection to Mediplast Oto-set with suction control
- Available in single and double bent versions
- Sterile, single packed in peel-pack

<b>Art.no.</b>	<b>mm</b>	<b>colour</b>	<b>ward.box/trsp.box</b>
----------------	-----------	---------------	--------------------------

### Double bent

6066101627	2,0 x 150	green	40/600
6066101635	3,0 x 150	blue	40/600
6066101643	4,0 x 150	white	40/600

### Short

6066101023	2,0 x 100	green	40/600
------------	-----------	-------	--------

### Single bent

6066100627	2,0 x 150	green	40/600
6066100635	3,0 x 150	blue	40/600
6066100643	4,0 x 150	white	40/600

ward.box = min.qty/order

Material: PVC/PP  
MDD Rule: Rule 5  
MDD Class: I sterile  
CE: Acc to regulation



## SUCTION FOR TONSILLECTOMY

- Single use, sterile, packed in peel-pack
- Medical grade stainless steel
- Transparent handle of PMMA with built in vacuum control

<b>Art. No.</b>	<b>mm</b>	<b>Ward box/trsp.box</b>
6066502040	4,0 x 190	20/300

ward box = minimum order qty

Material: PMMA/ Medical grade stainless steel  
 ("virgin steel")  
 MDD Rule: Rule 6  
 MDD Class: IIa  
 CE: Acc. to existing directive



## BRONCHOSCOPY SUCTION TIP

- Pediatric Bronchoscopy Suction Tip
- Single Use
- Sterile, single packed in peel pack

<u>Art.no.</u>	<u>mm</u>	<u>color</u>	<u>ward.box.</u>
6061626280	2,0 x 300	green	10

Material: PP/Medical Grade Stainless  
Steel (virgin steel)  
MDD Rule: Rule 6  
MDD Class: IIa  
CE: Acc to regulation



## **BRONCHO- / ESOPHAGO- SCOPIC SUCTION TIPS**

- Bronchoscopic- and Esophagoscopy- suction tips for adults and children
- For straight / rigid scopes
- Single use suction devices in different lengths and diameters
- New transparent handle for better control
- Sterile single packed in peel pack
- ETO sterilised

<b>Art No</b>	<b>Size mm</b>	<b>Ward box / trsp box.</b>
6061616242	2,0 x 400	10/90
6061616243	3,0 x 400	10/90
6061616244	4,0 x 400	10/90
6061616232	2,0 x 300	10/90
6061616233	3,0 x 300	10/90
6061616234	4,0 x 300	10/90
6061616252	2,0 x 500	10/90
6061616253	3,0 x 500	10/90
6061616254	4,0 x 500	10/90
6061616262	2,0 x 550	10/90
6061616263	3,0 x 550	10/90
6061616264	4,0 x 550	10/90

ward box = min qty / order

Material: PMMA/ Medical grade stainless steel (virgin steel)  
MDD Rule: Rule 6  
MDD Class: IIa  
EC: Acc. excisting directive





















## SUCTION TUBING

- A wide assortment of PVC and PVC free tubes
- Supple, resistant to folding and durable under vacuum pressure
- Available in 3 different packs:
  - Suction Tube Sterile - double packed in peel pack
  - Suction Tube Hygienical – hygienically packed in peel pack
  - Suction Tube on Roll – hygienically packed on roll











Material: PVC  
PVC free: XS119  
MDD Rule: Rule 2  
MDD Class: I, Is  
CE: Acc. to existing directives











## Suction tube, standard, serrated surface

Art. No.	Connector	Connector	Length	ID x OD	Hardness in Shore	Ch	Sterile	Box qty
6099900			200 cm	5,6 x 8,2	75	25	EO	50
6099903			250 cm					50
6099905			350 cm					40
6099910			200 cm	5,6 x 8,2	75	25	EO	50
6099915			350 cm					40
6099920			200 cm	5,6 x 8,2	75	25	EO	50
6099925			350 cm					40
6099930			200 cm	5,6 x 8,2	75	25	EO	50
6099935			350 cm					40
6099940			200 cm	5,6 x 8,2	75	25	EO	50
6099943			200 cm				Non-sterile	50
6099945			350 cm				EO	40
6099950			200 cm	7 x 10	75	30	EO	50
6099952			350 cm					40
6099959			200 cm	5,6 x 8,2	75	25	EO	50
6099957			200 cm	5,6 x 8,2	75	25	EO	50
6099958							Non-sterile	

## Suction tube, standard, serrated surface - PVC-Free

6088800			200 cm	5,6 x 8,2	75	25	EO	50
6088805			350 cm					40
6088810			200 cm	5,6 x 8,2	75	25	EO	50
6088840			200 cm	5,6 x 8,2	75	25	EO	50
6088850			200 cm	6 x 8	75	25	EO	50
6088860			200 cm	6 x 8	75	25	EO	50


## Suction tube, special, smooth surface

Art. No.	Connector	Connector	Length	ID x OD	Hardness in Shore	Ch	Sterile	Box qty
6099977			200 cm	5,6 x 8,2	75	25	EO	50
6099979			350 cm					40
6099980			200 cm	5,6 x 8,2	75	25	Non-sterile	1/50
6099982			200 cm	5,6 x 8,2	75	25	Non-sterile	1/50
6099983			350 cm					
6099985			200 cm	5,6 x 8,2	75	25	Non-sterile	1/50









## Bubble tube on roll, smooth tube

Art. No.	Connector	Connector	Length	ID x OD	Colour	Ch	Sterile	Box qty
6099970	-	-	30 m	3 x 5	Green	15	Non-sterile	1/5
6099973	-	-	30 m	5 x 7	Transparent	21	Non-sterile	1/5
6099975	-	-	30 m	3 x 5	Transparent	15	Non-sterile	1/5
6099976			30 m	7 x 9	Transparent	30	Non-sterile	1/5

## Suction tube, high flow aspiration, 360° rotating handle


Art. No.	Connector	Connector	Length	ID x OD	Hardness in Shore	Ch	Sterile	Box qty
6070002			200 cm	9,6 x 14,6	75	45	EO	15

## Suction tube, ENT, smooth tube, extra supple

Art. No.	Connector	Connector	Length	ID x OD	Hardness in Shore	Ch	Sterile	Box qty
6099990			200 cm	6 x 8	60 - 63	25	EO	50
6099992			350 cm					40
6099994			200 cm	6 x 8	60 - 63	25	EO	50
6099996			350 cm					40
6099987			200 cm	6 x 8	60 - 63	25	EO	50
6099998			200 cm					50
6099999			350 cm					40

## Suction tube special, with out connectors, transparent PVC, smooth tube

Art. No.	Connector	Connector	Length	ID x OD	Hardness in Shore	Ch	Sterile	Box qty
6099960	-	-	30 cm	7 x 10	68 - 70	30	EO	50
6099963	-	-	110 cm	7 x 10	68 - 70	30	EO	50
6099965	-	-	150 cm	7 x 10	68 - 70	30	EO	50
6099967	-	-	180 cm	7 x 10	68 - 70	30	EO	50

Art. No.	Connector	Length	ID x OD	-	-	Sterile	Box qty
6099981						Non-sterile	10



<b>Art.no.</b>		<b>Ward,/Trsp box</b>
----------------	--	-----------------------

Suction tube on roll

6091620043	4,0 x 6,0 mm	50/50
6091620050	5,0 x 8,0 mm	50/50
6091620068	6,0 x 9,0 mm	50/50
6091620076	7,0 x 10,0 mm	50/50
6091620084	8,0 x 11,0 mm	40/40
6091620090	8,0 x 12,0 mm	40/40
6091620092	9,0 x 13,0 mm	30/30
6091620100	10,0 x 14,0 mm	30/30

trsp box = min qty / order



## APPLICATORS

- Available in Aluminium or Steel
- Pure cotton end

<b>Art.no.</b>	<b>mm</b>	<b>ward.box/trsp.box</b>
601801183	alum. 0,9 x 150	100/10000
60201189	steel 0,9 x 150	100/10000
hygienically packed		

602011898	steel 0,9 x 150	20/100
sterile with 5 pcs in each pack		

ward.box = min.qty/order



## ANTRUM IRRIGATION SET

- 1,4 x 80 mm needle, 200 mm extension tube with luer lock connector
- Single use, sterile and single packed in peel pouch

<b>Art.no.</b>	<b>ward.box/trsp.box</b>
6066140146	10/300

ward.box = min.qty/order

Material: PP/ PVC/ Medical Grade Stainless  
Steel (virgin steel)  
MDD Rule: Rule 6  
MDD Class: IIa  
CE: Acc to regulation



## EAR FLUSHING TIP

- Fits on catheter coned syringes, has a safety wing
- Non-sterile, single packed in peel pouch

<b>Art.no.</b>	<b>ward.box/trsp.box</b>
6061220000	25/750

ward.box = min.qty/order

Material: ABS plastic  
MDD Rule: Rule 1  
MDD Class: I non-sterile  
CE: Acc to regulation



## SJECT FLUSHER SYRINGE

- Multi purpose flusher 100 ml
- Sterile, single packed in peel-pack

<b>Art.no.</b>	<b>ward.box/trsp.box</b>
6061990011	20/40

ward.box = min.qty/order

Material: PP  
MDD Rule: Rule 1  
MDD Class: I sterile  
CE: Acc to regulation



# MEDIPLAST®

## STANDARD AND CUSTOMIZED SETS, BALLS AND SWABS





## GAUZE SWABS HIGH-GRADE CLEAN

**MEDIPLAST®** can offer 8 layer gauze swabs in high-grade clean quality.

Our high-grade clean gauze swabs are packed in so called "after eight" boxes" or in paper wrapping.

Art No	Description	Ward box	Trsp box
<b>601055G</b>	Gauze swab 5 x 5 cm, "after eight" box	200	16000
<b>601075G</b>	Gauze swab 7,5 x 7,5 cm, "after eight" box	150	9000
<b>601010G</b>	Gauze swab 10 x 10 cm, "after eight" box	100	5000
<b>601020G</b>	Gauze swab 10 x 20 cm, "after eight" box	75	3000
<b>601220G</b>	Gauze swab 10 x 20 cm, <b>12 layer</b> , "after eight" box	50	2000
<b>6010255G</b>	Gauze swab 5 x 5cm, paper wrapping	100	40000
<b>6010275G</b>	Gauze swab 7,5 x 7,5cm, paper wrapping	100	10000
<b>6010210G</b>	Gauze swab 10 x 10cm, paper wrapping	100	10000
<b>60241021G</b>	Gauze swab 10 x 20 cm, paper wrapping, <b>24 layer</b>	50	1000



## GAUZE SWABS, STERILE, 8 LAYER

**MEDIPLAST®** can offer 8 layer gauze swabs sterile.

The sterile swabs are packed in 2 or 5-packs / pouch.

Art No	Description	Ward box	Trsp box
<b>6010552G</b>	Gauze swab 5 x 5 cm, 2-pack, sterile	120	4800
<b>6010555G</b>	Gauze swab 5 x 5 cm, 5-pack, sterile	125	5000
<b>6010752G</b>	Gauze swab 7,5 x 7,5 cm, 2-pack, sterile	120	4800
<b>6010755G</b>	Gauze swab 7,5 x 7,5 cm, 5-pack, sterile	125	5000
<b>6010102G</b>	Gauze swab 10 x 10 cm, 2-pack, sterile	120	4800
<b>6010105G</b>	Gauze swab 10 x 10 cm, 5-pack, sterile	125	3000
<b>601120G</b>	Gauze swab 7,5 x 7,5 cm, 20-pack, sterile	140	5600
<b>601127G</b>	Gauze finger bandage 3 x 38 cm, 2-pack	40	1920



## GAUZE SWABS STERILE, 12 AND 24 LAYER

**MEDIPLAST®** can offer 12 layer gauze swabs sterile.

The sterile swabs are packed in 2 or 5-packs / pouch.

Art No	Description	Ward box	Trsp box
<b>6012552G</b>	Gauze swab 5 x 5 cm, 2-pack, sterile	70	2800
<b>6012555G</b>	Gauze swab 5 x 5 cm, 5-pack, sterile	75	3000
<b>6012752G</b>	Gauze swab 7,5 x 7,5 cm, 2-pack, sterile	70	2800
<b>6012755G</b>	Gauze swab 7,5 x 7,5 cm, 5-pack, sterile	75	3000
<b>6012102G</b>	Gauze swab 10 x 10 cm, 2-pack, sterile	70	1680
<b>6012105G</b>	Gauze swab 10 x 10 cm, 5-pack, sterile	75	1800
<b>6012202G</b>	Gauze swab 10 x 20 cm, 2-pack, sterile	40	960
<b>6012205G</b>	Gauze swab 10 x 20 cm, 5-pack, sterile	40	960
<b>60241020G</b>	Gauze swab 10 x 20 cm, 2-pack, <b>24 layer</b>	20	480



## GAUZE BALLS

**MEDIPLAST®** has a wide assortment of gauze ball swabs. Non sterile gauze ball swabs are packed in bags. Sterile gauze ball swabs are packed in 5-pack / pouch. We also have a vaginal gauze ball swab on a wooden stick.

Art No	Description	Non Sterile / Sterile	Ward box	Trsp box
<b>601003G</b>	Gauze ball No 13, Ø 30 mm	Non sterile	100	4800
<b>601035G</b>	Gauze ball No 14, Ø 35 mm	Non sterile	100	4000
<b>601005G</b>	Gauze ball No 3, Ø 50 mm	Non sterile	100	1000
<b>601006G</b>	Gauze ball No 5, Ø 60 mm	Non sterile	100	800
<b>6010035G</b>	Gauze ball No 13, Ø 30 mm, 5-pack	Sterile	125	5000
<b>6010355G</b>	Gauze ball No 14, Ø 35 mm, 5-pack	Sterile	100	4800
<b>6010055G</b>	Gauze ball No 3, Ø 50 mm, 5-pack	Sterile	75	1800
<b>6010065G</b>	Gauze ball No 5, Ø 60 mm, 5-pack	Sterile	60	1440
<b>601080G</b>	Gauze flat No 3, 25 x 25 mm	Non sterile	500	12000
<b>601085G</b>	Gauze flat No 3, 25 x 25 mm, 5-pack	Sterile	300	12000
<b>601090G</b>	Gauze flat No 4, 30 x 30 mm, 5-pack	Sterile	300	12000
<b>609005</b>	Gauze flat 16 x 16 mm	Non sterile	1000	10000
<b>609004</b>	Gauze flat 8 x 8 mm	Non sterile	1000	20000
<b>608000V</b>	Vaginal gauze ball swab on a wooden stick	Non sterile	25	3500



## PRE-CUT GAUZE SWABS

**MEDIPLAST®** can offer non sterile, pre-cut gauze swabs, packed in paper wrapping.  
500 pcs/pack.

Art No	Description	Ward box	Trsp box
<b>601515G</b>	Pre-cut gauze swab 15 x 15 cm	500	20000
<b>602020G</b>	Pre-cut gauze swab 20 x 20 cm	500	10000
<b>603030G</b>	Pre-cut gauze swab 30 x 30 cm	500	5000
<b>604040G</b>	Pre-cut gauze swab 40 x 40 cm	500	2500
<b>605050G</b>	Pre-cut gauze swab 50 x 50 cm	500	2500



## NON WOVEN SWABS HIGH GRADE CLEAN

**MEDIPLAST®** can offer You 4 layer Non-Woven swabs, high grade clean. Our high grade clean Non-Woven swabs are packed in so called "after eight" boxes or in paper wrapping (30 g).

Art No	Description	Ward box	Trsp box
<b>601055NW</b>	Non-Woven swab 5 x 5 cm, "after eight" box	200	16000
<b>601075NW</b>	Non-Woven swab 7,5 x 7,5 cm, "after eight" box	150	9000
<b>601010NW</b>	Non-Woven swab 10 x 10 cm, "after eight" box	100	5000
<b>6010255NW</b>	Non-Woven swab 5 x 5 cm, Paper wrapping	100	40000
<b>6010275NW</b>	Non-Woven swab 7,5 x 7,5 cm, Paper wrapping	100	10000
<b>6010210NW</b>	Non-Woven swab 10 x 10 cm, Paper wrapping	100	10000
<b>601128NW</b>	Eye pad 5,6 x 7 cm	30	2100
<b>601130NW</b>	Eye pad 7,5 x 8 cm	30	1500



## NON WOVEN SWABS STERILE

**MEDIPLAST®** can offer You 4 layer Non-Woven swabs, sterile.

Our sterile Non-Woven swabs are packed in 2 or 5-packs / pouch (30 g).

Art No	Description	Ward box	Trsp box
<b>6010552NW</b>	Non-Woven swab 5 x 5 cm, 2-pack, sterile	120	4800
<b>6010555NW</b>	Non-Woven swab 5 x 5 cm, 5-pack, sterile	125	5000
<b>6010752NW</b>	Non-Woven swab 7,5 x 7,5 cm, 2-pack, sterile	120	4800
<b>6010755NW</b>	Non-Woven swab 7,5 x 7,5 cm, 5-pack, sterile	125	5000
<b>6010102NW</b>	Non-Woven swab 10 x 10 cm, 2-pack, sterile	120	2880
<b>6010105NW</b>	Non-Woven swab 10 x 10 cm, 5-pack, sterile	125	3000
<b>601129NW</b>	Eye pad 5,6 x 7 cm	25	2000
<b>601131NW</b>	Eye pad 7,5 x 8 cm	25	1500





## NON-WOVEN BALLS

**MEDIPLAST®** has a wide assortment of Non-Woven balls. Non sterile Non-Woven balls are packed in bags. Sterile Non-Woven balls are packed in 5-pack / pouch.

Art No	Description	Non Sterile / Sterile	Ward box	Trsp box
<b>601003NW</b>	Non-Woven ball No 13, Ø 30 mm	Non sterile	100	4800
<b>601035NW</b>	Non-Woven ball No 14, Ø 35 mm	Non sterile	100	4000
<b>601045NW</b>	Non-Woven ball, 45 mm	Non sterile	100	1000
<b>601005NW</b>	Non-Woven ball No 3, Ø 50 mm	Non sterile	100	1000
<b>601006NW</b>	Non-Woven ball No 5, Ø 60 mm	Non sterile	100	800
<b>6010035NW</b>	Non-Woven ball No 13, Ø 30 mm, 5-pack	Sterile	100	4000
<b>6010355NW</b>	Non-Woven ball No 14, Ø 35 mm, 5-pack	Sterile	75	3000
<b>6010455NW</b>	Non-Woven ball Ø 45 mm, 5-pack	Sterile	60	1440
<b>6010055NW</b>	Non-Woven ball No 3, Ø 50 mm, 5-pack	Sterile	60	1440
<b>6010065NW</b>	Non-Woven ball No 5, Ø 60 mm, 5-pack	Sterile	50	1200





## CUSTOMIZED SETS WITH BALLS

**MEDIPLAST®** can offer You sterile, set packed balls in gauze or in non woven, in bowls or trays. Both gauze and non woven balls are latex free.

Art No	Description	Ward box	Trsp box
<b>601122G</b>	Tray with 20 pcs gauze balls No 5, Ø 60 mm	10	40
<b>601123NW</b>	10 pcs non-woven balls No 5, Ø 60 mm, double packed	24	144
<b>601124NW</b>	Bowl with 15 pcs non-woven balls No 14, Ø 35 mm	16	64
<b>601125NW</b>	Tray with 10 pcs non-woven balls No 5, Ø 60 mm	28	168
<b>601126NW</b>	Bowl 250 ml with 5 pcs non-woven balls No 3, Ø 50 mm	30	180
<b>601133NW</b>	Bowl 500 ml with 10 pcs non-woven balls No 3, Ø 50 mm	16	64
<b>601134NW</b>	Bowl 1000 ml with 15 pcs non-woven balls No 3, Ø 50 mm	6	36
<b>601135NW</b>	Bowl 1000 ml with 20 pcs non-woven balls No 3, Ø 50 mm	6	36
<b>601136G</b>	Bowl 500 ml with 15 pcs gauze balls No 14, Ø 35 mm	16	64
<b>601137G</b>	Bowl 1000 ml with 20 pcs gauze balls No 5, Ø 60 mm	6	36
<b>601137</b>	2 kidney dishes, 1 bowl 250 ml, 1 bowl 500 ml, 2 bowl 1000 ml, 10 non-woven balls nr 3, Ø 50 mm	4	16
<b>601134</b>	Bowl 250 ml	30	180
<b>601135</b>	Bowl 500 ml	16	64
<b>601136</b>	Bowl 1000 ml	6	36
<b>601138</b>	kidney dish	12	48



## CATHETERIZATION SETS

**MEDIPLAST®** can offer You a wide assortment of sets for common, everyday procedures in healthcare. We have standardized series of sets but we can also customize sets according to specifications required (against a minimum quantity of 1000 pcs / year.) Here are some samples of standard catheterization sets:

Art No	Description	Ward box	Trsp box
<b>605011N</b>	2 trays 2 forceps 2 balls Ø 50 mm and 5 balls Ø 30 mm, gauze 2 swabs 10 x 10 cm, gauze 1 towel 50 x 80 cm and 1 towel 50 x 70 cm 1 towel with slit 50 x 80 cm 1 pair of nitrile gloves	10	40
<b>605013N</b>	2 kidney dishes och 1 bowl 160 ml 2 forceps 5 balls Ø 30 mm, and 2 balls Ø 35 mm, gauze 3 swabs 7,5 x 7,5 cm, gauze 1 towel 50 x 80 cm and 1 towel 50 x 70 cm 1 towel with slit 50 x 80 cm 1 pair of nitrile gloves	10	40
<b>605014N</b>	1 kidney dish 1 forcep 2 balls Ø 50 mm, gauze 2 swabs 7,5 x 7,5 cm, gauze 1 towel 50 x 70 cm 1 pair of nitrile gloves	10	40
<b>605016N</b>	1 kidney dish and 1 bowl 160 ml 2 forceps 5 balls Ø 30 mm, and 2 balls Ø 35 mm, gauze 3 swabs 7,5 x 7,5 cm, gauze 1 waste bag 1 towel 50 x 80 cm and 1 towel 50 x 70 cm 1 towel with slit 50 x 80 cm 1 pair of nitrile gloves	10	40



Art No	Description	Ward box	Trsp box
<b>605020N CLEAN PACK</b>	4 swabs 7,5 x 7,5 cm, non-woven 2 towels 50 x 70 cm 1 pair of nitrile gloves	60	360
<b>605021N</b>	1 kidney dish 1 bowl 120 ml 1 forcep 7 swabs 10 x 10 cm, gauze 1 towel 50 x 70 cm 1 towel with slit 50 x 80 cm	10	40
<b>605022N</b>	1 kidney dish 1 bowl 160 ml 1 forcep 3 swabs 7,5 x 7,5 cm, gauze in a bowl 9 swabs 7,5 x 7,5 cm, gauze in a kidney dish 2 towels 50 x 80 cm	10	40
<b>605023N</b>	1 kidney dish 1 bowl 160 ml 1 forcep 2 swabs 10 x 10 cm, non-woven 5 balls Ø 35 mm, non-woven 1 towel 50 x 80 cm 1 towel with slit 50 x 80 cm 1 pair of vinyl gloves	10	40
<b>605024N</b>	1 tray 1 forcep 2 balls Ø 35 mm, gauze 4 swabs 7,5 x 7,5 cm, gauze 1 towel 50 x 60 cm 1 towel with slit 50 x 80 cm 1 pair of nitrile gloves	28	168
<b>605025N</b>	1 tray 1 forcep 2 balls Ø 35 mm, non-woven 4 swabs 7,5 x 7,5 cm, non-woven 1 towel 50 x 60 cm 1 towel with slit 50 x 80 cm 1 pair of nitrile gloves	28	168



Art No	Description	Ward box	Trsp box
<b>605026N</b>	1 kidney dish 1 bowl 2 forceps 3 balls Ø 50 mm, non-woven 5 balls Ø 35 mm, non-woven 5 swabs 7,5 x 7,5 cm, non-woven 1 towel 60 x 70 cm 1 towel with slit 50 x 80 cm 1 pair of nitrile gloves	10	40
<b>605027N</b>	1 kidney dish 1 bowl 2 forceps 2 swabs 10 x 10 cm, gauze 7 balls Ø 35 mm, gauze 1 pair of nitrile gloves 1 towel	10	40
<b>605028N</b>	2 tray 2 forceps 3 swabs 7,5 x 7,5 cm, non-woven 3 balls Ø 50 mm non-woven 6 balls Ø 30mm non-woven 1 towel with slit 50 x 80 cm 1 towel 42 x 50 cm 1 towel 50 x 50 cm 1 pair of nitrile gloves	10	40



## WASHING SETS - STERILE

**MEDIPLAST®** can offer You a wide assortment of sets for common, everyday procedures in healthcare. We have standardized series of sets but we can also customize sets according to specifications required (against a minimum quantity of 1000 pcs / year.)

Art No	Description	Ward box	Trsp box
<b>605001T</b>	Washing set with 8 gauze balls, No 3, 1 forcep	18	108
<b>605002T</b>	Washing set with 8 gauze balls, No 3, 1 forcep in tray, towel 70x70 cm	18	108
<b>605004T</b>	Washing set with 8 gauze balls, No 3, 1 forcep and 1 towel	18	108
<b>605007T</b>	Washing set with 5 Non-Woven swabs 7,5 x 7,5 cm and a bowl 250 ml	30	180
<b>605008T</b>	Washing set with 10 Non-Woven swabs 10x10 cm and a bowl 250 ml, 2 towels 50x50 cm	24	144
<b>605012T</b>	Washing set with 1 kidney bowl, 1 bowl 250 ml, 1 bowl 500 ml, and 1 bowl 1000 ml with 10 non-woven balls Ø 50 mm	4	16
<b>605013T</b>	Washing set with 10 non-woven balls Ø 50 mm, forcep, tray, towel 60 x 90 cm	18	108
<b>605014T</b>	Washing set with 10 non-woven balls Ø 50 mm, bowl 1000 ml, towel 60 x 70 cm	6	36
<b>605136T</b>	Washing set with 8 gauze balls, No 3, 1 forcep/special and 1 towel	18	108
<b>605156T</b>	Washing set with 8 gauze balls, No 3, 1 pincette and towel	18	108



## WASHING SETS - UNSTERILE

**MEDIPLAST®** can offer You a wide assortment of sets for common, everyday procedures in healthcare. We have standardized series of sets but we can also customize sets according to specifications required (against a minimum quantity of 1000 pcs / year.)

### Standard Unsterile:

Art No	Description	Ward box	Trsp box
<b>605005T</b>	Washing set with 4 Non-Woven balls, No 14, 6 swabs 10 x 10 cm and a bowl 250 ml	90	90
<b>605006T</b>	Washing set with 4 Non-Woven balls, No 3 and a bowl 500 ml	50	50
<b>605009T</b>	Washing set with 12 Non-Woven balls 50 mm Ø in a bowl	40	40
<b>605010T</b>	Washing set with 6 Non-Woven balls 50 mm Ø in a bowl	150	150
<b>605011T</b>	Washing set with 16 Non-Woven balls 50 mm Ø in a bowl including pincette	40	40



## PROCEDURE SETS

**MEDIPLAST®** can offer You a wide assortment of sets for common, everyday procedures in healthcare. We have standardized series of sets but we can also customize sets according to specifications required (against a minimum quantity of 1000 pcs / year.) Here are some samples of standard wound care and dressing sets:

Art No	Description	Ward box	Trsp box
<b>605002G</b>	<b>Wound care set with:</b> 5 gauze balls, No 13 in tray with forcep	50	200
<b>605002NW</b>	<b>Wound care set with:</b> 5 Non-Woven balls, No 13 in tray with forcep	50	200
<b>605003G</b>	<b>Wound care set with:</b> 3 gauze balls, No 13, and 2 gauze swabs 5 x 5 cm in tray with forcep	50	200
<b>605003NW</b>	<b>Wound care set with:</b> 3 Non-Woven balls, No 13, and 2 gauze swabs 5 x 5 cm in tray with forcep	50	200
<b>601121F</b>	<b>Dressing set</b> 4 cotton applicators, 12 gauze swabs 10 x 10 cm and 5 pcs balls 60 mm Ø in kidney dish	10	40
<b>605019</b>	<b>Suture removal set</b> 1 tray, 1 suture knife, 3 gauze swabs 25 x 25mm, 1 pincette	20	320





## DIALYSIS SETS

**MEDIPLAST®** can offer You a wide assortment of sets for common, everyday procedures in healthcare. We have standardized series of sets but we can also customize sets according to specifications required (against a minimum quantity of 1000 pcs / year.) Here are just some samples of standard dialysis sets:

Art No	Description	Ward box	Trsp box
<b>606191D</b>	ON 1 Tray 2 Non-woven swabs 7,5 x 7,5 cm 5 Non-woven swabs 10 x 10 cm 2 Towels 50 x 50 cm OFF 4 Gauze balls Ø 50 mm 1 Towel 50 x 50 cm	<b>15</b>	<b>90</b>
<b>601010</b>	1 Tray 1 Bowl 60 ml 10 Gauze swabs 7,5 x 7,5 cm 1 Towel with slit 50 x 80 cm 1 Peang 1 tape/plaster 10 x 2,5 cm	<b>15</b>	<b>90</b>
<b>606193D</b>	ON 1 Tray 8 Non-woven swabs 7,5 x 7,5 cm 2 Syringes 5 ml luer 2 Syringes 10 ml luer 2 Towels 50 x 50 cm OFF 1 Tray 2 Syringes 5 ml luer 2 Syringes 10 ml luer 2 Stop cocks 1 Canula 1,2 x 50 mm 2 Towels 50 x 50 cm	<b>12</b>	<b>72</b>





## PRE-INJECTION SWAB

**MEDIPLAST®** can offer You Injection Swab Std and Pre-Injection Swab. The Injection Swab Std is made from unbleached cellulose. We can also offer a Pre-Injection Swab of a Non-Woven material, this version is single packed and prepared with 70% isopropyl alcohol.

Art no	Description	Ward box	Trsp box
<b>60615003</b>	Injection Swab Std 4 x 6 cm	250	18000
<b>60615004</b>	Pre-Injection Swab 3 x 6 cm, 70% IPA	100	2000



## Material Specification Gauze Products

Material specification:	100% cotton
No of Threads:	17 threads/cm <sup>2</sup>
Weight:	24 g/m <sup>2</sup>
Bleach method:	Peroxide bleach
Elastic band/ring in ball:	Poly-isoprene
Latex free:	Yes
Packing material Non-sterile:	Boxes - medical grade paper Ball in plastic bag – plastic bag in PE Paper pouch/paper bag – medical grade paper
Packing material Sterile:	Peel-pack - medical grade paper (PA) +
Swabs and balls	80 µm PE (inner pouch) Carton/boxes - Medical grade paper Non-recycled paper in all boxes. All boxes are non-bleached
Sterilisation method:	EO (ethylene oxide)



## Material Specification Non-Woven Products

Material specification:	70% rayon and 30% polyester
Weight:	30 g/m <sup>2</sup>
Bleach method:	Peroxide bleach
Elastic band/ring in ball:	Poly-isoprene
Latex free:	Yes
Packing material Non-sterile:	Boxes - medical grade paper Ball in plastic bag – plastic bag in PE Paper pouch/paper bag – medical grade paper
Packing material Sterile:	Peel-pack - medical grade paper (PA) +
Swabs and balls	80 µm PE (inner pouch) Carton/boxes - Medical grade paper Non- recycled paper in all boxes. All boxes are non-bleached
Sterilisation method:	EO (ethylene oxide)



## Material Specification Other Products

Tray:	PS (Polystyrene)
Pincette:	PS (Polystyrene)
Peang/Forcep:	PP (Polypropylene)
Folded towels/Wrapping towels:	20 g PE (Polyethylene + cellulose/paper)
Dialysis sets:	Sterilisation method: Gamma radiation





## ABDOMINAL SWABS AND SWABS X-RAY DETECTABLE, GAUZE

**MEDIPLAST®** can offer a high quality assortment of Abdominal swabs and swabs of gauze in various sizes. The gauze is 100% cotton, peroxide bleached. The X-ray proof thread is welded onto the Abdominal swabs and woven into the swabs.

The products are colour coded according to existing EN standards for OR's. The colour coding is found both on the checklist label and the pack label.

Art. No.	Description	Size cm	Layers/ply	Colour	Ward box / Trsp. box
<b>601301</b>	5-Pack G Sterile	30 x 45	4	Brown	45/540
<b>601302</b>	1-pack G Sterile	45 x 70	3	Grey	25/300
<b>601303</b>	5-pack G Sterile	45 x 70	3	Black	30/240
<b>601304</b>	5-pack G Sterile	10 x 60	4	Lilac	100/1200
<b>601305</b>	5-pack G Sterile	10 x 60	8	Green	45/540
<b>601306</b>	5-pack G Sterile	10 x 60	12	Red	30/360
<b>601307</b>	5-pack G Sterile	10 x 10	12	Orange	100/1200

ETO sterilized

EN standards: EN 550, EN 14079



## SURGICAL SPONGES X-RAY DETECTABLE, GAUZE

**MEDIPLAST®** can offer a high quality assortment of Surgical sponges of gauze in various sizes. The gauze is 100% cotton, peroxide bleached. The X-ray proof thread is welded inside the Surgical sponges.

The products are colour coded according to existing EN standards for OR's. The colour coding is found both on the checklist label and the pack label.

Art. No	Description	Size	Ø mm	Colour	Ward box / Trsp. box
<b>601311</b>	5-Pack G Sterile	0	13	Orange	250/2000
<b>601312</b>	5-pack G Sterile	1	25	Yellow	250/2000
<b>601313</b>	5-pack G Sterile	3	35	Blue	150/750
<b>601314</b>	5-pack G Sterile	4	40	Grey	150/750
<b>601315</b>	5-pack G Sterile	5	50	Black	75/600

ETO sterilized

EN standards: EN 550, EN 14079



## ABDOMINAL SWABS

### X-RAY DETECTABLE, NON WOVEN

**MEDIPLAST®** can offer a high quality assortment of Abdominal swabs of non woven material in various sizes. The non woven material is made of 70% rayon and 30% polyester. The X-ray proof thread is welded onto the Abdominal swabs. The products are colour coded according to existing EN standards for OR's. The colour coding is found both on the checklist label and the pack label.

Art. No	Description	Size cm	Layers/ply	Colour	Ward box / Trsp. box
<b>601320</b>	5-pack NW Sterile	10 x 40	4	Orange	125/1500
<b>601321</b>	5-pack NW Sterile	20 x 40	8	Green	65/780
<b>601322</b>	5-pack NW Sterile	30 x 40	12	Brown	45/540
<b>601323</b>	1-pack NW Sterile	40 x 65	16	Grey	30/360
<b>601324</b>	5-pack NW Sterile	40 x 65	16	Black	30/360

ETO sterilized

EN standards: EN 550, EN 1644, EN-1644-2



# MEDIPLAST®

## FLAT AND TUBULAR BANDAGES

**Hospidana** Made in Denmark

Medi plast AB  
Box 9504  
200 39 MALMÖ  
Sweden

[mediplast.info@mediplast.com](mailto:mediplast.info@mediplast.com)  
[www.mediplast.com](http://www.mediplast.com)



Phone: +46 40-671 23 00  
Fax: +46 40-671 23 10





## ELASTIC FIXATION BANDAGES

**MEDIPLAST®** can offer **Danamull** and **Danalast** elastic fixation bandages. Supporting fit over limbs and wound bandaged areas. The bandages are blood and salves resistant, and have good aeration. Autoclaveable. Latex free.

### DANAMULL:

Double weaving with woven edges  
70% Viscose, 30% Polyamide  
Elasticity: 100%

### PACKING:

Ward box: coated card board of recycled/new fibre mix.  
Trsp box: unbleached card board of recycled/new fibre mix.

Article No	Description	Pcs/ ward box	Pcs/ trsp.box
<b>6050601</b>	4 cm x 4 meters	20	1440
<b>6050602</b>	6 cm x 4 meters	20	960
<b>6050603</b>	8 cm x 4 meters	20	720
<b>6050604</b>	10 cm x 4 meters	20	600
<b>6050605</b>	12 cm x 4 meters	20	480

### DANALAST:

Single weaving  
60% Polyamide, 40% Viscose  
Elasticity: 100%

Article No	Description	Pcs/ ward box	Pcs/ trsp.box
<b>6050501</b>	4 cm x 4 meters	20	1440
<b>6050502</b>	6 cm x 4 meters	20	960
<b>6050503</b>	8 cm x 4 meters	20	720
<b>6050504</b>	10 cm x 4 meters	20	600
<b>6050505</b>	12 cm x 4 meters	20	480



## ELASTIC SUPPORT BANDAGES

**MEDIPLAST®** can offer **Danalastic** universal bandage to be used for support or light fixation by procedures like fixation of plaster strip.

The bandage is blood and salves resistant. Hand washable without loss of elasticity. Autoclaveable. Latex free.

### DANALASTIC:

63% Cotton  
37% Polyamide  
Elasticity: 70%

### PACKING:

Ward box: coated card board of recycled/new fibre mix.  
Trsp box: unbleached card board recycled/new fibre mix.

Article No	Description	Pcs/ ward box	Pcs/ trsp.box
<b>6050300</b>	4 cm x 5 meters	10	720
<b>6050301</b>	6 cm x 5 meters	10	480
<b>6050302</b>	8 cm x 5 meters	10	320
<b>6050303</b>	10 cm x 5 meters	10	280
<b>6050304</b>	12 cm x 5 meters	10	240
<b>6050305</b>	15 cm x 5 meters	10	200





## ELASTIC UNIVERSAL BANDAGE

**MEDIPLAST®** can offer **Danauniversal** to be used for support, relieve after injury and for ex. lymph oedema treatment. Skin friendly with high aeration. Crêpe type bandage that can also be used for light compression. Autoclaveable. Latex free.

### DANAUNIVERSAL:

Crêpe type bandage  
69% Cotton  
29% Polyamide  
2% Polyurethane  
Elasticity: 100%

### PACKING:

Ward box: coated card board of recycled/new fibre mix.  
Trsp box: unbleached card board recycled/new fibre mix.

Article No	Description	Pcs/ ward box	Pcs/ trsp.box
<b>6050401</b>	4 cm x 5 meters	10	720
<b>6050402</b>	6 cm x 5 meters	10	480
<b>6050403</b>	8 cm x 5 meters	10	320
<b>6050404</b>	10 cm x 5 meters	10	280
<b>6050405</b>	12 cm x 5 meters	10	240
<b>6050406</b>	15 cm x 5 meters	10	200



## COHESIVE BANDAGES

**MEDIPLAST®** can offer cohesive bandages for fixation, **Danamull Haft** in white, red or blue and, **Danalast Haft** in white. The bandage attaches itself to its own surface and not to the skin. This fixation pulp contains latex. The bandages have good aeration.

### **DANAMULL HAFT, white**

Double weaving  
46% Cotton  
27% Polyamide  
27% Viscose  
Elasticity: 80 – 100%

### **PACKING:**

Ward box: coated card board of recycle/new fibre mix.  
Trsp box: unbleached card board of recycled/new fibre mix.

Article No	Description	Colour	Pcs/ ward box	Pcs/ trsp.box
<b>6051111</b>	4 cm x 4 m	White	1	84
<b>6051112</b>	6 cm x 4 m	White	1	54
<b>6051113</b>	8 cm x 4 m	White	1	48
<b>6051114</b>	10 cm x 4 m	White	1	36

<b>6051101</b>	4 cm x 20 m	White	1	22
<b>6051102</b>	6 cm x 20 m	White	1	16
<b>6051103</b>	8 cm x 20 m	White	1	12
<b>6051104</b>	10 cm x 20 m	White	1	10

**DANAMULL HAFT, red and blue**

Double weaving  
46% Cotton  
27% Polyamide  
27% Viscose  
Elasticity: 80 – 100%

**PACKING:**

Ward box: coated card board of recycled/new fibre mix.  
Trsp box: unbleached card board of recycled/new fibre mix.

Article No	Description	Colour	Pcs/ ward box	Pcs/ trsp.box
<b>6051102R</b>	6 cm x 20 m	Red	1	16
<b>6051103R</b>	8 cm x 20 m	Red	1	12
<b>6051104R</b>	10 cm x 20 m	Red	1	10

<b>6051102B</b>	6 cm x 20 m	Blue	1	16
<b>6051103B</b>	8 cm x 20 m	Blue	1	12
<b>6051104B</b>	10 cm x 20 m	Blue	1	10

**DANALAST HAFT:**

Smooth woven bandage  
46% Cotton  
54% Polyamide  
Elasticity: 80 – 100%

**PACKING:**

Ward box: coated card board of recycled/new fibre mix.  
Trsp box: unbleached card board of recycled/new fibre mix.

Article No	Description	Pcs/ ward box	Pcs/ trsp.box
<b>6051211</b>	4 cm x 4 meters	1	84
<b>6051212</b>	6 cm x 4 meters	1	54
<b>6051213</b>	8 cm x 4 meters	1	48
<b>6051214</b>	10 cm x 4 meters	1	36

<b>6051201</b>	4 cm x 20 meters	1	24
<b>6051202</b>	6 cm x 20 meters	1	18
<b>6051203</b>	8 cm x 20 meters	1	12
<b>6051204</b>	10 cm x 20 meters	1	12



## ELASTIC TUBULAR BANDAGE

**MEDIPLAST®** can offer **Danatube**, an elastic tubular fixation bandage. Can also be used underneath plaster and cohesive bandages. Good aeration. Latex free. Length 20 meters.

### DANATUBE:

100% Unbleached cotton

### PACKING:

Ward box: coated card board of recycled/new fibre mix.

Trsp box: unbleached card board of recycled/new fibre mix.

Article No	Description	Pcs/ ward box	Pcs/ trsp.box
<b>6050721</b>	Size 0,1 = 1,5 cm, fingers, toes	1	64
<b>6050722</b>	Size 1,2 = 2,5 cm, fingers, toes with cotton wool, wadding	1	33
<b>6050723</b>	Size 3,4 = 4,0 cm, hand, small arms, childrens legs	1	29
<b>6050724</b>	Size 5,6 = 6,0 cm, large arms, normal legs	1	51
<b>6050725</b>	Size 7,8 = 8,0 cm, head, knee, small thighs	1	36
<b>6050726</b>	Size T1 = 10,0 cm, trunk bandage for children or large head	1	33
<b>6050727</b>	Size T2 = 16 cm, trunk bandage for adult	1	18



## TUBULAR BANDAGE FOR COMPRESSION

**MEDIPLAST®** can offer **Danagrip**, an elastic support and fixation bandage with compression effect. Can be used as single or double layer to achieve desired compression. Secures lasting compression even after long time use. Autoclavable. Holes and gashes can be cut out without risk for fringes. Hand washable 40°C. Can also be used underneath plaster and cohesive bandages. Good aeration. Very good elasticity, can be stretched 3 times its original width. Latex free, unbleached. Length 10 meters.

### **DANAGRIP, unbleached / skin colour:**

82% Unbleached cotton  
9% Lycra  
9% Polyamide

### **PACKING:**

Ward box: coated card board of recycled/new fibre mix.  
Trsp box: unbleached card board of recycled/new fibre mix.

Article No	Description	Pcs/ ward box	Pcs/ trsp.box
<b>6050901</b>	Size A = 4 cm, small arms, legs	1	36
<b>6050902</b>	Size B = 6,25 cm, small hands	1	24
<b>6050903</b>	Size C = 6,75 cm, hands, small arms, children's legs	1	24
<b>6050904</b>	Size D = 7,5 cm, normal legs	1	20
<b>6050905</b>	Size E = 8,75 cm, legs, knees	1	18
<b>6050906</b>	Size F = 10 cm, thighs	1	14
<b>6050907</b>	Size G = 12 cm, large thighs	1	14
<b>6050908</b>	Size J = 18 cm, small bodies	1	8
<b>6050909</b>	Size K = 22 cm, normal bodies	1	8
<b>6050910</b>	Size L = 32,5 cm, large bodies	1	6



## TUBULAR BANDAGES, WITHOUT COMPRESSION

**MEDIPLAST®** can offer **Danafast**, a tubular elastic bandage without compression. Can be used as a fixation bandage underneath plaster and also for dermatology applications.

**Danafast** is elastic width wise and is therefore easy to put on. Delivered on roll, for easy cut of desired length, can be used as single or double layer, no extra fixation is needed. Good aeration. Colour coded with a woven line which makes choice of size easy.

**Danafast** bandages are produced according to Drug Tariff 46:1989 standards. Latex free. Length 10 meters.

### DANAFAST:

92% Unbleached viscose  
5% Unbleached lycra  
3% Polyamide

### PACKING:

Ward box: coated card board of recycled/new fibre mix.  
Trsp box: unbleached card board of recycled/new fibre mix.

Article No	Description	Pcs/ ward box	Pcs/ trsp.box
<b>6050751</b>	RED = 3,0-3,8 cm, small arms and legs	1	35
<b>6050752</b>	GREEN = 4,5-5,5 cm, feet, hands, normal arms, small legs	1	24
<b>6050753</b>	BLUE = 7,5-8,2 cm, normal legs, large arms, thighs	1	24
<b>6050754</b>	YELLOW = 10,2-12,5 cm, small trunks, large thighs	1	20
<b>6050755</b>	BEIGE = 17-19 cm, adult trunks	1	18



## ELASTIC TUBULAR NET BANDAGE

**MEDIPLAST®** can offer **Dananet**, a large mesh, latex free tubular fixation net bandage. Can be used for multiple fixation tasks where secure fixation is sought after. Delivered on roll for easy cut to desired length. Supple fit and good elasticity. Hand washable 40°C. Length 25 meters.

### DANANET:

75% Polyamide  
25% Lycra

### PACKING:

Ward box: coated card board  
of recycled/new fibre mix  
Trsp box: unbleached card board of  
recycled/new fibre mix.

Article No	Description	Pcs/ ward box	Pcs/ trsp.box
<b>6050821</b>	Size 0,5 = 8 mm, small fingers	1	66
<b>6050822</b>	Size 1 = 14 mm, fingers, toes	1	66
<b>6050823</b>	Size 2 = 17 mm, hands, multiple fingers	1	40
<b>6050824</b>	Size 3 = 20 mm, hands, feet	1	40
<b>6050825</b>	Size 4 = 23 mm, hands, feet, knees	1	40
<b>6050826</b>	Size 5 = 26 mm, feet, elbows, knees, head	1	40
<b>6050827</b>	Size 5,5 = 36 mm, head, small shoulders	1	38
<b>6050828</b>	Size 6 = 48 mm, trunks, children	1	38
<b>6050829</b>	Size 7 = 52 mm, shoulders, chest, hip	1	34
<b>6050830</b>	Size 8 = 57 mm, shoulders, chest, hip	1	34





## ELASTIC BANDAGE PRIMALAST

Elastic "Prima" bandage for support and relieve and for light compression.

Material: Cotton 100%, peroxide bleached.

Elasticity: 90%.

Characteristics: Manufactured according to EN ISO 13485:2000. Special woven edges to avoid fringing.

Unstretched length approx 2 m. No shrinking.

Non sterile: Autoclavable in steam autoclave at 134 degrees C.

Sterile: ETO sterilized and delivered in SBS bag

<b>Non sterile:</b>	<b>Art No</b>	<b>Size</b>	<b>Ward box/Trsp box</b>
	607028	6 cm x 4 m	20/400
	607029	8 cm x 4 m	20/320
	607030	10 cm x 4 m	20/240
	607034	12 cm x 4 m	20/240
	607032	15 cm x 4 m	20/160
<b>Sterile:</b>	<b>Art No</b>	<b>Size</b>	<b>Ward box/Trsp box</b>
	607028S	6 cm x 4 m	24/144
	607029S	8 cm x 4 m	24/144
	607030S	10 cm x 4 m	24/96
	607034S	12 cm x 4 m	24/96
	607032S	15 cm x 4 m	24/96

Ward box. = min order qty





## LOW ELASTIC BANDAGE GASELAST

Low elastic bandage for fixation of dressing. Non sterile.

Material  
100% Cotton. Autoclavable.

<b>Art no</b>	<b>Size</b>	<b>Ward box / Trsp box</b>
607011	6 cm x 4 m	20/520
607012	8 cm x 4 m	20/560
607013	10 cm x 4 m	20/480

Ward box = min order qty



## SURGICAL SPONGES

### X-RAY DETECTABLE, NON WOVEN

**MEDIPLAST®** can offer a high quality assortment of Surgical sponges of non woven material in various sizes. The non woven material is made of 70% rayon and 30% polyester. The X-ray proof thread is welded inside the Surgical sponges. The products are colour coded according to existing EN standards for OR's. The colour coding is found both on the checklist label and the pack label.

Art. No	Description	Size	Ø mm	Colour	Ward box / Trsp. box
<b>601327</b>	5-pack NW Sterile	1	25	Yellow	175/1050
<b>601328</b>	5-pack NW Sterile	3	30	Blue	150/900
<b>601329</b>	5-pack NW Sterile	4	35	Grey	75/450
<b>601330</b>	5-pack NW Sterile	5	40	Black	75/450

ETO sterilized

EN standards: EN 550, EN 1644, EN-1644-2



# INFUSION AND TRANSFUSION





## INFUSION SET STD

- 3 std tube lengths
- Alternatives with detachable 3-way stopcock
- PVC or PVC-free material
- Roll clamp with maximum precision for steady drip flow rate
- Spike without silicone coating for secure connection to infusion bags/bottles
- Fluid filter 15my. Hydrophobic filter 0,2my
- Protection cap with hydrophobic filter simplifies priming of sets and reduces contamination risk
- All sets with rotating luer-lock



Art no	Description	Ward box/Trsp box 10/100
60MR60500150	non vented 180 cm inf set	
60MR60500180	non vented 220 cm inf set	
60MR60500151	non vented 320 cm inf set	

60MR60500154	vented 180 cm inf set
60MR60500184	vented 220 cm inf set
60MR60500155	vented 320 cm inf set

Art no	Description PVC-free	Ward box/Trsp box 10/100
60MR60500152	non vented 180 cm inf set	
60MR60500182	non vented 220 cm inf set	
60MR60500160	non vented 320 cm inf set	

60MR60500156	vented 180 cm inf set
60MR60500186	vented 220 cm inf set
60MR60500161	vented 320 cm inf set

All products are single packed in peel-pack  
Trsp box = min qty/order



### Infusion sets including 3-way stopcock and extension tube

Vented or non vented, PVC or PVC-Free material.

Detachable 3-way stopcock of lipid resistant Polycarbonate

Infusion set according to specification above.

Art no	Description	Ward box/Trsp box 10/100
60MR60500141	non vented 180 cm inf set, 3-way stopcock + 15 cm	
60MR60500139	non vented 220 cm inf set, 3-way stopcock + 30 cm	
60MR60500140	vented 180 cm inf set, 3-way stopcock + 15 cm	

Art no	Description PVC-free	Ward box/Trsp box 10/100
60MR60500158	non vented 180 cm inf set, 3-way stopcock +15 cm	
60MR60500138	vented 180 cm inf set, 3-way stopcock +15 cm	
60MR60500183	non vented 220 cm inf set, 3-way stopcock +15 cm	
60MR60500187	vented 220 cm inf set, 3-way stopcock +15 cm	

All products are single packed in peel-pack

Trsp box = min qty/order

<u>Infusion set 20 drps aqua dest = 1 ml</u> Protection cap: PE and Polyamide Hydrophobic, fluid filter: Polyamide Spike, roller clamp and luer: ABS Drip chamber: PP	Tube: Transp. PVC or Polyefine alloy (XS118) ID/OD 3,0/4,1 mm <u>3-way stopcocks:</u> Lipid resistant Polycarbonate, 72 Hrs, pressure resistant 2,75 bar	Sterilization method: ETO ISO 8536-4 EN 980 Luer EN 20594-1 MDD Rule: Rule 2 MDD Class: IIa CE: According to existing directives
---	--	---



## MIXING SETS

- A closed system for preparation and administration of antibiotics and chemotherapy drugs
- Simple and safe to use
- Will make the drug handling safer. The mixing sets prevent spilling and unnecessary exposure
- Hydrophobic filter on the protection cap eliminates risk for leakage and contamination of the luer-lock connection
- Rotating luer-lock on all models

*For mixing of dry substances and transfer of solution from ampoule to infusion bag with a closed procedure*



*60MR03508357C  
For preparation of stem solutions, will also give the possibility to adapt dose volume*



Art no	Description	Ward box/Trsp box
60MR03508357	180 cm, 1-way	10/100
60MR03508357C	180 cm, 1-way for stem solutions	10/100
60MR03508358	180 cm, 2-way	10/100
60MR03508359	180 cm, 4 way	10/50
60MR03508057	1 way, 180 cm PVC-free tube	10/100
60MR03508058	2-way, 180 cm PVC-free tube	10/50
60MR03508059	4 way, 180 cm PVC-free tube	10/50

All products are single packed in peel-pack  
Trsp box = min qty/order

Material:  
Protection caps: PE  
Drip chambers: PP  
Spike, roller clamp, luer: ABS  
Hydrophobic filter, infusion filter: Polyamide  
Clave: Polyester, Silicon and ABS

Tube: PVC or Polyfine Alloy (XS118)

Sterilization method: ETO  
ISO 8536-4 and 8536-12 EN 980  
Luer EN 20594-1  
MDD Rule: Rule 2  
MDD Class: IIa  
CE: According to existing directives





# INFUSION SET - SPECIAL

## Infusion set with integrated back check valve in the luer-lock connector

- Used in parallell infusions
- Helps to make infusion therapy safer and easier
- Prohibits back flow. Closes automatically by back pressure
- No risk for infusion solutions/drugs to be mixed
- Gives steady infusion flow rate
- Reduces risk for stop in cannulas/catheters



Art no	Ward box/Trsp box
60MP60500170 non vented inf set 180 cm	10/100

## Infusion set with Clave

- Vented infusion set with Y-connector and MicroClave, injection membrane
- Clave is a needle free connector with a small dead space: 0,04 ml. Fits all luer-connectors.
- Increases safety for staff and patient
- Minimizes risk for contamination
- No risk for leakage/spilling
- Simplifies injections during ongoing infusion



Art no	Ward box/Trsp box
60MP60500154C vented inf set 180 cm	10/100
60MP60500154Cv vented inf set 180 cm	10/100
2 back check valves	

## Infusion set with integrated infusion filter (0,2 my) för Taxol infusion, PVC-free

Art no	Ward box/Trsp box
60MP60500168 non vented inf set 180 cm	10/50



All products are single packed in peel-pack  
Trsp box = min qty/order

<p>Infusion set 20 drps aqua dest = 1 ml</p> <p>Protection cap: PE</p> <p>Hydrophobic air filter: Polyamide</p> <p>Spike, roller clamp and luer: ABS</p> <p>Drip chamber: PP</p> <p>Infusion filter: 15my Polyamide</p> <p>Tube: PVC eller Polyefine Alloy (XS118)</p>	<p>3-way stopcocks:</p> <p>Lipid resistant Polycarbonate, 72 Hrs,</p> <p>Back check valve: Silicon</p> <p>Clave: Polyester, Silicon and ABS</p>	<p>Sterilization method: ETO</p> <p>ISO 8536-4 and 8536-12 EN 980</p> <p>Luer EN 20594-1</p> <p>MDD Rule: Rule 2</p> <p>MDD Class: IIa</p> <p>CE: According to existing directives</p>
--	---	--



# INFUSION SETS - TIVA

- Infusion sets with integrated manifolds, back check valves and patient connections
- Suitable for wards where parallel infusions is common, also adapted for Total IntraVenous Anaesthesia (TIVA)
- The sets are available with 2 or 3 stopcocks, made of a lipid resistant material
- Available in PVC free material



Colour coded manifolds gives easy and secure overview of the different infusions

## Non vented, 180 cm inf set

Art no	Description	Ward box/Trsp box 10/50
60MR035021	2-way manifold, 1 back check valve	
60MR035031	3-way manifold, 1 back check valve	

Art no	Description PVC-free	Ward box/Trsp box 10/50
60MR03502	2-way manifold, 3 back check valve	
60MR03503	3-way manifold, 4 back check valve	

## TIVA-sets complete, incl manifolds, back check valves and extension tube for syringe pump, 150 cm inf set

Art no	Description PVC-free	Ward box/Trsp box 5/50
60MR03502150	2-way manifold, 3 back check valve	
60MR03503150	3-way manifold, 4 back check valve	



All products are single packed in peel-pack  
Trsp box = min qty/order

<u>Infusion set 20 drps aqua dest = 1 ml</u> Protection cap: PE and Polyamid Hydrophobic, fluid filter: Polyamide Spike, roller clamp, luer: ABS Drip chamber: PP Tube: Transp PVC or PVC-free Polyefine alloy (XS118), ID/OD 3,0/4,1 mm	<u>3-way stopcocks:</u> Lipid resistant Polycarbonate, 72 Hrs, pressure resistant 2,75 bar <u>Back check valve:</u> Silicon and ABS, dead space: 0,18 ml, Pressure resistant up to 4 bar (=58 psi)	Sterilization method: ETO  ISO 8536-4 and 8536-12 EN 980 Luer EN 20594-1 MDD Rule: Rule 2 MDD Class: IIa CE: According to existing directives
---	---	---





## 3-WAY STOPCOCKS, MANIFOLDS AND EXTENSION TUBING

### 3-Way stopcock PVC-free

Polycarbonate, lipid resistant 72 hrs, Pressure tested 2,75bar.  
Rotating luer-lock

#### Art no

60MP0436

Ward box/40 pcs Trsp box/600 pcs



### 3-Way stopcock with Extension tube

Manifold of lipid resistant Polycarbonate, 72 hrs, pressure tested 2,75 bar. Tube of pvc with rotating luer-lock. Protection cap with hydrophobic filter

#### Art no

Rotating luer-lock

60MR04100401 tube 15 cm

60MR04100430 tube 30 cm

60MR04100405 tube 50 cm

60MR04100410 tube 100 cm

Ward box/10 pcs Trsp box/100 pcs



#### With low volume tube:

60MR04100501, rotating luer-lock

DEHP-free PVC tube 15 cm. ID 1,0 mm

Ward box/10 pcs Trsp.box/100 pcs

### 3-Way stopcock with Extension tube PVC-free

With rotating luer-lock. Manifold of lipid resistant Polycarbonate, 72 hrs, pressure tested 2,75 bar. Protection cap with hydrophobic filter

#### Art no

60MR04100610 tube 10 cm

60MR04100650 tube 50 cm

60MR04100601 tube 100 cm

Ward box/10 pcs Trsp box/100 pcs

### Manifolds on plate PVC-free

Lipid resistant Polysulfone. 72 hrs. Pre glued tube, 180 cm with rotating luer-lock. Protection cap with hydrophobic filter

#### Art no

60MR04103180 3 pcs 3-way stopcocks

60MR04105180 5 pcs 3-way stopcocks

Ward box/10 pcs Trsp box/50 pcs





## Manifolds PVC-free

Lipid resistant Polycarbonate. 72 hrs.

Glued, rotating luer-lock. Protection cap with hydrophobic filter

### Art no

60MR04102000 2 pcs 3-way stopcocks

60MR04104000 4 pcs 3-way stopcocks

Ward box/40 pcs Trsp box/600 pcs



### Art no

60MR04104015 4 pcs 3-way stopcocks 15 cm tube

60MR04104015V 4 pcs 3-way stopcocks 15 cm tube, back check valve

Ward box/10 pcs Trsp box/100 pcs

## Extension tube for infusion

Fixed luer-lock, PVC or PVC-free. ID/OD 3/4, 1 mm

Protection cap with hydrophobic filter

### Art no

60MR4108303 length 30 cm PVC

60MR04108310 length 100 cm PVC

60MR04106310 length 100 cm PVC-free

60MR04108317 length 200 cm PVC

Ward box/10 pcs Trsp box/100 pcs

## Extension tube for syringe pump

Material: DEHP-free PVC. ID/OD 1,0/2,0 mm

### Art no

60MP1501 150 cm, priming volume 1,2ml, M/F

60MP2001 200 cm, priming volume 1,6 ml, M/F

60MP3001 300 cm, priming volume 2,4 ml, M/F

Ward box/10 pcs Trsp box/100 pcs

### Art no

Material: PE. ID/OD 1,0/2,0 mm

60MP150PE 150cm, priming volume 1,2ml, MF

60MP200PE 200cm, priming volume 1,6ml, MF

Ward box/10 pcs Trsp box/50 pcs

All products are single packed in peel-pack

Trsp box = min qty/order

Sterilization method: ETO

ISO 8536-4 and 8536-12 EN 980  
Luer EN 20594-1  
MDD Rule: Rule 2  
MDD Class: IIa  
CE: according to existing directive



# INFUSION ACCESSORIES

## Transfer spike

- Transfer spike for all types of Infusion/additive bottles
- Single use, sterile single packed in peel pack



Art no	Description	Ward box / trsp box
60MP04000380	Transfer spike	30/450

Material: ABS

## Infusion filter

Applicable for Taxol- and/or Remicade-infusions

- PVC-free, 0,2 my with 20 cm tube
- Back check valve
- Single use, sterile single packed in peel pack



Art no	Description	Ward box / trsp box
60MP04100168	Infusion filter with back check valve	10/100

60MR04100168	Infusion filter	10/100
--------------	-----------------	--------

Materials:

Filter membrane: Polysulphone, Filter house: acrylic, LuerLock: ABS, tubing: Polyefine, Back check valve: Silicone

Trsp box = min order qty

Sterilization: EO  
ISO 8536-4 and 8536-12 EN 980

MDD Rule: Rule 2  
MDD Class: IIa  
CE: Acc to excisting directives



# TRANSFUSION SETS

- Extra large filter area, 16,3 cm<sup>2</sup>  
Gives an effective filtration and allows a high transfusion-flow rate
- Available in different tube lengths
- PVC or PVC-Free material
- Protection cap with hydrophobic filter simplifies filling of sets and reduces contamination risk
- All models with rotating luer-lock
- Available with 3-way stopcock



## Transfusion sets

Art no	Description	Ward box/Trsp box 10/100
60MR60500200	180 cm inf set	
60MR60500300	180 cm inf set, two spikes	
60MR60500210	220 cm inf set	
60MR60500201	320 cm inf set	

Art no	Description PVC-free	Ward box/Trsp box 10/100
60MR60500202	180 cm inf set	
60MR60500212	220 cm inf set	
60MR60500203	320 cm inf set	

## Transfusion sets incl 3-way stopcock and extension tube

3-way stopcock in lipid resistant Polycarbonate

Art no	Description	Ward box/Trsp box 10/100
60MR60500214	220 cm inf set + 3-way stopcock, 15 cm	
60MR60500230	180 cm inf set + 3-way stopcock, 30 cm	

Art no	Description PVC-free	Ward box/Trsp box 10/100
60MR60500231	180 cm inf set + 3-way stopcock, 15 cm	

All products are single packed in peel-pack

Trsp box = min qty/order

Transfusion set. 20 drps aqua dest = 1 ml Protection cap: PE Hydrophobic filter: Polyamide Spike, roller clamp and luer: ABS Drip chamber: PP	Transfusion filter: 200 my. Polyamide Tube: PVC eller Polyfine Alloy (XS118)  3-way stopcock: Lipid resistant Polycarbonate, 72 Hrs	Sterilization method: ETO ISO 1135-4 EN 980 Luer EN 20594-1 MDD Rule: Rule 2 MDD Class: IIa CE: According to existing directives
---	---	---

Medi plast AB  
P.O Box 9504  
SE-200 39 MALMÖ  
SWEDEN

mediplast.info@mediplast.com  
www.mediplast.com



Phone: +46 40-671 23 00  
Fax: +46 40-671 23 10



## MEDIFUSOR - PRESSURE CUFF

For infusion and transfusion therapy where it is important with constant high flow rate, I e in the ER, ICU and OR's typically.

- Available in Nylon coated PU (Polyurethane) which fits both 500 ml and 1000 ml infusion bags.
- Also available in transparent PU which also fits for 500ml resp. 1000 ml infusion bags. The transparent version simplifies control of fluid level in the bag.

Art. No	Description	Ward box/Trsp. box
604805	Grey nylon 500-1000 ml	1/1
604807	Transparent 500 ml	1/1
604810	Transparent 1000 ml	1/1
Spare part		
604811	Pump with 3-v valve fits all pressure cuffs above	1/1

Material:  
Art No 604805 Nylon coated PU  
Art No 604807 Transparent PU  
Art No 604810 Transparent PU  
Tubes and pump in PVC



MDD Class: I  
CE: Acc to existing directives

Medi plast AB  
P.O. Box 9504  
SE-200 39 MALMÖ  
Sweden

mediplast.info@mediplast.com  
www.mediplast.com



Phone: +46 40-671 23 00  
Fax: +46 40-671 23 10



















## SUCTION TUBING











- A wide assortment of PVC and PVC free tubes
- Supple, resistant to folding and durable under vacuum pressure
- Available in 3 different packs:
  - Suction Tube Sterile - double packed in peel pack
  - Suction Tube Hygienical – hygienically packed in peel pack
  - Suction Tube on Roll – hygienically packed on roll

Material: PVC  
PVC free: XS119  
MDD Rule: Rule 2  
MDD Class: I, Is  
CE: Acc. to existing directives









## Suction tube, standard, serrated surface

Art. No.	Connector	Connector	Length	ID x OD	Hardness in Shore	Ch	Sterile	Box qty
6099900			200 cm	5,6 x 8,2	75	25	EO	50
6099903			250 cm					50
6099905			350 cm					40
6099910			200 cm	5,6 x 8,2	75	25	EO	50
6099915			350 cm					40
6099920			200 cm	5,6 x 8,2	75	25	EO	50
6099925			350 cm					40
6099930			200 cm	5,6 x 8,2	75	25	EO	50
6099935			350 cm					40
6099940			200 cm	5,6 x 8,2	75	25	EO	50
6099943			200 cm				Non-sterile	50
6099945			350 cm				EO	40
6099950			200 cm	7 x 10	75	30	EO	50
6099952			350 cm					40
6099959			200 cm	5,6 x 8,2	75	25	EO	50
6099957			200 cm	5,6 x 8,2	75	25	EO	50
6099958							Non-sterile	

## Suction tube, standard, serrated surface - PVC-Free

6088800			200 cm	5,6 x 8,2	75	25	EO	50
6088805			350 cm					40
6088810			200 cm	5,6 x 8,2	75	25	EO	50
6088840			200 cm	5,6 x 8,2	75	25	EO	50
6088850			200 cm	6 x 8	75	25	EO	50
6088860			200 cm	6 x 8	75	25	EO	50


## Suction tube, special, smooth surface

Art. No.	Connector	Connector	Length	ID x OD	Hardness in Shore	Ch	Sterile	Box qty
6099977			200 cm	5,6 x 8,2	75	25	EO	50
6099979			350 cm					40
6099980			200 cm	5,6 x 8,2	75	25	Non-sterile	1/50
6099982			200 cm	5,6 x 8,2	75	25	Non-sterile	1/50
6099983			350 cm					
6099985			200 cm	5,6 x 8,2	75	25	Non-sterile	1/50









## Bubble tube on roll, smooth tube

Art. No.	Connector	Connector	Length	ID x OD	Colour	Ch	Sterile	Box qty
6099970	-	-	30 m	3 x 5	Green	15	Non-sterile	1/5
6099973	-	-	30 m	5 x 7	Transparent	21	Non-sterile	1/5
6099975	-	-	30 m	3 x 5	Transparent	15	Non-sterile	1/5
6099976			30 m	7 x 9	Transparent	30	Non-sterile	1/5

## Suction tube, high flow aspiration, 360° rotating handle


Art. No.	Connector	Connector	Length	ID x OD	Hardness in Shore	Ch	Sterile	Box qty
6070002			200 cm	9,6 x 14,6	75	45	EO	15

## Suction tube, ENT, smooth tube, extra supple

Art. No.	Connector	Connector	Length	ID x OD	Hardness in Shore	Ch	Sterile	Box qty
6099990			200 cm	6 x 8	60 - 63	25	EO	50
6099992			350 cm					40
6099994			200 cm	6 x 8	60 - 63	25	EO	50
6099996			350 cm					40
6099987			200 cm	6 x 8	60 - 63	25	EO	50
6099998			200 cm					50
6099999			200 cm	6 x 8	60 - 63	25	EO	50
			350 cm					40

## Suction tube special, with out connectors, transparent PVC, smooth tube

Art. No.	Connector	Connector	Length	ID x OD	Hardness in Shore	Ch	Sterile	Box qty
6099960	-	-	30 cm	7 x 10	68 - 70	30	EO	50
6099963	-	-	110 cm	7 x 10	68 - 70	30	EO	50
6099965	-	-	150 cm	7 x 10	68 - 70	30	EO	50
6099967	-	-	180 cm	7 x 10	68 - 70	30	EO	50

Art. No.	Connector	Length	ID x OD	-	-	Sterile	Box qty
6099981						Non-sterile	10





<b>Art.no.</b>		<b>Ward,/Trsp box</b>
----------------	--	-----------------------

Suction tube on roll

6091620043	4,0 x 6,0 mm	50/50
6091620050	5,0 x 8,0 mm	50/50
6091620068	6,0 x 9,0 mm	50/50
6091620076	7,0 x 10,0 mm	50/50
6091620084	8,0 x 11,0 mm	40/40
6091620090	8,0 x 12,0 mm	40/40
6091620092	9,0 x 13,0 mm	30/30
6091620100	10,0 x 14,0 mm	30/30

trsp box = min qty / order




# SUCTION CATHETERS

## TYPE MÜLLY

**STRAIGHT, COUDÉ (ANGLED TIP)**

**AND WITH VACUUM CONTROL**



- Colour coded connections for easy identification CH size
- Coudé catheter enables choice of direction for left or right bronchi. Used for drainage of bronchi
- Suction catheter with vacuum control gives an atraumatic suction treatment and enables finger tip control and intermittent suction technique
- Catheter lengths straight – 25 cm, 40 cm, 53 cm
- Catheter length coudé (angled) – 53 cm
- Catheter length for suction catheter with vacuum control – 35 cm for models with cm graduation and 53 cm for models without cm graduation
- Material: Medical Grade PVC / PVC-free
- ETO sterilised

Art No	Size (CH)	Ward box / Trsp box
--------	-----------	---------------------

### Suction catheter straight, 25 cm

60130525	No 05	grey	100/600
60130625	No 06	green	100/600
60130825	No 08	blue	100/600
60131025	No 10	black	100/600
60131225	No 12	white	100/600
60131425	No 14	light green	100/600
60131625	No 16	orange	100/600

### Suction catheter straight, 40 cm

60130540	No 05	grey	100/600
60130640	No 06	green	100/600
60130840	No 08	blue	100/600
60131040	No 10	black	100/600
60131240	No 12	white	100/600
60131440	No 14	light green	100/600
60131640	No 16	orange	100/600

### Suction catheter straight, 53 cm

60130653	No 06	green	100/600
60130853	No 08	blue	100/600
60131053	No 10	black	100/600
60131253	No 12	white	100/600
60131453	No 14	light green	100/600
60131653	No 16	orange	100/600
60131853	No 18	red	100/600

trsp box = min qty / order

Material: PVC, SH A 78  
MDD Rule: Rule 1  
MDD Class: I sterile  
CE: Acc. existing directive

Art No	Size (CH)	Ward box / Trsp box
--------	-----------	---------------------

**Suction catheter straight, 25 cm, PVC-Free**

60880525	No 05	grey	100/600
60880625	No 06	green	100/600
60880825	No 08	blue	100/600
60881025	No 10	black	100/600
60881225	No 12	white	100/600
60881425	No 14	light green	100/600
60881625	No 16	orange	100/600

**Suction catheter straight, 40 cm, PVC-Free**

60880540	No 05	grey	100/600
60880640	No 06	green	100/600
60880840	No 08	blue	100/600
60881040	No 10	black	100/600
60881240	No 12	white	100/600
60881440	No 14	light green	100/600
60881640	No 16	orange	100/600

**Suction catheter straight, 53 cm, PVC-Free**

60880653	No 06	green	100/600
60880853	No 08	blue	100/600
60881053	No 10	black	100/600
60881253	No 12	white	100/600
60881453	No 14	light green	100/600
60881653	No 16	orange	100/600

trsp box = min qty / order

Material: PVC-free XS118  
MDD Rule: Rule 1  
MDD Class: I sterile  
CE: Acc. existing directive

Art No	Size (CH)	Ward box / Trsp box
--------	-----------	---------------------

### Suction catheter type Coudé, 53 cm

60170853	No 08	blue	100/600
60171053	No 10	black	100/600
60171253	No 12	white	100/600
60171453	No 14	light green	100/600
60161653	No 16	orange	100/600
60161853	No 18	red	100/600

### Suction catheter with vacuum control, cm-graduation, 35 cm


60151805	No 05	grey	100/600
60151806	No 06	green	100/600
60151808	No 08	blue	100/600

### Suction catheter with vacuum control, 53 cm

60150653	No 06	green	100/600
60150853	No 08	blue	100/600
60151053	No 10	black	100/600
60151253	No 12	white	100/600
60151453	No 14	light green	100/600
60151653	No 16	orange	100/600
60151853	No 18	red	100/600

trsp box = min qty / order

Material: PVC, SH A 78  
MDD Rule: Rule 1  
MDD Class: I sterile  
CE: Acc. existing directive



## DUODENAL-, STOMACH- AND RECTAL TUBES

- Tubes in Medical Grade PVC with colour coded connectors for easy identification of the CH size of the tube
- Duodenal tubes type Ryle and Levin with distance markings at 45/55/65/75 cm from the distal end
- Stomach tubes with distance markings at 45/55/65 cm from the distal end
- ETO sterilised
- Sterile, single packed in peel-pack

Material: PVC, SH A 75  
MDD Rule: Rule 1  
MDD Class: I sterile  
CE: Acc. existing directives

<b>Art No</b>	<b>Size (Ch)</b>	<b>Colour</b>	<b>Ward,/Trsp box</b>
---------------	------------------	---------------	-----------------------

**Duodenalsonder**

Levin - with x-ray line, 125 cm

60612108	Nr 08	blue	25/300
60612110	Nr 10	black	25/300
60612112	Nr 12	white	25/300
60612114	Nr 14	green	25/300
60612116	Nr 16	orange	25/300
60612118	Nr 18	red	25/300

Ryle - with barium sulphate tip

With x-ray line, 125 cm

60612712	Nr 12	white	25/300
60612714	Nr 14	green	25/300
60612716	Nr 16	orange	25/300

**Stomach tubes**

Standard - with x-ray line, 80 cm

60612616	Nr 16	orange	25/300
60612618	Nr 18	red	25/300
60612620	Nr 20	yellow	25/300
60612622	Nr 22	pink	25/300
60612624	Nr 24	blue	25/300
60612626	Nr 26	maroon	25/300
60612628	Nr 28	grey	25/300
60612630	Nr 30	green	25/300

**Rectal tubes, 40 cm**

60614122	Nr 22	pink	25/300
60614124	Nr 24	blue	25/300
60614126	Nr 26	maroon	25/300
60614128	Nr 28	grey	25/300

trsp box = min qty / order



## CATHETERS NELATON & TIEMANN

- For intermittent bladder drainage
- Single use, Medical Grade PVC
- Length, 40 cm male, and 23 cm female
- Sterile, single packed in peel-pack
- Nelaton catheter 40 cm available in 6 ch - 20 ch
- Nelaton catheter 23 cm available in 6 ch - 16 ch
- 2 opposite side eyes for maximum drainage
- Tiemann catheters (bent tip) 40 cm available in 8 ch - 18 ch

Material: PVC, SH A 75  
MDD Rule: Rule 1  
MDD Class: I sterile  
CE: Acc. existing directives



<b>Art No</b>	<b>Size (CH)</b>	<b>Conn. Color</b>	<b>Ward/Trsp box</b>
---------------	------------------	--------------------	----------------------

**Nelaton Catheters 40 - cm (male)**

60613106	No 06	green	100/600
60613108	No 08	blue	100/600
60613110	No 10	black	100/600
60613112	No 12	white	100/600
60613114	No 14	light green	100/600
60613116	No 16	orange	100/600
60613118	No 18	red	100/600
60613120	No 20	yellow	100/600

**Nelaton Catheters 23 - cm (female)**

60613406	No 06	green	100/600
60613408	No 08	blue	100/600
60613410	No 10	black	100/600
60613412	No 12	white	100/600
60613414	No 14	light green	100/600
60613416	No 16	orange	100/600

**Tiemann Catheters - 40 – cm**

60613608	No 08	blue	100/600
60613610	No 10	black	100/600
60613612	No 12	white	100/600
60613614	No 14	light green	100/600
60613616	No 16	orange	100/600
60613618	No 18	red	100/600

ward box = minimum order qty



## GUEDEL AIRWAYS

- Colour coded Guedel Airways for easy size identification
- Hygienically packed in bulk, set or in single packed peel pack
- Made of Polyethylene

Art. No.	No	Length mm	Color	Ward box/Trsp. box
60605104	000	40 mm	pink	100/1000
60605105	00	50 mm	blue	100/1000
60605106	0	60 mm	black	100/1000
60605107	1	70 mm	white	100/1000
60605108	2	80 mm	green	100/1000
60605109	3	90 mm	yellow	100/1000
60605110	4	100 mm	red	100/1000
60605111	5	110 mm	orange	100/1000
60605104-1	000	40 mm	pink	1/100/1000
60605105-1	00	50 mm	blue	1/100/1000
60605106-1	0	60 mm	black	1/100/1000
60605107-1	1	70 mm	white	1/100/800
60605108-1	2	80 mm	green	1/100/700
60605109-1	3	90 mm	yellow	1/100/700
60605110-1	4	100 mm	red	1/100/600
60605111-1	5	110 mm	orange	1/100/500
60605112	000-5		set in box	1/1

Ward box = minimum order qty

Material: Polyethylene MDD Rule: Rule 1  
MDD Class: I  
CE: Acc. existing directive



# TUBE GUARD

## BITE BLOCK FIXATION OF ET-TUBES

- Safe, hygienic and easy to use
- Easily squeezed onto the ET tube
- Occupies minimum space in the mouth
- No force on the tongue back and down
- Patient friendly
- Simplifies the mouth hygiene routine
- No risk for the cuff tube to be squeezed thanks to the unique design
- Made of Polypropylene
- Latexfree



*TubeGuard – unique fixation and bite block for ET-tubes*

<b>Art no</b>	<b>Size</b>	<b>Ward box/Trsp box</b>
601060	Fits ET-tube 6 and 6,5	50/1200
601070	Fits ET-tube 7, 7,5 and 8	50/1200

Ward box = minimum order qty



## MUCUS TRAP SET


Mucus trap set for mucus sample collection in a closed system, which makes handling easy and safe.

- Complete set with suction catheter with softly rounded, CH 14, tip and female cone for connection with suction tube.
- Transparent 25 ml-graduated trap/container with extra lid.
- “All in one”, timesaving pack, where only a suction tube needs to be connected.
- Single packed in peel-pack, ETO sterilized.

Art. No.	Description	Ward box/Trsp. box
6066560	Mucus trap set	50/200
6066561	Mucus trap set, PVC-free	50/200

Ward box = min order qty

Material: Suction catheter/tube: PVC / Polyefin Alloy Trap/container: Polypropylene Lid: ABS Filter: Hydrofobic		MDD Rule: 5 MDD Class: IIa CE: Acc to excisting directives
---	--	--



## MUCUS TRAP SET TRACHEAL

Mucus trap set for mucus sample collection in a closed system, which makes handling easy and safe, specialized for trachea specimen collection.

- “All in one” timesaving set for connection in between the normal suction tube and suction device.
- Transparent 10 ml-graduated trap/container with extra lid.
- Complete set with female cone for connection with suction device with male connector ( f. ex. bronchoscopy suction etc etc), and one with suction control handle for connection to the suction tube.
- Single packed in peel-pack, ETO sterilized.
- Self adhesive label for patient data is enclosed.

<b>Art. No</b>	<b>Description</b>	<b>Ward box/Trsp. box</b>
6066565	Mucus trap set Tracheal	50/200
6066566	Mucus trap set Tracheal PVC-free	50/200

Ward box = min order qty

Material: Suction catheter/tube: PVC / Polyefin Alloy Trap/container: AS Copolymer Lid: Polypropylene Filter: Hydrofobic		MDD Rule: 5 MDD Class: IIa CE: Acc to existing directives
--	--	---



## PVC TUBING ON ROLL

- Available in different lengths and diameters
- Hygienically packed, non-sterile
- Soft and Supple, Sh-A of 73

<b>Art.no.</b>	<b>mm</b>	<b>Roll Length</b>
6091620043	4,0 x 6,0	50 m
6091620050	5,0 x 8,0	50 m
6091620068	6,0 x 9,0	50 m
6091620076	7,0 x 10,0	50 m
6091620084	8,0 x 11,0	40 m
6091620090	8,0 x 12,0	40 m
6091620092	9,0 x 13,0	30 m
6091620100	10,0 x 14,0	30 m

1 roll = min. order qty

Material: PVC  
MDD Rule: Rule 1  
MDD Class: I non-sterile  
CE: Acc to regulation



## SUCTION TUBE FOR LAPORASCOPY

- For direct connection to Smoke Evac-200L
- Single use, sterile, single packed in peel pack (or bulk: 60M1010011)
- Made of PVC, 7 x 10 mm (ID/OD)

<u>Art.no.</u>		<u>Ward box/trsp.box</u>
60TB1010011	4,5 m	100/100
60M1010011	4,5 m	100/100

ward box = minimum qty / order

Material: PVC  
MDD Rule: Rule 3  
MDD Class: I sterile, or none sterile  
CE: Acc. to regulation



Material: PE/PVC  
/SAN  
MDD Rule: Rule 7  
MDD Class: IIa  
CE: Acc to regulation

## SINGLE USE, ORTHOPEDIC SUCTION SET

60QP09061 Containing:

- Suction tube fem – fem, 3,5m, PVC
- Transparent handle including filter, SAN (StyreneAcrylNitril)
- Extra filter, SAN
- Suction tip, two ways, 23 cm, PVC
- Suction tip, two ways, 15cm, PVC

<u>Art.no</u>	<u>description</u>	<u>ward.box/trsp.box</u>
60QP09061	Orthopedic Suction Set	20/20
60QP090611	As above, no Suction Tube	20/20
Accessories, single packed (sterile)		
60QP090612	Suction tip 23 cm	20/80
60QP090613	Suction tip, 15 cm	20/80
60QP090614	Filter	20/80





## PLEURA SET FOR THORACENTESIS

- A closed system which is reliable and easy to use
- Drainable 2 liter bag
- The volume of the punctured fluids can be read of the 60ml syringe or the bag
- The 60 ml syringe connects with the 3 way stop cock to the tube
- The large inner diameter of the tubes makes it possible for drainage of large quantities of fluids and secretion
- Sterile, single packed in peel-pac
- The set contains two separate canulas one 1,6 x 80 mm and one 2,0 x 80 mm

**Art.no.**


601-80

**ward.box./trsp.box**

20/20

Min qty = 20pcs

Material: PVC/ Medical Grade Stainless Steel  
(virgin steel)  
MDD Rule: Rule 7  
MDD Class: IIa  
CE: Acc to regulation



## GASTRIC LAVAGE SET VENTRICULAR SET

- The set consist's of two transparent bags, tubing, and a connection tube specially adapted to the gastric probe
- The volume of the two bags are 3,5 respectively 5 liters
- Total tubing length of 3 m
- The set is delivered hygienically packed

**Art.no.**  
6054010101

**ward.box /trsp box**  
10/10


Material: PVC  
MDD Rule: Rule 1  
MDD Class: I non-sterile  
CE: Acc to regulation



**NEW!**  
**Small dispenser**

## VOMIT BAG AND DISPENSER

- Simple handling
- 1,5 litre bag with graduation
- Easy to close with a simple handgrip
- Small contamination risk
- Single use, packed in bags of 50 pcs
- Wall mounted metal dispenser system



Article No		Ward box/Trsp box
601132K	Vomit bag	50/400
601133K	Dispenser for 50 bags	1/1
601131K	Dispenser for 25 bags	1/1

Trsp box = minimum order qty

Material: Polyethylene PE  
MDD Class: I  
CE: Acc. existing directive



## CANULAS, SPECIAL AND LUMBAL/SPINAL

- Hypodermic needles with specially bevelled points
- For numerous applications.
- Single use, sterile and single packed
- With metal hub and luer-lock
- With mandrin
- Protective caps made of polypropylene

**Special needles**

<b>Art.no.</b>	<b>mm</b>	<b>G</b>	<b>Bevel</b>	<b>ward.box/trsp.box</b>
6062590715	0,7 x 150	22	14	20/600
6062590720	0,7 x 200	22	14	20/600
6062590811	0,8 x 110	21	14	20/600
6062600093	0,9 x 110	20	14	20/600
6062620091	0,9 x 120	20	14	20/600
6062630090	0,9 x 150	20	14	20/600
6062630192	0,9 x 200	20	14	20/600
6062600135	1,3 x 110	18	20	20/600
6062610126	1,2 x 115	18	14	20/600
6062600184	1,8 x 110	15	20	20/600
6062610209	2,0 x 115	14	20	20/600
6062600259	2,5 x 110	13	20	20/600

**Lumbal/Spinal needles**

<b>Art.no.</b>	<b>mm</b>	<b>G</b>	<b>Bevel</b>	<b>ward.box/trsp.box</b>
6062570478	0,7 x 40	22	24	20/600
6062570775	0,7 x 70	22	24	20/600
6062570973	0,7 x 90	22	24	20/600
6062570980	0,7 x 120	22	24	20/600
6062570999	0,9 x 90	20	24	20/600
6062571000	0,9 x 120	20	24	20/600
6062571020	0,9 x 150	20	24	20/600
6062571120	1,2 x 100	18	24	20/600
6062572020	1,2 x 150	18	24	20/600
6062510169	0,6 x 25	23	24	20/600

**Medioplast needles**

<b>Art.no.</b>	<b>mm</b>	<b>G</b>	<b>Color</b>	<b>ward.box/trsp.box</b>
6063560080	0,8 x 80	21	green	100/10.000
6063560205	2,1 x 80	14	green	100/10.000
6062640440 (dental)	0,4 x 40	27	grey	100/10.000

Material: PP/PE/Medical Grade Stainless Steel  
(virgin steel)  
MDD Rule: Rule 6  
MDD Class: IIa  
CE: Acc to regulation



# OPTIDYNAMIC


- Used to measure lumbal pressure
- Improved dynamic pressure sensitivity
- Improved possibility to detect minor spinal abstractions
- Easy reading
- The fluid manometer holds only 0,05 ml / 100 mm
- Can be used as a siphon drawing liquid into collecting tubes
- Luer-lock fittings

**Art.no.**                      **ward.box/trsp.box**

6061650008

25/225

ward.box = min.qty/order




Material: PE/PVC  
MDD Rule: Rule 1  
MDD Class: I sterile  
CE: Acc to regulation




## FORCEPS

- Autoclavable
- Radiopaque
- Lockable



**NOTE! To sterilize and store with jaws open.**

**A handy instrument for a wide field of applications.**



- The Mediplast® universal, disposable, lockable forceps for numerous applications.
- The instrument is a combination of forceps and clamps (max tube diameter 7 mm).
- Autoclavable and radiopaque.
- Forceps are suitable for inclusion in the hospital's own sets for different purposes.
- Supplied hygienically packed.

Product No.	Product	Ward/ Transport pack
6070000013	Forceps	250/1000

Material: PP  
MDD Rule: Rule 1  
MDD Class: I non-sterile  
CE: Acc to regulation





## PINCETTES

**MEDIPLAST®** can offer You two different Pincettes. A green, simpler pincette that is well suited for ex. holding a swab/ball to clean a wound. The white one has a finer tip and gives a more stable grip, suitable for both wound cleaning and for instance suture removal where the pincette needs to be more rigid. Both models are made of Polystyrene.

Art no	Description	Ward box	Trsp box
<b>601170</b>	Pincette green, non sterile. 50 pcs / plastic pouch.	50	300
<b>601171</b>	Pincette green, sterile. Single packed.	50	300
<b>601172</b>	Pincette white, non sterile. 50 pcs / plastic pouch.	50	300
<b>601173</b>	Pincette white, sterile. Single packed.	50	300



## 50/60ML SYRINGES

- Available with catheter cone or luer lock connector
- Single use, sterile and single packed in peel pouch

<u>Art.no.</u>		<u>ward.box/trsp.box</u>
6063156000	luer lock	60/360
60MP20002660	catheter cone	60/360

ward.box = min.qty/order

Material: PP  
MDD Rule: Rule 2  
MDD Class: IIa  
CE: Acc to regulation




## SJECT FLUSHER SYRINGE

- Multi purpose flusher 100 ml
- Sterile, single packed in peel-pack

<b>Art.no.</b>	<b>ward.box/trsp.box</b>
6061990011	20/40

ward.box = min.qty/order

Material: PP  
MDD Rule: Rule 1  
MDD Class: I sterile  
CE: Acc to regulation



## MEDIPLAST® AVLOX SYSTEM

**The Avlox® System is used to close all standard incisions in the chest, abdomen, neck and limbs**


- The Mediplast® Avlox System is a "gold standard" to use in thyroidectomy operations
- The system is non-traumatic and time saving
- The speed of application reduces operating time
- The sterile disposable format requires a minimum of preparation both for application and removal
- Scarring is reduced to a minimum
- Easy clip removal with minimal patient discomfort
- The clips are made of medical grade stainless steel and kept sterile in the applicator
- The material is polycarbonate and stainless steel medical grade

### **Mediplast Remover®**

Instrument for removal of the clips.

Material: Blue polycarbonate – harmless to the environment

<b>Art.no.</b>	<b>description</b>	<b>ward.box./trsp.box</b>
6070010130	Avlox	15/450
6070010510	Remover	20/300



Avlox, 12

Remover

Material: PC/Stainless Steel  
MDD Rule: Rule 4, Subrule 1  
MDD Class: I Sterile  
CE: Acc to regulation

Material: Nylon  
MDD Rule: Rule 1  
MDD Class: I Sterile  
CE: Acc to regulation



**LATEX-FREE**

## FLUSHING SETS

- **TUR-P SET MEDIPLAST®** urologic flushing set. Single use, sterile packed in inner PE pouch + outer peel pack. Y-connector with 2 spikes for flushing drip bags, roll clamp and large drip chamber with air guard enabling large flows without risk of air bubbles in the tube. For connection with universal connector in silicone to cystoscope or catheter.

Art.No: 60Y-147 20pcs/box  
 2-way Y-connector, 2 spikes for flushing drip bags,  
 150 cm PVC-tube (ID 4,8 x OD 6,8 mm)

- **FLUSHING SET MEDIPLAST®**  
 For flush drip in for instances laporascopies.  
 Single use, sterile packed in inner PE pouch + outer peel pack. With spike and roll clamp, universal connector of silicone for connection to troachard.

Art.No: 60Y-280 20pcs/box  
 1-way, 280 cm PVC-tube (ID 4,8 x OD 6,8 mm)

Material: PVC/ABS  
 MDD Rule: Rule 3  
 MDD Class: IIa  
 CE: Acc. to regulation




## CONNECTION TUBE FOR NEPHROSTOMY

- Connection Tube PVC 2,9 x 4,1 mm
- Green, soft, female end, fits all types of urine bags
- Rotating luer-lock in other end for easy connection
- Single use product. Sterile, packed in peel-pack

<b>Art.no</b>	<b>Description</b>	<b>ward.box/trsp.box</b>
60MR223115	15 cm, LL, female universal connector	10/300
60MR222915	15 cm, LL, female universal connector, 3 way stopcock, at LL end	10/300
60MR222815	15 cm, LL, female universal connector, 3 way stopcock, at universal end	10/300
60MR222830	30 cm, LL, female universal connector, 3 way stopcock, at universal end	10/300
6060300001	LL-25, 25 cm, LL female universal connector	10/300

Material: PVC, Polycarbonate, ABS  
MDD Rule: Rule 2  
MDD Class: I s  
CE: Acc to regulation



## CYSTOSCOPY TUBING

- For direct connection to cystoscope or troachard
- Single use, sterile, single packed in peel pack
- Made of PVC, 7 x 10 mm (ID/OD)

<b>Art.no.</b>		<b>Ward box/trsp.box</b>
6061626262	1,1 m	25/50
6061626263	1,5 m	25/50

ward box = minimum qty / order

Material: PVC  
MDD Rule: Rule 3  
MDD Class: IIa sterile  
CE: Acc. to regulation





# GLOVES





Medi plast is part of Procurator AB, a Swedish group of companies with focus on personal protection equipment, medical technical products and office and paper supplies. The group has around 300 employees and a total turnover of approx. 110 MEUR

Gloves for personal protection is the biggest product group in the Procurator Group, our yearly turnover in this product group is approx. 25 MEUR/Year and growing.

In the assortment today we have approx. 470 different models of gloves, of which 86 are of our own brand – 14 of these are single use gloves.

Medi plast sells these types of gloves to the healthcare and on export, Procurator sells the gloves suited for industry and trade in Scandinavia.



**SWEHAND** 



**SWEHAND**  **GRIPEX**

## Examination glove in LATEX - POWDER FREE

A tanned, yellow glove that is strong and supple. Rolled collar and textured surface. For use in hospitals, care centres, laboratories, home care and some medicine and food processing industry. Double chlorinated for minimal allergic reaction.

**AQL:** 1,5

### **Chemical Data/Material**

- Natural rubber, latex

### **Pull Tenacity/Elasticity:**

>7,5 N and > 700%

### **Classification / Standard**

- Product standard EN 455 part 1, 2 & 3
- Quality system ISO 9002, EN 46002
- CE marked, approving authority TÜV
- Protein content < 50 micrograms/g

**Length:** Sizes XS - XL - 240 mm

### **Thickness:**

Collar: 0,11 mm ( $\pm 0,02$  mm)

Inner hand: 0,16 ( $\pm 0,02$  mm)

<b>Circumference:</b>	Size XS	- 160 mm ( $\pm 5$ mm)
	Size S	- 170 mm ( $\pm 5$ mm)
	Size M	- 190 mm ( $\pm 5$ mm)
	Size L	- 210 mm ( $\pm 5$ mm)
	Size XL	- 230 mm ( $\pm 5$ mm)

**Colour:** Tanned yellow

## **Sizes and packaging information**

Art No	Size	Inner box	Transport box
40550453060	XS (6)	100	1000
40550453070	S (7)	100	1000
40550453080	M (8)	100	1000
40550453090	L (9)	100	1000
40550453100	XL (10)	100	1000

ver EN 070531





**SWEHAND**  **NIMEX**



## Examination glove in NITRILE - POWDER FREE

An examination glove that is extra durable and has a high threshold against puncture. Minimal risk for allergical reaction. Rolled collar and textured surface for better grip. For use in hospitals, care centres, laboratories, home care and some medicine and food processing industry. Free from Latex!

**AQL:** 1,5

**Pull tenacity/Elasticity:**

>7,5 N and > 700%

**Length:** Sizes XS - XL - 240 mm

**Thickness:**

Inner hand; 0,14 mm ( $\pm 0,02$  mm)

Collar; 0,11 ( $\pm 0,02$  mm)

### **Chemical Data/Material**

- Nitrile

### **Classification / Standard**

- Product standard EN 455 part 1, 2 & 3
- Quality system ISO 9002, EN 46002
- CE marked, approving authority TÜV

<b>Circumference:</b>	Size XS	- 150 mm ( $\pm 5$ mm)
	Size S	- 170 mm ( $\pm 5$ mm)
	Size M	- 190 mm ( $\pm 5$ mm)
	Size L	- 210 mm ( $\pm 5$ mm)
	Size XL	- 230 mm ( $\pm 5$ mm)

### **Sizes and packaging information**

Art No	Colour	Size	Inner box	Transport box
40531521060	Blue	XS (6)	100	1000
40531521070	Blue	S (7)	100	1000
40531521080	Blue	M (8)	100	1000
40531521090	Blue	L (9)	100	1000
40531521100	Blue	XL (10)	100	1000
40531523060	White	XS (6)	100	1000
40531523070	White	S (7)	100	1000
40531523080	White	M (8)	100	1000
40531523090	White	L (9)	100	1000
40531523100	White	XL (10)	100	1000

ver EN 070531





**SWEHAND**  **PEVEX**

## Examination glove in VINYL - POWDER FREE

A transparent examination glove, suitable for work in hospitals, care centres, laboratories, home care and some medicine and food processing industry. Free from latex!

**AQL:** 1,5

### **Chemical Data/Material**

- Poly-Vinyl-Chloride (PVC)

### **Pull tenacity / Elasticity:**

>9 MPA and > 350%

### **Classification / Standard**

- Product standard EN 455 part 1, 2 & 3
- Quality system ISO 9002, EN 46002
- CE marked, approving authority: TÜV

**Length:** Size XS - 245 mm ( $\pm 5$  mm)  
Size S - XL - 250 mm ( $\pm 10$  mm)

### **Thickness:**

Fingers; 0,06 – 0,15 mm ( $\pm 5$  mm)  
Inner hand; 0,08 – 0,15 mm ( $\pm 5$  mm)

**Circumference:**

Size XS	- 170 mm ( $\pm 5$ mm)
Size S	- 175 mm ( $\pm 5$ mm)
Size M	- 195 mm ( $\pm 5$ mm)
Size L	- 215 mm ( $\pm 5$ mm)
Size XL	- 235 mm ( $\pm 5$ mm)

**Colour:** Transparent

### **Sizes and packaging information**

Art No	Size	Inner box	Transport box
40560453060	XS (6)	100	1000
40560453070	S (7)	100	1000
40560453080	M (8)	100	1000
40560453090	L (9)	100	1000
40560453100	XL (10)	100	1000

ver EN 070531



**SWEHAND**  **PEVEX**

## **Examination glove VINYL - POWDER FREE - PHTHALATE FREE**

A transparent examination glove, internal coating and PU covered for better comfort. Minimal risk for allergic reaction. For use in hospitals, care centres, laboratories, home care and some medicine and food processing industry. Free from latex!

**AQL:** 1,5

**Pull tenacity/Elasticity:**

>9 MPA and > 350%

**Length:** Sizes XS - XL - 250 mm ( $\pm 10$  mm)

**Thickness:**

Fingers; 0,10 mm ( $\pm 0,03$  mm)

Collar; 0,10 mm ( $\pm 0,03$  mm)

Inner hand; 0,11 mm ( $\pm 0,03$  mm)

**Circumference:**

Size XS/6	- 150 mm ( $\pm 10$ mm)
Size S/7	- 170 mm ( $\pm 10$ mm)
Size M/8	- 190 mm ( $\pm 10$ mm)
Size L/9	- 210 mm ( $\pm 10$ mm)
Size XL/10	- 230 mm ( $\pm 10$ mm)

**Chemical Data/Material**

- Poly-Vinyl-Chloride (PVC)

**Classification / Standard**

- Product standard EN 455 part 1, 2 & 3
- Quality system ISO 9002 and EN 46002
- CE marked, approving authority: TÜV

**Colour:** Transparent

### **Sizes and packaging information**

Art No	Size	Inner box	Transport box
40560462060	XS (6)	100	1000
40560462070	S (7)	100	1000
40560462080	M (8)	100	1000
40560462090	L (9)	100	1000
40560462100	XL (10)	100	1000

ver EN 070531

Medioplast AB  
P.O Box 9504  
SE-200 39 MALMÖ  
SWEDEN

mediplast.info@mediplast.com  
www.mediplast.com



Phone: +46 40-671 23 00  
Fax: +46 40-671 23 10





**SWEHAND**  **PEVEX**

## Examination glove in STRETCH-VINYL - POWDER FREE

A white, very soft and supple examination glove, it has near latex feeling to it and offers 30% more elasticity than other vinyl gloves. The Stretch-quality enables better finger top feeling. For use in hospitals, care centres, laboratories, home care and some medicine and food processing industry. Free from latex!

**AQL:** 1,5

**Pull Tenacity/Elasticity:**

>11 MPA and > 450%

**Length:** Size XS - XL - 250 mm ( $\pm 10$  mm)

**Thickness:**

Fingers; 0,10 mm ( $\pm 0,03$  mm)

Collar; 0,10 mm ( $\pm 0,03$  mm)

Inner hand; 0,11 mm ( $\pm 0,03$  mm)

**Circumference:**

Size XS	- 150 mm ( $\pm 10$ mm)
Size S	- 170 mm ( $\pm 10$ mm)
Size M	- 190 mm ( $\pm 10$ mm)
Size L	- 210 mm ( $\pm 10$ mm)
Size XL	- 230 mm ( $\pm 10$ mm)

**Chemical Data/Material**

- Poly-Vinyl-Chloride (PVC)

**Classification / Standard**

- Product standard EN 455 part 1, 2 & 3
- Quality system ISO 9002 and EN 46002
- CE marked, approving authority: TÜV

**Colour:** Transparent White

### Sizes and packaging information

Art No	Size	Inner box	Transport box
40560460060	XS (6)	100	1000
40560460070	S (7)	100	1000
40560460080	M (8)	100	1000
40560460090	L (9)	100	1000
40560460100	XL (10)	100	1000

ver EN 070531





**SWEHAND**  **POLLY**



## Glove in POLYETHYLENE - POWDER FREE

A transparent glove with coinage surface for better comfort and grip. Suitable for hospitals, care centres, laboratories, home care and certain medicine and food processing industry. Latex free!

**Length:** Size 8 - 30 cm  
Size 10 - 30 cm  
Size 10 - 53 cm (40531112100)

**Thickness:** 0,025 mm

**Colour:** Transparent

### **Chemical Data/Material**

- Polyethylene

### **Classification / Standard**

- Product standard EN 420 Category 1
- CE marked
- ISO 9002

### **Sizes and packing information**

Art No	Size	Length	Inner box	Transport box
40531115080	8	30 cm	100	5000
40531115100	10	30 cm	100	5000
40531112100	10	53 cm	100	5000

ver EN 070531







**SWEHAND**  **POLLYSOFT**



## Glove in SAM-POLYMER - POWDER FREE

A transparent, pink glove with coinage surface for better comfort and grip. Minimum risk for allergical reaction. SAM-Polymer material makes the glove extra supple. For use in hospitals, care centres, laboratories, home care and some medicine and food processing industry. Latex free!

**Length:** Size 8 - 35 cm  
Size 9 - 35 cm  
Size 10 - 35 cm

**Thickness:** 0,035 mm

**Colour:** Transparent, Pink

### **Chemical Data/Material**

- SAM-Polymer; LDPE, LLPE, EVA

### **Classification / Standard**

- Product standard EN 420 Category 1
- CE marked
- ISO 9002

### **Sizes and packaging information**

Art No	Size	Length	Inner box	Transport box
40531113080	8	35 cm	100	5000
40531113090	9	35 cm	100	5000
40531113100	10	35 cm	100	5000

ver EN 070531





**SWEHAND**  **POLEX**

## Glove in POLYETHYLENE - POWDER FREE

A transparent glove with coinage surface. An extremely strong and durable glove with good comfort and grip. Minimum risk for allergical reaction. Suited very well for "wet work". For use in hospitals, care centres, laboratories, home care and for the food processing industry. An environmental alternative without Phthalate softeners.

**Length:** Size 7 - 10 - 28 cm

### **Chemical Data/Material**

- Polyethylene

**Thickness:** 0,08 mm

### **Classification / Standard**

- Product standard EN 420 Category 1
- CE marked
- ISO 9002

**Colour:** Transparent

### **Sizes and packaging information**

Art No	Size	Inner box	Transport box
40531114070	7	150	1500
40531114080	8	150	1500
40531114090	9	150	1500
40531114100	10	150	1500

ver EN 070531





## EXAMINATION GLOVES

### Examination Gloves in COPOLYMER - POWDER FREE

A transparent glove in copolymer fastened on a paper sheet. A soft glove with excellent finger tip feeling. Applicable for daily tasks in health care. Available sterile or hygienically packed.

**Length:** S - L - 28 cm  $\pm$  2%

**Thickness:** 0,03 mm

**Sterilization method:** Gamma radiation

**Colour:** Transparent

#### **Chemical data / Material**

Copolymer of polyethylene, LDPE

#### **Classification / Standard**

- EN ISO 13485:2003
- EN 552      • EN 455-1
- EN 556      • EN 455-2

#### **Sizes and Packings**

Art No	Size / Shape	Ward box	Trsp box (min PO qty)
	<b>Non-sterile glove on paper sheet</b>		
	<b>Single pc packed</b>		
604500	Small, nr 6 - 6 1/2	100 pcs	1000 pcs
604501	Medium, No 7 - 7 1/2	100 pcs	1000 pcs
604502	Large, No 8 - 8 1/2	100 pcs	1000 pcs
	<b>Sterile glove on paper sheet</b>		
	<b>Single pc packed</b>		
604510	Small, No 6 - 6 1/2	100 pcs	2400 pcs
604511	Medium, No 7 - 7 1/2	100 pcs	2400 pcs
604512	Large, No 8 - 8 1/2	100 pcs	2400 pcs
	<b>Sterile glove on paper sheet</b>		
	<b>Pair packed</b>		
604520	Small, No 6 - 6 1/2	50 pairs	1200 pairs
604521	Medium, No 7 - 7 1/2	50 pairs	1200 pairs
604522	Large, No 8 - 8 1/2	50 pairs	1200 pairs

ver EN 070531



## KNITTED GLOVE COTTON

### Glove of unbleached cotton with Latex Burls

An environmental friendly glove in unbleached cotton, with latex burls on the palm surface, which can be washed in maximum 40°C. Suitable for goods picking and transportation work, very good fit and very good grip. Colour coded sizes.

#### Technical data

**Length:** Size XS - 210 mm ( $\pm 5$ )  
 Size S - 220 mm ( $\pm 5$ )  
 Size M - 230mm ( $\pm 5$ )  
 Size L - 240 mm ( $\pm 5$ )  
 Size XL - 250 mm ( $\pm 5$ )

#### **Chemical Data/Material**

- Cotton
- Latex burls

**Colour:** Nature cotton, nature latex

#### **Classification / Standard**

- Product standard EN 420 Category 1
- CE marked
- ISO 9002

#### **Sizes and packaging information**

Art No	Size	Colour code	Inner box	Transport box
40530724060	XS (6)	Green	12 pairs	300 pairs
40530724070	S (7)	Lilac	12 pairs	300 pairs
40530724080	M (8)	Orange	12 pairs	300 pairs
40530724090	L (9)	Black	12 pairs	300 pairs
40530724100	XL (10)	Grey	12 pairs	300 pairs

ver EN 070531



## KNITTED GLOVE

### Glove of Cotton / Polyester

A thin knitted unbleached glove that can be washed in maximum 40°C. Suitable as inner glove, under protective glove.

#### **Chemical Data/Material**

- Cotton / Polyester

**Colour:** Nature cotton, unbleached

#### **Classification / Standard**

- Product standard EN 420 Category 1
- CE marked
- ISO 9002

Art No	Size		Inner box	Transp. box
40530692080	8 (27 cm)	Thin Cotton / polyester	12 pairs	600 pairs
40530692100	10 (28 cm)	Thin Cotton / polyester	12 pairs	600 pairs

### Glove of Cotton

A knitted unbleached glove that can be washed in maximum 40°C. Suitable as inner glove, under protective glove.

#### **Chemical Data/Material**

- Cotton

**Colour:** Nature cotton, unbleached

#### **Classification / Standard**

- Product standard EN 420 Category 1
- CE marked
- ISO 9002

Art No	Size		Inner box	Transp. box
40530698100	10 (27 cm)	Thin Cotton	12 pairs	600 pairs
40530697080	8 (23 cm)	Medium dense Cotton	12 pairs	600 pairs
40530697100	10 (27 cm)	Medium dense Cotton	12 pairs	600 pairs
40530615080	8 (23 cm)	Dense Cotton	12 pairs	600 pairs
40530615100	10 (25 cm)	Dense Cotton	12 pairs	600 pairs
40530618080	8 (24 cm)	Dense cotton, with wristlet	12 pairs	600 pairs
40530618100	10 (25 cm)	Dense cotton, with wristlet	12 pairs	600 pairs

ver EN 070531