



SPIGGLE & THEIS
Medizintechnik

MICRO SUCTION TUBES

- For Otological Surgery
- For Nasal Surgery

Disposable Suction Tubes for Otological Surgery

Art. No.	Length	Diameter	Colour	Connection
301004	60 mm	0,4 mm	grey	Luer
301007	60 mm	0,7 mm	black	Luer
301010	60 mm	1,0 mm	white	Luer
301014	60 mm	1,4 mm	red	Luer
302004	70 mm	0,4 mm	grey	Luer
302006	70 mm	0,6 mm	blue	Luer
302007	70 mm	0,7 mm	black	Luer
302008	70 mm	0,8 mm	green	Luer
302010	70 mm	1,0 mm	white	Luer
302012	70 mm	1,2 mm	pink	Luer
302014	70 mm	1,4 mm	red	Luer
302016	70 mm	1,6 mm	white	Luer
301020	70 mm	2,0 mm	green	Luer
303004	80 mm	0,4 mm	grey	Luer
303007	80 mm	0,7 mm	black	Luer
303014	80 mm	1,4 mm	red	Luer
303020	80 mm	2,0 mm	green	Luer
303025	80 mm	2,5 mm	transparent	Luer-Lock
303030	80 mm	3,0 mm	transparent	Luer-Lock
303040	80 mm	4,0 mm	transparent	Luer-Lock
304010	100 mm	1,0 mm	white	Luer
304014	100 mm	1,4 mm	red	Luer
304020	100 mm	2,0 mm	green	Luer
304025	100 mm	2,5 mm	transparent	Luer-Lock
304030	100 mm	3,0 mm	transparent	Luer-Lock

Disposable Suction Tubes for Nasal Surgery

Art. No.	Length	Diameter
3012020	120 mm	2,0 mm
3012025	120 mm	2,5 mm
3012030	120 mm	3,0 mm
3012035	120 mm	3,5 mm
3012040	120 mm	4,0 mm
3012050	120 mm	5,0 mm
3014020	140 mm	2,0 mm
3014025	140 mm	2,5 mm
3014030	140 mm	3,0 mm
3014035	140 mm	3,5 mm
3014040	140 mm	4,0 mm
3014050	140 mm	5,0 mm
3016020	160 mm	2,0 mm
3016025	160 mm	2,5 mm
3016030	160 mm	3,0 mm
3016035	160 mm	3,5 mm
3016040	160 mm	4,0 mm
3016050	160 mm	5,0 mm

Suction Handle

Art. No.	Description
301000	With tip control, Luer connection, sterile

Package
50 pcs./box



- Angled cannula made of stainless steel
- Luer-/ Luer-Lock-connection
- Atraumatic rounded cannula point
- Single sterile packaging
- Wide range of different sizes
- Packaging of 50 pieces per box



- Tip control
- Blue tube connection
- Angled cannula made of stainless steel
- Atraumatic rounded cannula point
- Single sterile packaging
- Packaging of 50 pieces per box

Another quality product selected and stocked by:

WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489
e enquiries@bamford.co.nz www.bamford.co.nz

When performance counts

