




SPIGGLE & THEIS
Medizintechnik



NASAL SPLINTS

- External Nasal Splints
- Internal Nasal Splints
- Septal Buttons

Another quality product selected and stocked by:

WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489
e enquiries@bamford.co.nz w www.bamford.co.nz

When performance counts





Internal Silicone Nasal Splint with Airway

- Integral airway allows post-operative nasal breathing
- Silicone prevents adhesion
- Sterile, ready to use
- Soft material for easy insertion and removal
- Punched suture holes
- Surface with lotus effect

Art.No.	Packaging
2030100	5 pairs / box
2030100-50	50 pairs / box



Internal Silicone Nasal Splint without Airway

- Sterile, ready to use
- Punched suture holes
- Soft material for easy insertion and removal
- Surface with lotus effect

Art.No.	Packaging
2030110	5 pairs / box



Internal Nasal Splint, Fluoroplastic

- Bivalve splint, provides septal support
- Fluoroplastic prevents adhesion
- Easy insertion and removal
- Punched suture holes
- Sterile, ready to use
- Available in two sizes and thicknesses

Art.No.	Size	Thickness	Packaging
2041025	standard	0.25 mm	5 pairs / box
2041025-S	standard	0.25 mm	10 pieces / box
2041050	standard	0.50 mm	5 pairs / box
2041050-S	standard	0.50 mm	10 pieces / box
2041075	standard	0.75 mm	5 pairs / box
2041075-S	standard	0.75 mm	10 pieces / box
2042025	large	0.25 mm	5 pairs / box
2042025-S	large	0.25 mm	10 pieces / box
2042050	large	0.50 mm	5 pairs / box
2042050-S	large	0.50 mm	10 pieces / box
2042075	large	0.75 mm	5 pairs / box
2042075-S	large	0.75 mm	10 pieces / box
2043050	64 x 30 mm	0.50 mm	5 pairs / box
2043025	64 x 38 mm	0.25 mm	5 pairs / box



Septal Button, Silicone

- For non-surgical closure of septal perforation
- Soft material for easy insertion and removal
- Silicone prevents adhesion
- No scabbing, epistaxis or other septal perforation symptoms
- Sterile, ready to use
- Two designs available

Art.-No.	Ø	Feature	Packaging
2031100	3.2 cm	one part	1 piece
2031110	3.2 cm	two parts	1 piece




External Aluminium Nasal Splints

External aluminium nasal splints belong to the basic equipment found in ENT-Surgery, Casualty and Accident and Emergency Units. The material properties offer user and patient many advantages over standard plaster of paris bandaging. They can be applied quickly and efficiently, provide stable protection and are optically much less conspicuous.

- Padded Aluminium Nasal Splint
- Easily malleable
- Self-adhesive back
- Four sizes available
- Self-adhesive Dorsum foam pad included



Art.No.	A	x	B	x	C	Packaging
2060100	50	x	33	x	28 mm	10 pieces / box
2060200	56	x	38	x	30 mm	10 pieces / box
2060300	64	x	40	x	36 mm	10 pieces / box
2060400	66	x	42	x	45 mm	10 pieces / box



External Aluminium Nasal Splints – Sets


Each box includes:

- 5 x nasal splints of the indicated dimensions
- 5 x foam protective bases
- 5 x alcohol swabs
- 5 x Skin Prep™
- 5 x nose strips
- 5 x Remove™ swabs

Equipment box (2060500), contents:

- 5 x alcohol swabs
- 5 x Skin Prep™
- 5 x nose strips
- 5 x Remove™ swabs

Art.No.	A	x	B	x	C	Packaging
2060100-K	50	x	33	x	28 mm	1 box
2060200-K	56	x	38	x	30 mm	1 box
2060300-K	64	x	40	x	36 mm	1 box
2060400-K	66	x	42	x	45 mm	1 box
2060500						1 box



External Thermoplastic Nasal Splints

Here, thermoplastic is used to describe a substance that changes its' consistency and stability under temperature influence. When heated in a water bath, the material changes its' colour from white to transparent and becomes malleable under the influence of warmth. The process is reversed by cooling i.e. the product sets in the modelled shape. Using such material for an external nasal splint allows individual fitting according to the patients anatomy.

- Much less time and work required than with plaster of paris
- Includes three-part foam bridge
- Five sizes in two different thicknesses
- Self-adhesive back (NOT 2051010 and 2050010)

Art.No.	A	x	B	x	C	Thickness	Packaging
2050100	56	x	38	x	30 mm	2.4 mm	10 pieces / box
2050200	64	x	40	x	36 mm	2.4 mm	10 pieces / box
2050300	66	x	42	x	45 mm	2.4 mm	10 pieces / box
2051100	56	x	38	x	30 mm	1.6 mm	10 pieces / box
2051200	64	x	40	x	36 mm	1.6 mm	10 pieces / box
2051300	66	x	42	x	45 mm	1.6 mm	10 pieces / box

Art.No.	A	x	B	Thickness	Feature	Packaging
2050000	80	x	50 mm	2.4 mm	cut to size	10 pieces / box
2050010	100	x	75 mm	2.4 mm	cut to size	10 pieces / box
2051000	80	x	50 mm	1.6 mm	cut to size	10 pieces / box
2051010	100	x	75 mm	1.6 mm	cut to size	10 pieces / box

