



NON CONFORMANCE module

MEDITRAX
SMART TRACKING SYSTEM



features : non conformance tracking




accountability & traceability

The Meditrax Non Conformance Module is a systematic, controlled manual documentation system for the accurate recording and control of all medical devices that have been rejected from the sterilisation process.

The module ensures that items rejected from the sterilization process are segregated, accounted for and visually identified as non-conforming, therefore removing them from potential use.

Clear identification of non-conforming items via the non conformance label ensures inadvertent use of potentially harmful products.



manual recording

components

The MEDITRAX® Non Conformance Module starter pack includes :

Qty	Description	Code
1	MEDITRAX® label applicator	MHG1800
1	non gassing ink roller	IRT/NG05
2 pads	non conformance sheets	TTS203
5 rolls	non conformance labels	TTS/201
1	non conformance area sign	TTS3005

implementation

MEDITRAX® Non Conformance Module includes at no extra cost ongoing :-

- staff training
- auditing of system
- quality assurance consultation
- maintenance program

summary

The MEDITRAX® Non Conformance Module benefits

- Tracks devices rejected from the sterilisation process
- Batch record data captured
- Instrument data captured
- Staff data linked to the above
- All materials supplied
- Easy to learn

quick reference

MEDITRAX® Non Conformance Module includes ;



BATCH RECORD



ITEM RECORD



STERILISER
DATA



NON CONFORMANCE



STAFF DATA
RECORDED

related MEDITRAX® products :









- Total Traceability System
- Batch Recall Plus II module
- Sterile Release module
- Patient record traceability module





NON CONFORMANCE module

consumables for MEDITRAX® Non Conformance Module :

Code		Description	Packaging	Order Qty
MHG1800		MEDITRAX® label applicator	1 each	
IRT/NG05		non gassing ink roller	1 each	
873170		standard ink roller	1 each	
TTS203		non conformance sheets	50/pad 1 pad	
MH105A/30		non conformance labels	700/roll	
TTS300S		non conformance area sign	1 each	

accountability & traceability

training

Austmel provides MEDITRAX® System Inservice at the point of implementation. Courses are run on a person-to-person basis ensuring rapid absorption. Advanced courses are designed to train the trainer where required. MEDITRAX® training systems comply with the AS/NZS 4187:2003 standards and documentation.

one-on-one MEDITRAX® customer support

Austmel Pty Ltd is able to offer one-on-one support during office hours on our 1800 074 085 service. This is a valuable service to healthcare personnel who require immediate technical support to use the system. This service is a complete one-on-one person to person level. Additionally IT Remote Access Tech-Support provides quick on-line diagnosis and repair of systems to clients requiring rapid response under maintenance contracts.

documentation

MEDITRAX® documentation is built on the philosophy of paper waste minimisation, taking into account policies and procedures relating to the processing of items being sterilised in healthcare facilities.

All aspects of the complex requirements are clearly and easily followed.

background

Austmel Pty. Ltd is an Australian owned company established in 1974. Austmel have been directly involved in the supply and development of Sterility Assurance products for over 30 years. Our unique products and systems extend into both the Health Care and Industrials markets, including :

- MEDITRAX® Accountability & Traceability Systems
- Biological and Chemical Sterilization Indicators
- Bowie-Dick Tests
- Sterilization Packaging
- Quality Assurance Consultation
- TMI Electronic Data Loggers
- Surgical Disposables

T 1800 074 085 customer support

summary

The MEDITRAX® Non Conformance Module benefits

- Tracks devices rejected from the sterilisation process
- Batch record data captured

- Instrument data captured
- Staff data linked to the above
- All materials supplied
- Easy to learn

quick reference

MEDITRAX® Non Conformance Module includes ;



related MEDITRAX® products :

- Total Traceability System
- Cold Sterilisation module
- Flash Steriliser module
- Cold Sterilisation module

