



**SPIGGLE & THEIS**  
Medizintechnik

## Non-Stick Bipolar Forceps

*Quality you can rely on.*



### **NON-STICK BIPOLAR FORCEPS**

High-tech forceps for a low-tech price

An adhesion when coagulating with bipolar forceps is a big problem for the surgeon. Our **NON-STICK bipolar forceps** solve this problem easily and reliably.

#### **Advantages:**

- No adhesion due to gold plated solid silver tips
- Precise working without interruptions
- Perfect parallel motion through guide pin
- Parallel tips giving optimal coagulation
- Autoclavable
- Standard adapter grants compatibility to all current HF systems
- Quality made in Germany

#### **Function:**

- No punctual overheating – Tissue detaches easily from forceps.
- With correct usage these characteristics remain over the entire life span of the instrument.

Another quality product selected and stocked by:

WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489  
e enquiries@bamford.co.nz www.bamford.co.nz

When performance counts





**Bipolar forceps, straight, guide pin**

Art. Nr./No.	Länge / Length	Spitze / Tip
80-100-16	16 cm	0,5 mm
80-100-19	19,5 cm	0,5 mm
80-101-16	16 cm	1,0 mm
80-101-19	19,5 cm	1,0 mm
80-101-22	22 cm	1,0 mm
80-103-16	16 cm	2,0 mm
80-103-19	19,5 cm	2,0 mm



**Bipolar forceps, straight, angled tip, guide pin**

Art. Nr./No.	Länge / Length	Spitze / Tip
80-105-16	16 cm	0,5 mm
80-105-19	19,5 cm	0,5 mm
80-106-16	16 cm	1,0 mm
80-106-19	19,5 cm	1,0 mm
80-107-16	16 cm	2,0 mm
80-107-19	19,5 cm	2,0 mm



**Bipolar forceps, bayonet, guide pin**

Art. Nr./No.	Länge / Length	Spitze / Tip
80-110-16	16 cm	0,5 mm
80-110-19	19,5 cm	0,5 mm
80-111-16	16 cm	1,0 mm
80-111-19	19,5 cm	1,0 mm
80-112-16	16 cm	2,0 mm
80-112-19	19,5 cm	2,0 mm



**Bipolar forceps Meuser, 60° curved, guide pin**

Art. Nr./No.	Länge / Length	Spitze / Tip
80-121-18	18 cm	2,0 mm



**Bipolar forceps, bayonet, tip upwards angled, guide pin**

Art. Nr./No.	Länge / Length	Spitze / Tip
80-118-19	19,5 cm	1,0 mm
80-119-19	19,5 cm	2,0 mm



Bipolar forceps, angled, very fine,  
0.3 mm, guide pin

Art. Nr./No.	Länge / Length	Spitze / Tip
80-123-16	16,5 cm	0,3 mm



Bipolar forceps, bayonet, guide pin

Art. Nr./No.	Länge / Length	Spitze / Tip
80-125-03	25 cm	0,3 mm
80-125-07	25 cm	0,7 mm
80-125-12	25 cm	1,2 mm
80-125-20	25 cm	2,0 mm



Bipolar forceps, bayonet, tip downwards angled,  
guide pin

Art. Nr./No.	Länge / Length	Spitze / Tip
80-128-03	25 cm	0,3 mm
80-128-07	25 cm	0,7 mm
80-128-12	25 cm	1,2 mm
80-128-20	25 cm	2,0 mm



Bipolar forceps, bayonet, tip upwards angled,  
guide pin

Art. Nr./No.	Länge / Length	Spitze / Tip
80-127-03	25 cm	0,3 mm
80-127-07	25 cm	0,7 mm
80-127-12	25 cm	1,2 mm
80-127-20	25 cm	2,0 mm



Bipolar forceps, bayonet, tip upwards angled,  
long tip, guide pin

Art. Nr./No.	Länge / Length	Spitze / Tip
80-129-12	25 cm	1,2 mm