

# P<sup>3</sup>

[www.p3-medical.com](http://www.p3-medical.com)

## SINGLE USE VENTED TRACHEAL TUBE INTRODUCER LATEX FREE

VTT Oct 10

- With lumen to allow ventilation.
- Single use - reduces the risk of cross infection.
- Length 750mm, 14ch. Supplied with 15mm connector.
- Extra tip graduations allow accurate judgement of introducer depth within trachea.
- Individually packed sterile, 10 units within highly visible outer tube for ease of identification in an emergency.
- Robust outer tube prevents damage in transit and storage.

- Twin distal eyes ensure oxygen flow.



Code	Description	Dimensions
BV14750	Vented tracheal tube introducer	14ch x 750mm



Another quality product selected and stocked by:

WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489  
e [enquiries@bamford.co.nz](mailto:enquiries@bamford.co.nz) [www.bamford.co.nz](http://www.bamford.co.nz)

When performance counts **BAMFORD** delivers

