

InfuseIT™

Disposable Pressure Infusor

Latex
Free



500 ml

1000 ml

3000 ml

Applications

- Surgery
- Intensive Care
- Radiology, Special Procedures
- Cardiac Catheterization Laboratory
- Emergency/Trauma Department
- Labor & Delivery
- One Day Surgery
- Pre-Operative Care
- Post-Operative Care

CE 0086



5350 SW 1st Lane, Ocala, FL 34474 USA
Tel: +1 352 854 8080; Fax: +1 352 854 7480
www.zefon.com

MM00009 Rev 1 2007-10

EC REP

Blue Box Medical Limited
Unit 29, New Forest Enterprise Centre,
Rushington Business Park, Chapel Lane
Totton, Southampton, SO40 9LA
Phone: +023 8066 9000; Fax: +023 8066 9002

Disposable Pressure Infusor

Features:

- Unique pressure gauge is easy to read and durable.
- Convenient IV Pole Loop Hanger is sturdy and easy to hang onto IV poles - even for short personnel.
- IV Fluid Bag Hook is durable and designed to fit all standard bags snugly.
- Durable Construction - will not rip or tear.
- Clear sleeve around bag provides a visible window to see what fluids are being pressurized.
- Palm fitted, efficient hand bulb makes inflation easy.
- Choice of stopcock valve or bulb thumbwheel valve.



Easy to Read Pressure Gauge

Unique gauge provides an easy to read 360 degree window, which can be read from far away.



Starting Position

Displays the pressure at 0mmHG before any pressure is applied.



Normal Pressure

Indicator marks are located at the 175mmHG & 300mmHG levels. Green color at 300mmHG indicates normal proper inflation.



Overpressurization

Red color indicates overpressurization. When pressure reaches this area an automatic bleed valve will engage reducing the pressure to a safe level.

Product Number	Description
ZIT-500	InfuseIT 500ml with Stopcock valve, 10/bx
ZIT-502	InfuseIT 500ml with Thumbwheel valve, 10/bx
ZIT-1000	InfuseIT 1000ml with Stopcock valve, 10/bx
ZIT-1002	InfuseIT 1000ml with Thumbwheel valve, 10/bx
ZIT-3000	InfuseIT 3000ml with Stopcock valve, 10/bx
ZIT-3002	InfuseIT 3000ml with Thumbwheel valve, 10/bx

Distributed by:



5350 SW 1st Lane
Ocala, FL 34474 USA
Tel: +1 3528548080
Fax: +1 3528547480
www.zefon.com