

VitalCheck™

Reusable Blood Pressure Cuffs

For over ten years Zefon International has been designing and manufacturing superior blood pressure cuffs for leading blood pressure measurement instrument manufacturers. With a complete range of sizes, materials and connectors, Zefon brand cuffs represent the best value to bring consistent performance, patient comfort and economy to your medical setting.



Features

- Reusable; Durable construction for long life and multiple patient use.
- Cost Effective; Bladder and sleeve design allows sleeves to be replaced as they become worn or contaminated.
- Made from Urethane material that is soft and comfortable for the patient.
- Meets AAMI and AHA Cuff Sizing Recommendations.
- Color coded cuff sleeve makes it easy to identify cuff size.
- Latex Free.
- Designed as a direct replacement for OEM cuffs provided by monitor manufacturers.
- Full range of connectors and adapters allow standardization of cuffs across the whole hospital.
- Easy connection to electronic NIBP monitors or manual blood pressure measurement systems.

Latex Free **CE**




5350 SW 1st Lane, Ocala, FL 34474 USA
Tel: +1 3528548080; Fax: +1 3528547480
www.zefon.com

MM00005 Rev 4 2008-06

EC REP Blue Box Medical Limited
Unit 29, New Forest Enterprise Centre,
Rushington Business Park; Chapel Lane
Totton, Southampton, SO40 9LA United Kingdom
Phone: +023 8066 9000; Fax: +023 8066 9002

Ordering Information


Single Tube Cuffs



| | | | CAS 740/750 Criticare 508 Datascopes Accutor Plus & Passport 2 Philips/HP Siemens Spacelabs | BCI Colin Datascopes Passport XG & Expert Nihon Kodan Lifescope 9 & 12 Welch Allyn Vital Signs/Atlas (Before July 2000) | Welch Allyn Vital Signs/Atlas (after July 2000) | Welch Allyn Protocol Propaq |
|-------------|------------|--------------------|--|---|---|------------------------------------|
| | | | | | | |
| Cuff Size | Color Code | Limb Circumference | 1-Tube Bayonet Connector | 1-Tube Female Luer Connector | 1-Tube Screw Type Connector | 1-Tube Female Twist Lock Connector |
| Infant | Brown | 9-15cm | ZVC4901HP | ZVC4901F | ZVC4901D | ZVC4901MQF |
| Child | Green | 14-21.5cm | ZVCC4401HP | ZVCC4401F | ZVCC4401D | ZVCC4401MQF |
| Small Adult | Light Blue | 20.5-28.5cm | ZVCC4301HP | ZVCC4301F | ZVCC4301D | ZVCC4301MQF |
| Adult | Navy Blue | 27-35cm | ZVCC4201HP | ZVCC4201F | ZVCC4201D | ZVCC4201MQF |
| Large Adult | Plum | 34-44cm | ZVCC4101HP | ZVCC4101F | ZVCC4101D | ZVCC4101MQF |
| Thigh | Brown | 42-52cm | ZVC4001HP | ZVC4001F | ZVC4001D | ZVC4001MQF |

All cuffs are packaged with 5 cuffs per box

Double Tube Cuffs



| | | | Critikon® Dinamap® Criticare 507 Invivo Research MDE Nihon Kodan Lifescope L Spacelabs | Datex-Ohmeda | GE Medical Eagle 4000 Tram 100/200/300 | Manual |
|-------------|------------|--------------------|---|---|--|--|
| | | | | | | |
| Cuff Size | Color Code | Limb Circumference | 2-Tube Screw Type Connector | 2-Tube Male/Female Twist Lock Connector | 2-Tube Female Twist Lock Connector | 2-Tube Manual Inflation with Male Luer |
| Infant | Brown | 9-15cm | ZVC4902D | ZVC4902MQMF | ZVC4902MQF | ZVC4912L |
| Child | Green | 14-21.5cm | ZVCC4402D | ZVCC4401MQMF | ZVCC4402MQF | ZVCC4412L |
| Small Adult | Light Blue | 20.5-28.5cm | ZVCC4302D | ZVCC4302MQMF | ZVCC4302MQF | ZVCC4312L |
| Adult | Navy Blue | 27-35cm | ZVCC4202D | ZVCC4202MQMF | ZVCC4202MQF | ZVCC4212L |
| Large Adult | Plum | 34-44cm | ZVCC4102D | ZVCC4102MQMF | ZVCC4102MQF | ZVCC4112L |
| Thigh | Brown | 42-52cm | ZVC4002D | ZVC4002MQMF | ZVC4002MQF | ZVC4012L |

All cuffs are packaged with 5 cuffs per box

Distributed By:



5350 SW 1st Lane,
Ocala, FL 34474 USA
Tel: +1 3528548080
www.zefon.com