

DuraCheck™

Reusable Bladderless Blood Pressure Cuffs



For over ten years Zefon International has been designing and manufacturing superior blood pressure cuffs for leading blood pressure measurement instrument manufacturers. With a complete range of sizes, materials and connectors, Zefon brand cuffs represent the best value to bring consistent performance, patient comfort and economy to your medical setting.

Features

- Reusable; Durable construction for long life and multiple patient use.
- Bladderless design is accurate and easy to clean.
- Made from durable nylon material that is comfortable and strong.
- Meets AAMI and AHA Cuff Sizing Recommendations.
- Color coding makes it easy to identify cuff size.
- Latex Free.
- Complete selection of connectors and adapters allow standardization of cuffs across the whole hospital.
- Easy connection to electronic NIBP monitors or manual blood pressure measurement systems.

Latex Free CE




5350 SW 1st Lane, Ocala, FL 34474 USA
Tel: +1 3528548080; Fax: +1 3528547480
www.zefon.com





MM00002 Rev 5 2008-06

EC REP Blue Box Medical Limited
Unit 29, New Forest Enterprise Centre,
Rushington Business Park; Chapel Lane
Totton, Southampton, SO40 9LA United Kingdom
Phone: +023 8066 9000; Fax: +023 8066 9002

Ordering Information

Single Tube Cuffs








Monitor			
CAS 740/750 Criticare 508 Datascope Accutor Plus & Passport 2 Philips/HP Siemens Spacelabs	Colin Datascope Passport XG & Expert Nihon Koden Lifescope 9 & 12 Welch Allyn Vital Signs/Atlas (Before July 2000)	Welch Allyn Vital Signs/Atlas (after July 2000)	Welch Allyn Protocol Propaq
			
1-Tube Bayonet Connector	1-Tube Female Luer Connector	1-Tube Screw Type Connector	1-Tube Female Twist Lock Connector
ZDC4401HP	ZDC4401F	ZDC4401D	ZDC4401MQF
ZDC4301HP	ZDC4301F	ZDC4301D	ZDC4301MQF
ZDC4201HP	ZDC4201F	ZDC4201D	ZDC4201MQF
ZDC4221HP	ZDC4221F	ZDC4221D	ZDC4221MQF
ZDC4101HP	ZDC4101F	ZDC4101D	ZDC4101MQF
ZDC4121HP	ZDC4121F	ZDC4121D	ZDC4121MQF
ZDC4001HP	ZDC4001F	ZDC4001D	ZDC4001MQF

* Available late summer 2008

All cuffs are packaged with 5 cuffs per box

Double Tube Cuffs



Monitor			
Critikon® Dinamap® Criticare 507 Invivo Research Nihon Koden (Lifescope I) Spacelabs	Datex-Ohmeda	GE Medical Eagle 4000 Tram 100/200/300	Manual
			
2-Tube Screw Type Connector	2-Tube Male/Female Twist Lock Connector	2-Tube Female Twist Lock Connector	2-Tube Manual Inflation with Male Luer
ZDC4402D	ZDC4402MQMF	ZDC4402MQF	ZDC4412L
ZDC4302D	ZDC4302MQMF	ZDC4302MQF	ZDC4312L
ZDC4202D	ZDC4202MQMF	ZDC4202MQF	ZDC4212L
ZDC4222D	ZDC4222MQMF	ZDC4222MQF	ZDC4232L
ZDC4102D	ZDC4102MQMF	ZDC4102MQF	ZDC4112L
ZDC4122D	ZDC4122MQMF	ZDC4122MQF	ZDC4132L
ZDC4002D	ZDC4002MQMF	ZDC4002MQF	ZDC4012L

* Available late summer 2008

All cuffs are packaged with 5 cuffs per box

Distributed By:



5350 SW 1st Lane,
Ocala, FL 34474 USA
Tel: +1 3528548080
www.zefon.com