

FLEXOPEEL® Seethrough Packaging Robust and Bacteriaproof.



Our Service

- Full service from advice on the packaging and punctual delivery to comprehensive customer support on site by our sales representatives
- Well-founded technical know-how
- High-grade flexographic printing in up to four colours

Customer-specific packaging solutions in individual sizes with print as required.

Please ask your sales representative!



A combination of Tyvek® and laminated film makes our **FLEXOPEEL®** packaging very robust and – in comparison to paper – they are also not sensitive to moisture: they remain absolutely impermeable to bacteria. Extremely high strength, excellent gas/air permeability and a user-friendly peel system are further **FLEXOPEEL®** advantages.

Suitable for sterilisation by

- ethylene oxide
- irradiation
- plasma (uncoated Tyvek®)

SENGEWALD 

clinic innovations

Sengewald Klinikprodukte GmbH
Adlerstrasse 2
83101 Rohrdorf-Thansau
Germany
Telephone: +49-80 31-72 24-0
Telefax: +49-80 31-72 24-95
www.sengewald.de
infosengewald@pregis.com

more quality
more safety
more flexibility

Our Range

- **FLEXOPEEL®** seethrough pouches and tubes made of Tyvek® and various types of laminated film (e.g. PET-O/PE or PA/PE film)
- available with coated or uncoated Tyvek®
- specially for products with sharp edges
- with and without indicator

Properties

- suitable for sterilisation with ethylene oxide and irradiation
- extremely high puncture resistance
- high Tyvek® permeability to gas and air enable a fast and safe sterilisation and ventilation
- insensitive to water in comparison with paper: bacteriaproof even in a moist environment
- long-term storage ability of sterile products due to the impermeability to bacteria
- sterile and user-friendly extraction by peeling open the packaging
- particularly with coated Tyvek® grades, absolutely fibre and splinter-free peeling open is achieved
- fast identification of the contents
- individually printable in up to four colours with customer-specific layout and marking as desired
- also suitable for plasma sterilisation as **FLEXOPEEL®** pouches with uncoated Tyvek®

FLEXOPEEL® packaging is in accordance with the DIN 58953, DIN EN 868 and DIN EN ISO 11607 standards

FLEXOPEEL® Seethrough Packaging

Standard Stock Range

Made of PET-O/PE and Tyvek® for sterilisation by ethylene oxide and irradiation

Code #	Outside Dimensions	Packaging Unit
FLEXOPEEL® Seethrough Pouches without Indicator		
11 08 20	80 x 200 mm	2000
11 10 25	100 x 250 mm	1000
11 15 25	150 x 250 mm	1000
11 20 30	200 x 300 mm	500
FLEXOPEEL® Seethrough Tubes with Ethylene Oxide Indicator		
21 10 15	100 mm x 150 m	3
21 15 15	150 mm x 150 m	2
21 20 15	200 mm x 150 m	2
21 30 15	300 mm x 150 m	1
21 38 15	380 mm x 150 m	1