

FLEXOPEEL® Header Pouches The Variety with the Flap.



Our Service

- Full service from advice on the packaging and punctual delivery to comprehensive customer support on site by our sales representatives
- Well-founded technical know-how
- High-grade flexographic printing in up to four colours

Customer-specific packaging solutions in individual sizes with print as required.

Please ask your sales representative!

Our **FLEXOPEEL®** packaging is available with an opening header made of Tyvek®. They make the packaging of voluminous goods in particular safe and user-friendly. The advantages of the **FLEXOPEEL®** materials such as great strength, insensitivity to moisture in comparison with paper or very good gas and air permeability also distinguish the header pouches.

Suitable for sterilisation by

- ethylene oxide
- irradiation

SENGEWALD 
clinic innovations

Sengewald Klinikprodukte GmbH
Adlerstrasse 2
83101 Rohrdorf-Thansau
Germany
Telephone: +49-80 31-72 24-0
Telefax: +49-80 31-72 24-95
www.sengewald.de
infosengewald@pregis.com

more quality
more safety
more flexibility

Our Range

- **FLEXOPEEL®** header pouches made of various laminated films (e.g. coextruded PE film or a combination of PE/PA/PE) with a peelable header made of coated Tyvek®
- Front side made of laminated film, divided rear side made of laminated film and a Tyvek® header
- Very well suited to voluminous products
- With and without gas indicator

Properties

- Suitable for sterilisation by gas and irradiation
- Extremely high puncture resistance
- High gas and air permeability through the Tyvek® enable fast and safe sterilisation and ventilation
- A user-friendly peel system at the header as well as an overlap at the filling aperture ensure easy opening and filling of the pouch
- Insensitivity to moisture in comparison to paper: bacteriproof, even in a moist environment
- Due to bacterial impermeability long-term storage of the sterile products is possible
- Sterile extraction due to absolutely fibre and splinter-free peeling open
- High seal strength for the packaging of voluminous products
- Fast identification of the contents
- Individually printable in up to four colours: customer-specific layout and markings as desired

FLEXOPEEL® header pouches are in accordance with the DIN 58953, DIN EN 868 and DIN EN ISO 11607 standards

FLEXOPEEL® Header Pouches

Standard Stock Range

PE/PA/PE film and coated Tyvek® for sterilisation by ethylene oxide and irradiation

Code #	Outside Dimensions	Packaging Unit
FLEXOPEEL® Header Pouches without Indicator		
16 11 01	230 x 400 mm	1000
16 11 02	330 x 400 mm	1000
16 11 03	370 x 580 mm	1000
16 11 04	400 x 620 mm	1000
16 11 05	350 x 430 mm	1000
16 11 06	450 x 640 mm	1000
16 11 07	580 x 700 mm	500
16 11 08	500 x 700 mm	500
16 11 09	650 x 900 mm	500

