

FLEXOTEAR® Pouches

The Robust Solution.



Our Service

- Full service from advice on the packaging and punctual delivery to comprehensive customer support on site by our sales representatives
- Well-founded technical know-how
- High-grade flexographic printing in up to four colours

These customer-specific packaging solutions are only manufactured as tailor-made custom products. Please ask your sales representative!



FLEXOTEAR® packaging is made out of transparent plastic film on both sides. That makes them particularly robust and their contents can be identified at a glance.

The linear tear open effect ensures very easy package opening.

FLEXOTEAR® pouches with Tyvek® valves are particularly suited to sterilisation by ethylene oxide. They are also available without Tyvek® valves and are then solely for sterilisation by irradiation. (**LINEARTEAR®**)

Suitable for sterilisation by

- ethylene oxide
- irradiation

SENGEWALD 
 clinic innovations

Sengewald Klinikprodukte GmbH
 Adlerstrasse 2
 83101 Rohrdorf-Thansau
 Germany
 Telephone: +49-80 31-72 24-0
 Telefax: +49-80 31-72 24-95
 www.sengewald.de
 infosengewald@pregis.com

more quality
 more safety
 more flexibility

Our Range

- **FLEXOTEAR®** pouches made of coextruded PE film and Tyvek® valves
- up to two valves are possible
- with opening aid and special **LINEARTEAR®** effect (i.e. easy linear tearing open at a freely definable tear notch)
- **FLEXOTEAR®** packaging is not available as a standard stock range.

FLEXOTEAR® pouches without Tyvek® valves can also be purchased under the **LINEARTEAR®** brand name. These packagings are also distinguished by opening aids and the special **LINEARTEAR®** effect and are solely suitable for sterilisation by irradiation.

Properties

- Suitable for sterilisation by irradiation and, with Tyvek® valves, also for ethylene oxide sterilisation.
- Long-term storage is possible due to the highest degree of bacterial impenetrability, also in moist environments
- Fast, easy opening of the packaging due to opening aids and the special **LINEARTEAR®** effect
- No need to stick indicator tapes or labels onto the packaging, as ethylene oxide indicators can be printed onto the Tyvek® valves
- Ideal for voluminous products due to sturdy seals
- Seethrough on both sides of the pouches: fast identification of the contents
- Format limits: min. 100 x 200 mm and max. 700 x 550 mm
- Individually printable with up to four colours: customer-specific layout and markings as desired

FLEXOTEAR® packaging is in accordance with the DIN 58953, DIN EN 868 and DIN EN ISO 11607 standards