

MEDIPEEL® Seethrough Packaging Economical and Reliable.



Our Service

- Full service from advice on the packaging and punctual delivery to comprehensive customer support on site by our sales representatives
- Well-founded technical know-how
- High-grade flexographic printing in up to four colours

Customer-specific packaging solutions in individual sizes with print as required.

Please ask your sales representative!

MEDIPEEL® seethrough pouches and tubes are not only a clear but also a safe affair: they exceed the requirements of the valid standards. Optical defect controls and indicators for the relevant methods of sterilisation provide additional safety. A targeted choice from our extensive range, as well as environmentally responsible materials, make **MEDIPEEL®** a sensible solution for the packaging of sterile medical devices.

Suitable for sterilisation with

- steam (MEDIPEEL®D)
- ethylene oxide (MEDIPEEL®D and R)
- formalene (MEDIPEEL®D and R)
- irradiation (MEDIPEEL®R)

SENGEWALD 

clinic innovations

Sengewald Klinikprodukte GmbH
Adlerstrasse 2
83101 Rohrdorf-Thansau
Germany
Telephone: +49-80 31-72 24-0
Telefax: +49-80 31-72 24-95
www.sengewald.de
infosengewald@pregis.com

more quality
more safety
more flexibility

Our Range

MEDIPEEL® Seethrough Pouches:

- ideal for large series, e.g. in the C.S.S.D. and industry
- for light products without sharp edges
- made of special medical paper and PET-O/PP film for sterilisation with steam, ethylene oxide and formalene (MEDIPEEL®D) or PET-O/PE film for sterilisation with ethylene oxide and irradiation (MEDIPEEL®R)
- with and without gusset

MEDIPEEL® Self-Adhesive Seethrough Pouches:

- ideal for practices and small sterilisation departments
- self-adhesive: saves the sealing machine
- for light products without sharp edges
- made of special medical paper and PET-O/PP film

MEDIPEEL® Seethrough Tubes:

- ideal for small series
- for light products without sharp edges
- made of special medical paper and PET-O/PP film for sterilisation with steam, ethylene oxide and formalene (MEDIPEEL®D) or PET-O/PE film for sterilisation with ethylene oxide and irradiation (MEDIPEEL®R)
- with and without gusset

Properties

- suitable for sterilisation with steam, ethylene oxide and formalene (MEDIPEEL®D) or PET-O/PE film for sterilisation with ethylene oxide and irradiation (MEDIPEEL®R)
- tinted film makes it easier to detect optical defects in the seals
- print is completely visible through the film
- simple and clean peeling thanks to advice on the peel direction
- colour indicators for steam, ethylene oxide or formalene sterilisation outside the fill area (MEDIPEEL®D), no migration onto the contents is possible, time-consuming sticking on of labels is unnecessary
- PVC-free film and oxygen-bleached paper
- individually printable with up to four colours: customer-specific layout and markings on request

MEDIPEEL® packaging is in accordance with the DIN 58953, DIN EN 868 and DIN EN ISO 11607 standards

MEDIPEEL® D Seethrough Packaging

Standard Stock Range

Made of special medical paper and PET-O/PP film for sterilisation with steam, ethylene oxide and formalene.

Code #	Outside Dimensions	Packaging Unit
MEDIPEEL® D Seethrough Pouches, without Gusset		
10 06 19	60 x 190 mm	3500
10 06 24	60 x 240 mm	3500
10 08 15	80 x 150 mm	2500
10 08 20	80 x 200 mm	3500
10 08 25	83 x 250 mm	2500
10 08 30	83 x 300 mm	2500
10 10 25	103 x 250 mm	2000
10 10 30	106 x 310 mm	2000
10 10 38	106 x 380 mm	1000
10 10 59	100 x 590 mm	1000
10 15 20	150 x 200 mm	2000
10 15 30	150 x 300 mm	2000
10 15 38	150 x 380 mm	1000
10 20 28	200 x 280 mm	1000
10 20 38	200 x 380 mm	1000
10 25 30	250 x 300 mm	500
10 25 40	250 x 400 mm	1000
10 25 48	250 x 480 mm	500
10 30 50	300 x 500 mm	500
MEDIPEEL® D Seethrough Pouches, with Gusset		
19 10 28	100 x 50 x 280 mm	1000
19 10 38	100 x 50 x 380 mm	1000
19 15 38	150 x 50 x 380 mm	500
19 20 38	200 x 50 x 380 mm	500
19 20 48	200 x 50 x 480 mm	500
19 25 48	250 x 65 x 480 mm	500
19 30 58	300 x 65 x 580 mm	500
MEDIPEEL® D Seethrough Tubes, without Gusset		
20 05 20	50 mm x 200 m	6
20 08 20	75 mm x 200 m	4
20 10 20	100 mm x 200 m	3
20 15 20	150 mm x 200 m	2
20 20 20	200 mm x 200 m	2
20 25 20	250 mm x 200 m	2
20 30 20	300 mm x 200 m	1
20 38 20	380 mm x 200 m	1
20 42 20	420 mm x 200 m	1
MEDIPEEL® D Seethrough Tubes, with Gusset		
30 08 10	75 x 25 mm 100 m	4
30 10 10	100 x 50 mm 100 m	3
30 15 10	150 x 50 mm 100 m	2
30 20 10	200 x 50 mm 100 m	2
30 25 10	250 x 65 mm 100 m	2
30 30 10	300 x 65 mm 100 m	1
30 38 10	380 x 80 mm 100 m	1
30 42 10	420 x 80 mm 100 m	1
MEDIPEEL® Self-Adhesive Seethrough Pouches, without Gusset		
17 08 18	80 x 180 mm	2000
17 08 27	80 x 270 mm	2000
17 10 30	100 x 300 mm	2000
17 10 60	100 x 600 mm	1000
17 15 30	150 x 300 mm	1000
17 15 40	150 x 400 mm	1000
17 20 34	200 x 340 mm	500
17 30 40	300 x 430 mm	500