



# Endotracheal Tube

100% Silicone for Great Comfort and Flexibility



**FORTUNE MEDICAL INSTRUMENT CORP.**

6Fl., No. 29, Sec. 2, Jungjeng E. Rd., Danshuei Jen, Taipei Hsien, TAIWAN

Email: [info@fortunemed.com](mailto:info@fortunemed.com) [sales1@fortunemed.com](mailto:sales1@fortunemed.com)

TEL: (886)-2-2624 2233 FAX: (886)-2-2624-2266

**[www.FortuneMed.com](http://www.FortuneMed.com)**



0123



**ISO 13485 & CE Certified**

Another quality product selected and stocked by:

WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489  
e [enquiries@bamford.co.nz](mailto:enquiries@bamford.co.nz) [www.bamford.co.nz](http://www.bamford.co.nz)


When performance counts



**Fortune®** Silicone Endotracheal Tube is made of 100% medical grade silicone. It provides superior biocompatibility and flexibility. Great transparency allows easy visual inspection. Graduation and X-ray opaque line help depth and location confirmation. A pilot balloon is equipped to monitor balloon status. Low pressure cuff models are available to reduce the pressure on trachea wall and to provide greater comfort to patient.

#### Silicone Endotracheal Tubes, Cuffless:

PED-Magill REF. No.	PED-Murphy REF. No.	Size/I.D.	O.D.	Length	Description
1510-0025	1515-0025	2.5 mm	4.0 mm	140 mm	Cuffless
1510-0030	1515-0030	3.0 mm	4.7 mm	160 mm	Wire reinforced
1510-0035	1515-0035	3.5 mm	5.3 mm	180 mm	X-ray opaque line
1510-0040	1515-0040	4.0 mm	6.0 mm	200 mm	Magill series: Magill tip
1510-0045	1515-0045	4.5 mm	6.7 mm	220 mm	Murphy series: Murphy eye
1510-0050	1515-0050	5.0 mm	7.3 mm	240 mm	



#### Silicone Endotracheal Tubes, OC Type (Ordinary cuff, Curved) :

OC-Magill REF. No.	OC-Murphy REF. No.	Size/I.D.	O.D.	Length	Description
1520-0045	1525-0045	4.5 mm	6.7 mm	220 mm	Ordinary cuff
1520-0050	1525-0050	5.0 mm	7.3 mm	240 mm	Curved
1520-0055	1525-0055	5.5 mm	8.0 mm	270 mm	Pilot balloon
1520-0060	1525-0060	6.0 mm	8.7 mm	280 mm	X-ray opaque line
1520-0065	1525-0065	6.5 mm	9.3 mm	290 mm	
1520-0070	1525-0070	7.0 mm	10.0 mm	300 mm	
1520-0075	1525-0075	7.5 mm	10.7 mm	310 mm	
1520-0080	1525-0080	8.0 mm	11.3 mm	320 mm	*OC-Magill series: Magill tip
1520-0085	1525-0085	8.5 mm	12.0 mm	330 mm	*OC-Murphy series: Murphy eye
1520-0090	1525-0090	9.0 mm	12.7 mm	340 mm	
1520-0095	1525-0095	9.5 mm	13.3 mm	350 mm	
1520-0100	1525-0100	10.0 mm	14.0 mm	360 mm	



#### Silicone Endotracheal Tubes, LC Type (Low pressure cuff, Curved) :

LC-Magill REF. No.	LC-Murphy REF. No.	Size/I.D.	O.D.	Length	Description
1530-0045	1535-0045	4.5 mm	6.7 mm	220 mm	Low pressure cuff
1530-0050	1535-0050	5.0 mm	7.3 mm	240 mm	Curved
1530-0055	1535-0055	5.5 mm	8.0 mm	270 mm	Pilot balloon
1530-0060	1535-0060	6.0 mm	8.7 mm	280 mm	X-ray opaque line
1530-0065	1535-0065	6.5 mm	9.3 mm	290 mm	
1530-0070	1535-0070	7.0 mm	10.0 mm	300 mm	
1530-0075	1535-0075	7.5 mm	10.7 mm	310 mm	
1530-0080	1535-0080	8.0 mm	11.3 mm	320 mm	*LC-Magill series: Magill tip
1530-0085	1535-0085	8.5 mm	12.0 mm	330 mm	*LC-Murphy series: Murphy eye
1530-0090	1535-0090	9.0 mm	12.7 mm	340 mm	
1530-0095	1535-0095	9.5 mm	13.3 mm	350 mm	
1530-0100	1535-0100	10.0 mm	14.0 mm	360 mm	



#### Silicone Endotracheal Tubes, OW Type (Ordinary cuff, Wire reinforced) :

OW-Magill REF. No.	OW-Murphy REF. No.	Size/I.D.	O.D.	Length	Description
1540-0045	1545-0045	4.5 mm	6.7 mm	220 mm	Ordinary cuff
1540-0050	1545-0050	5.0 mm	7.3 mm	240 mm	Wire reinforced
1540-0055	1545-0055	5.5 mm	8.0 mm	270 mm	Pilot balloon
1540-0060	1545-0060	6.0 mm	8.7 mm	280 mm	
1540-0065	1545-0065	6.5 mm	9.3 mm	290 mm	
1540-0070	1545-0070	7.0 mm	10.0 mm	300 mm	
1540-0075	1545-0075	7.5 mm	10.7 mm	310 mm	
1540-0080	1545-0080	8.0 mm	11.3 mm	320 mm	*OW-Magill series: Magill tip
1540-0085	1545-0085	8.5 mm	12.0 mm	330 mm	*OW-Murphy series: Murphy eye
1540-0090	1545-0090	9.0 mm	12.7 mm	340 mm	
1540-0095	1545-0095	9.5 mm	13.3 mm	350 mm	
1540-0100	1545-0100	10.0 mm	14.0 mm	360 mm	



#### Silicone Endotracheal Tubes, LW Type (Low pressure cuff, Wire reinforced) :

LW-Magill REF. No.	LW-Murphy REF. No.	Size/I.D.	O.D.	Length	Description
1550-0045	1555-0045	4.5 mm	6.7 mm	220 mm	Low pressure cuff
1550-0050	1555-0050	5.0 mm	7.3 mm	240 mm	Wire reinforced
1550-0055	1555-0055	5.5 mm	8.0 mm	270 mm	Pilot balloon
1550-0060	1555-0060	6.0 mm	8.7 mm	280 mm	
1550-0065	1555-0065	6.5 mm	9.3 mm	290 mm	
1550-0070	1555-0070	7.0 mm	10.0 mm	300 mm	
1550-0075	1555-0075	7.5 mm	10.7 mm	310 mm	
1550-0080	1555-0080	8.0 mm	11.3 mm	320 mm	*LW-Magill series: Magill tip
1550-0085	1555-0085	8.5 mm	12.0 mm	330 mm	*LW-Murphy series: Murphy eye
1550-0090	1555-0090	9.0 mm	12.7 mm	340 mm	
1550-0095	1555-0095	9.5 mm	13.3 mm	350 mm	
1550-0100	1555-0100	10.0 mm	14.0 mm	360 mm	



O: ordinary cuff, L: low pressure cuff, W: wire reinforced