



ENT
CONSUMABLES

Overview

CoreMed is an Australian based medical distribution company, specialising in the sales and marketing of ENT products.

We provide a professional, efficient and experienced approach towards the product and service requirements of hospital operating theatres, emergency departments and medical practices.

We are proud to be the Exclusive Distributor for Spiggle & Theis in Australia and New Zealand.



Contents	Page
1. Nasal Splints	1
2. Disposable Suction Tubes	2
3. Ventilation Tubes	3
4. Sterile Covers	4
5. Bandages	5-6

Internal Nasal Splints

- Sterile, ready to use
- Punched suture holes
- Soft material for easy insertion and removal
- Surface with Lotus leaf effect (Silicone splints only)

Code	Feature	Packaging	Material
2030100	standard / with airway	5 pairs / box	silicone
2030110	standard / without airway	5 pairs / box	silicone
2042025	large / 0.25 mm thickness	5 pairs / box	fluoroplastic
2042050	large / 0.50 mm thickness	5 pairs / box	fluoroplastic



Septal Button, Silicone

- Sterile, ready to use
- For non-surgical closure of septal perforation
- Soft material for easy insertion and removal
- Silicone prevents adhesion
- No scabbing, epistaxis or other septal perforation symptoms

Code	Size	Feature	Packaging
2031100	3.2 cm	one part	1 piece
2031110	3.2 cm	two parts	1 piece



External Aluminium Nasal Splints

- Easily malleable
- Self-adhesive back, self-adhesive Dorsum foam pad included
- 10 pieces per box

Code	Dimensions	Size
2060200	56 x 38 x 30 mm	small
2060300	64 x 40 x 36 mm	medium
2060400	66 x 42 x 45 mm	large



External Thermoplastic Nasal Splints

- Much less time and work required than with plaster of paris
- Includes three-part foam bridge
- Self-adhesive back
- 10 pieces per box

Code	A x B x C	Thickness	Size
2050100	56 x 38 x 30 mm	2.4 mm	small
2050200	64 x 40 x 36 mm	2.4 mm	medium
2050300	66 x 42 x 45 mm	2.4 mm	large
2051100	56 x 38 x 30 mm	1.6 mm	small
2051200	64 x 40 x 36 mm	1.6 mm	medium
2051300	66 x 42 x 45 mm	1.6 mm	large



Disposable Suction Tubes - Micro

Code	Length	Diameter	Colour	Connection
301004	60 mm	0.4 mm	grey	Luer
301007	60 mm	0.7 mm	black	Luer
301010	60 mm	1.0 mm	white	Luer
301014	60 mm	1.4 mm	red	Luer
302004	70 mm	0.4 mm	grey	Luer
302006	70 mm	0.6 mm	blue	Luer
302007	70 mm	0.7 mm	black	Luer
302008	70 mm	0.8 mm	green	Luer
302010	70 mm	1.0 mm	white	Luer
302012	70 mm	1.2 mm	pink	Luer
302014	70 mm	1.4 mm	red	Luer
302016	70 mm	1.6 mm	white	Luer
301020	70 mm	2.0 mm	green	Luer
303004	80 mm	0.4 mm	grey	Luer
303007	80 mm	0.7 mm	black	Luer
303014	80 mm	1.4 mm	red	Luer
303020	80 mm	2.0 mm	green	Luer
303025	80 mm	2.5 mm	transparent	Luer-Lock
303030	80 mm	3.0 mm	transparent	Luer-Lock
303040	80 mm	4.0 mm	transparent	Luer-Lock
304010	100 mm	1.0 mm	white	Luer
304014	100 mm	1.4 mm	red	Luer
304020	100 mm	2.0 mm	green	Luer
304025	100 mm	2.5 mm	transparent	Luer-Lock
304030	100 mm	3.0 mm	transparent	Luer-Lock



- Angled cannula made of stainless steel
- Luer / Luer-Lock connection
- Atraumatic rounded cannula point
- Single sterile packaging
- Wide range of different sizes
- Packaging of 50 pieces per box

Suction Handle

Code	Description	Package
301000	With tip control, Luer connection, sterile	50 pcs./box



Disposable Suction Tubes

Code	Length	Diameter
------	--------	----------

3012020	120 mm	2.0 mm
3012025	120 mm	2.5 mm
3012030	120 mm	3.0 mm
3012035	120 mm	3.5 mm
3012040	120 mm	4.0 mm
3012050	120 mm	5.0 mm

3014020	140 mm	2.0 mm
3014025	140 mm	2.5 mm
3014030	140 mm	3.0 mm
3014035	140 mm	3.5 mm
3014040	140 mm	4.0 mm
3014050	140 mm	5.0 mm

3016020	160 mm	2.0 mm
3016025	160 mm	2.5 mm
3016030	160 mm	3.0 mm
3016035	160 mm	3.5 mm
3016040	160 mm	4.0 mm
3016050	160 mm	5.0 mm



- Tip control
- Blue tube connection
- Angled cannula made of stainless steel
- Atraumatic rounded cannula point
- Single sterile packaging
- Packaging of 50 pieces per box

3010025	100mm	2.5 mm
3010030	100mm	3.0 mm
3010040	100mm	4.0 mm

- Rounded tip
- Curved
- Atraumatic rounded tip
- With suction control
- Sinus surgery
- Packaging of 50 pieces per box



3025020	250 mm	2.0 mm
3025030	250 mm	3.0 mm

For microlaryngoscopy (MLs):



- Distal opening (adherence protection)
- Packaging of 10 pieces per box

Vent tubes

Titanium

Code	Type	Inner Diameter	Feature	Packaging
10076	Collar Button	0.76 mm		10pcs / box
10100	Collar Button	1.00 mm		10pcs / box
10100F	Collar Button	1.00 mm	w/wire	10pcs / box
10125	Collar Button	1.25 mm		10pcs / box
10125F	Collar Button	1.25 mm	w/wire	10pcs / box
10150	Collar Button	1.50 mm		10pcs / box
10150F	Collar Button	1.50 mm	w/wire	10pcs / box
101100	Shepard	1.00 mm		10pcs / box
101100F	Shepard	1.00 mm	w/wire	10pcs / box
101125	Shepard	1.25 mm		10pcs / box
101125F	Shepard	1.25 mm	w/wire	10pcs / box



Silicone

Code	Type	Inner Diameter	Feature	Packaging
103100	T-tube, 12 mm	1.10 mm		10pcs / box
103150	T-tube, 5 mm	1.10 mm		10pcs / box
103200	Grommet T-tube	1.27 mm	w/tab	10pcs / box
103300	Paparella I	1.10 mm	w/tab	10pcs / box
103310	Paparella II	1.50 mm	w/tab	10pcs / box



Fluoroplastic

Code	Type	Inner Diameter	Feature	Packaging
104100	Collar Button	1.00 mm		10pcs / box
104100F	Collar Button	1.00 mm	w/wire	10pcs / box
104125	Collar Button	1.25 mm		10pcs / box
104125F	Collar Button	1.25 mm	w/wire	10pcs / box
104150	Collar Button	1.50 mm		10pcs / box
104150F	Collar Button	1.50 mm	w/wire	10pcs / box
104115	Shepard	1.15 mm		10pcs / box
104115F	Shepard	1.15 mm	w/wire	10pcs / box
104200	Shah	1.15 mm	w/tab	10pcs / box
104210	Shah	1.15 mm		10pcs / box



Sterile Covers


For Leica and Wild-Leitz microscopes:

Code	Dimensions	Shape	Lense	Feature	Packaging
2500110H	115 x 165 mm	D	70 mm	angled lense for Leica	10 pieces / box
2500120H	115 x 260 mm	D	70 mm	angled lense for Leica	10 pieces / box
2500130H	137 x 265 mm	D	70 mm	angled lense for Leica	5 pieces / box
2500140F	115 x 300 mm	D	70 mm	angled lense for Leica	10 pieces / box



For Zeiss, Kaps or Moller-Wedel microscopes:

Code	Dimensions	Shape	Lense	Feature	Packaging
2548050F	115 x 180 cm	D	48 mm		10 pieces / box
2548090H	115 x 210 cm	D	48 mm		10 pieces / box
2565070F	50 x 160 cm	A	65 mm		10 pieces / box
2565020H	115 x 180 cm	D	65 mm		10 pieces / box
2565030H	115 x 260 cm	D	65 mm		10 pieces / box
2565040H	115 x 300 cm	D	65 mm		10 pieces / box
2565050H	115 x 210 cm	D	65 mm		10 pieces / box
2565080H	125 x 67 x 155 cm	B	65mm		10 pieces/box
2565090H	137 x 260 cm	D	65 mm		5 pieces / box
2565120F	50 cm		65 mm	cap-shaped	10 pieces / box



Sundry Covers

Code	Dimensions	Type	Packaging
2500010	universal	Microscope handgrips	20 pairs / box
2800070	8 x 150 cm	Drill cable cover	30 pieces / box

External nasal bandage

- For post-operative nasal surgery
- Absorbs nasal secretion
- Cotton compress contains a stretch net bandage
- Fits the nasal contour
- Non-sterile

Code	Size	Packaging
10904	universal	40 pieces / box



Ear Bandages

- High quality material
- Ready for application
- Applied within seconds
- Optimal patient comfort
- Single, sterile packaging



Ear bandage one-sided

- Easy and quick application
- Post-operative bandage
- One-sided
- Compress included
- Flexible and comfortable
- Single, sterile packaging

Code	Size	Packaging
35000	small	10 pieces / box
35001	medium	10 pieces / box
35002	large	10 pieces / box



Ear bandage one-sided with Velcro

Code	Size	Packaging
35180	universal	10 pieces / box



Ear bandage two-sided

Code	Size	Packaging
35030	small	5 pieces / box
35029	medium	5 pieces / box
35028	large	5 pieces / box



Ear bandage for otoplasty

Code	Size	Packaging
35033	small	5 pieces / box
35032	medium	5 pieces / box
35031	large	5 pieces / box



Elastic head bandages

- For post-operative head surgery
- Elastic stretch bandage with face opening
- Simple application
- Non-sterile

Code	Size	Packaging
90006	universal	25 pieces / box



Elastic compression bandage

- Compression for incisions in the area below the ear
- With ear slit and fixation strap
- Velcro sealing strip
- Simple application
- Non-sterile

Code	Size	Packaging
90007	universal	5 pieces / box





Exclusive Distributor for Spiggle & Theis

PO Box 565, Balmain, NSW 2041, Australia P +61 02 9810 6947 F +61 02 8072 1556 admin@coremed.com.au
www.coremed.com.au