

Sterile Suction Tubing

- Single-use soft and flexible medical grade plastic transparent tubing
- Thick tubing walls offer safety and collapse resistance
- Smooth inner bore surface allowing for easy flow of suction
- Tubing is comfortable to handle, malleable and resistant to kinking and coiling
- Ridges on outer tubing surface minimizes friction at contact
- Universal female connectors on both ends fit all suction outlets and devices
- 10mm outer diameter and 7 mm inner diameter with tubing total length of 3.6m
- Individually pre-packed in peel down pouches ready to use
- Sterile and Latex free

Product No.	Description	Quantity
DEF3194	Suction Tubing 7mm x 3.6m Sterile	20 per box



Another quality product selected and stocked by:

WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489
e enquiries@bamford.co.nz www.bamford.co.nz

When performance counts

