



# STERIS SYSTEM I® module


**MEDITRAX**  
SMART TRACKING SYSTEM



## accountability & traceability

The MEDITRAX® Steris System I® module is a systematic controlled documentation and tracking system for the accurate recording and control of all medical devices processed and distributed from the proprietary Steris System I® process.

This module records data relating to processing the device, linking individual instruments to batch cycle information and patient records. Additionally the responsible person for each individual process signs a serialized log book linking these events in chronological sequence.



## features : instrument tracking

## standard recording

### components

The MEDITRAX® Steris System I® module includes :

Qty	Description	Code
1	label applicator	MHG1800
1	label applicator CD instructions	CDI800
5 rolls	Steris System I® batch labels	M/ST/003
1000	Steris System I® patient labels	M/ST/002
5 pads	patient record sheets	TTS/201
1	patient label dispenser	ASS/80R
1	anti bacterial pen	AB/210
1	Steris System I® log book	M/ST/001
1	Steris System I® ID label	M/ST/004
1	policies & procedures manual	M/ST/PP005

### implementation

MEDITRAX® Steris System I® module includes ongoing :-

- staff training
- review of system
- quality assurance consultation
- maintenance program at no extra cost.

### summary

The MEDITRAX® Steris System I® module benefits

- Customised for Steris System I® Process
- Batch record data captured
- Patient record data captured
- Instrument data captured
- Staff data linked to the above
- All materials supplied
- Easy to learn

The MEDITRAX® Steris System I® Log Book is manufactured from Sterisyn waterproof synthetic paper for essential record protection in harsh wet environments.

### quick reference

MEDITRAX® Steris System I® module includes ;



ENDOSCOPY



BATCH RELEASE



PATIENT RECORD



STAFF DATA RECORDED

## related MEDITRAX® products :

- Total Traceability System
- Batch Recall Plus II module
- Flash Sterilization module
- TrayBax Computerised Tracking & Traceability System





# STERIS SYSTEM I® module

## accountability & traceability

consumables for MEDITRAX® Steris System I® module :

Code		Description	Packaging	Order Qty
MHG1800		MEDITRAX® label applicator	1 each	
IRT/ING05		MEDITRAX® non gassing ink roller	1 each	
873170		MEDITRAX® standard ink roller	1 each	
TTS/201		MEDITRAX® patient record sheets	5 pads/packet 50 sheets/pad	
M/ST/003		MEDITRAX® Steris System I® batch labels	5 rolls/packet 1000 labels/roll	
M/ST/002		MEDITRAX® Steris System I® patient labels	1000/roll	
AB/210		MEDITRAX® anti bacterial pen	1 each	
M/ST/001		MEDITRAX® Steris System I® log book	1 each	

### training

Austmel provides MEDITRAX® System Inservice at the point of implementation. Courses are run on a person-to-person basis ensuring rapid absorption. Advanced courses are designed to train the trainer where required. MEDITRAX® training systems comply with the AS/NZS 4187 standards and documentation.

### one-on-one MEDITRAX® customer support

Austmel Pty Ltd is able to offer one-on-one support during office hours on our toll free phone service. This is a valuable service to healthcare personnel who require immediate technical support to use the system. This service is a complete one-on-one person to person level. Additionally IT Remote Access Tech-Support provides quick on-line diagnosis and repair of systems to clients requiring rapid response under maintenance contracts.


### documentation

MEDITRAX® documentation is built on the philosophy of paper waste minimisation, taking into account policies and procedures relating to the processing of items being sterilised in healthcare facilities. All aspects of the complex requirements are clearly and easily followed.

### background

Austmel Pty Ltd is an Australian owned company established in 1974. Austmel have been directly involved in the supply and development of Sterility Assurance products for over 30 years. Our unique products and systems extend into both the Health Care and Industrials markets, including :

- MEDITRAX® Accountability & Traceability Systems
- Biological and Chemical Sterilization Indicators
- Bowie-Dick Tests
- Sterilization Packaging
- Quality Assurance Consultation
- TMI Electronic Data Loggers
- Surgical Disposables



**T** 1800 074 085 Aus or 0800 400 151 NZ **customer support**

### summary

The MEDITRAX® Steris System I® module starter pack features :

- Customised for Steris System I® Process
- Batch record data captured
- Patient record data captured
- Instrument data captured
- Staff data linked to the above

- All materials supplied
- Easy to learn

The MEDITRAX® Steris System I® Log Book is manufactured from Sterisyn waterproof synthetic paper for essential record protection in harsh wet environments.

### quick reference

MEDITRAX® Steris System I® module includes :



ENDOSCOPY



BATCH RELEASE



PATIENT RECORD



STAFF DATA RECORDED

### related MEDITRAX® products :

- Total Traceability System
- Batch Recall Plus II module
- Flash Sterilization module
- TrayBax Computerised Tracking & Traceability System

