

Suprasorb® Liquacel

Hydroactive fibre dressing, sterile

Product Description

Suprasorb® Liquacel is a needle bonded fabric which consists of cellulose fibres and cellulose-ethyl-sulphonate fibres. In contact with wound exudate Suprasorb® Liquacel forms a cohesive gel which traps exudate, debris and bacteria and transfers them to a secondary dressing. This keeps the exudate away from the wound edge, reduces the risk of maceration and promotes wound healing.

This product is presented sterile.



Use

- assists in protecting against maceration of the wound edge and the surrounding skin
- effective exudate management assists moist wound healing and the healing process
- uncomplicated application under moist conditions, as it can be removed in one piece
- when dry can be cut to wound shape/size and/or folded if required
- can be used as a wound filler

Features

The cellulose-ethyl-sulphonate fibres interact with the wound exudate or blood to form a gel which creates a moist wound healing environment.

Available in a dressing or rope.

Storage

Stored appropriately - not above 25 °C, in a dry place protected from light Suprasorb® Liquacel has a shelf life of 3 years.

Ordering Information

Product	Size	Code
Suprasorb® Liquacel	5cm x 5cm	33435 **
Suprasorb® Liquacel	10cm x 10cm	33436 **
Suprasorb® Liquacel	15cm x 15cm	33437 *
Suprasorb® Liquacel Rope	2cm x 45cm	33438 *

* supplied in packs of 5
** supplied in packs of 10

Indications

Suitable for the treatment of acute and chronic wounds including:

- Lower leg ulcers
- Pressure ulcers
- Diabetic ulcers
- Lacerations, incisions and abrasions
- Second degree burns
- Traumatic wounds and wounds that are prone to bleeding, for example following surgical or mechanical debridement
- Post surgical wounds, such as skin graft donor sites
- Wounds healing by secondary intention if they are exudating and superficial to deep

Contra-Indications

- 3rd degree burns
- Known sensitivity to the product or its components.