

Product Data T-mat™




Cairn Technology's T-mat is a disposable absorbent mat that collects and binds liquid hygienically – promoting a cleaner working environment.

Some surgical and investigative procedures can generate considerable amounts of fluid. The T-mat minimises the impact of these fluids:

T-mat

- ✓ Polymer gelling technology rapidly absorbs and gels liquids, controlling and locking liquid ready for simple, safe disposal
- ✓ Traps liquids on contact
- ✓ Simple, inexpensive and effective way of collecting and holding liquid hygienically and safely
- ✓ Only a few millimetres thick
- ✓ No handling of liquid Clinical Waste
- ✓ Quick and effective waste disposal
- ✓ Gelled waste is disposed of as Clinical Waste



Dimensions: 750 x 360mm

Weight: Approx. 100g

Use: Single use
 Do not reuse
 Do not cut outer sheath
 Do not use if outer sheath damaged

Absorbency
 : 8 litres of water or
 3.2 litres of 0.9% saline solution

Applications: Absorbs aqueous liquids including blood & urine

Contains: Inner: cellulose & sodium polyacrylate
 Outer sheath: polypropylene fibre

Sterile: No

Disposal: After use, treat as Clinical Waste

Packaging: 100 mats per box

Supplier: Cairn Technology Ltd, Sheffield

