

TRACOE *vario*

January 2012

Transparent and metal-free TRACOE *vario* tracheostomy tube is now also available in **extra-long** as REF 461!



NEW!

Advantages:

- The extra-long tube is especially convenient for obese patients.
- The extra-long tube with adjustable neck flange and X-ray contrast line is metal-free and therefore applicable in MRI. Only in the Luer valve of the cuff control balloon is a small metal spring which needs to be taped away during examination.
- Through the variably adjustable neck flange with rotating wings a universal adjustment to any patient's neck anatomy is possible (e.g. thick neck, postoperative oedema, asymmetrical necks).

REF 461 TRACOE *vario* with adjustable neck flange and low pressure cuff.
The tracheostomy tube is available in sizes 7, 8, 9 and 10 mm.

Dimensions REF 461

Size	ID (mm)	OD (mm)	Length TL* (mm)	Length AB** (mm)
07	7.0	9.7	104	110
08	8.0	11.2	123	128
09	9.0	12.3	126	134
10	10.0	13.7	135	146

*TL: length along centre line - ** AB: length along outer bend

For further information please visit our homepage: www.tracoe.com