




TOTAL TRACEABILITY SYSTEM


MEDITRAX
SMART TRACKING SYSTEM



total traceability system

The supply of sterile instruments to the hospital is the main objective of the CSD and success is critical. To prove sterility, policy and procedures must be strictly adhered to. The MEDITRAX® Total Traceability System (TTS) was designed to provide documentation with tight control, to dramatically lessen the potential of process errors and highlight any deficiencies, reducing the probability of patient harm.

The MEDITRAX® TTS provides control by giving a unique identity to instruments and instrument sets, in the form of a specially designed TTS label. The TTS label provides the ability to trace items through the sterilization process to and from the patient. The information collected by the system is stored in a simple effective filing system, extending to a comprehensive archival system.



features : device tracking to patient

standard recording

components

The MEDITRAX® Total Traceability System includes :

Qty	Description	Code
1	label applicator	MHG1800
1	label applicator CD instructions	CD1800
2	non gassing ink roller	IRTNG05
5 rolls	batch labels	MH102A
1000	MEDITRAX® labels	TTS2
250	batch cycle cards	TTS202/B
1	sterilizer ID labels	MH113L
1	applicator ID labels	MH111S
50	patient record sheet	TTS201
50	daily audit sheets	TTSAM11
50	weekly audit sheets	TTSAM05
50	batch clips	TTSS20
1	clipboard for batch cards	TTSAM007
10	sterilizer pen	TTS209P
1	authorised transfer pad	TTSATD01
50	non conformance report	TTSNC/50
50	non conformance sheets	TTS203
1	non conformance log	TTSNC/60
1 roll	non conformance labels	MH105A
6	binders	TTSS01
1	medimarc bullet tip marker	MM02T

implementation

MEDITRAX® Total Traceability System includes ongoing :-

- staff training
- review of system
- quality assurance consultation
- maintenance program at no extra cost.

summary

The MEDITRAX® Total Traceability System benefits

- Full audit trail traceability
- Batch record data captured
- Patient record data captured
- Instrument data captured
- Sterilizer data captured
- Non conformance monitored
- Staff data linked to the above
- Closed loop accountability
- All materials supplied
- Easy to learn

quick reference

MEDITRAX® Total Traceability System module includes ;



BATCH RECORD



6849 3421 9877 3419
SERIAL # RECORDED



PATIENT RECORD



STAFF DATA
RECORDED

related MEDITRAX® products :

- TTS Non Conformance module
- TTS Multi Batching module
- Archive module
- TrayBax Computerised Sterility Management System



Austmel Pty Ltd A Head Office Claude House,16 Jade Drive P PO Box 2345, Nerang QLD 4211, AUSTRALIA T 1800 074 085 F 07 5597 5244 W www.austmel.com E info@austmel.com
New Zealand T 0800 400 151

Austmel Pty Ltd A Head Office Claude House,16 Jade Drive P PO Box 2345, Nerang QLD 4211, AUSTRALIA T 1800 074 085 F 07 5597 5244 W www.austmel.com E info@austmel.com
New Zealand T 0800 400 151



- TTS Non Conformance module
- TTS Multi Batching module
- Archive module
- TrayBax Computerised Sterility Management System

related MEDITRAX® products :



TOTAL TRACEABILITY SYSTEM



total traceability system

consumables for MEDITRAX® Total Traceability System :

Code		Description	Packaging	Order Qty
MHG1800		MEDITRAX® label applicator	1 each	
IRTNG05		MEDITRAX® non gassing ink roller	1 each	
873170		MEDITRAX® standard ink roller	1 each	
MHI02A		MEDITRAX® batch labels	850/roll <div>green red yellow cerise blue</div>	
TTS2		MEDITRAX® labels	1000/packet	
TTS202/B or TTS202/5		MEDITRAX® batch cycle cards	250/packet 500/box	
TTS201		MEDITRAX® patient record sheets	5 pads/packet 50 sheets/pad	
TTSAM11		MEDITRAX® daily audit sheets	50/pad	
TTSAM05		MEDITRAX® weekly audit sheets	50/pad	
TTSS20		MEDITRAX® batch clips	50/packet	
TTSAM007		MEDITRAX® clipboard for batch cards	1 each	
TTS209P or TTS209P25		MEDITRAX® sterilizer pens	1 each 25/packet	
TTSATD01		MEDITRAX® authorised transfer pad	50/pad	
TTSNC/50		MEDITRAX® non conformance reports	50/pad	
TTS203		MEDITRAX® non conformance sheets	50/pad	
TTSNC/60		MEDITRAX® non conformance log	1 each	
MHI05A		MEDITRAX® non conformance label	1000/roll 5 rolls/packet	
TTSS0		MEDITRAX® binders (specify sterilizer numbers)	6/box <div>1 2 3 4 5 6</div>	
MM02T		MEDITRAX® medimarc bullet tip marker	1 each <div>black blue green red</div>	

training

Austmel provides MEDITRAX® System Inservice at the point of implementation. Courses are run on a person-to-person basis ensuring rapid absorption. Advanced courses are designed to train the trainer where required. MEDITRAX® training systems comply with the AS/NZS 4187 standards and documentation.

one-on-one MEDITRAX® customer support

Austmel Pty Ltd is able to offer one-on-one support during office hours on our toll free phone service. This is a valuable service to healthcare personnel who require immediate technical support to use the system. This service is a complete one-on-one person to person level. Additionally IT Remote Access Tech-Support provides quick on-line diagnosis and repair of systems to clients requiring rapid response under maintenance contracts.

documentation

MEDITRAX® documentation is built on the philosophy of paper waste minimisation, taking into account policies and procedures relating to the processing of items being sterilized in healthcare facilities. All aspects of the complex requirements are clearly and easily followed.

background

Austmel Pty Ltd is an Australian owned company established in 1974. Austmel have been directly involved in the supply and development of Sterility Assurance products for over 30 years. Our unique products and systems extend into both the Health Care and Industrials markets, including:

- MEDITRAX® Accountability Traceability Systems
- Biological and Chemical Sterilization Indicators
- Bowie-Dick Tests
- Sterilization Packaging
- Quality Assurance Consultation
- TMI Electronic Data Loggers
- Surgical Disposables

summary

The MEDITRAX® Total Traceability System benefits

- Full audit trail traceability
- Batch record data captured
- Patient record data captured
 - Instrument data captured
 - Sterilizer data captured
 - Non conformance monitored
 - Closed loop accountability
 - Staff data linked to the above
 - All materials supplied
 - Easy to learn

