

# Tracoe Cuff Pressure Monitor



## TRACOE cuff pressure monitor

**Ref 720** 1x TRACOE cuff pressure monitor and 1 connection tube Ref 702

The TRACOE cuff pressure monitor is used to fill and monitor the pressure of high-volume, low pressure cuffs of tracheostomy and endotracheal tubes.

Its outstanding features and benefits:

- Easy use
- Precise, reproducible pressure indication
- Extra large, easy to read scale from 0 to 100 mbar
- Widely spaced scale markings with colour coded pressure range

The TRACOE cuff pressure monitor comes with a connection tube with male and female Luer lock connectors on its ends (Ref 702).

## Connection Tube

**Ref 702** 10 units supplied in individual sterile packages

Tube for connecting a tracheostomy or endotracheal tube to the cuff pressure monitor, with Luer lock connectors.

Length: 1m

## Connection Tube

**Ref 701** 10 units supplied in individual sterile packages

Tube for connecting a tracheostomy or endotracheal tube to the cuff pressure monitor, with Luer lock connectors.

Length: 2m

Another quality product selected and stocked by:

WM Bamford & Co Ltd p 0800-BAMFORD f 04-569-6489  
e enquiries@bamford.co.nz www.bamford.co.nz

When performance counts

