

Two-Stage Breast Expander Quick Reference Poster

ENHANCED SURFACE

All Expanders are available with an enhanced surface (textured). Simply add an “es” after the first 4 numbers.

(e.g. 3612es-14)

TEARDROP BREAST EXPANDERS

<u>VOLUME (CC)</u>	<u>CATALOG #</u>	<u>WIDTH (CM)</u>	<u>HEIGHT (CM)</u>	<u>PROJECTION (CM)</u>
500	3612es-22-2i	10	13	8.5
600	3612es-23-2i	11	14	10
850	3612es-24-2i	12	16	11
1200	3612es-26-2i	14	18	13

ROUNDED TEARDROP BREAST EXPANDERS

<u>VOLUME (CC)</u>	<u>CATALOG #</u>	<u>WIDTH (CM)</u>	<u>HEIGHT (CM)</u>	<u>PROJECTION (CM)</u>
400	3612ES-60-2i	10	11	8.5
550	3612ES-61-2i	11	12	9.5
650	3612ES-62-2i	11.5	12	10
850	3612ES-63-2i	12.5	13.5	11
1200	3612ES-64-2i	14	14	12
1300	3612ES-65-2i	15	15	13

PORT TYPE

Remote or Internal;
Internal port is designated with an “i” after the catalog number.

(e.g. 3612es-22-2i)

WIDE BASE - OVAL BREAST EXPANDERS

<u>VOLUME (CC)</u>	<u>CATALOG #</u>	<u>WIDTH (CM)</u>	<u>HEIGHT (CM)</u>	<u>PROJECTION (CM)</u>
450	3611ES-204-2i	12	10	9
650	3611ES-200-2i	13	11	10
1000	3611ES-199-2i	15	13	11.5

CUSTOM EXPANDERS

Special modifications can be made for special needs, including suture tabs, modified port exits, or a pre-attached port.

*Modifications should be noted when ordering

*Prices may vary

WIDE BASE - ROUND BREAST EXPANDERS

<u>VOLUME (CC)</u>	<u>CATALOG #</u>	<u>WIDTH (CM)</u>	<u>HEIGHT (CM)</u>	<u>PROJECTION (CM)</u>
500	3611ES-38-2i	11	11	9
900	3611ES-27-2i	13	13	12.5