

The Poppa™

Universal Connection System

The Poppa™ NIBP Cuff Connection System

A true universal connection system for NIBP Cuffs that allows standardization of connections on NIBP monitors and blood pressure cuffs.



Features

- Unique design makes connecting and disconnecting easy, yet still provides a leak-free seal.
- Safety—The POPPA™ has purposely been designed to be completely incompatible with luer connections to avoid any possible misconnection with intravenous lines and catheters.
- No more splitters, Y-pieces or other adapting type systems.
- Follows the Joint Commission recommendations to avoid purchasing non-intravenous devices that can physically mate with a female luer IV line adapter.*
- Improve inventory management by stocking one style cuff for each size range that can connect to any NIBP monitor.
- In-line connection keeps The Poppa™ safely away from patient tampering and cross-contamination.
- Simple, easy to clean and disinfect design with NO moving parts.

*Tubing misconnections—a persistent and potentially deadly occurrence. (Joint Commission Sentinel Event Alert, Issue 36, April 3rd, 2006.) www.jointcommission.org

LateX Free CE



5350 SW 1st Lane, Ocala, FL 34474 USA
Tel: +1 3528548080; Fax: +1 3528547480
www.zefon.com

MM00010 Rev 1 2008-07



EC REP Blue Box Medical Limited
Unit 29, New Forest Enterprise Centre,
Rushington Business Park; Chapel Lane
Totton, Southampton, SO40 9LA, United Kingdom
Phone: +023 8066 9000; Fax: +023 8066 9002

Ordering Information

Getting started with The Poppa™ is as easy as 1, 2, 3

1. Outfit all monitors in the facility with The Poppa™.
2. Purchase cuffs with The Poppa™ connector installed.
3. Easily connect to any monitor.

Poppa Monitor Connectors

Connector	Illustration	Pack Size	Product Number
Poppa 1-Tube Monitor Connector (green)		10/pk	BPC-P001SF
		100/pk	BPC-P001SFB
Poppa 2-Tube Monitor Connector (green)		10/pk	BPC-P002DF
		100/pk	BPC-P002DFB

Poppa Connectors Pre-installed on NIBP Cuffs

All cuffs with Poppa connectors are configured as 1 Tube cuffs and come with a white Poppa connector pre-installed.



Cuff Size	Color Code	Size Range	CUFF STYLE		
			ComfortCheck™ Soft Fabric Disposable BP Cuffs	StatCheck™ Vinyl Disposable BP Cuffs	DuraCheck™ Reusable BP Cuffs
Child	Green	11-22cm	ZCC4401P	ZSC4401P	ZDC4401P
Small Adult	Light Blue	16-28cm	ZCC4301P	ZSC4301P	ZDC4301P
Adult	Dark Blue	26-33cm	ZCC4201P	ZSC4201P	ZDC4201P
Adult Long*	Dark Blue	26-33cm	ZCC4221P	ZSC4221P	ZDC4221P
Large Adult	Red	33-41cm	ZCC4101P	ZSC4101P	ZDC4101P
Large Adult Long*	Red	33-41cm	ZCC4121P	ZSC4121P	ZDC4121P
Thigh	Brown or Dark Blue	39-55cm	ZCC4001P	ZSC4001P	ZDC4001P

*Long size cuffs available summer 2008.

Distributed By:



5350 SW 1st Lane,
Ocala, FL 34474 USA
Tel: +1 3528548080
www.zefon.com